
## Choosing the right magnetic beads





## Three types of magnetic particles available from Cytiva

	Sera-Mag™ magnetic beads	SeraSil-Mag™ silica beads	Mag Sepharose™ magnetic beads
Core material	Carboxyl	Silica	Agarose
Surface chemistry options	<ul> <li>Streptavidin coated/ Streptavidin blocked</li> <li>Amine-blocked</li> <li>Oligo (dT) coated</li> <li>NeutrAvidin™</li> <li>Protein A/G</li> </ul>	Silica (OH)	<ul> <li>Streptavidin</li> <li>Protein A or G</li> <li>NHS</li> <li>His</li> </ul>
Size	1 μm	400 nm and 700 nm	37–100 μm
Size distribution	Monodisperse	Monodisperse	Polydisperse
Surface	Solid but irregular (cauliflower like)	Uniform and spherical	Porous
Capacity	Medium	High	High
Application focus	<ul> <li>Nucleic acid extraction and cleanup</li> <li>Direct conjugation of ligands such as antibodies, enzymes, or oligos</li> <li>mRNA purification</li> <li>Enrichment of target nucleic acid sequences for NGS</li> <li>Size selection of nucleic acids</li> <li>Antibody purification with combined A/G affinities</li> </ul>	<ul> <li>Traditional nucleic acid extraction using chaotropic salts</li> <li>Existing Sera-Xtracta™ extraction kits use of SeraSil-Mag™ silica beads to deliver high quality performance</li> </ul>	<ul> <li>Antibody purification, screening, immunoprecipitation, and pull-down</li> <li>Direct conjugation of ligands such as antibodies, enzymes, and oligos</li> <li>Affinity purification of histidine tagged recombinant proteins</li> </ul>

## Comparison of magnetic bead surface chemistries and applications

Туре	Properties	Applications	Variations	Product	Pack size	Product code
Sera-Mag™ magnetic beads						
Sera-Mag™ carboxylate-		Conjugation or direct binding	High-speed version available:  • Sera-Mag™ SpeedBeads carboxylate-modified	Sera-Mag™ Carboxylate- Modified [E3] Magnetic Beads	15 mL	44152105050250
modified magnetic beads	for direct capture	<ul><li>applications:</li><li>Covalent attachment</li></ul>			100 mL	44152105050350
	Surface suitable for conjugation	<ul> <li>Affinity purification and</li> </ul>	•		1000 mL	44152105050450
		pull-down		Sera-Mag™ SpeedBead	15 mL	65152105050250
	containing amino groups	<ul> <li>Nucleic acid isolation and purification</li> </ul>		Carboxylate-Modified [E3] Magnetic Beads	100 mL	65152105050350
		<ul> <li>NGS size selection</li> </ul>			1000 mL	65152105050450
				Sera-Mag™ Carboxylate-	15 mL	24152105050250
				Modified [E7] Magnetic Beads	100 mL	24152105050350
					1000 mL	24152105050450
				Sera-Mag™ SpeedBead Carboxylate-Modified [E7] Magnetic Beads	15 mL	45152105050250
					100 mL	45152105050350
					1000 mL	45152105050450
Sera-Mag™ amine-blocked	Surface suitable for conjugation	Conjugation applications, similar to carboxylate-modified beads	<ul> <li>High-speed version available:</li> <li>Sera-Mag™ SpeedBeads amine-blocked</li> </ul>	Sera-Mag™ SpeedBeads 1 mL Amine-Blocked Magnetic Beads 5 mL 100 ml	1 mL	19152104011150
magnetic beads	through covalent bonding  Non-surfactant, non-protein-				5 mL	19152104010150
	blocked surface				100 mL	19152104010350
	Low non-specific binding					
Sera-Mag™ Oligo(dT)-coated	Hybridizes with mRNA poly-A	<ul> <li>mRNA binding applications:</li> <li>mRNA extraction and purification</li> <li>RT-PCR</li> <li>cDNA library construction</li> <li>Subtractive hybridization</li> </ul>		Sera-Mag™ Oligo (dT)	1 mL	38152103011150
magnetic beads	tails				5 mL	38152103010150
	High colloidal stability				100 mL	38152103010350
		NGS (RNA sequencing)				

Туре	Properties	Applications	Variations	Product	Pack size	Product code
Sera-Mag™ streptavidin-coated magnetic beads	Binds biotinylated ligands such	Immunoassay and molecular	<ul> <li>High-speed version available:</li> <li>Sera-Mag™ SpeedBeads streptavidin-coated</li> <li>Biotin binding ranges:</li> <li>2500 to 3500 pmol/mg</li> <li>3500 to 4500 pmol/mg</li> <li>4500 to 5500 pmol/mg</li> </ul>	Sera-Mag™ Streptavidin 2500 to 3500 (Low) pmol per mg Magnetic Beads	1 mL	30152103011150
	as proteins, nucleic acids, and peptides	<ul><li>biology applications:</li><li>Sample preparation and</li></ul>			5 mL	30152103010150
	Covalently bound streptavidin	assay development for genomics and proteomics			100 mL	30152103010350
	coating  Fact reaction kinetics	genomics and proceomics		Sera-Mag™ Streptavidin	1 mL	30152104011150
	Fast reaction kinetics			3500 to 4500 (Med.) pmol per mg Magnetic Beads	5 mL	30152104010150
	Low non-specific binding				100 mL	30152104010350
	High throughput and precision			Sera-Mag™ Streptavidin	1 mL	30152105011150
				4500 to 5500 (High) pmol per mg Magnetic Beads	5 mL	30152105010150
					100 mL	30152105010350
				Sera-Mag™ SpeedBeads Streptavidin 3500 to 4500 (Med.) pmol per mg	1 mL	66152104011150
					5 mL	66152104010150
					100 mL	66152104010350
Sera-Mag™ streptavidin-	Binds biotinylated ligands such	High-specificity biotin binding	High-speed version available:  • Sera-Mag™ SpeedBeads streptavidin-blocked	Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Beads	1 mL	21152104011150
blocked magnetic beads	as proteins, nucleic acids, and peptides	applications: molecular and immunodiagnostics			5 mL	21152104010150
	Non-surfactant, non-protein- blocked surface	NGS library preparation			100 mL	21152104010350
	Lower non-specific binding than streptavidin-coated beads via additional blocking of non-specific binding sites					
Sera-Mag™ NeutrAvidin™- coated magnetic beads	Binds biotinylated ligands such as proteins, nucleic acids, and peptides	Alternative to Streptavidin in immunoassay and molecular biology applications:  • Sample preparation and assay development for genomics and proteomics	<ul> <li>High-speed version available:</li> <li>Sera-Mag™ SpeedBeads     NeutrAvidin™-coated</li> <li>Biotin binding range:</li> <li>3500 to 4500 pmol/mg</li> </ul>	Sera-Mag SpeedBeads NeutrAvidin™-Coated Magnetic Beads	1 mL 5 mL	78152104011150 78152104010150
	Fast reaction kinetics				100 mL	78152104010350
	Low non-specific binding					
	High throughput and precision					

Туре	Properties	Applications	Variations	Product	Pack size	Product code
Sera-Mag™ protein A/G magnetic beads	Binds IgA and IgG proteins	<ul><li>Antibody isolation applications:</li><li>Affinity purification and pull-down</li></ul>		Sera-Mag SpeedBeads Protein A/G Magnetic Beads	1 mL	17152104011150
	Coating based on IgA/IgG fusion protein				5 mL	17152104010150
	Broad binding capabilities	<ul> <li>Immunoprecipitation</li> </ul>			100 mL	17152104010350
SeraSil-Mag™ silica beads						
SeraSil-Mag™ silica-coated	Reversibly binds nucleic acids	Applications with low sample	Available in 400 μm or 700 μm	Beads  Beads  SeraSil-Mag™ 700 Silica  Beads  Beads	5 mL	29357369
magnetic beads	based on salt concentration  Monodisperse particles with	amounts  Nucleic acid extraction	particle sizes		60 mL	29357371
	narrow size ranges	for molecular diagnostics applications such as qPCR			450 mL	29357372
					1000 mL	29705862
					5 mL	29357373
					60 mL	29357374
					450 mL	29357375
					1000 mL	29705861
Mag Sepharose™ magnetic beads						
His-Mag Sepharose™ Ni	Highly cross-linked spherical agarose (Sepharose), including magnetite IMAC immobilized with nickel	Small-scale purification and screening of histidine-tagged proteins from different sources			2 × 1 mL	28967388
magnetic beads					5 × 1 mL	28967390
					100 mL	29104065
His Mag Sepharose™ excel	Strip resistant ligand with strongly bound nickel for immobilized metal ion affinity chromatography (IMAC)	Small-scale capture and purification of histidine-tagged proteins secreted into eukaryotic cell culture supernatants		His Mag Sepharose™ excel Beads	2 × 1 mL	17371220
magnetic beads					5 × 1 mL	17371221
					10 × 1 mL	17371222
NHS Mag Sepharose™ magnetic	Coupling of antibodies, aptamers, and proteins through primary amino groups on the molecules to the NHS ligand on NHS Mag Sepharose™	Enrichment of target protein for further downstream analyses such as mass spectometry (MS) and electrophoresis techniques		NHS Mag Sepharose™	500 μL	28944009
beads				Beads	4 × 500 μL	28951380

Туре	Properties	Applications	Variations	Product	Pack size	Product code
Protein A Mag Sepharose™ magnetic beads	Maximum binding capacity due to dense coating of Protein A  Optimized capacity for enrichment or immunoprecipitation requiring only low amounts of antibody needed	Enrichment of target proteins via immunoprecipitation or pull- down assays  Optimised for downstream analyses such as mass spectometry (MS) and electrophoresis techniques		Protein A Mag Sepharose™ Beads	500 μL 4 × 500 μL	28944006 28951378
Protein A Mag Sepharose™ Xtra magnetic beads	Maximum binding capacity due to dense coating of Protein A	High capacity small-scale purification and screening of monoclonal and polyclonal antibodies from various species		Protein A Mag Sepharose™ Xtra Beads	2 × 1 mL 5 × 1 mL	28967056 28967062
Protein G Mag Sepharose™ magnetic beads	Maximum binding capacity due to dense coating of Protein G  Optimized capacity for enrichment or immunoprecipitation requiring only low amounts of antibody needed	Enrichment of target proteins via immunoprecipitation or pull- down assays  Optimised for downstream analyses such as mass spectometry (MS) and electrophoresis techniques		Protein G Mag Sepharose™ Beads	500 μL 4 × 500 μL	28944008 28951379
Protein G Mag Sepharose™ Xtra magnetic beads	Maximum binding capacity due to dense coating of Protein G	High capacity small-scale purification/screening of monoclonal and polyclonal antibodies from various species		Protein G Mag Sepharose™ Xtra Beads	2 × 1 mL 5 × 1 mL	28967066 28967070

## cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva. Sera-Mag, SeraSil-Mag, Sepharose, and Sera-Xtracta are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

NeutrAvidin is a trademark of Thermo Fisher Scientific. Any other third-party trademarks are the property of their respective owners.

© 2022 Cytiva

For local office contact information, visit cytiva.com/contact

CY31839-22Sep22-SG

