Palltronic® Flowstar V

INTEGRITY TEST INSTRUMENT

Fast, accurate integrity testing with plug-and-play network integration

The Palltronic® Flowstar V instrument ensures accurate filter integrity testing (FIT) with up to a 20% reduction in test duration compared to previous generations. It has been designed to comply entirely with the FDA's 21 CFR Part 11 guidance on electronic records and signatures and provides the user with automation and network integration features that can save them time and effort.

Designed to be used in any cleanroom environment, the Palltronic Flowstar V instrument has been developed and validated in accordance with the current GAMP guidelines. It also comes complete with the requisite controls and other features for operation in a 21 CFR Part 11 compliant environment. Additional data integrity controls as well as electronic records and signatures further support quality risk management.

The instrument has been designed such that test fluid contact areas are manufactured from FDA-listed material.

The human machine interface (HMI) is a 307 mm (12.1 in.) color touchscreen that can be adjusted for ease-of-use, whilst the machine is also fitted with LED strips to support test status and result visibility whilst in use.



Fig 1. Palltronic Flowstar V integrity test instrument.



Features and benefits

	Features	Benefits
Technology	Direct flow measurement	Faster testing High accuracy of test results Calibration based on traceable flow references
Ease-of-use	Automatic control of the test system	Reduced test time under stable test conditions
	Storage of test parameters and recipes	Reduced setup times Avoids input errors
	Light instrument with a small footprint	Portable Easy to locate close to point of use
	Intuitive user interface with 12 in. (1280 × 800 pixels) color touchscreen and adjustable screen angle	Enhanced readability Intuitive design improving ease of use for data entry and navigation Adjustable positioning of display angle to improve ease of interaction and aid comfort when in use
	LED alarm strip	Enables quick check of instrument status whilst in operation, even from a distance
Quality and regulatory	Developed in accordance with GAMP 5	Assurance instrument complies with GAMP standard
	IP54 rated design	Robust and easy to clean
	Compliance to 21 CFR Part 11	Meets regulatory requirements
	Electronic records and electronic signature	Requirement for electronical batch documentation
	Automated data entry through barcode reader	Mitigates risks of manual entry errors of test parameters (quality risk management)
	Test repeat control feature	Ensures compliance with the work instructions (SOP) Mitigates risks of entry errors Improves ease of use
	Data security and audit trail	Facilitates the detection of unauthorized record editing
	Test parameters and records can be stored electronically	Ensures permanent records
	Fade-resistant printouts	Ensures permanent records
Network integration and automation	Possibility to connect to OPC-UA, Profibus, DeviceNet, Profinet, Serial, PAS-X MES, WiFi 2.4 GHz	Faster and simplified integration to local network Improved flexibility for remote connection
	Lightweight Direct Access Protocol (LDAP)	Increased safety margins with lower risks in forgetting user ID and password
	USB and LAN networking connectivity	Connectivity ease of use for flash drives, printers and LAN networking

Secures compliance to data integrity and 21 CFR Part 11 requirements

Electronic Records and Electronic Signature (ERES) features

The instrument has all features and requirements required to be operated in accordance with 21 CFR Part 11 and other guidelines like the Medicines and Healthcare Products Regulatory Agency (MHRA) related to data integrity and electronic signature safety.

Access management helps prevent uncontrolled access to critical functions on the Palltronic Flowstar V integrity test instrument.

Access management and password management allows the electronic signing of test records. The Palltronic Flowstar V integrity test instrument allows:

- Configurable access control management (four access levels Super Administrator, Administrator, Supervisor, and Operator)
- · Configurable password management

Audit trail features

The instrument captures all data needed for audit purposes, facilitating the review process regardless of frequency, while complying with all data integrity requirements issued by the FDA, MHRA, and Pharmaceutical Inspection Cooperation Scheme (PIC/S). This ease of review plays a major role for regulated companies that are expected to regularly review data with a direct impact on batch release and patient safety. The Palltronic Flowstar V integrity test instrument fulfils ALCOA PLUS principles which require that data be attributable, legible, contemporaneous, original, accurate, complete, consistent, enduring, and available. ALCOA has served as a framework for guaranteeing data integrity and good documentation practice (GDP) by regulated industries.

In addition, the Palltronic Flowstar V integrity test instrument has an event audit trail which logs activities on the instrument and a record audit trail which includes test program records, access management records, user records, date/time change records and configuration records.

The primary records of the Palltronic Flowstar V integrity test instrument, including the test results, cannot be altered after generation and storage. Only signatures and any comments with the signature can be added. With this approach, the test result is a static record.

Finally, with the different data backup possibilities, including manual and automatic backup options, and the expanded data storage capacity on the instrument, the Palltronic Flowstar V integrity test instrument ensures compliance with all data backup, retention and retrieval requirements.



Fig 2. Integrity testing a filter capsule.

State-of-the-art automation protocols and remote control

Supplementary software

Data management system (DMS) Pro software

The DMS Pro software is designed to synchronize and back up data across one or more Palltronic Flowstar V instruments simultaneously. Features include:

- Synchronization of date and time across instruments
- · Synchronization of test programs across instruments
- · Automated, continuous backup of audit trails and test records
- · Use of domain log in for Active Directory passwords



Fig 3. DMS Pro software provides users with information on filter integrity test instruments from across their facility.

mPath™ operational insights software (OIS)

The mPath OIS is a browser-based analytics solution that monitors filter integrity test activity across multiple Palltronic Flowstar V instruments. It can improve productivity by providing insights into:

- Test results by instrument, filter product and individual operators
- · Troubleshooting false-failure integrity test results
- · Preventive maintenance and calibration scheduling




Fig 4. mPath operational insights software monitors filter integrity test activity across multiple instruments and provides insights on test results, troubleshooting, and other key analytics.

Automation connectivity

The Palltronic Flowstar V instrument has been designed with manufacturing execution system compatibility in mind. Through the use of automation adaptors and software licenses, the instrument can work with a variety of protocols to:

- · Access real-time FIT progress and results
- Maintain traceability of filter products across batch and lot numbers
- Enforce individual FIT protocol rules
- Automated retrieval of test records through plug and play PAS-X MES integration

Supplementary hardware

Palltronic Flowstar V MUX instrument

The Palltronic Flowstar V MUX instrument connects to the Palltronic Flowstar V integrity test instrument⁽¹⁾ and allows a user to test up to 16 separate filter cartridges or capsules sequentially, reducing the need for a user to be present at frequent intervals to connect and disconnect filters to be tested.



Fig 5. Palltronic Flowstar V MUX instrument.

(1) The Palltronic Flowstar V instrument must be operating on at least firmware version 2.1.

Technical specifications

Physical dimensions	Weight: 11 kg (24.3 lb) width × height × depth: 350 × 263 × 450 mm (13.8 × 10.4 × 17.7 in.)	
Filter tests	Forward flow test (FF) Bubble point test (BP) Combined forward flow / bubble point test (FF/BP) Water intrusion test (WIT)	
Other tests	Leak test (LT) (flow measurement based for volume < 50 L) Pressure decay test (PD) (leak test for volumes up to 200 L)	
Function tests	Self-test Flow check test Network test	
Other functions	Test program transfer function Configuration export/import function Access management export/import function User list export/import function Test result export function Full back up function (including all above functions)	
Operating system	Linux ⁽¹⁾	
Data capacity	Test program storage: > 500 programs Test result storage: > 10 000 tests Active user storage: > 250 users	
Language options ⁽²⁾	Chinese Dutch English French German Italian Japanese Portuguese Russian Spanish US English	
Communication ports	1× CAN 1× ethernet port (RJ45) 2× USB (standard) 1× USB (under screen on the front panel) 1× custom printer interface	
Automation protocols	OPC UA Industry standard Profibus, Profinet, Serial, and DeviceNet (via external USB modules)	
Measuring range	Forward flow test: 0.1 to 1000 mL/min Water intrusion test: 0.03 to 50 mL/min water Bubble point test: 400 to 6500 mbar (5.8 to 94.3 psi)	
Accuracy	Forward flow test: \pm 3% of measurement or \pm 0.05 mL/min whichever is greater Water intrusion test: \pm 3% of measurement or \pm 0.02 mL/min, whichever is greater	
Resolution	Forward flow test 0.1 mL/min Water intrusion test 0.01 mL/min Bubble point test 50 mbar (0.7 psi) Pressure decay test: 1 mbar	
Thermal printer	See detailed description and specifications on page 7	
Calibration limits	Flow measurement: ± 3% of measurement Pressure measurement: ± 30 mbar (10 bar full scale output, calibrated up to 7 bar)	
Electrical data	Voltage: automatically adjusted between 100 and 240 V AC Input frequency: 50 Hz / 60 Hz Power: typically 75 W, peak 150 W Fuse: 3.15 A, slow blow External vent valve: 24 V DC External pressure transducer: relative (to gauge) pressure measurement, 4 mA to 20 mA signal	
Touchscreen	Size: 12.1 in. LCD (1280 × 800) Features: stylus pen, capacitive touch, 30° adjustable angle, color, illuminated background, adjustable contrast, enclosed multicolour LED strip	

Pneumatic specifications	Maximum gas supply pressure: 8000 mbar (116 psi) Minimum gas supply: – self-test: 3000 mbar (43.5 psi) – flow range 0.01 to 149 mL/min 1000 mbar (14.5 psi) above test pressure – flow range 150 to 1000 mL/min 2000 mbar (29.0 psi) above test pressure Test pressure range: 50 to 6500 mbar (0.7 to 94.3 psi)
Pneumatic connections	1× Staubli IN male 1× Staubli OUT female 1× Staubli male, 1× external vent 1× internal vent
Accessory printer	Palltronic Flowstar printer accessory connected via custom printer interface Lifetime of the printout > 10 years depending on storage conditions Paper width: 112 mm (4.4 in.) Paper roll diameter: 50 mm (1.9 in.)
Peripherals	One digital out for external vent valve Power and control via CAN bus connector Automation connectivity to network via USB adapter box
Environmental conditions	Splash-proof: IP54 ⁽³⁾ Operating temperature: +5 to +40°C Storage temperature: -20 to +50°C Humidity: 80% up to 31°C, 50% at 40°C Cleaning fluid compatibility – 70%/30% ethanol/water, 70%/30% isopropyl alcohol/water Altitudes up to 3000 m ⁽⁴⁾ (9800 ft)

The source code which falls under the GNU General Public License (GPL) is stored on the USB flash drive supplied with the instrument.

Contact your local Cytiva representative for other languages.

To meet the full requirements of IP54, splash-proof electrical connections to the instrument are necessary. These are available as accessories; please contact us for details.

For altitudes higher than 2000 m, we recommend an on-site calibration for best accuracy of measurement.

Thermal printer

The thermal printer for the Palltronic Flowstar V integrity test instrument delivers accurate printouts on long-lifetime Mitsubishi HiTec thermal paper.

Printer installation

The printer can be used on a table next to the Palltronic Flowstar V instrument or be set in the indentation on the top of the Palltronic Flowstar V instrument.



Fig 6. Thermal printer as delivered.



Fig 7. Thermal printer installed on the top of the Palltronic Flowstar V instrument.

Technical specifications

Environmental conditions (printer only)

Operating	5 to 40°C, 20% to 85% relative humidity (RH) non-condensing -40 to +70°C, 5% to 95% RH non-condensing	
Storage		
IP rating (when transparent lid closed)	IP 20	
Overall maximum dimensions (without cable)	Width × height × depth: 184 × 77 × 165 mm	
Printer language	Printer should be capable to handle the printer language PCL5e or PostScript Printer language is proprietary and compatible with Palltronic Flowstar V instrument	

Ordering information

Product	Product code
Palltronic Flowstar V filter integrity tester	FFS05
Palltronic Flowstar V accessory printer	FFS-PRNT
Palltronic Flowstar V barcode reader	FFS-BARCODE
Palltronic Flowstar V DMS Pro software (1)	FFS-DMS
Palltronic Flowstar V filter integrity tester with DMS Pro software	FFS-FFS05-DMS
Palltronic Flowstar V with automation adapter (Profibus)	FFS-AUTADAPTBUS
Palltronic Flowstar V with automation adapter (DeviceNet)	FFS-AUTADAPTDNET
Palltronic Flowstar V with automation adapter (Profinet)	FFS-AUTADAPTPNET
Palltronic Flowstar V with automation adapter (Serial)	FFS-AUTADAPTSBUS
Palltronic Flowstar V protective transportation case	FFS-TC
Palltronic Flowstar V thermal printer paper rolls	S44.5.05
Palltronic Flowstar V power cable EU	FFS-PCEU
Palltronic Flowstar V power cable UK	FFS-PCUK
Palltronic Flowstar V power cable US	FFS-PCUS
Palltronic Flowstar V power cable Japan	FFS-PCJP
Palltronic Flowstar V power cable Switzerland	FFS-PCCH
Palltronic Flowstar V power cable Australia	FFS-PCAUS
Palltronic Flowstar V power cable China	FFS-PCZH
Palltronic Flowstar V USB stick	FFS-USB
Palltronic Flowstar V touchscreen pen	FFS-PEN
Palltronic Flowstar V external vent valve	FFS-BHVALVE
Palltronic Flowstar V MUX sequential filter integrity test accessory	FFS05MUX
Palltronic Flowstar V PAS/X software (1)	FFS05-PASX
mPath operational insights software (1) installation package	FFS-OIS
License key to activate software, expires after 1 year	FFS-OISKEY1
License key to activate software, expires after 3 years	FFS-OISKEY3
License key to renew subscription, expires after 1 year	FFS-OISKEY1N
License key to renew subscription, expires after 3 years	FFS-OISKEY3N

⁽¹⁾ For all marked software products, orders are required to mention the serial number of the Palltronic Flowstar V instrument for Cytiva to generate a complimentary license key for the software.



Fig 8. Bar code reader simplifies data entry and reduces operator input errors.

Protect your equipment and optimize productivity

We have more than 70 years of experience in providing maintenance, calibration and repair, ensuring seamless operation with minimal downtime and regulatory compliance throughout the lifecycle of your critical equipment.

Our range of service packages keeps your equipment protected and well-maintained, and includes itemized pay-as-you-go services, start-up care and training packages and a variety of post-warranty service plans that include priority response times, discounts for emergency repairs and flexible payment options. Our service plans provide total peace of mind and worry-free support throughout the coverage period.

cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva. mPath is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Active Directory is a trademark of Microsoft Corporation. DeviceNet is a trademark of Odva, Inc. GAMP is a trademark of International Society for Pharmaceutical Engineering, Inc . Linux is a registered trademark of Linux Torvalds in the United States, Modbus is a trademark of Modicon, PostScript is a trademark of Adobe Systems, Inc. Palltronic is a registered trademark of Pall Corporation. Profibus and Profinet are trademarks of Profibus Nutzerorganisation E.v. Staubli is a trademark of Staubli International AG. Any other trademarks are the property of their respective owners. © 2024 Cytiva

For local office contact information, visit cytiva.com/contact

CY40365-17Dec24-DF

