Cytiva™ Protein Select™ resin

AFFINITY CHROMATOGRAPHY

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
Deconvoluted mass spectrum

Protein mass with tag	22 141 Da
Expected protein mass after tag removal	16 876 Da
Obtained protein mass after purification with Cytiva Protein Select resin	16 876.7 Da

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Matrix	Highly cross-linked agarose
Particle size, d50V ¹	~ 60 µm
Ligand	Protein ligand derived from <i>E. coli</i>
Coupling chemistry	Epoxy
Yield ²	Up to 20 mg protein/mL resin
Chemical stability and compatibility	Stable in commonly used aqueous buffers
pH stability	Operational: 2 to 10 Regeneration/CIP: 4 M urea in 0.1 M NaOH
Recommended flow velocity ³	Recommended: 600 cm/h Maximum: 600 cm/h
Storage	2°C to 8°C, 20% ethanol
Delivery conditions	20% ethanol

¹ Median particle size of the cumulative volume distribution.

² Yield will depend on protein design and cleavage-time (hold time). Yield is the amount of cleaved target protein in the eluate.

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Experimental conditions with Cytiva Protein Select tag methodology

Sample: IL-1 β (2nd Pro amino acid changed to Phe amino acid to improve cleavage kinetic – see Practical considerations)
Column: HiTrap Protein Select column 1 mL
Flow rate: 4 mL/min
Cleavage time: 4 h
Buffer: 2 × PBS, pH 7.4
Load: *E. coli* lysate containing IL-1 β tagged with Cytiva Protein Select tag
System: ÄKTA pure™

Experimental conditions with histidine-tag methodology

1. IMAC bind-elute affinity purification of His-tagged IL-1 β

Sample: *E. coli* lysate containing His-tagged IL-1 β
Column: HisTrap HP 1 mL
Flow rate: 1 mL/min
Binding buffer: 2 × PBS, pH 7.4
Washing buffer: 2 × PBS, 20 mM imidazole
Elution buffer: 2 × PBS, 500 mM imidazole
System: ÄKTA pure

2. Buffer exchange before tag cleavage


Sample: Eluted His-tagged IL-1 β in IMAC elution buffer:
2 × PBS, 500 mM imidazole
PD-10 desalting column with Sephadex™ G-25 resin
Column: PD-10 desalting column with Sephadex™ G-25 resin
Sample buffer after buffer exchange: 30 mM Tris, 150 mM NaCl, pH 7.5

3. Tag cleavage

Protease: His-tagged TEV protease
Tag cleavage: Sample mixture of His-tagged IL-1 β and TEV protease on end-over-end rotation
Tag cleavage buffer: 30 mM Tris, 150 mM NaCl, pH 7.5
Tag cleavage temperature: Cold room
Tag cleavage time: Overnight

4. IMAC flow-through purification to recover the tag-free protein

Sample: Mixture of His-tagged IL-1 β and TEV protease after cleavage
Column: HisTrap HP 1 mL
Sample loading flow rate: 0.5 mL/min
Binding and washing buffer: 30 mM Tris, 150 mM NaCl, pH 7.5
Column washing flow rate: 1 mL/min
Elution buffer: 2 × PBS, 500 mM imidazole
Elution flow rate: 1 mL/min

- Totally cleaved IL-1 β protein from HiTrap Protein Select column 1 mL
- His-IL-1 β purified on HisTrap HP 1 mL followed by desalting
- TEV-protease cleavage of His-IL-1 β
- IL-1 β (flow-through) from HisTrap HP containing residual amino acids

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	His-tag	Cytiva Protein Select tag
Purity level	98% (contains tag residual amino acids)	98% (no remaining residual tag amino acids)
Time	~ 3 days	< 1 day
Number of steps	4	1
Comparison	- Extra time and effort to cleave and remove the tag - Some tag amino acids stay on the protein after cleavage	- Self cleavage of the tag - Completely traceless tag cleavage

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	AIEC + HIC	Cytiva Protein Select resin + AIEC
Convenience	- Developing a HIC chromatography step required extensive screening and optimization. - Required the addition of a non-chromatography step (heat treatment).	- This affinity chromatography does not require optimization. - No need for buffer screening.
Time spent for developing the process	Developing the full process was two times faster when Cytiva Protein Select resin was used in the first step	
Purity level	95%	Close to 100%

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
Flow rate: 5 mL/min

Sample load: 150 mL

Hydrophobic interaction chromatography

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Sample: Cleaved Interleukin-1 β purified with Cytiva Protein Select resin
Column: HiScale 16/40 Superdex 75 Increase
System: ÄKTA pure 25
Buffer: 30 mM Tris, 300 mM NaCl, pH 7.5
Purity analysis: Analytical SEC and SDS-PAGE

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- Top challenges in recombinant protein purification process development
- Tag affinity chromatography challenges faced by researchers today

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Format	Volume	Product code
Bulk resin	25 mL	17542101
	100 mL	17542102
	500 mL	17542103
HiTrap column	1 × 1 mL	17542151
	5 × 1 mL	17542152
	1 × 5 mL	17542153
	5 × 5 mL	17542154

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cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva. ÄKTA, ÄKTA avant, ÄKTA pure, AxiChrom, Capto, HiScale, HisTrap, HiTrap, ImageQuant, Protein Select, Sephadex, Sepharose, Superdex, and Tricorn are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Coomassie is a trademark of Thermo Fisher Scientific. Tween is a trademark of Croda Group of Companies. Any other third-party trademarks are the property of their respective owners.

Memo 30070145

© 2023 Cytiva

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY38320-05Dec23-DF