

## Procedure

# Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles

Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles (Table 1) use a non-surfactant, non-protein blocking reagent and provide low nonspecific binding (NSB), and high binding capacity for biotinylated target molecules.

Typical applications include:

- improvement and simplification of ligand binding
- affinity purifications
- immunoprecipitation
- protein interaction studies
- DNA: protein pulldowns
- purification of biotin-labeled proteins and nucleic acids
- other molecular biology applications

Biotinylated molecules are bound to the magnetic particles which are removed from the solution using a magnetic field. For manual processing, a simple magnetic stand can be used. An automated platform can also be used for magnetic purification.

Streptavidin is a  $M_r$  60 000 protein from *Streptomyces avidinii* (1). The protein is a tetramer containing four biotin-binding sites and is covalently coupled to the surface of blocked magnetic particles. There are two to three biotin binding sites available for each streptavidin molecule bound to the microparticle surface. Unlike avidin, streptavidin has a low isoelectric point (pI=5) and no carbohydrate groups, resulting in low non-specific binding. Furthermore, the protein is coupled to magnetic particles that exhibit very low non-specific binding in the presence of complex biological samples. The affinity between streptavidin and biotin is very high, requiring harsh conditions for disruption such as SDS-PAGE reducing sample buffer. Therefore, it is possible to elute binding partners in an interaction complex without co-eluting the biotinylated component.

Sera-Mag™ SpeedBeads combine the advantages of a high surface area, high affinity and high specific activity. They are colloiddally stable in the absence of a magnetic field. However, the particles can be separated rapidly and completely from suspension when a magnetic field is applied. Binding of biotinylated ligands to streptavidin groups on the surface is easily accomplished using standard avidin-biotin technology.

In the past, achieving high activity and stable binding of solid phase ligands has been a major difficulty. Compounds that are difficult to attach to microparticle surfaces by conventional means may be amenable to biotinylation. Due to the high affinity of the avidin-biotin reaction, binding biotinylated compounds to Sera-Mag™ SpeedBeads may improve specific activity. In such cases, biotinylation may be carried out in aqueous or organic solvent. Then, the biotin derivative can be bound to Sera-Mag™ SpeedBeads simply by mixing in appropriate buffer conditions. For example, nucleic acids which adsorb poorly to microparticle surfaces are readily bound to Sera-Mag™ SpeedBeads particles after biotinylation.

The use of magnetic particles as a solid phase support in immunoassays and molecular biology applications is well documented. Standard protocols are available to biotinylate a wide range of ligands including proteins, nucleic acids, haptens, peptides, and other molecules

Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles are uniform, colloiddally stable, monodispersed, non-porous superparamagnetic spheres made by a proprietary core shell process. The core is a carboxylate-modified particle made by free radical emulsion polymerization of styrene and acid monomer. Magnetite ( $Fe_3O_4$ ) is coated onto this core particle and then encapsulated with propriety polymers. Finally, the surface is blocked with a proprietary method to help prevent the nonspecific binding of proteins.

Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles are nominal 1  $\mu m$  particles with highly active streptavidin covalently bound to the surface. They are supplied at approximately 1% solids (10 mg/mL) in 0.05 % sodium azide.

**Table 1.** Characteristics of Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles

Composition	Streptavidin monolayer covalently coupled to bead surface
Magnetization	Superparamagnetic (no magnetic memory)
Mean diameter	1 $\mu m$ (nominal)
Bead concentration	10 mg/mL (bead weight/volume); 1% solids
Binding capacity (per mg of bead)	~ 3500 pmol biotinylated fluorescein
Particle density	~ 2.0 g/cm <sup>3</sup>

**Note:** Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles are not supplied in RNase-free solutions

## Important information before using Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles

- Do not freeze or dry Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles. This causes the particles to aggregate and lose binding activity.
- After labeling proteins or nucleic acids with biotin, remove unincorporated biotin with a desalting column. Free biotin will reduce the binding capacity of the particles.
- To minimize protein degradation, include protease inhibitors in the preparation of cell lysate.
- A low pH elution may be used for single-use applications. To limit leaching of streptavidin, do not exceed 10 minutes for the elution step in either manual or automated protocols.
- Boiling the magnetic particles in SDS-PAGE reducing sample buffer is acceptable for single-use applications. Boiling causes microparticle aggregation and loss of binding activity.
- Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles can be used successfully with mass spectrometry because the non-specific binding is very low.

## Procedure for manual immunoprecipitation using a biotinylated antibody

### Additional materials required

- Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles
- 1.5 mL microcentrifuge tubes
- Buffer set:
  - Hybridization buffer: 3.5x SSC, 5mM EDTA, 0.2% SDS
  - Wash buffer: 0.5x SSC
  - Elution buffer: Water
- Alternate elution buffer: SDS-PAGE reducing sample buffer
- Biotinylated antibody
- Antigen sample
- Cell lysis buffer (used to prepare antigen sample)
- Magnetic stand for 1.5 mL tube (e.g., MagRack 6, Cytiva code number 28948964)

## Prewashing Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles

**Note:** To ensure homogeneity, mix the particles thoroughly before use by repeated inversion, gentle vortexing or using a rotating platform.

1. Aliquot 50 µL (0.5 mg) of Sera-Mag™ SpeedBeads Streptavidin-Blocked Magnetic Particles into a 1.5 mL microcentrifuge tube.
2. Place the tube into a magnetic stand to collect the particles against the side of the tube. Remove and discard the supernatant.
3. Add 1 mL of binding/wash buffer to the tube. Invert the tube several times or vortex gently to mix. Collect the particles with a magnetic stand, then remove and discard the supernatant.

**Note:** Do not allow the particles to dry. If necessary, store them in binding/wash buffer prior to proceeding with purification protocol.

## Immunoprecipitation

**Note:** This protocol is a general guideline for immunoprecipitation and requires optimization for each application.

4. Combine antigen sample with 10 µg of biotinylated antibody. Incubate 1 to 2 hours at room temperature or overnight at 4°C with mixing.

**Note:** Dilute each sample to a minimum volume of 300 µL with cell lysis buffer or binding/wash buffer.

5. Add the antigen sample/biotinylated antibody mixture to a 1.5 mL microcentrifuge tube containing prewashed magnetic beads (see above) and incubate at room temperature for 1 hour with mixing.
6. Collect the particles with a magnetic stand, remove the supernatant and save for analysis.
7. Add 300 µL of binding/wash buffer to the tube and gently mix. Collect the particles and then discard the supernatant. Repeat twice.
8. Elution buffer recovery of antigen: Add 100 µL of elution buffer to the tube. Incubate the tube at room temperature with mixing for 5 minutes. Magnetically separate the particles and save the supernatant containing target antigen.

**Note:** If a low pH elution buffer is selected for elution, streptavidin may leach from the particles. Low pH elution buffers are effective for most antibody-antigen interactions. However, to ensure efficient release of target antigen from the antibody, prerinse the particles with 300 µL 0.1% Tween™-20 in water (no buffering capacity) before adding low pH elution buffer.

Alternate Elution: SDS-PAGE reducing sample buffer recovery of antigen: Add 100 µL of SDS-PAGE reducing sample buffer to the tube and heat the samples at 96°C to 100°C in a heating block for 5 minutes. Magnetically separate the particles and save the supernatant containing the target antigen.

**Note:** If SDS-PAGE buffer is selected for elution, the eluate will contain streptavidin monomers and dimers and biotinylated antibody along with target antigen.

## Troubleshooting

Problem	Possible cause	Solution
Low protein recovery	Proteolysis of sample	Add protease inhibitors
	Not enough magnetic beads used for capture	Increase the amount of magnetic beads used for capture
	Insufficient target protein present in sample	Increase amount of antigen sample
	Free biotin present in sample	Dialyze biotinylated antibody or pass it through a desalting column prior to binding to the magnetic beads
Protein does not elute	Elution conditions are too mild	Increase incubation time with elution buffer or use more stringent elution buffer
Multiple, nonspecific bands appear in eluted sample	Non-specific protein binding to the magnetic beads	Add 50–200 mM NaCl to the binding/wash and /or elution buffers
Magnetic beads aggregate	Magnetic beads were frozen or centrifuged	Handle the beads as directed in the instructions
	Buffer used is incompatible with magnetic beads	

## References

1. Chalet, I. and Wolf, F.J. The properties of streptavidin, a biotin-binding protein produced by *Streptomyces*. *Arch Biochem Biophys.* **106**, 1–5 (1964).

## cytiva.com/sera-mag

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva. Sera-Mag is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Tween is a trademark of the Croda Group of Companies. Any other third-party trademarks are the property of their respective owners.

© 2022 Cytiva

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

CY28490-11May22-PD

