



Rapid purification of GST-tagged proteins from large sample volumes

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

Rapid purification of GST-tagged proteins from large sample volumes

Lars Haneskog, Anna Heijbel, and Hans J. Johansson, GE Healthcare, Uppsala, Sweden

Johan Öhman and Eva Rupp-Thuresson, Bivitrum, Stockholm, Sweden

Introduction

The expression of recombinant proteins or domains fused to Glutathione S-transferase (GST) using the pGEX vectors allows for easy purification by affinity chromatography on Glutathione Sepharose™ under mild conditions. The GST-tag can be removed by a site-specific protease either on the column before elution, or after elution. We have developed a chromatography medium, Glutathione Sepharose 4 Fast Flow, for purification of GST-tagged proteins. It has high binding capacity for GST (10 mg/ml medium), and high rigidity makes the medium suitable for large-scale chromatography.

A GST-tagged protein was purified and cleaved by protease on a GSTrap™ FF 5-ml prepacked column to obtain a pure protein without the GST-tag. A second GST-tagged protein, expressed at a low level, was rapidly purified on a 34-ml column by application of a large sample volume at a high flow rate. The high volumetric flow rate allowed preparation to be done four times faster than before.

Summary and conclusions

- One-step purification on GSTrap FF was used to obtain a highly pure SH2 domain. SH2 was expressed as a GST-tagged protein and purified with on-column proteolytic removal of the GST-tag.
- Glutathione Sepharose 4 Fast Flow was used for rapid purification of a GST-tagged glycoprotein (expressed at a low level) by the application of a large sample volume to the column. This procedure offers a realistic alternative to the need to obtain higher expression levels by the optimization of the fermentation process or the preparation of a new gene construct.



Materials and methods

Samples

1. Clarified homogenate of *E. coli* containing the SH2 domain of a phosphatase fused with GST (SH2-GST) (M_r 37 000).
2. Clarified medium from a culture of Human Embryo Kidney cells (HEK293) expressing a M_r 120 000 glycosylated and secreted protein. Preliminary experiments revealed expression levels of 0.5 to 1.5 μ g GST-tagged protein per ml of culture medium.

Chromatography

Chromatography was done at room temperature on Glutathione Sepharose 4 Fast Flow using ÄKTAexplorer™ 10 and ÄKTAexplorer 100. GSTrap FF 1-ml and GSTrap FF 5-ml prepacked columns were used for purification and on-column cleavage of SH2-GST.

A 34-ml Glutathione Sepharose 4 Fast Flow medium bed (1.6 x 17 cm) was packed in a XK 16/20 column and used for purification of the low-expressed eukaryotic protein.

On-column cleavage with protease

After sample application and wash, the GSTrap FF 5-ml column was disconnected from the ÄKTAexplorer and filled using a syringe with 7 ml of 20 U/ml thrombin protease from GE Healthcare in binding buffer. The column was reconnected and left overnight (14 hours) at room temperature before elution.

Mass spectrometry analysis

A desalted sample of purified SH2 domain was adjusted to 50% acetonitrile/0.5% trifluoroacetic acid, before mixing with saturated sinapinic acid in 50% acetonitrile/0.5% trifluoroacetic acid, and air drying. Mass spectrometry analyses were performed in linear m/z window on a MALDI-TOF-MS instrument.

Results

Purification and on-column cleavage of an SH2-GST fusion protein

Chromatography of 2 ml of clarified *Escherichia coli* (*E. coli*) homogenate on GSTrap FF 1-ml, for estimation of expression level, yielded about 2 mg of SH2-GST fusion protein (Fig 1A). The eluted material mostly contained SH2-GST, although a small amount of GST was detected (Fig 1B). No detectable contaminants were recovered.

Column: GSTrap FF, 1 ml
Sample: 2 ml clarified *E. coli* homogenate expressing a M_r 37 000 SH2-GST fusion protein
Binding buffer: 150 mM NaCl, 20 mM phosphate buffer, pH 7.3
Elution buffer: 20 mM reduced glutathione, 50 mM Tris-HCl, pH 8.0
Flow rate: 2 ml/min (sample application and washing) and 0.5 ml/min (elution)
System: ÄKTAexplorer 10




Fig 1A. Purification of SH2 phosphatase domain-GST fusion protein.




Fig 1B. SDS-PAGE analysis on PhastGel™ Gradient 8–25, silver staining.

To obtain the pure SH2 domain without its GST-tag, chromatography was scaled-up on GSTrap FF 5 ml, and on-column cleavage was done overnight with thrombin before elution of released SH2 domain with binding buffer (Fig 2A). The eluted SH2 domain fraction contained 2 mg of protein and the GST fraction eluted by glutathione contained 4 mg. SDS-PAGE indicated that the prepared SH2 domain was pure and that the protease cleavage was complete (Fig 2B). Mass spectrometry revealed essentially two peaks corresponding to the single-charged (m/z 12 472) and double-charged protein (m/z 6 241), which agrees with the expected M_r of the SH2 domain (Fig 3). The spectra contained no other peaks in the m/z window used (insert in Fig 3).

Column: GSTrap FF, 5 ml
Sample: 100 ml clarified *E. coli* homogenate expressing a M_r 37 000 SH2-GST fusion protein
Binding buffer: 150 mM NaCl, 20 mM phosphate buffer, pH 7.3
Elution buffer: 10 mM reduced glutathione, 50 mM Tris-HCl, pH 8.0
Protease treatment: 20 U/ml thrombin protease (GE Healthcare) for 14 hours at room temperature
Flow rate: 10 ml/min (sample application and washing) and 2.5 ml/min (elution)
System: ÄKTApexplorer 10





Fig 2A. Purification of SH2 domain with concomitant removal of the GST-tag.




Fig 2B. SDS-PAGE analysis on PhastGel Gradient 8–25, silver staining.




Fig 3. MALDI-TOF-MS analysis of the SH2 domain.

Purification of a eukaryotic GST-tagged protein for crystallography

A large volume (1.5 liters) of clarified cell culture medium from HEK293 cells expressing small amounts of a M_r 120 000 glycosylated and secreted protein was applied to a 34-ml Glutathione Sepharose 4 Fast Flow column. The column was washed with binding buffer and 1 mg of pure protein was eluted by a step-gradient of glutathione (Fig 4), all within five hours. The protein was concentrated by ultrafiltration (cut-off M_r 10 000) and used for successful crystallization trials (Fig 5).

Column: 1.6 x 17 cm Glutathione Sepharose 4 Fast Flow (34 ml) packed in XK 16/20 column
Sample: 1500 ml clarified cell culture medium HEK293 expressing a M_r 120 000 glycosylated and secreted protein
Binding buffer: 140 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, 1.8 mM KH₂PO₄, pH 7.3
Elution buffer: 20 mM reduced glutathione, 50 mM NaCl, 50 mM Tris-HCl, pH 8.0
Flow rate: 10 ml/min (sample application and washing) and 1 ml/min (elution)
System: ÄKTApexplorer 100




Fig 4A. Purification of low-expressed eukaryotic protein.




Fig 4B. SDS-PAGE analysis on 4% to 12% gel followed by silver staining.



Fig 5. Crystals of the eukaryotic GST-tagged protein obtained in the initial crystallization trial.

www.gehealthcare.com/hitrap

GE Healthcare Bio-Sciences AB
 Björkgatan 30
 751 84 Uppsala
 Sweden

ÅKTAexplorer, Drop Design, Sepharose, PhastGel, and GSTRap are trademarks of GE Healthcare companies. GE, imagination at work, and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare that supplies them. General Electric Company reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2006 General Electric Company - All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Europe GmbH
 Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
 Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
 Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Pacific Tel +85 65 62751830 Fax +85 65 62751829 • Australasia Tel +61 2 8820 8299 Fax +61 2 8820 8200 • Austria Tel 01 57606 1613 Fax 01 57606 1614 • Belgium Tel 0800 73 890 Fax 02 416 8206 • Canada Tel 1 800 463 5800 Fax 1 800 567 1008 • Central & East Europe Tel +43 1 972 720 Fax +43 1 972 722 750 • Denmark Tel +45 70 25 24 50 Fax +45 45 16 2424 • Eire Tel 1 800 709992 Fax +44 1494 542010 • Finland & Baltics Tel +358 9 512 3940 Fax +358 9 512 39439 • France Tel 01 69 35 67 00 Fax 01 69 41 98 77 • Germany Tel 0800 9080 711 Fax 0800 9080 712 • Greater China Tel +852 2100 6300 Fax +852 2100 6338 • Italy Tel 02 26001 320 Fax 02 26001 399 • Japan Tel 81 3 5331 9336 Fax 81 3 5331 9370 • Korea Tel 82 2 6201 3700 Fax 82 2 6201 3803 • Latin America Tel +55 11 3933 7300 Fax +55 11 3933 7304 • Middle East & Africa Tel +30 210 96 00 687 Fax +30 210 96 00 693 • Netherlands Tel 0800-82 82 82 1 Fax 0800-82 82 82 4 • Norway Tel +47 815 65 777 Fax +47 815 65 666 • Portugal Tel 21 417 7035 Fax 21 417 3184 • Russia, CIS & NIS Tel +7 495 956 5177 Fax +7 495 956 5176 • Spain Tel 902 11 72 65 Fax 935 94 49 65 • Sweden Tel 018 612 1900 Fax 018 612 1910 • Switzerland Tel 0848 8028 10 Fax 0848 8028 11 • UK Tel 0800 515 313 Fax 0800 616 927 • USA Tel +1 800 526 3593 Fax +1 877 295 8102

