

Purification and chromatographic characterisation of an integral membrane protein

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit cytiva.com/contact

Purification and chromatographic characterisation of an integral membrane protein

- One-step purification of an integral membrane protein
- Chromatography at low temperature and in the presence of non-ionic detergents
- Chromatographic characterisation of the purified protein on the same system

Summary

ÄKTA_{FPLC}[™] controlled by UNICORN[™] software was used at cold room temperature for one-step purification and for characterisation of an integral membrane protein in the presence of non-ionic detergents. The target protein was a recombinant histidine-tagged cytochrome *bo*₃ ubiquinol oxidase expressed in *Escherichia coli*. Detergent extracts of bacterial membranes were applied to HiTrap[™] Chelating. Adsorbed proteins were eluted by an imidazol gradient in two well separated peaks, which both contained cytochrome *bo*₃.

The charge and size homogeneity of the purified protein was assessed under non-denaturing conditions, in the presence of non-ionic detergent by anion exchange chromatography on Mono Q[™] HR 5/5 and by gel filtration on Superdex[™] HR 10/30. A single sharp peak was observed for both fractions upon salt gradient elution of proteins adsorbed to the ion-exchanger, but minor peaks were detected for the low-imidazol fraction. Similarly, gel filtration gave essentially one single peak corresponding to monomeric protein. These results indicate that the preparations were homogenous in size and charge.

The application can be performed with the more recent ÄKTA_{purifier}[™] instead of the ÄKTA_{FPLC}. The Mono Q HR 5/5 and Superdex HR 10/30 are replaced by Mono Q 5/50 GL and Superdex 200 10/300 GL. HisTrap[™] HP or HisTrap FF crude can be used instead of HiTrap Chelating.




Fig. 1. Overview of the procedure.

Introduction

Detergent extracted membrane proteins can often be purified by the fractionation methods used for water-soluble proteins, although all steps of the purification must be performed in the presence of detergent. Membrane proteins in detergent solution frequently tend to aggregate or become denatured. Therefore, the concentration and type of detergent that is suitable for a particular fractionation must be tested for each situation. Due to the instability and tendency to associate there is often a need for fast purification protocols applied at low temperature.



Structure determination of proteins by crystallography and other high-resolution studies on membrane proteins usually demand that the protein can be obtained at high purity and homogeneity. Purity is usually taken to mean absence of contaminants, whereas homogeneity means similarity of the molecules in terms of conformation, association state, charge, etc. Purification of proteins should therefore not only be followed by purity analysis such as SDS-PAGE, but also by characterisation of the size and charge of the protein under non-denaturing conditions.

This application note shows the use of ÄKTA_{FPLC} with UNICORN control for chromatographic preparation and characterisation of a recombinant histidine-tagged integral membrane protein on prepacked 1 ml and 5 ml HiTrap Chelating columns at low temperature and in the presence of non-ionic detergent. The homogeneity of the preparations was characterised by anion exchange and gel filtration on Mono Q HR 5/5 (now replaced by Mono Q 5/50 GL) and Superdex 200 HR 10/30 (now replaced by Superdex 200 10/300 GL) columns.


Results and discussion

One-step purification

Membrane protein extract in dodecyl- β -D-maltoside solution was applied to a 1 ml prepacked HiTrap Chelating column loaded with nickel ions (Fig 2a). Elution was done in the presence of detergent by applying a gradient of imidazol concentration. Two distinct peaks were obtained and the material in fractions corresponding to both peaks had the characteristic colour of cytochrome bo_3 . SDS-PAGE showed that both fractions contained the four subunits expected for cytochrome bo_3 (Fig 3). Fraction 1 gave additional bands indicating contaminations. Fraction 2 showed only traces of contaminants. Scaling-up five times using 5 ml HiTrap Chelating gave the same result, consisting of 25 mg and 18 mg protein in Fraction 1 and 2, respectively (Fig 2b). Both cytochrome bo_3 preparations were kept in the refrigerator.

Sample: Detergent extracts of Escherichia coli membranes (see Materials and methods)
Binding buffer: 20 mM Tris-HCl, 5 mM imidazol, 0.03% dodecyl- β -D-maltoside and 300 mM NaCl, pH 7.5
Elution buffer: 20 mM Tris-HCl, 500 mM imidazol, 0.03% dodecyl- β -D-maltoside and 300 mM NaCl, pH 7.5
Gradient: 0 to 60% in 20 column volumes
Temperature: 5 °C
System: ÄKTA_{FPLC} with UNICORN control
Sample loop: Superloop 10 ml or 50 ml

a)
Sample volume: 8.5 ml
Column: HiTrap Chelating 1 ml, Ni²⁺-loaded
Flow: 1 ml/min



b)
Sample volume: 48.5 ml
Column: HiTrap Chelating 5 ml, Ni²⁺-loaded
Flow: 5 ml/min




Fig. 2. Purification on HiTrap Chelating.




Fig. 3. SDS electrophoresis on PhastSystem using PhastGel 8–25%, silver staining.


Charge characterisation

The cytochrome bo_3 preparations were assessed for purity and charge homogeneity under nondenaturing conditions by anion exchange chromatography on Mono Q HR 5/5 in the presence of the non-ionic detergent octyl β -D-glucoside. Fraction 1 eluted in one main peak (Fig 4a), but several small peaks were also observed that may be caused by trace contaminants. This agrees with the SDS-PAGE analysis (Fig 3). Fraction 2 gave a single sharp peak on Mono Q HR 5/5 (Fig 4b). These results confirm the high purity observed by SDS-PAGE and show that the material in Fraction 2 was highly homogenous by charge.

Column: Mono Q HR 5/5
 Binding buffer: 20 mM Tris-HCl, 1% octyl- β -D-glucoside, pH 7.5
 Elution buffer: 20 mM Tris-HCl, 1% octyl- β -D-glucoside, 1 M NaCl, pH 7.5
 Gradient: 0 to 50%, in 20 column volumes
 Flow: 1 ml/min
 Temperature: 5 °C
 System: ÄKTA_{APL}C with UNICORN control
 Sample loop: Superloop 10 ml

a)

Sample: 2.5 ml fraction 1-material diluted 10 times with binding buffer



Size characterisation

The cytochrome bo_3 preparations were analysed under non-denaturing conditions for size characterisation by gel filtration on Superdex HR 10/30 in the presence of octyl- β -D-glucoside (Fig 5). The main peak positions for both fractions were identical. Fraction 2 eluted in a narrower zone than Fraction 1, but there was no indication of large aggregates in the samples. The origin of the small unresolved peaks forming shoulders on the main peaks is unknown, but may correspond to optical effects caused by protein-free detergent micelles.

The column was calibrated by application of standard proteins from the LMW and HMW Gel filtration kits. Apparent molecular weights of the detergent-protein complexes were calculated from their elution positions by the Molecular weight calibration function. A molecular weight of $160\,000 \pm 5\,000$ was obtained for the detergent-protein complexes eluted in the main peak. This value is probably underestimated due to the dynamic structure and flexibility of the detergent moiety of the complexes, but indicated that cytochrome bo_3 is inserted in the detergent micelles as monomer. Chromatography in 1% octyl- β -D-glucoside revealed detergent-protein complexes of a molecular weight of $410\,000 \pm 3\,000$ (data not shown), indicating that cytochrome bo_3 formed dimers under these conditions.

b)

Sample: 3 ml fraction 2-material diluted 10 times with binding buffer




Fig. 4. Chromatographic characterisation on Mono Q HR 5/5.

Column: Superdex 200 HR 10/30

Sample application: 300 μ l applied via 0.1 ml sample loop

Buffer: 20 mM Tris-HCl, 1.5 % octyl- β -D-glucoside, 150 mM NaCl, pH 7.5


Flow: 0.25 ml/min

Temperature: 5 $^{\circ}$ C

System: ÄKTA_{FPLC} with UNICORN control


a)

Sample: Fraction 1 from HiTrap Chelating 1 ml



b)

Sample: Fraction 2 from HiTrap Chelating 1 ml



Materials and methods

Materials

All columns and sample loops were from GE Healthcare.

Sample preparation

Recombinant histidine-tagged cytochrome bo_3 from *Escherichia coli* was expressed in *E. coli* (1). This 130 000 Dalton enzyme is composed by four subunits (molecular weights 58 000, 33 000, 22 000, 17 000) of which the subunit II contain the histidine-tag on the carboxy-terminus. Bacterial membranes were prepared essentially as described by Rumbley *et al.* (1) and were solubilised by addition of 10 ml 1% dodecyl- β -D-maltoside, 300 mM NaCl in 20 mM Tris-HCl buffer, pH 7.5 per 1 g sedimented membranes followed by stirring on ice for 30 min. Non-solubilised material was removed by centrifugation at $40\,000 \times g$ for 20 min.

Chromatography

ÄKTA_{FPLC} with UNICORN was used for chromatography. The system was equipped with Pump P-920, Monitor UPC 900 with a 2 mm pathlength and 280 nm filter, Monitor pH/C-900, Fraction Collector Frac-900, mixer M-900 with 0.6-ml mixing chamber, Buffer A inlet valve (IV-908), Column position valves (PV-908), and an outlet valve (PV-908). Buffers were prepared at 22 $^{\circ}$ C to given pH. Chromatography was performed at 5 $^{\circ}$ C. Experimental details are given in the figures.

Acknowledgments

The purification protocol was developed in corporation with Dr. So Iwata and Mr. Jeff Abramson at Department of Biochemistry, Uppsala University, Sweden. This work was supported by EU grant BIO4-CT97-2357 (DG12-SSMI) "New approaches in the crystallisation and crystallography of membrane proteins".

Fig. 5. Chromatographic characterisation on Superdex 200 HR 10/30.

References

1. Rumbley, J. N., Furlong Nickels, E. and Gennis, R. B. (1997) One-step purification of histidine-tagged cytochrome bo_3 from *Escherichia coli* and demonstration that associated quinone is not required for the structural integrity of the oxidase. *Biochim. Biophys. Acta* 1340, 131–142.

Ordering information

Products	Quantity	Code no.
ÄKTA _{FFLC}		18-1118-67
HiTrap Chelating 1 ml	5	17-0408-01
HiTrap Chelating 5 ml	1	17-0409-01
Mono Q 5/50 GL	1	17-5166-01
Superdex 200 10/300 GL	1	17-5175-01
Superloop 10 ml	1	19-7585-01
Superloop 5 ml	1	19-7850-01

Related products

ÄKTApurifier 10		28-4062-64
Superdex 200 5/150 GL	1	28-9065-61
HiTrap HP 1 ml	5	17-5247-01
HiTrap HP 5 ml	1	17-5248-01
HiTrap IMAC HP 1 ml	5	17-0920-03
HiTrap IMAC HP 5 ml	5	17-0920-05

Accessories

Gel Filtration LMW Calibration Kit	1	28-4038-41
Gel Filtration HMW Calibration Kit	1	28-4038-42

Related product literature

Literature	Code no.
Application note: Purification of a labile, oxygen-sensitive enzyme for crystallisation and 3D structure determination	18-1128-91
Rapid optimisation and development of an automated two-step purification procedure for monoclonal IgG antibodies	18-1128-93
Data File: ÄKTAexplorer™ System Series	18-1124-09
Data File: ÄKTApurifier	18-1119-48
Data File: Monitor UPC-900	18-1128-40
Data File: UNICORN control system	18-1111-20
Purifying Challenging Proteins Handbook	28-9095-31
Gel Filtration Handbook	18-1022-18
Recombinant Protein Purification Handbook	18-1142-75

For contact information for your local office,
please visit, www.gelifesciences.com/contact

www.gelifesciences.com/trap

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden



imagination at work

GE, imagination at work and GE monogram are trademarks of General Electric Company.

ÅKTA_{FLC}, ÅKTApurifier, ÅKTAexplorer, HiTrap, HisTrap, Mono Q, PhastGel, PhastSystem, Superdex, Superloop, and UNICORN are trademarks of GE Healthcare companies.

Histidine-tagged protein purification. Purification and preparation of fusion proteins and affinity peptides comprising at least two adjacent histidine residues may require a license under US patent numbers 5,284,933 and 5,310,663 and equivalent patents and patent applications in other countries (assignee: Hoffman La Roche, Inc).

All third party trademarks are the property of their respective owners.

© 1998-2007 General Electric Company – All rights reserved.

First published March 1998

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information

GE Healthcare Limited, Amersham Place, Little Chalfont, Buckinghamshire,
HP7 9NA, UK

GE Healthcare Bio-Sciences Corp., 800 Centennial Avenue, P.O. Box 1327, Piscataway,
NJ 08855-1327, USA

GE Healthcare Europe GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences KK, Sanken Bldg., 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo, 169-0073 Japan