



# One-step purification of monoclonal IgM from cell culture supernatants

**Intellectual Property Notice:** The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

## cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.  
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.  
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.  
All other third-party trademarks are the property of their respective owners.  
© 2020 Cytiva  
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.  
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY13098-01 May20-AN

# One-step purification of monoclonal IgM from cell culture supernatants

*Martin Hall, Annika Bergenstråhle, and Anna Heijbel*

*GE Healthcare, 751 84 Uppsala, Sweden*

## Introduction

The development of a general purification protocol for immunoglobulin M (IgM) has met with problems; the unique variation in the amino acid sequence of the variable regions gives different monoclonal antibodies and different chromatographic characteristics. This problem can be circumvented by using thiophilic interaction chromatography for the purification of IgM. The exact mechanism of thiophilic interaction between protein and ligand is not fully understood, but adsorption is promoted by water-structuring salts (1).

HiTrap™ IgM Purification HP columns are packed with a thiophilic adsorption medium consisting of the ligand 2-mercaptopyridine coupled to Sepharose™ High Performance. The column is well suited for purification of IgM, and the main application area for HiTrap IgM Purification HP is purification of monoclonal IgM from hybridoma cell culture supernatants.

## Summary

- Monoclonal immunoglobulin M (IgM) from hybridoma cell culture supernatants was purified in a single chromatography step by using the HiTrap IgM Purification HP column.
- The binding of IgM to the HiTrap column was accomplished by adding ammonium or potassium sulphate (0.8 to 1.0 M or 0.5 M final concentration, respectively) to the IgM-containing cell culture supernatant.
- Bound IgM was eluted with 20 mM sodium phosphate, pH 7.5.
- The eluted material was estimated to contain more than 80% IgM, as judged by gel filtration.

## Acknowledgments

We gratefully acknowledge Dr. Nils Linde, Medicago AB, Sweden, Imma Brogren, Dept. of Genetics and Pathology, Uppsala University, Sweden, and Anette Holtz, Phadia, Uppsala, Sweden, for their contributions to the work.



Methods and results

IgM purification

All chromatographic separations were performed on ÄKTAexplorer™ and ÄKTApurifier™ chromatography systems. Binding of IgM to the HiTrap IgM Purification HP column was obtained in all samples tested (purified human polyclonal IgM, purified murine monoclonal IgM, and cell culture supernatants from hybridoma cell lines expressing monoclonal murine IgM). Samples up to 100 ml of IgM-containing (10 µg/ml) cell culture supernatant were applied to the column, without detecting any IgM in the flow-through (Table 1). This gives a dynamic capacity of at least 1 mg IgM/ml medium. The IgM elutes when applying a salt-free buffer. Any remaining material on the column can be eluted with 30% isopropanol. See Figure 1 for a typical chromatogram.

Up to 100% of the applied IgM was recovered in the eluate (Table 1) from samples containing 0.5 M potassium sulphate or 0.8 M ammonium sulphate. The purity of the eluted IgM was 80% or higher, as judged by analytical gel filtration (Figures 3 and 4b).

Column: HiTrap IgM Purification HP, 1 ml (0.7 × 2.5 cm)  
Sample: 75 ml of cell culture supernatant containing α-S IgM, 0.5 M potassium sulphate  
Binding buffer: 20 mM sodium phosphate, 0.5 M potassium sulphate, pH 7.5  
Elution buffer: 20 mM sodium phosphate, pH 7.5  
Cleaning buffer: 20 mM sodium phosphate, 30% isopropanol, pH 7.5  
Flow rate: 1 ml/min




Fig 1. Purification of murine monoclonal α-Shigella (α-S) IgM on HiTrap IgM Purification HP.




Fig 2. Non-reducing SDS-PAGE analysis on PhastGel™ Gradient 4–15, silver stained.




Fig 3. Gel filtration analysis of eluted α-S IgM.

Table 1. Recovery of anti-Shigella (α-S) IgM measured by ELISA

Sample	α-S IgM conc. (mg/ml)	Volume (ml)	Recovery (%)
Cell culture supernatant	0.01	100	100
Flowthrough pool, 0 to 90 ml	not detectable	90	0
Flowthrough pool, 90 to 100 ml	not detectable	10	0
Eluate	0.1 to 0.3	9	> 100
Cleaning eluate	approx. 0.1	1	approx. 10

## Effects of increasing salt concentrations

If the recovery of IgM is not satisfactory, the ammonium sulphate concentration can be raised to 1.0 M in the sample. This improves the recovery of IgM at the expense of the purity, as more contaminants in the sample bind and co-elute. At ammonium sulphate concentrations of 1.2 M or higher, there is a strong increase in binding of contaminants from the cell culture medium (Dulbecco's Modified Eagle Medium with 5% Fetal Calf Serum) to the column (Fig 4).

## References

1. Porath, J. and Belew, M. *Trends in Biotechnology*, **5**, 225–229 (1987).

Column: Superdex 200, 10 × 300 mm

Samples: **A.** 200 µl cell culture supernatant containing a murine monoclonal  $\alpha$ -sheep red blood cell IgM was applied to the HiTrap IgM Purification HP column at different ammonium sulphate concentrations (0.8 to 1.2 M)


**B.** 200 µl HiTrap IgM Purification HP eluate, 0.8 M ammonium sulphate in sample and binding buffer

**C.** 200 µl HiTrap IgM Purification HP eluate, 1.0 M ammonium sulphate in sample and binding buffer

**D.** 200 µl HiTrap IgM Purification HP eluate, 1.2 M ammonium sulphate in sample and binding buffer

Buffer: 20 mM sodium phosphate, 0.15 M NaCl, pH 7.2

Flow rate: 0.4 ml/min



**Fig 4.** Gel filtration analysis of cell culture supernatant (A) and eluate (B-D) from HiTrap IgM Purification HP.

[www.gehealthcare.com/hitrap](http://www.gehealthcare.com/hitrap)

GE Healthcare Bio-Sciences AB

Björkgatan 30

751 84 Uppsala

Sweden

Drop Design, Sepharose, Superdex, HiTrap, PhastGel, ÄKTAexplorer, and ÄKTApurifier are trademarks of GE Healthcare companies. GE, imagination at work, and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare that supplies them. General Electric Company reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2006 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd  
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp  
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK  
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

Asia Pacific Tel +85 65 62751830 Fax +85 65 62751829 • Australasia Tel +61 2 8820 8299 Fax +61 2 8820 8200 • Austria Tel 01/57606 1613 Fax 01/57606 1614 • Belgium Tel 0800 73 890 Fax 02 416 8206 • Canada Tel 1 800 463 5800 Fax 1 800 567 1008 • Central & East Europe Tel +43 1 972 720 Fax +43 1 972 722 750 • Denmark Tel +45 70 25 24 50 Fax +45 45 16 2424 • Eire Tel 1 800 709992 Fax +44 1494 542010 • Finland & Baltics Tel +358 9 512 3940 Fax +358 9 512 39439 • France Tel 01 69 35 67 00 Fax 01 69 41 98 77 • Germany Tel 0800 9080 711 Fax 0800 9080 712 • Greater China Tel +852 2100 6300 Fax +852 2100 6338 • Italy Tel 02 26001 320 Fax 02 26001 399 • Japan Tel 81 3 5331 9336 Fax 81 3 5331 9370 • Korea Tel 82 2 6201 3700 Fax 82 2 6201 3803 • Latin America Tel +55 11 3933 7300 Fax +55 11 3933 7304 • Middle East & Africa Tel +30 210 96 00 687 Fax +30 210 96 00 693 • Netherlands Tel 0800-82 82 82 1 Fax 0800-82 82 82 4 • Norway Tel +47 815 65 777 Fax +47 815 65 666 • Portugal Tel 21 417 7035 Fax 21 417 3184 • Russia, CIS & NIS Tel +7 495 956 5177 Fax +7 495 956 5176 • Spain Tel 902 11 72 65 Fax 935 94 49 65 • Sweden Tel 018 612 1900 Fax 018 612 1910 • Switzerland Tel 0848 8028 10 Fax 0848 8028 11 • UK Tel 0800 515 313 Fax 0800 616 927 • USA Tel +1 800 526 3593 Fax +1 877 295 8102



imagination at work