

Benzamidine Sepharose 4 Fast Flow (high sub)

HiTrap Benzamidine FF (high sub)

AFFINITY CHROMATOGRAPHY

HiTrap™ Benzamidine FF (high sub) are prepacked, ready to use 1 ml and 5 ml columns containing Benzamidine Sepharose™ 4 Fast Flow (high sub). These products are excellent tools for removal of proteolytic activity from a protein or peptide preparation or purification of trypsin and trypsin-like serine proteases such as thrombin and enterokinase (Table 1).

HiTrap Benzamidine FF (high sub) is also convenient to use for removal of enzymes (e.g., thrombin) after cleavage of the tag from a recombinant fusion protein.

Fast, simple and easy separations are provided by the combination of the prepacked column with a high binding capacity affinity medium that has good flow properties. HiTrap Benzamidine FF (high sub) columns can be operated with a syringe, pump or liquid chromatography system.

- Fast and convenient to use
- Very high binding capacity
- Simple operation with a syringe, a pump, an ÄKTA™ system, or other chromatography systems
- Scaling-up is easily achieved by connecting several 1 ml or 5 ml columns in series or using lab packages

Medium characteristics

Benzamidine Sepharose 4 Fast Flow (high sub) is based on highly cross-linked 4% agarose enabling high flow rates without losing binding capacity.

p-Aminobenzamidine (pABA), a synthetic inhibitor of trypsin and trypsin-like serine proteases, is covalently coupled via an amide bond to a long spacer arm attached to Sepharose 4 Fast Flow via a stable ether linkage (Fig 2).

Table 2 shows the main characteristics of HiTrap Benzamidine FF (high sub) and Benzamidine Sepharose 4 Fast Flow (high sub).

Figure 3 shows the pH stability of the medium.



Fig 1. Benzamidine Sepharose 4 Fast Flow (high sub) and prepacked HiTrap Benzamidine FF (high sub) are designed for removal and/or purification of serine proteases.




Fig 2. Partial structure of Benzamidine Sepharose 4 Fast Flow (high sub).

Table 1. Examples of different serine proteases

Enzyme	Source	M _r	pI
Thrombin	Bovine pancreas	23 345	10.5
Trypsin	Human plasma:		
	chain A	5 700	
	chain B	31 000	7.1
Urokinase	Human urine	54 000	8.9
Enterokinase	Porcine intestine:		
	heavy chain	134 000	4.2
	light chain	62 000	
Plasminogen	Human plasma	90 000	6.4–8.5
Prekallikrein	Human plasma	nd	nd
Kallikrein	Human plasma	86 000	nd
Kallikrein	Human saliva	nd	4.0

Column characteristics

HiTrap Benzamidine FF (high sub) is prepacked 1 ml and 5 ml columns made of polypropylene, a biocompatible material that does not interact with biomolecules. Top and bottom frits are manufactured from porous polyethylene. The column is delivered with a stopper on the inlet and a twist-off end on the outlet. Several columns can be connected in series for easy scale-up. HiTrap columns can not be opened or refilled.


Table 2. Main characteristics of HiTrap Benzamidine 4 FF (high sub) and Benzamidine Sepharose 4 FF (high sub)

Matrix	highly cross-linked 4% spherical agarose	
Average particle size	90 µm	
Ligand	p-aminobenzamidine (pABA)	
Ligand concentration	≥ 12 µmol p-aminobenzamidine/ml medium	
Spacer	14-atom	
Binding capacity	≥ 35 mg trypsin/ml medium	
Column dimensions (i.d.×h)	0.7 × 2.5 cm (1 ml) 1.6 × 2.5 cm (5 ml)	
Column volumes	1 ml and 5 ml	
Recommended flow rates	1 ml/min and 5 ml/min for 1 and 5 ml columns respectively	
Column hardware pressure limit	5 bar (0.5 MPa, 70 psi)	
Maximum flow rates*	4 ml/min and 20 ml/min for 1 and 5 ml columns respectively	
Chemical stability	All commonly used aqueous buffers	
pH stability**		
	cleaning	pH 1 to 9
	working	pH 2 to 8
Storage temperature	4°C to 8°C	
Storage buffer	20% ethanol in 0.05 M acetate buffer, pH 4	

* H₂O at room temperature


** The ranges given are estimates based on our knowledge and experience. Please note the following: pH stability, cleaning refers to the pH interval for regeneration, cleaning-in-place and sanitization procedures.

pH stability, working refers to the pH interval where the medium is stable over a long period of time without adverse effects on its subsequent chromatographic performance. See Fig 3.

**Fig 3.** Benzamidine Sepharose 4 Fast Flow (high sub) has been tested for trypsin capacity after being stored at pH 1–11 at ambient temperature for 24 h, 1 week and 4 weeks.

Operation

Like all HiTrap columns, HiTrap Benzamidine FF (high sub) is quick and convenient to use. Instructions and connectors are included with each pack of columns. In general, the separation can be easily achieved with a syringe (using the luer adapter provided). Figure 4 A–C illustrates this technique. Alternatively, the column can be operated using a laboratory pump or a chromatography system.

**Fig 4.** Using HiTrap Benzamidine FF (high sub) with a syringe. (A) Prepare buffers and sample. Remove the column's top cap and twist off the end. Wash and equilibrate. (B) Load the sample and begin collecting fractions. (C) Wash, elute and continue collecting fractions.

Scaling-up

HiTrap Benzamidine FF (high sub) is available as 1 ml and 5 ml columns. Two or more columns can easily be connected in series by screwing the end of one into the top of the next for fast scaling-up. Note that the backpressure will increase. Benzamidine Sepharose 4 Fast Flow (high sub) is also available as lab packages for further scaling-up.

Applications

Removal of trypsin-like serine proteases from human plasma

Proteases included in human plasma can damage the sample if not removed. This application shows the use of HiTrap Benzamidine FF (high sub) for removal of trypsin-like serine proteases in one step. The 1 ml column was equilibrated with 10 ml binding buffer, 1 ml human plasma was loaded and the column was washed with 10 ml binding buffer before a one step elution using 5 ml elution buffer.

The results show that according to the activity test, almost all trypsin-like serine protease activity is removed from the human plasma sample and bound to the column.

Sample: 1 ml human plasma filtered through a 0.45 µm filter
Column: HiTrap Benzamidine FF (high sub), 1 ml
Binding buffer: 20 mM Tris-HCl, 0.5 M NaCl, pH 7.4
Elution buffer: 50 mM glycine, pH 3.0
0-100% elution buffer in one step
Flow rate: 1.0 ml/min
Instrumentation: ÄKTAE explorer 10
Protease activity: S-2288 from Haemochrom Diagnostica AB A_{405} measurement. The activity is presented as the proteolytic activity/mg protein




Fig 5. Removal of trypsin-like serine proteases from human plasma using HiTrap Benzamidine FF (high sub), 1 ml.

Removal of thrombin after on-column cleavage of a GST-tagged protein

1. Isolation of GST-tagged protein on GSTRap FF

The SH2-GST fusion protein sample (2 ml, clarified *E. coli* homogenate) was loaded on a GSTRap™ FF 1 ml column equilibrated in binding buffer. Non-binding proteins were washed out using 10 ml binding buffer.

2. On-column cleavage using thrombin

Between the SH2 protein and the GST-tag, a thrombin cleavage site was incorporated. To separate the protein from its GST-tag, the serine protease, thrombin, was dissolved in binding buffer (1 ml, 20 units/ml) and applied to the column using a syringe. The column was sealed with the supplied connectors and left to incubate for 2 hours at room temperature.

Note: Cleavage conditions have to be optimized for each case. For a sensitive protein fast cleavage may be preferred over complete cleavage.

3. Sample purification

HiTrap Benzamidine FF (high sub) 1 ml column was equilibrated with water and binding buffer before being placed after the GSTRap FF column. Placed in series, the columns were washed with 7 ml binding buffer and later with 5 ml high salt buffer. This procedure allowed the cleaved SH2 protein and thrombin to be washed out from the GSTRap FF column and thrombin to bind when passing the HiTrap Benzamidine FF (high sub) column. Collected fractions contained pure cleaved SH2 protein (Fig 7).

4. Elution of HiTrap Benzamidine FF (high sub)

HiTrap Benzamidine FF (high sub) 1 ml column was eluted competitively using 10 ml of 20 mM p-aminobenzamidine in binding buffer. Since p-aminobenzamidine by itself gives a high absorbance at 280 nm an enzymatic activity assay for detection of thrombin was used (Fig 6).

5. Elution of GSTRap FF

GSTRap FF 1 ml column was eluted using 10 ml of 20 mM reduced glutathione in 50 mM Tris, pH 8.0 (reduced glutathione has a slight absorbance at 280 nm).

Thrombin activity assay

A small amount of thrombin was used in this on-column cleavage of SH2-GST. The presence of thrombin could neither be determined by SDS-PAGE analysis nor by absorbance (280 nm), due to the small amount and the usage of p-aminobenzamidine for elution. Therefore the chromogenic substrate S-2238 (Haemochrom Diagnostica AB) was used to determine the thrombin activity and thereby verify the removal of protease.

Results

GSTRap FF efficiently bound and purified the SH2-GST protein. On-column cleavage using thrombin was almost complete in only 2 hours. (To achieve full cleavage the incubation time or amount of protease can be increased). HiTrap Benzamidine FF (high sub) 1 ml bound the total added thrombin, as shown by the activity measurement in Figure 6, resulting in pure cleaved SH2 protein (2 mg) without contamination of used protease (Fig 7). This whole procedure is completed in less than one day.

Sample: 2 ml clarified *E. coli* homogenate expressing SH2-GST (M_r 37 000) with a thrombin cleavage site
Columns: GSTrap FF, 1 ml
 HiTrap Benzamidine FF (high sub), 1 ml
Binding buffer: 20 mM Na phosphate, 0.15 M NaCl, pH 7.5
High salt wash buffer: 20 mM Na phosphate, 1.0 M NaCl, pH 7.5
HiTrap Benzamidine FF (high sub) elution buffer: 20 mM p-aminobenzamidine in binding buffer
GSTrap FF elution buffer: 50 mM Tris, 20 mM reduced glutathione, pH 8.0
Flow rate: 0.5 ml/min
System: ÄKTAprime
Protease treatment: 1 ml, 20 units/ml thrombin (Cytiva) for 2 hours at room temperature
Thrombin activity: S-2238 from Haemochrom Diagnostica AB, A₄₀₅ measurement




Fig 6. Removal of thrombin after on-column cleavage of a GST-tagged protein.

Cleaning

The cleaning protocol has to be designed for each protein. General recommendation is to use a solution of guanidine hydrochloride to remove precipitated or denatured substances. For hydrophobically bound substances a solution of nonionic detergent or ethanol is recommended. Short term use of pH 1 to 9 is possible. For example, more than half of the trypsin capacity is available after 1 week (168 h) contact time at pH 9 (Fig 3). However, prolonged exposure to pH greater than 8 and lower than 2 should be avoided due to hydrolysis of the ligand at high pH and decomposition of the matrix at low pH.

Storage

HiTrap Benzamidine FF (high sub) and Benzamidine Sepharose 4 Fast Flow (high sub) are supplied in 0.05 M acetate buffer, pH 4 containing 20% ethanol as a bacteriostat. Storage temperature is 4°C to 8°C.




Fig 7. SDS-PAGE on ExcelGel™ SDS Gradient 8% to 18%, Coomassie™ staining.

Ordering information

Product	Quantity	Code number
HiTrap Benzamidine FF (high sub)	2 × 1 ml	17-5143-02
HiTrap Benzamidine FF (high sub)	5 × 1 ml	17-5143-01
HiTrap Benzamidine FF (high sub)	1 × 5 ml	17-5144-01
Benzamidine Sepharose 4 Fast Flow (high sub)	25 ml	17-5123-10

Related products	Quantity	Code number
HiTrap Desalting	1 × 5 ml	29-0486-84
HiTrap Desalting	5 × 5 ml	17-1408-01
HiPrep™ 26/10 Desalting	1 × 53 ml	17-5087-01
HiPrep 26/10 Desalting	4 × 53 ml	17-5087-02
GSTrap FF	2 × 1 ml	17-5130-02
GSTrap FF	5 × 1 ml	17-5130-01
GSTrap FF	1 × 5 ml	17-5131-01
GSTrap FF	5 × 5 ml	17-5131-02
Thrombin	500 units	27-0846-01
Factor Xa	400 units	27-0849-01

Accessories	Quantity	Code number
1/16" male/Luer female*	2	18-1112-51
Tubing connector flangeless/M6 female	2	18-1003-68
Tubing connector flangeless/M6 male	2	18-1017-98
Union 1/16" female/M6 male	6	18-1112-57
Union M6 female/1/16" male	5	18-3858-01
Union luerlock female/M6 female	2	18-1027-12
HiTrap/HiPrep, 1/16" male connector for ÄKTA design	8	28-4010-81
Stop plug female, 1/16" [†]	5	11-0004-64
Fingertight stop plug, 1/16" [‡]	5	11-0003-55

* One connector included in each HiTrap package

[†] Two, five, or seven stop plugs female included in HiTrap packages depending on products.

[‡] One fingertight stop plug is connected to the top of each HiTrap column at delivery.

Related literature	Quantity	Code number
The Recombinant Protein Purification Handbook, Principles and Methods	1	18-1142-75
Affinity Chromatography Handbook, Principles and Methods	1	18-1022-29
Affinity Columns and Media, Selection Guide	1	18-1121-86
Convenient Protein Purification, HiTrap Column Guide	1	18-1129-81

cytiva.com/hitrap

cytiva.com/protein-purification

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. ÄKTA, ExcelGel, GSTrap, HiPrep, HiTrap, and Sepharose are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Coomassie is a trademark of Thermo Fisher Scientific LLC. All other third party trademarks are the property of their respective owner.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

CY14297-11Sep20-DF

