



Highly sensitive phosphopeptide analysis using Ettan MDLC and a linear ion trap mass spectrometer

Intellectual Property Notice: Part of GE Healthcare's Life Sciences business was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

Highly sensitive phosphopeptide analysis using Ettan MDLC and a linear ion trap mass spectrometer

Key words: *phosphopeptide • liquid chromatography • Ettan MDLC • biocompatible flowpath • reversed-phase chromatography (RPC) • dual-mode RPC • high-throughput method • neutral loss • MS/MS • MS³ • linear ion trap*

Ettan™ MDLC coupled to a linear ion trap mass spectrometer was used to characterize phosphopeptides in the low femtomole range. Phosphopeptides in a complex biological sample were sequenced, and the phosphorylation sites were identified.

Introduction

Because of the importance of phosphorylation and dephosphorylation to the function of many metabolically active proteins (1, 2), there is a great interest in improving the methods for reliable analysis of phosphopeptides. Unfortunately, phosphoproteins often exist in low abundance in cells, and the corresponding phosphopeptides resulting from digestion exhibit low ionization efficiency in positive electrospray ionization (ESI). Furthermore, the frequently low stoichiometry of phosphorylation in proteins (3) makes it important to increase the sensitivity of the analysis.

One way to increase the sensitivity is to ensure that no peptides are lost or degraded during the analysis. Because Ettan MDLC is a biocompatible system, there are no metal ions in the flowpath that can complex with phosphate groups (4). In addition, the biocompatible materials that make up the pumpheads and flowpath promote the longevity of the system, especially when used with salt-containing buffers.

We developed a method for femtomole-scale phosphopeptide analysis using Ettan MDLC and a Finnigan™ LTQ™ linear ion trap mass spectrometer. The phosphopeptide analysis procedure consisted of the following steps:

- Desalting the peptide sample on an RPC trap column using Ettan MDLC
- Separation of the peptide sample on a nanoscale analytical RPC column using Ettan MDLC
- Data-dependent neutral loss scan in the Finnigan LTQ mass spectrometer
- Phosphopeptide identification by database searching using TurboSEQUENT™ protein identification software

The method was developed using β-casein trypsin digest, either alone or mixed with a six-protein digest. To test the method, phosphorylation sites in a complex biological sample—mouse brain tissue—were identified.

Materials

Products used

Ettan MDLC	18-1176-44
NAP™ 10 Columns	17-0854-01
PlusOne™ DTT	17-1318-01
PlusOne Iodoacetamide	RPN6302
PlusOne Tris	17-1321-01
Trypsin, sequencing grade	17-6002-75



Other products required

Finnigan LTQ mass spectrometer (Thermo Electron)

TurboSEQUENT protein identification software (Thermo Electron)

PepMap™ C18 RPC trap column, 300 µm i.d. x 5 mm, 5 µm, 100 Å (LC Packings)

PepMap C18 RPC analytical column, 75 µm i.d. x 150 mm, 5 µm, 100 Å (LC Packings)

Zorbax™ 300-SB C18 trap column, 300 µm i.d. x 5 mm, 3 µm (Agilent)

Zorbax 300-SB C18 analytical column, 75 µm i.d. x 150 mm, 3 µm (Agilent)

β-Casein, bovine (Sigma-Aldrich)

Protein Mixture Digest (LC Packings)

Formic acid, ultrapure (Merck Suprapur™)

Water, HPLC grade

Acetonitrile, HPLC grade

Ammonium bicarbonate (Merck)

Methods

Sample preparation

Approximately 0.5 mg of mouse brain tissue was prepared according to the following procedure. The proteins were dissolved by adding 1 ml of 9 M urea with 50 mM DDT to the tissue and allowing it to incubate for 60 min at 20 °C. A 1-ml aliquot consisting of 8 M urea, 250 mM TrisHCl, and 125 mM iodoacetamide, pH 8.8, was added. The mixture was allowed to incubate for 60 min at 20 °C. A 1-ml aliquot of this mixture was buffer exchanged with 20 mM ammonium bicarbonate, pH 7.8, on a NAP-10 desalting column. The protein sample was digested with trypsin (concentration ratio of 50:1) for 4 h. The trypsin was then inactivated by adding formic acid to the sample.

β-Casein was dissolved in 50 mM ammonium bicarbonate, and then digested with trypsin.

A 1-fmol/µl solution of a six-protein digest (cytochrome C, lysozyme, alcohol dehydrogenase, bovine serum albumin, apotransferrin, and β-galactosidase; LC Packings) was mixed with a 1-fmol/µl solution of β-casein digest (equimolar ratio).

High-throughput nanoscale liquid chromatography

Ettan MDLC was used in the high-throughput configuration using two RPC trap columns 300 µm i.d. x 5 mm (PepMap, C18, 5 µm, 100 Å or Zorbax 300-SB C18) for online desalting and sample cleanup, followed by two nanoscale RPC analytical columns 75 µm i.d. x 150 mm (PepMap, C18, 5 µm, 100 Å or Zorbax 300-SB C18) for high-resolution separation.

One set of RPC trap/analytical columns was equilibrated while the second set separated the sample. The mobile phases were A: HPLC-grade water with 0.1% formic acid, and B: 84% HPLC-grade acetonitrile with 0.1% formic acid. A gradient of 0–56% B was run at a flow rate of 250 nl/min for 48 min.

Mass spectrometry

A Finnigan LTQ linear ion trap mass spectrometer equipped with a nanospray interface was used as the detector for the peptides that eluted from the RPC column. The MS method consisted of a cycle combining one full MS scan with three MS/MS events (25% collision energy) followed by an MS³ event (35% collision energy) that was triggered upon detection of -98, -49, or -32.7 Da from the precursor (neutral loss of phosphoric acid). Dynamic exclusion duration was set to 30 s.

Peptide identification

The MS/MS and MS³ spectra from all the runs were searched using TurboSEQUENT protein identification software (Thermo Electron). Modifications were set to allow for the detection of oxidized Met (+16), carboxyamidomethylated Cys (+57), phosphorylated Ser, Thr, and Tyr (+80), and dehydrated Ser and Thr (-18). The peptide matches were filtered based on cross-correlation scores (Xcorr) of 1.5, 2, and 2.5 for charge states 1⁺, 2⁺, and 3⁺, respectively. All spectra in the result lists were manually confirmed. The spectra were considered good when 50% of the b and y ions were assigned, including amino acids with phosphate groups.

Results and discussion

A 200-fmol trypsin digest of β-casein was analyzed by LC-MS/MS for instrument tuning and method development. The averaged spectra at a retention time of 28 min are shown in Figure 1 (top). The signal observed at *m/z* 982.5 is due to the loss of the neutral fragment from the precursor ion at *m/z* 1032.4 (Fig 1, middle). The neutral loss in MS/MS automatically triggered an MS³ experiment, which gives a very rich and clear peptide fragmentation (Fig 1, bottom). Interpretation of the data identified with very high confidence a β-casein phosphopeptide (FQpSEEQQTEDELQD, *m/z* 1031).

Based on the findings from method development, we decided to lower the collision energy in MS/MS to 15–25% to enhance neutral loss and avoid complete fragmentation of the peptides. To obtain amino acid sequence information, the neutral loss ion was further fragmented in MS³ mode. The collision energy for MS³ was set to 35%, which led to good fragmentation of the selected phosphopeptides (5).

To determine whether phosphopeptides could be detected at low concentrations in mixtures, a mixture of β-casein trypsin digest and six-protein digest was analyzed by the described method. The sample was injected in different

concentrations corresponding to 20, 2, and 1 fmol of β -casein. The resulting spectra were subjected to database searches using TurboSEQUEST software. The phosphopeptide was easily identified ($X_{\text{corr}} > 3.2$) even at 1 fmol, with as many as 15 out of 16 b ions and 15 out of 16 y ions assigned (Fig 2). The S/N of 1 fmol of phosphopeptide in the extracted ion chromatogram was from 5 to 10, which is above the detection limit. The phosphopeptide ion at m/z 1031 was easily detected in the averaged spectrum (Fig 3).

Table 1 illustrates the linearity of the method in the low femtomole range, indicating that there were no losses of the phosphopeptide in the system. This demonstrates the biocompatibility of Ettan MDLC: there are no metal ions in the system that can complex with phosphate groups. There was a slight deviation at 1 fmol, which was not surprising considering that the phosphopeptide concentration was close to the detection limit.

To test the method further, it was applied to a complex biological sample—mouse brain tissue—to search for phosphoproteins. We had earlier identified 2786 nonphosphorylated proteins in this sample by 2-D LC-MS/MS and TurboSEQUEST database searches. The base peak ion chromatogram and all MS^3 events are shown in Figure 4. The number of MS^3 events was typically around 50 per run. Notably, some of the neutral losses from precursors that eluted late in the LC run did not originate from a phosphopeptide. These false positives were clearly found during the manual spectra investigation but not always in the automatic database searches. In total, seven putative phosphorylation sites were assigned from the mouse database search of the MS^3 spectra (Table 2). At least two scans and four consecutive runs confirmed the sequences for the first two phosphopeptides shown in the list.

Table 1. LC retention times and area of the β -casein phosphopeptide for various amounts of β -casein/six-protein digest injected.

Run	Amount of protein injected (fmol)	Average peak area	Area %
1-4	1	39553	8.12
5-8	2	46627	9.57
9-10	20	487020	100.00




Fig 1. MS spectra from 200 fmol of β -casein digest. Averaged MS spectra from the peak at 28 min retention time (top). MS^2 neutral loss spectrum of the doubly charged phosphopeptide at m/z 1032.4 (middle). MS^3 spectrum of the MS^2 peak at m/z 982.5 (bottom).




Fig 2. TurboSEQUEST search results from 1 fmol of β -casein/six-protein digest. The phosphopeptide was found and sequenced, and the phosphorylation site, serine 3, is denoted with #.




Fig 3. Averaged spectra at 24 min retention time of β -casein phosphopeptide from a blank, 1-fmol, 2-fmol, and 20-fmol injections of β -casein tryptic digest.




Fig 4. Total ion chromatogram (top) and MS³ ion trace (bottom) of mouse brain tissue.




Fig 5. A) Full MS spectrum at 28 min retention time from a mouse brain sample. B) Resulting loss of phosphoric acid from the triply charged ion, m/z 420. C) MS³ spectrum of the m/z 367 ion.

Where a phosphopeptide was identified in several runs, an average of the Xcorr values is shown in Table 2. Dihydropyrimidinase-like protein 3 (collapsing response mediator), myristoylated alanine-rich protein kinase C, and MARCKS-like protein (6) are all known phosphoproteins involved in cell signaling in the brain. Drebrin A is a neuron-specific protein that is downregulated in the brains of patients with Alzheimer's disease or Down's syndrome (7). Examples of the MS, MS/MS, and MS³ spectra from the dihydropyrimidinase-like protein 3 phosphopeptide are shown in Figure 5. In addition to the precursor ion, a b or y ion can also show a loss of phosphoric acid. This might account for some of the unknown intense ions in the spectra.

A TurboSEQUEST search was also performed on all MS/MS spectra. The phosphopeptides in Table 2 were confirmed in the MS/MS spectra. Many nonphosphorylated peptides from the phosphoproteins were also found. MS/MS spectra from phosphopeptides are often dominated by the neutral loss ion, as mentioned earlier, which results in less sequence information than from the nonphosphorylated peptides. Even so, more than 40 possible sequences with one or more phosphate groups—either indicated by dehydrated serine or threonine, or by the loss of phosphate—were identified apart from those in Table 2. Additional spectrum interpretation is needed to confirm this data.

Table 2. TurboSEQUEST result list of the putative phosphoproteins and phosphorylation sites from the MS³ spectra of mouse brain tissue.

Protein	Sequence	Mass of MH ⁺ (Da)	Average Xcorr value
dihydropyrimidinase-like 3; unc-33-like phosphoprotein	GpSPTRPNPPVR	1159.6	2.63
T-cell receptor alpha V region t2c6	DSALYYCAL(p)SGD(p)SNNR	1787.7	2.83
myristoylated alanine rich protein kinase C subst	LSGFpSFKK	895.5	2.67
glutamyl-prolyl ⁻ tRNA	SQGSGLPSSGGAGEGQGPK	1542.7	2.52
MARCKS-like protein	LSGLpSFKR	889.5	2.45
Similar to microtubule-associated protein 2	VAIIRpTPPKSPATPK	1557.9	3.06
Drebrin A	LSpVLHR	890.5	2.10

Conclusions

Ettan MDLC in combination with the Finnigan LTQ mass spectrometer is an ideal platform to study phosphopeptides because of the biocompatibility of Ettan MDLC and the speed and MSⁿ capabilities of the Finnigan LTQ system. Phosphopeptides could be sequenced and identified with a high level of confidence at low femtomole concentrations. Brain phosphoproteins were found in complex mixtures with a high level of confidence.

To further increase the sensitivity of phosphopeptide analysis, an enrichment strategy can be implemented, such as immobilized metal affinity chromatography (IMAC). The resolution of the analysis might also be increased by adding a cation-exchange column to the first-dimension separation (8).

References

1. Hunter, T. Signaling—2000 and beyond. *Cell* **100**, 113–127 (2000).
2. Cohen, P. Protein kinases—the major drug targets of the twenty-first century? *Nat. Rev. Drug Discovery* **1**, 309–315 (2002).
3. Ficarro, S. B. *et al.* Phosphoproteome analysis by mass spectrometry and its application to *Saccharomyces cerevisiae*. *Nat. Biotechnol.* **20**, 301–305 (2002).
4. Andersson, L. and Porath, J. Isolation of phosphoproteins by immobilized metal (Fe³⁺) affinity chromatography. *Anal. Biochem.* **154**, 250–254 (1986).
5. Ruse, C. *et al.* Integrated analysis of the human cardiac transcriptome, proteome, and phosphoproteome. *Proteomics* **4**, 1505–1516 (2004).
6. McNamara, R. K. and Lenox, R. H. Distribution of the protein kinase C substrates MARCKS and MRP in the postnatal developing rat brain. *J. Comparative Neurology* **397**, 337–356 (1998).
7. Kobayashi, R. *et al.* Antisense knockdown of drebrin A, a dendritic spine protein, causes stronger preference, impaired pre-pulse inhibition, and an increased sensitivity to psychostimulant. *Neurosci. Res.* **49**, 205–17 (2004).
8. Beausoleil, S. A. *et al.* Large-scale characterization of HeLa cell nuclear phosphoproteins. *Proc. Natl. Acad. Sci. USA* **101**, 12130–12135 (2004).

Related products

Offline Fractionation Kit

11-0007-20

Asia Pacific	Tel: +852 2811 8693	Fax: +852 2811 5251
Australasia	Tel: + 61 2 9899 0999	Fax: +61 2 9899 7511
Austria	Tel: 01/57606-1619	Fax: 01/57606-1627
Belgium	Tel: 0800 73 888	Fax: 03 272 1637
Canada	Tel: 1 800 463 5800	Fax: 1 800 567 1008
Central, East, & South East Europe	Tel: +43 1 982 3826	Fax: +43 1 985 8327
Denmark	Tel: 45 16 2400	Fax: 45 16 2424
Finland & Baltics	Tel: +358-(0)9-512 39 40	Fax: +358 (0)9 512 39 439
France	Tel: 01 69 35 67 00	Fax: 01 69 41 96 77
Germany	Tel: 0761/4903-490	Fax: 0761/4903-405
Italy	Tel: 02 27322 1	Fax: 02 27302 212
Japan	Tel: +81 3 5331 9336	Fax: +81 3 5331 9370

Latin America	Tel: +55 11 3933 7300	Fax: +55 11 3933 7304
Middle East & Africa	Tel: +30 210 9600 687	Fax: +30 210 9600 693
Netherlands	Tel: 0165 580 410	Fax: 0165 580 401
Norway	Tel: 815 65 555	Fax: 815 65 666
Portugal	Tel: 21 417 7035	Fax: 21 417 3184
Russia & other C.I.S. & N.I.S	Tel: +7 (095) 232 0250, 956 1137	Fax: +7 (095) 230 6377
South East Asia	Tel: 60 3 8024 2080	Fax: 60 3 8024 2090
Spain	Tel: 93 594 49 50	Fax: 93 594 49 55
Sweden	Tel: 018 612 1900	Fax: 018 612 1910
Switzerland	Tel: 0848 8028 12	Fax: 0848 8028 13
UK	Tel: 0800 616928	Fax: 0800 616927
USA	Tel: +1 800 526 3593	Fax: +1 877 295 8102

www.amershambiosciences.com/mdlc

GE Healthcare
Amersham Place
Little Chalfont
Buckinghamshire
HP7 9NA
UK

General Electric Company reserves the right, subject to any regulatory approval if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information. © 2005 General Electric Company—All rights reserved. GE and GE Monogram are trademarks of General Electric Company. Amersham Biosciences, Ettan, NAP, and PlusOne are trademarks of GE Healthcare Limited. Finnigan, LTQ, and TurboSEQUEST are trademarks of Thermo Electron Corporation. PepMap is a trademark of LC Packings, a Dionex Company. Suprapur is a trademark of Merck KGaA. Zorbax is a trademark of Agilent Technologies. Amersham Biosciences Corp., a General Electric Company going to market as GE Healthcare.



GE imagination at work