



# Ni Sepharose 6 Fast Flow—a new high capacity nickel-charged medium for batch purification of histidine-tagged proteins

**Intellectual Property Notice:** The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

## cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.  
© 2020 Cytiva  
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.  
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

# Ni Sepharose 6 Fast Flow—a new high capacity nickel-charged medium for batch purification of histidine-tagged proteins

M. Carlsson, A. Bergh, A. Heijbel, J. Lundqvist, U. Saplakoglu\*, K. Torstenson, and K. Öberg

GE Healthcare, Björkgatan 30, 751 84 Uppsala, Sweden

\*GE Healthcare, 800 Centennial Ave, PO Box 1327, NJ 08855-1327, USA

## Introduction

Immobilized Metal Affinity Chromatography (IMAC) is an excellent chromatography technique for purification of histidine-tagged proteins. Pure protein can be obtained from a pretreated clarified cell extract in a single step. Ni Sepharose™ 6 Fast Flow is a new medium precharged with nickel for optimized purification of histidine-tagged proteins. The medium is useful for batch/gravity-flow. Nickel leakage from the medium is negligible, which allows repeated purification of histidine-tagged proteins without recharging the medium.

## Ni Sepharose 6 Fast Flow for batch purifications

- Compatible with commonly used reducing agents, such as DTT, DTE, TCEP, and  $\beta$ -mercaptoethanol.
- Stable in a wide range of denaturants, detergents, and buffer systems.
- Protein binding capacity higher than that found with comparative media from other manufacturers.

## Compatibility with reducing and denaturing agents, detergents, additives, and buffer substances

Ni Sepharose 6 Fast Flow is stable in a wide range of reducing and denaturing agents, detergents, and additives. The medium is, for example, stable in reducing agents such as DTT at concentrations up to 5 mM (Table 1).



**Table 1.** Ni Sepharose 6 Fast Flow is stable in the following reducing agents, denaturing agents, detergents, additives, and buffer substances at the concentrations given

Reducing agents*	5 mM DTE
	5 mM DTT
	20 mM $\beta$ -mercaptoethanol
	5 mM TCEP
	10 mM reduced glutathione
Denaturing agents	8 M urea <sup>†</sup>
	6 M Gua-HCl <sup>†</sup>
Detergents	2% Triton™ X-100 (nonionic)
	2% Tween™ 20 (nonionic)
	2% NP-40 (nonionic)
	2% cholate (anionic)
	1% CHAPS (zwitterionic)
Other additives	500 mM imidazole
	20% ethanol
	50% glycerol
	100 mM Na <sub>2</sub> SO <sub>4</sub>
	1.5 M NaCl
	1 mM EDTA <sup>‡</sup>
60 mM citrate <sup>‡</sup>	
Buffer substances	50 mM sodium phosphate, pH 7.4
	100 mM Tris-HCl, pH 7.4
	100 mM Tris-acetate, pH 7.4
	100 mM HEPES, pH 7.4
	100 mM MOPS, pH 7.4
	100 mM sodium acetate, pH 4 <sup>‡</sup>

\* For best results, perform a blank run before including reducing agents in the sample/buffers.

<sup>†</sup> Tested for one week at 40 °C.

<sup>‡</sup> The strong chelator EDTA has been used successfully in some cases at 1 mM. Generally, chelating agents should be used with caution (and only in the sample, not in the buffers). Any metal-ion stripping may be counteracted by addition of a small excess of MgCl<sub>2</sub> before centrifugation/filtration of the sample. Note that stripping effects may vary with applied sample volume.

## Repeated purifications of a histidine-tagged protein on the same column with buffers and sample containing 2 mM or 5 mM DTT

Column: HisTrap™ FF 1-ml

Sample: (Histidine)<sub>6</sub>-tagged maltose binding protein in *E. coli* extract (conc. ~1 mg/ml, M<sub>r</sub> ~43 000)


Sample volume: 7 ml

Binding buffer: 2 or 5 mM DTT, 20 mM sodium phosphate, 25 mM imidazole, 500 mM NaCl, pH 7.4

Elution buffer: 2 or 5 mM DTT, 20 mM sodium phosphate, 500 mM NaCl, pH 7.4

Flow rate: 1 ml/min

A blank run without DTT was performed before the purifications.



## Results:

- Purity and recovery of the protein were unaffected by DTT.

## High binding capacity of Ni Sepharose 6 Fast Flow


The batch purification performance of Ni Sepharose 6 Fast Flow was compared to the performance of two Qiagen™ products—Ni-NTA Superflow™ and Ni-NTA Agarose. Batch purification of histidine-tagged green fluorescent protein (GFP-[His]<sub>6</sub>) was performed using empty disposable PD-10 columns packed with the different media.

### Materials and methods

*E. coli* BL21 cells were lysed by sonication in binding buffer (20 mM sodium phosphate, 500 mM NaCl, 20 mM imidazole, pH 7.3). After clarification, pure GFP-[His]<sub>6</sub> was added to give 2.25 mg GFP-[His]<sub>6</sub>/ml of sample.



The media were prepared to give 50% slurries in binding buffer. A sample volume of 4 ml (approx. 9 mg GFP-[His]<sub>6</sub>) was added to 0.2 ml or 1 ml of the 50% slurries (corresponding to 0.1 ml and 0.5 ml of the media respectively) and incubation at room temperature was performed on a shaker for 1 h. The media with bound GFP-[His]<sub>6</sub> were packed into empty disposable PD-10 columns and washed with four column volumes (CV) of binding buffer. Elution was performed using five and seven CVs of elution buffer (20 mM sodium phosphate, 500 mM NaCl, 500 mM imidazole, pH 7.3) for 0.5 and 0.1 ml medium respectively. The amounts of eluted GFP-[His]<sub>6</sub> per ml medium were estimated by 280 nm measurements. The purity of pooled fractions was analyzed by SDS-PAGE.

The protocol described is summarized schematically below.






## Results

- As estimated from absorbance readings, the amount of eluted GFP-(His)<sub>6</sub> and the % yield was higher for Ni Sepharose 6 Fast Flow than for the other media.



*(Due to the kinetics of batch purification, the amount of GFP-(His)<sub>6</sub> was not directly proportional to the quantity of the medium used.)*

- Purity of GFP-(His)<sub>6</sub> obtained with 0.1 ml medium was similar for the three media, as determined by SDS-PAGE.




Lane 1: LMW markers  
 Lane 2: Flow through  
 Lane 3: Wash 1  
 Lane 4: Wash 2  
 Lane 5: Elution fraction 1  
 Lane 6: Elution fraction 2  
 Lane 7: Elution fraction 3  
 Lane 8: Elution fraction 4


*(Similar purities were observed also with 0.5 ml of the media.)*

## Results

- A higher binding capacity results in reduced costs as a consequence of the lower volume requirement of Ni Sepharose 6 Fast Flow.



(Data based on packed bed breakthrough experiments.)



## Conclusions

- High binding capacity of Ni Sepharose 6 Fast Flow reduces overall costs of histidine-tagged protein purification.
- Repeated exposure of the medium to reducing agents such as DTT did not affect binding capacity or purity of the histidine-tagged protein.
- Higher yield of batch purified GFP-(His)<sub>6</sub> was obtained using Ni Sepharose 6 Fast Flow compared with Ni-NTA Superflow and Ni-NTA Agarose.

## Acknowledgement

We thank the owners of the clones used in this work: MBP-(His)<sub>6</sub> provided by Pharmacia Diagnostics, Uppsala, Sweden, GFP-(His)<sub>6</sub> provided by Dr. David Drew, Dept. of Biochemistry and Biophysics, Stockholm University.

[www.chromatography.amershambiosciences.com](http://www.chromatography.amershambiosciences.com)

GE Healthcare  
Amersham Biosciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

HisTrap and Sepharose are trademarks of GE Healthcare Ltd, a General Electric Company. GE and GE monogram are trademarks of General Electric Company. Triton is a trademark of Union Carbide Chemicals and Plastics Co. Tween is a trademark of ICI Americas Inc. Qiagen is a trademark of the Qiagen group. Superflow is a trademark of Sterogene Bioseparations Inc.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

Protocols of the comparative studies are found at [www.amershambiosciences.com/protocol-his](http://www.amershambiosciences.com/protocol-his). All experiments were performed at GE Healthcare Bio-Sciences, Protein Separations laboratories.

Purification and preparation of fusion proteins and affinity peptides comprising at least two adjacent histidine residues may require a license under US pat 5,284,933 and US pat 5,310,663, including corresponding foreign patents (assignee: Hoffman La Roche, Inc).

© 2005 General Electric Company – All rights reserved.

Amersham Biosciences AB, a General Electric company going to market as GE Healthcare.

