

## High-throughput screening for elution conditions on Capto MMC using PreDictor plates

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva<sup>™</sup> brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

#### cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.

Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners. © 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit <a href="https://contact.com/contact">cytiva.com/contact</a>

CY13156-04May20-AN

### **GE Healthcare**

# High-throughput screening for elution conditions on Capto MMC using PreDictor plates

Carina Engstrand, Tryggve Bergander and Karol Lacki GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden

#### Introduction

The increasing demands on the Biopharma industry with shorter time lines and cost constraints have challenged the use of conventional column chromatography for early process development. High-throughput screening of chromatographic conditions can substantially decrease the time and amount of sample needed during this process. The work presented here describes the use of PreDictor<sup>™</sup> plates, 96-well filter plates filled with chromatography media, for parallel screening of elution conditions on the multimodal cation exchange medium Capto<sup>™</sup> MMC using Design of Experiments (DoE). The results are compared with data generated in chromatography columns. Comparisons of time and sample required to perform the study in the two different formats are made.

#### **Experimental approach**

In this elution study, buffer ionic strength (BIS), pH, and concentration of two eluting salts were varied using DoE (MODDE<sup>TM</sup>, version 8.0) within a range of conditions listed in Table 1.

Theoretically, the batch adsorption/desorption occurring in the wells of PreDictor plates involves the same mass transfer mechanisms as in column chromatography (1). To verify that data derived from PreDictor plates give the same qualitative information as data generated by column chromatography, a corresponding study was performed in a Tricorn<sup>™</sup> 5/100 column (column volume (CV) 2 ml). 
 Table 1. The factors, associated ranges, and responses measured in this study

Factors	Range	Responses
Salt concentration	0.25, 1.00, 1.75 M	Recovery in 1st elution (for plates) or in 3 CV (for column)
Salt type	NaCl or NH <sub>4</sub> Cl	
Buffer Ionic Strength (BIS) <sup>1</sup>	0.026, 0.163, 0.300 M	
рН²	5.75, 6.25, 6.75	

<sup>1</sup> Calculated from the Henderson-Hasselbach equation without correction for pKa shifts, and at different ionic strengths. Obtained by addition of sodium hydroxide.

 $^{\rm 2}$   $\,$  The buffers used were sodium malonate (pH 5.75 and 6.25) and sodium phosphate (pH 6.75).

#### **Materials and Methods**

Prototype PreDictor Capto MMC plates filled with 6 µl chromatography medium per well were used to study elution conditions for bovine serum albumin (BSA). Equilibration buffer was 50 mM sodium acetate, 250 mM NaCl, pH 4.75. Removal of liquid from the plates was done using a vacuum manifold.

#### **Design of Experiments**

For the experimental design the factors and responses listed in Table 1 were chosen. The experimental layout in the PreDictor plate was done according to the RED-MUP design in MODDE. The design was a full factorial design in multiple levels in order to resolve main, interactive, and curvature effects. The software used for the design and evaluation was MODDE version 7.0 for column experiments and MODDE version 8.0 for plate experiments (2).



#### Prestudy: Binding step

To identify experimental settings (e.g., incubation time and protein concentration) for the elution study, a pre-study of the binding step was performed. Batch uptake curves at four different protein concentrations were monitored by adding sample at different times, followed by simultaneous removal and measurement of unbound protein. The dynamic binding capacity at 10% breakthrough (QB10%) for BSA had previously been determined in columns to be approximately 35 mg/ml chromatography medium at 1 minute residence time. A suitable load is often considered 70% of QB10%, in this case 25 mg/ml. Thus, for the experimental procedure, conditions giving a binding capacity of 25 mg/ml should be chosen.



**Fig 1.** Batch uptake curves for BSA at four different concentrations at various incubation times.

#### Elution study Chromatography in 96-well format

After equilibration of the chromatography media in the wells, 200 µl of 1 mg BSA/ml equilibration buffer was loaded onto each well in the PreDictor plate and incubated on an orbital shaker for 10 minutes. The wells were washed three times with equilibration buffer, whereafter elution was done using the different conditions described in the experimental set-up. The elution step was repeated three times. Eluted material was collected in a UV-readable microplate (Costar).

#### **Evaluation of data**

The absorbance at 280 nm was used to calculate the concentrations of the elution pools according to Lambert-Beer's law. The recovery of protein was then calculated using Equation 1. The recovery in the first elution was used as the response in MODDE.

Equation 1: m = C × V, recovery (%) = 100 ×  $m_{recovered}$ /  $m_{loaded}$ 

#### Column chromatography

Tricorn 5/100 columns with a CV of 2 ml were used for the column experiments (3). Equilibration buffer was 50 mM sodium acetate, 250 mM NaCl, pH 4.75. The sample was BSA (4 mg/ml) in equilibration buffer. The dynamic binding capacity at 10% breakthrough (QB10%), 1 minute residence time, had previously been determined using frontal analysis. For the recovery runs, 70% of QB10% was loaded to the column.

#### Results

From the batch uptake in Figure 1, it was apparent that an initial concentration of 1 mg/ml and an incubation time of 10 min resulted in a robust response with a capacity of approximately 25 mg/ml. Thus, this was used as the sample concentration and the incubation time in the elution study in PreDictor plates.

Evaluation of the recoveries in MODDE showed that the results in PreDictor plates were similar to those obtained in column chromatography (Fig 2).



**Fig 2 A.** Contour plots for the recovery of BSA in A) PreDictor plates and B) Tricorn column (data from reference 3). Recovery is plotted as a function of salt concentration and buffer ionic strength (BIS) at three different pH values for the two salt types NaCl and NH<sub>4</sub>Cl. Experimental data points are shown as black dots.

Both studies define the best conditions for high recoveries at pH 6.75 at buffer ionic strength around 0.23 M and a concentration of  $NH_4Cl$  above 0.5 M. When NaCl is the preferred eluting salt, higher salt concentrations are needed (i.e., above 1 M). The difference in the two models could be due to a significantly larger elution volume relative the amount of chromatography media in the plate compared to the column experiment.

A comparison of the experimental time and sample consumption for this experiment performed in plates and in columns is shown in Figure 3.



**Fig 3.** Estimate of experimental time and sample consumption for this screening using PreDictor plates and conventional column chromatography (Tricorn 5/100, CV 2 ml).

#### References

- 1. Application Note. Screening of loading conditions on Capto S using a new highthroughput format, PreDictor plates. 28-9258w-40 AA.
- 2. MODDE software vesion 7.0 and version 8.0, Umetrics AB, www.umetrics.com.
- 3. Application Note. Optimizing elution conditions on Capto MMC using Design of Experiments. 11-0035-48 AA.

#### Conclusions

The optimal conditions for high recoveries of BSA on Capto MMC obtained on PreDictor plates show good correlation with conditions identified using traditional packed bed chromatography. At the same time PreDictor plates allow for shorter experimental time and lower sample consumption, making them a time-saving and accurate tool for screening of chromatographic conditions during process development.

For contact information for your local office, please visit, www.gelifesciences.com/contact

www.gelifesciences.com/bioprocess

GE Healthcare Bio-Sciences AB Björkgatan 30 751 84 Uppsala Sweden GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Capto, PreDictor, and Tricorn are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

The Tricorn column and components are protected by US design patents USD500856, USD506261, USD500555, USD495060 and their equivalents in other countries.

© 2008 General Electric Company – All rights reserved. First published Mar 2008.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information

GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden GE Healthcare Bio-Sciences Corp., 800 Centennial Avenue, P.O. Box 1327, Piscataway,

NJ 08855-1327, USA GE Healthcare Europe GmbH, Munzinger Strasse 5,

D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences KK, Sanken Bldg., 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo, 169-0073 Japan

![](_page_4_Picture_13.jpeg)