



ReadyToProcess increases facility capacity and shortens change-over time

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit cytiva.com/contact

CY13170-05May20-AN

ReadyToProcess increases facility capacity and shortens change-over time

Ingemar Daniels, Product Manager, GE Healthcare, Uppsala, Sweden.

Kristoffer Hansson, Manager of Cell Culture Purification and

*Sandeep Kristiansson, Principal Scientist, Cell Culture Purification,
CMC Biologics, Copenhagen, Denmark.*

Poster presented at ASME: 2nd Annual Bioprocess Technology Seminars & Exhibition — Europe
Stockholm, Sweden, 23–25 June, 2008



Introduction

Contract manufacturers work to tight schedules, producing multi-products on several production lines. Cleaning is a critical element that impinges on valuable production time. CMC Biologics in Denmark has evaluated ReadyToProcess™ as a means of increasing uptime through eliminating some cleaning steps.

The cleaning evaluation study considered the man-hours needed to clean equipment in a production campaign, comparing conventional chromatography systems and columns with products from the ReadyToProcess platform. The evaluation is based on modelling and data from normal routines within the CMC facility. Columns and systems were of comparable dimensions.

Initial preparation and cleaning procedures using conventional chromatography solutions

When using conventional chromatography solutions, a campaign equipment rationale must be written for all equipment before production can begin. This document provides an evaluation of the predefined acceptance limits to determine whether they are adequate with respect to the total equipment surface area in the campaign, the estimated number of doses to be produced, and the daily dose size for both previous batches and for the planned batch. The risk of product carry-over is calculated, evaluated and also documented in the campaign equipment rationale.

Cleaning processes are implemented as part of the initial preparations for a campaign when equipment is set-up, system inlets and outlets installed, column tube cleaned and packed with media, and HETP testing carried out. It is also included during a campaign, between batches. Cleaning and cleaning validation follow approved protocols and provide proof that equipment is properly clean.

At the end of a campaign a cleaning verification is performed. This involves following a prescribed cleaning protocol that defines the parameters used to verify cleaning. At CMC, these parameters are:

- Rinse samples
 - TOC (Total Organic Carbon)
 - Endotoxin
 - Bioburden
- Swab samples
 - TOC
- Visual inspection




Figure 1. Overview of sampling and response times for cleaning results at CMC.

Table 1. Steps and time allocated to preparing and cleaning a chromatography system. Note that in this study specific time and costs for preparation of CIP solutions are considered negligible in relation to facility costs.

Timing	Activity	Estimated time needed	Comment	Step eliminated in a single use process
Before use of new equipment	Produce cleaning protocol	1 day		✓
Before start of a campaign	Campaign equipment rationale	0.5 day	Includes QA review	✗
Before start of a campaign	Equipment set up	0.5 day		✗
Before start of a campaign	Cleaning before use	0.5 day	Includes preparation of CIP solutions	✓
During a campaign	Cleaning in between lots	1 day		✗
After a campaign	Preparation of sampling	1–2 hours	Training in protocol, wash of sampling glasses	✓
After a campaign	Cleaning after use	0.5 day		✓
After a campaign	Rinse sampling, drain of system and swab sampling	1 day	Includes preparation of CIP solutions	✓
After a campaign	QC assays	0.5–1 day	Two persons half a day	✓
After a campaign	Cleaning report	1 hour	Includes review	✓

The procedure also includes the analytical assays, reviewing results and approving each assay. The bottleneck in the cleaning verification process is the time taken for assay results to be ready. Preliminary bioburden results can be read 3 days after sampling, but a final result is not available until after 5 days (Fig 1).

Conventional vs single-use

For a production campaign using conventional chromatography equipment, the total time needed to clean a system is 6–7 man days (Table 1) and for a column 5–6 man days (Table 2). However, the equipment is not




Figure 2. Estimated increase in the number of batches in a 14-day and 7-day production campaign schedule with single-use products.

Table 2. Steps and time allocated to preparing and cleaning a chromatography column. Note that in this study specific time and costs for preparation of CIP solutions are considered negligible in relation to facility costs.

Timing	Activity	Estimated time needed	Comment	Step eliminated in a single use process
Before use of new equipment	Produce cleaning protocol	1 day	Only for new equipment. The cleaning protocol is reused.	✓
Before start of a campaign	Campaign equipment rationale	0.5 day	Includes QA review	✗
Before start of a campaign	Equipment set up	0.5 day		✓
Before start of a campaign	Cleaning before use	0.5 day	Includes preparation of CIP solutions	✓
Before start of a campaign	Column packing and HETP test	1 day	Includes preparation of packing and HETP solutions	✓
During a campaign	Cleaning in between lots	0.5 day		✗
After a campaign	Preparation of sampling	1–2 hours		✓
After a campaign	Column unpacking, exchange of spare parts etc.	1 day	Two persons half a day	✓
After a campaign	Cleaning after use and sampling after use	1 day	Two persons half a day, includes preparation of CIP solutions	✗ (Column can be cleaned and stored until batch release)
After a campaign	QC assays	0.5–1 day	Actual depends on if assays for several instruments are run at the same time	✓
After a campaign	Cleaning report	1 hour	Includes review	✓

ready for use in another production campaign until the results from the post-campaign cleaning are approved, which is at least 5 facility days after CIP and sampling was done (Fig 1).

The tables do not take into consideration time for handling failed cleanings and deviations. The consequence of a failed cleaning is that time has to be spent on evaluation and follow up of errors. In a worst case scenario, part of the cleaning or QC work may have to be repeated. Extra time for failure is not included in production schedules and can thus affect production capacity. Within the industry, the average failure rate for cleaning is about 10%.

Employing single-use columns and systems offers a number of distinct advantages that save time within a facility. ReadyToProcess columns are pre-packed and pre-sanitized, so there is no need for column packing or cleaning procedures pre- and post-campaign. As a consequence, the risk of a cleaning failure and any subsequent rescheduling of production are avoided, and the amount of documentation and analysis in the cleaning process is reduced. By their nature, single-use products also remove the risk of product cross-contamination.

Increases capacity by 30%

Removing time-consuming steps with single-use products provides an opportunity to increase production capacity Tables 1&2. For a facility with a 49-week production schedule (3 weeks allocated to maintenance and calibration) and production campaigns of 14 days, eliminating cleaning steps can mean a reduction in change-over from 7 to 4 days, which would give an increase in the number of batches produced per year from 16 to 19. For 7-day production campaigns, the same reduction in change-over time would increase the maximum number of batches per year from 24 to 31, or by 30% (Fig 2).

Single-use products will also benefit multi-product facilities that demand a fast change-over between production campaigns. This is particularly important for facilities where the upstream capacity is larger than the downstream capacity, for example where both microbial fermentation and mammalian cell culture bioreactors deliver harvest to the same downstream production line.

Faster route to clinic

By eliminating some of the time-consuming steps from a campaign – within initial set-up, cleaning, analysis and documentation – downtime can be turned into uptime and production capacity within a facility increased. In cases where rapid facility change-over is demanded, single-use products will be especially useful. The ReadyToProcess platform will be a good choice for early clinical phases or for potent products. A standard design with fixed column sizes makes the concept valuable for many applications and for a faster route to market.

ReadyToProcess

The ReadyToProcess platform simplifies and accelerates bioprocessing, reducing upfront investment and shortening development time. The platform has features that effectively eliminate the need to clean, sterilize, or validate multiple-use systems in the manufacturing process. Products are designed to enable lean and responsive biopharmaceutical development and production with assured safety and cost-efficiency, from cell culture and fermentation to purification. The ReadyToProcess platform comprises fluid handling, filtration and chromatography solutions, as well as the WAVE Bioreactor™ and WAVE Mixer™.



CMC Biologics

CMC is a contract manufacturing organization, with manufacturing facilities in Europe (CMC Biologics A/S) and the USA (CMC ICOS Biologics Inc) offering a wide range of integrated cGMP manufacturing services using microbial fermentation and mammalian cell culture processes. CMC's headquarters are located in Copenhagen, Denmark, part of the Medicon Valley region and its US facility is based in Seattle, WA. The range of services offered for both mammalian cell culture and microbial fermentation include the following:

- Process development and scale-up of protein processes
- Production of biopharmaceuticals for: pre-clinical, clinical trial phase I, II, III, and market supply
- Consultancy in facility design, quality and regulatory issues
- New technology evaluation
- Analytical Development & Quality Control



For contact information for your local office,
please visit, www.gelifesciences.com/contact

www.gelifesciences.com/readytoprocess

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

ReadyToProcess, WAVE Bioreactor and WAVE Mixer are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners.

© 2008 General Electric Company – All rights reserved.

First published Oct. 2008.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information

GE Healthcare Bio-Sciences AB,
Björkgatan 30,
751 84 Uppsala, Sweden

GE Healthcare Bio-Sciences Corp.,
800 Centennial Avenue,
P.O. Box 1327, Piscataway,
NJ 08855-1327, USA

GE Healthcare Europe GmbH,
Munzinger Strasse 5,
D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences KK,
Sanken Bldg., 3-25-1, Hyakunincho,
Shinjuku-ku, Tokyo,
169-0073 Japan



imagination at work