



Impact of sterile media filtration on the performance of cells cultured in WAVE Bioreactors

Intellectual Property Notice: Part of GE Healthcare's Life Sciences business was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY16289-04Sep20-PT

Impact of Sterile Media Filtration on the Performance of Cells Cultured in WAVE Bioreactors

Seul A Bae¹, Jean-François Hamel¹ & Jonathan Royce²

¹Massachusetts Institute of Technology, 66 Ames St., Cambridge, MA 02142 ²GE Healthcare Bio-Sciences Corp. 14 Walkup Drive, Westborough, MA 01581

Introduction

Disposable bioreactors have gained widespread acceptance in cell culture applications because they provide a flexible resource for multi-product facilities and speed the production of bio-molecules to the market. Disposable bioreactors offer a number of advantages for researchers including reduction of preparation time, elimination of cleaning and sterilization time, and ease of use. (Fries et al.) Of the many models of disposable bioreactors, WAVE bioreactors are the most popular—with over 2000 units in operation worldwide.

The operation of WAVE bioreactors is very straightforward. However, the feeding of cell culture media to the bioreactor continues to be a manual process which is associated with a heightened risk of contamination. Usually, the use of tube-fusers or welders ensures a sterile environment inside the reactors, inoculation bottles and tubes, but misuse or lack of this equipment may result in contamination without any simple method to test for sterility.


To minimize the risk of contamination, we propose that an ULTA™ Pure HC capsule filter be attached to the disposable Cellbag™ so that media can be fed directly through the filter after the initial inoculation. If successful, the pre-installation of the filter would eliminate the possibility of contamination during media feeding and the time spent performing a separate media sterilization step. The following experiments test the comparability of mammalian cell growth and viability in ULTA Pure HC filtered media versus media sterilized by other means.

Materials

The following equipment and materials were used in these experiments.

Consumables

Standard Cellbag 10	GE Healthcare	CB0010L10-02
Custom Cellbag 10 w/ULTA	GE Healthcare	CB0010L10GE01-1
ULTA Pure HC 2-inch capsule	GE Healthcare	KMPGH9202HH



Equipment

WAVE EH20/50	GE Healthcare	28-4115-03
Biomedical Analyzer	Nova Biomedical	Bioprofile 400

Cell Culture Media & Supplements

SAFC EX-CELL	SAFC biosciences	C4567
DMEM (Meth-A)	Cellgro	50-014-PB
DMEM (Hybridoma)	Gibco	10313
FBS (Meth-A)	Gibco	103B
FBS (Hybridoma)	Atlanta Biologicals	S11550
Penicillin/streptomycin	Cellgro	30-002-CI
Sodium pyruvate	Cellgro	25-000-CI
L-glutamine	MP Biomedicals	101806
GlutaMax	Gibco	35050
Glucose	EMD	DX0145-3
Sodium bicarbonate	Mallinckrodt	7412-12

Table 1 – Filter Description

Type:	GE Healthcare ULTA Pure HC
Format:	2-inch capsule w/pleated membrane
Prefilter membrane:	0.6 µm polyethersulfone prefilter
Final membrane:	0.2 µm polyethersulfone sterilizing-grade
Surface area:	0.047 m ²

Methods

Media analysis

Three common cell culture media were prepared:

Base Media	Supplements (per liter)
Ex-cell ACF CHO medium	100mL FBS 20mL L-glutamine 10mL glucose 10mL penicillin/streptomycin 10mL sodium pyruvate
DMEM medium for hybridoma WN cells	100 mL FBS 10 mL GlutaMAX 10 mL penicillin/streptomycin
DMEM medium for Meth-A cells	100 mL FBS 10 mL penicillin/streptomycin 10 mL sodium pyruvate 0.87 g L-glutamine 2.0 g sodium bicarbonate 2.5 g glucose


5L of each medium was warmed to 37°C and then filtered using an ULTA Pure HC filter. Filtered medium was collected to one 5L bottle and four 2mL samples were taken before filtration and after filtering every 0.5L of medium. Metabolite concentrations in the samples were analyzed using a NOVA Biomedical analyzer (Bioprofile 400).

Cell growth performance

	WAVE 1A	WAVE 1B	WAVE 2A	WAVE 2B
Objective	Compare mammalian cell growth in filtered and unfiltered media.			
Cell Type	Meth-A			
Run Mode	Batch		Fed-Batch	
Inoculation volume	1.25L cells and media were inoculated, 3.75L of sterile medium was added for final volume of 5L.	1.25L cells and media were inoculated, 3.75L of unsterile medium was added through ULTA Pure HC filter.	1.25L cells and medium were inoculated sterilely.	1.25L cells and medium were inoculated sterilely.
Media addition	No extra medium was added.		1.25L of medium was added sterilely when VCD reached 4e5 cells/mL. 2.5L of medium was added the next day for a final working volume of 5L.	1.25L of unsterile medium was added through the ULTA Pure HC filter when VCD reached 4e5 cells/mL. 2.5L of unsterile medium was added the next day for a final working volume of 5L.
Temperature	37°C			
Rock angle/speed	Initial rock angle was set at 6° and rock speed at 17 rpm. Rock speed was later increased to 21rpm to increase dissolved oxygen level.		Initial rock angle was set at 5° and speed at 14rpm. Rock angle and speed were adjusted accordingly after each media addition.	
Sampling	Samples of 3-4mL were taken from each bag everyday for CEDEX and NOVA analysis			


Results

Media Analysis



Media composition was not affected by filtration. The small changes in pH and osmolality observed are likely due to the media's carbonate buffering systems and the fact that the filtration was not performed in a carbon dioxide enriched environment. Note that this test result represents a worst-case scenario since the filter loading (~100 L/m²) is far lower than typical filter capacities for cell culture media (>1000 L/m²).

Cell Growth Performance



No statistically significant differences in cell growth were observed. Maximum achieved viable cell density was within ±11% for both experiments. The eventual decline in viability observed in both experiments represents the natural life cycle of cells grown in batch or fed-batch reactors and is not attributable to the use of ULTA filters or Cellbag bioreactor chambers.

Conclusions

Analysis of cell culture media composition before and after filtration indicates only minor changes in the media which are not associated with changes in cell growth performance.

The data presented here indicate that the integration of ULTA Pure HC capsule filters does not affect the growth of mammalian cells cultured in WAVE Bioreactors.

The use of ULTA Pure HC capsule filters improves the usability of the Cellbag by eliminating aseptic handling steps which are not only time consuming but also are associated with a heightened risk of contamination.

Future Work

As a next step, we plan to extend this work to a cell line which produces an extracellularly expressed protein. By doing so, we can evaluate not only cell density but also protein production.

Literature Cited

Fries, Serena et al, "Evaluation of Disposable Bioreactors – Rapid Production of Recombinant Proteins By Several Animal Cell Lines", *BioProcess International*, October 2005.

Acknowledgments

We thank Carrie Miller (GE Healthcare) for her assistance in designing the custom Cellbag. We would also like to thank Laura D'Aoust, David Whitehead, and Sunitha Kilaru (MIT) who helped with the cell culture work and reactor setups.

For Further Information

Please contact jonathan.royce@ge.com or sabae@mit.edu

© 2009 General Electric Company – All rights reserved. All goods and services are sold subjects to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information. GE, imagination at work, and GE monogram are trademarks of General Electric Company. WAVE Bioreactor, Cellbag and ULTA are trademarks of GE Healthcare companies. All third party trademarks are property of their respective owners. First published 2008.