

## Gradient Kit P-50



**Fig 1.** Gradient Kit P-50 mounted on the right-hand side of HiLoad Pump P-50.

The Gradient Kit P-50 is designed for use with HiLoad™ Pump P-50 for gradient formation at high flow rates. The kit includes the Solenoid Valve PSV-50, bracket for mounting the valve on Pump P-50, inlet tubing and tubing for connecting the valve to the pump. The valve operates as a switch valve controlled from Gradient Programmer GP-250 Plus or Liquid Chromatography Controller LCC-501 Plus. A second valve can be mounted to allow programmed control of sample application. For optimal results with gradient elution, a separate mixer (size depending on the flow rate used) should be included in the system.

These Instructions refer to the mounting, connection and use of the Solenoid Valve PSV-50. Refer to the respective User Manuals for further details of Pump P-50, Gradient Programmers, Liquid Chromatography Controllers and Mixers.

### Unpacking

Check the contents of the transport package against the enclosed Packing List. Notify your Pharmacia Biotech representative immediately if any parts are missing or damaged.


### Installation

The mounting bracket for holding Solenoid Valve PSV-50 is fitted on the front of Pump P-50, on either the right or left of the pump head. Figure 1 shows the valve mounted on the right.

#### Fitting the mounting bracket

Push the black rubber holder into place in the mounting bracket. If you are only going to use one valve, fit the holder on the side of the bracket which will be nearest the





**Fig 2.** Tubing connections for valve(s). Note the different orientation of the valves according to whether they are mounted on the right or left of the pump head. Placing and connections for the optional sample application valve are shaded.

Note that the valve ports marked IN are used in this arrangement as outlet ports from the sample application valve to the switch valve and from the switch valve to the pump.

pump head (see Fig 2). If you intend to use a second valve for sample application, fit both holders into the mounting bracket.

The orientation of the holders in the mounting bracket is determined by whether you intend to fit the valve on the right or left of the pump head. Make sure the notches in the holders can accommodate the inlet and outlet tubing arrangement as shown in Figure 2.

Remove the screw from underneath the front edge of the pump on the side where you intend to fit the mounting bracket. Fit the mounting bracket using the longer screw provided with the Gradient Kit.

### Tubing connections

Fit inlet and outlet tubing to the valve(s) before placing the valve(s) in the holders. Some of the tubing connections are awkward to get at when the valves are in place in the holders.

Fit inlet and outlet tubing to the valve(s) according to Figure 2. Tighten all tubing connections by hand or gently with a key.

Use the short (88 mm) tubing provided with the Gradient Kit to connect the switch valve to the pump. A second short tubing is required for connecting between the two valves when a sample application valve is used.

Always use I.D. 1.9 mm tubing on the inlet side of the pumphead and I.D. 1.2 mm tubing on the outlet side. Keep all tubing lengths as short as possible to minimize gradient deterioration.

Place the valve(s) in the holder(s) and connect the pump tubing to the inlet port on the pump head.

### Electrical connections


Connect the REMOTE contact on Pump P-50 to the PUMP A contact on GP-250 Plus or LCC-501 Plus using the Communication Cable (Code No. 19-6005-02). Make sure the PUMP P contact on the controller is disconnected.

Fit the 3-pole connector to the valve cable, connecting the cable leads to the outermost contacts on the connector. Connect the switch valve to the SWITCH VALVE contact on the controller. Connect the sample application valve, if used, to one of the SOLENOID VALVE contacts on the controller using a 2-pole connector.

The default wash program on the controller causes the switch valve to switch between equal volumes of buffer A and B with a cycle time of 3 s for about 80 s at 30 ml/min.

## Gradient Mixers

For gradient work, a Mixer should be fitted between the outlet from the pump and the chromatography column. Choose the Mixer size according to the flow rate in the system, and set the mixing volume calibration value in the controller according to the Mixer size.



**Fig 3.** Electrical connections for Pump P-50 and switch valve PSV-50 used in gradient work. Connections for the optional sample application valve are shaded.

**Table 1.** Recommended Mixer sizes.

<b>Flow rate</b>	<b>Mixer</b>	<b>Mixing volume</b>
0–5 ml/min	0.6 ml	60 µl
3–50 ml/min	6.0 ml	600 µl

## Ordering information

<b>Designation</b>	<b>Code No.</b>
Gradient Kit P-50	19-1992-50
<b>Replacement parts</b>	
Solenoid Valve PSV-50	19-1994-01
Holder PSV-50	19-1992-55
Connector, 3-pole female (for switch valve)	18-0012-52
Connector, 2-pole female (for sample valve)	19-8629-01
Remote communication cable 1.5 m (controller–pump)	19-6005-02
Inlet tubing I.D. 1.9 mm, 1.5 m	19-1992-08
<b>Accessories</b>	
Mixer 0.6 ml, 24 V	19-6700-01
Mixer 6 ml, 110 V AC	18-3500-01
Mixer 6 ml, 220 V AC	18-3501-01

www.gehealthcare.com

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

HiLoad is a trademark of GE Healthcare companies. GE, imagination at work and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

© 2006 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
SE-751 84 Uppsala  
Sweden

GE Healthcare Europe GmbH  
Munzinger Strasse 5  
D-79111 Freiburg  
Germany

GE Healthcare UK Ltd  
Amersham Place  
Little Chalfont  
Buckinghamshire, HP7 9NA  
UK

GE Healthcare Bio-Sciences Corp  
800 Centennial Avenue  
P.O. Box 1327  
Piscataway, NJ 08855-1327  
USA

GE Healthcare Bio-Sciences KK  
Sanken Bldg, 3-25-1, Hyakunincho  
Shinjuku-ku, Tokyo 169-0073  
Japan

Asia Pacific Tel: +85 65 62751830 Fax: +85 65 62751829 • Australasia Tel: +61 2 8820 8299 Fax: +61 2 8820 8200 • Austria Tel: 01 /57606 1613 Fax: 01 /57606 1614 • Belgium Tel: 0800 73 890 Fax: 02 416 8206 • Canada Tel: 1 800 463 5800 Fax: 1 800 567 1008 • Central, East, & South East Europe Tel: +43 1 972 720 Fax: +43 1 972 722 750 • Denmark Tel: +45 70 25 24 50 Fax: +45 45 16 2424 • Eire Tel: 1 800 709992 Fax: +44 1494 542010 • Finland & Baltics Tel: +358 9 512 3940 Fax: +358 9 512 39439 • France Tel: 01 69 35 67 00 Fax: 01 69 41 98 77 • Germany Tel: 0800 9080 711 Fax: 0800 9080 712 • Greater China Tel: +852 2100 6300 Fax: +852 2100 6338 • Italy Tel: 02 26001 320 Fax: 02 26001 399 • Japan Tel: 81 3 5331 9336 Fax: 81 3 5331 9370 • Korea Tel: 82 2 6201 3700 Fax: 82 2 6201 3803 • Latin America Tel: +55 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 210 96 00 687 Fax: +30 210 96 00 693 • Netherlands Tel: 0800-82 82 82 1 Fax: 0800-82 82 82 4 • Norway Tel: +47 815 65 777 Fax: +47 815 65 666 • Portugal Tel: 21 417 7035 Fax: 21 417 3184 • Russia & other C.I.S. & N.I.S Tel: +7 495 956 5177 Fax: +7 495 956 5176 • Spain Tel: 902 11 72 65 Fax: 935 94 49 65 • Sweden Tel: 018 612 1900 Fax: 018 612 1910 • Switzerland Tel: 0848 8028 10 Fax: 0848 8028 11 • UK Tel: 0800 515 313 Fax: 0800 616 927 • USA Tel: +1 800 526 3593 Fax: +1 877 295 8102



imagination at work