

ÄKTExpress

Cue Cards




Contents

System overview	2
UNICORN for ÄKTExpress, overview	4
ÄKTExpress workflow, overview	5
1 Creating a method plan for purification	6
2 Preparing solutions, system and columns	9
3 Starting a run	12
4 Monitoring a run	14
5 Column and system procedures	15
6 Evaluating the results	18
Getting help	19
Appendix A	20
Appendix B Flow charts	24

System overview

Column block - for purification runs


- 1 - AC or IEX column for sample 1
- 2 - AX or IEX column for sample 2
- 3 - AC or IEX column for sample 3
- 4 - AC or IEX column for sample 4 or
- AC or IEX column for 2nd AC/IEX step or
- DS column for 1st DS step out of 2 DS-steps
- 5 - Desalting column for first or last DS step or
- Gel filtration column



See also "A2 Column positions" on page 20

Column valve

Column connection




Left holder

- Gel filtration column
- HiPrep™ Desalting column

On bench


- Pump rinse solution

Control panel




Sled for microplate

Deep well microplate, 24 or 96 square shaped wells




Loop valve

Superloop™ connection



Right holder


- Superloop



Tubing overview

Detailed information on inlets and outlets for a specific run need to be checked in the "2.1 Print out a summary" on page 9. The below image is an example, showing inlets and outlets for a 4-step protocol with 2 samples.

Tubing overview, example 4-step protocol



Detailed flow charts

See "Appendix B Flow charts" on page 24



Safety


Read the ÄKTAexpress Operating Instructions before using the system.

UNICORN for ÄKTExpress, overview




Method Editor wizard

The Method Editor wizard is used to create Method Plans for Purify or Prepare and Maintain. Input information: protocol type, columns, system and column procedures.



Default parameters

View or modify default parameters





System Control wizard

The System Control wizard is used to start runs based on Method Plans. Input information: Systems to use, number of samples per system and sample information.



Prints Summary check list.


System	system 1	
Total run	7.0 hour (total)	
As sensor		
Inlet position	Solution/Buffer	Volume (ml)
A1	AC/EX Binding Buffer	500 ml
A2	Stop or DIP AC/EX, Step 1	0 ml
A3	AC/EX Wash Buffer before Elution	0 ml
A4	D5/RT Buffer	1940 ml
A5	Water	200 ml
A6	0.5M NaOH	0 ml
A7	Tag Clearance Buffer	54 ml
A8	20% EtOH	200 ml
B1	AC/EX Elution Buffer	53 ml
B2	Recharge	0 ml
S1		30 ml
S2		30 ml
S3		30 ml
S4		30 ml
Column position	Type of column	Column
1	Affinity	HiTrap_HP_1_wf0
2	Affinity	HiTrap_HP_1_wf0
3	Affinity	HiTrap_HP_1_wf0



ÄKTexpress workflow, overview

"1 Creating a method plan for purification" on page 6

- Select Prepare and Maintain or Purify
- Select Columns
- Include system and column procedures (page 7 and page 15)
 - Preparation
 - During Purification run
 - Post run




"2 Preparing solutions, system and columns" on page 9

- Print out a summary
- Prepare buffers and inlets
- Connect columns and Superloop

"3 Starting a run" on page 12

- Run setup
 - Select method plan
 - Select number of systems and samples
 - Add sample IDs and sample information
 - Prepare microplate
 - Prepare samples
- Start the run
- Be prepared for any manual actions (Load samples into loops, protease or sample into Superloop, move sample inlets)




System system 1		
Total run	7.8 hour (max)	
Alt sensor		
Inlet position	Solution/Buffer	Volume (ml)
A1	AC/REX Binding Buffer	501
A2	Strip or CP AC/REX Step 1	0 ml
A3	AC/REX Wash Buffer before Elution	0 ml
A4	DS/RF Buffer	1940 ml
A5	Water	268 ml
A6	0.5M NaOH	0 ml
A7	Tag Cleavage Buffer	54 ml
A8	20% Ethanol	268 ml
B1	AC/REX Elution Buffer	53 ml
B2	Recharge	0 ml
S1		30 ml
S2		30 ml
S3		30 ml
S4		30 ml
Column position	Type of column	Column
1	Alfrity	HiTrap_HPL_1.ml ()
2	Alfrity	HiTrap_HPL_1.ml ()
3	Alfrity	HiTrap_HPL_1.ml ()

"4 Monitoring a run" on page 14

"5 Column and system procedures" on page 15

- Column procedures (page 15)
 - Equilibrate/re-equilibrate
 - Affinity blank run
 - Ion exchange blank run
 - Metal Ion Stripping (chelating columns)
 - Metal Ion charging/re-charging
 - CIP columns
 - Remove/fill columns with ethanol
- System procedures (page 16)
 - Manual filling of the inlet tubing
 - Purge pump with methanol
 - Clean sample inlet tubing
 - Clean pump
 - Wash frac tubing
 - Rinse all outlet tubing
 - Remove ethanol from the system
 - Fill system with ethanol
 - Clean system

"6 Evaluating the results" on page 18



1 Creating a method plan for purification

1.1 Overview

Use the Method Editor wizard to create a new method plan or to change an existing method plan.




- Click the Method Wizard icon in the Method Editor module.
- Select New or open an existing method plan to edit.
- Select Purify and click Next to proceed.

In addition to the basic method steps, it is possible to include most system and/or column procedures in the purification method, (e.g., system and column CIP).

System and column procedures can also be run as separate prepare and maintain methods.

For more guidance on protocol selection, see the User Manual.

For guidance on run times, see "A3 Typical run times" on page 20.



1.2 Protocols supported

For purification of	use either of these protocols
tagged proteins ¹	AC AC-DS/GF AC-DS ² -IEX AC-DS ² -IEX-DS/GF
double-tagged proteins ¹	AC-DS ² -AC AC-DS ² -AC-DS/GF
tagged proteins in need of desalting or ion exchange first ¹	DS-AC DS-AC-DS/GF IEX-DS ² -AC IEX-DS ² -AC-DS/GF
proteins with removed tag	AC ³ /DS/IEX/GF DS-IEX IEX-DS/GF DS-IEX-DS/GF IEX-DS ² -IEX-DS/GF

¹ Automated on-column tag cleavage can be performed in protocols starting with AC or IEX.

² Optional desalting step.

³ Collect the flowthrough fraction as the target protein will not bind to the AC-column.

1.3 Automated on-column tag cleavage

Automated on-column tag cleavage can be performed in protocols starting with AC or IEX.

- In the Purify method plan, check the On-column tag cleavage box on the Purify 1 page to include the cleavage in the method plan.
- Choose Incubation time and select Fill Columns with Cleavage buffer if preferred. To condition the column before the tag cleavage, it can be filled with cleavage buffer before protease is injected.

1.4 Column selection

For guidance on column selection and position, see “A1 Columns supported (examples)” on page 20, “A2 Column positions” on page 20 and the User Manual.

Affinity columns

1 or 5 ml column

- Depends on required yield.
- If many purification steps are used, start with a good amount of protein and probably with a 5 ml column.
- For large sample volumes, use a 5 ml column to save time.
- To increase protein purity, overload the column.

Ion exchange columns

Required resolution

The smaller bead size the higher resolution:

- HiTrap™ Q HP and SP HP: 34 µm
- RESOURCE™ Q and S: 15 µm
- Mono Q™ and Mono S™ 5/50 GL: 10 µm

The longer column the higher resolution.

Anion exchange, e.g., RESOURCE Q

- The pH of the buffer should be at least 1 pH unit above the pI of the protein(s).

Cation exchange, e.g., RESOURCE S

- The pH of the buffer should be at least 1 pH unit below the pI of the protein(s).

Desalting columns

HiTrap or HiPrep

- Depends on loading volumes:
 - Up to 3 ml: Use 2 × HiTrap Desalting.
 - Up to 13 ml¹: Use HiPrep Desalting.

Gel filtration columns

Superdex™ 75 or 200

- Separation range:
 - 3–70 kDa: Use Superdex 75 prep grade.
 - 10–600 kDa: Use Superdex 200 prep grade.

Sephacryl™ S-100, S-200 or S-300

- Separation range:
 - 1–100 kDa: Use Sephacryl S-100.
 - 5–250 kDa: Use Sephacryl S-200.
 - 10–1500 kDa: Use Sephacryl S-300.

HiLoad™ and HiPrep 16/60² or 26/60³

- Loading volume:
 - up to 5 ml: Use HiLoad or HiPrep 16/60.
 - up to 13 ml³: Use HiLoad or HiPrep 26/60.

¹ If you are expecting peaks larger than 7.5 ml from the affinity step, it is recommended to use a double loop, see “Appendix B Flow charts” on page 24.

² An optional column holder is available.

³ An optional column holder is available.

Tip

Grouping samples

Considerations when running several separation systems in parallel:

To maximize the success rate, run similar proteins on the same system, for example:

- To maximize GF performance, run small proteins <70 kD on one system and large proteins >70kD on another system.
- To maximize IEX performance, run proteins with low pI on one system and proteins with high pI on another system.

1.5 Buffer selection

For guidance on buffer selection, see “A4 Buffers” on page 22

1.6 Procedures possible to include

It is possible to include many system and column procedures in the purification run, see the list below.

Further explanations can be found in the Method Wizard help texts and in the User Manual.

Preparations:

- Fill sample inlet tubing with buffer
- Remove ethanol from system/columns
- Column equilibration
- Column blank run (affinity and ion exchange columns)
- Guided loading of Superloop

During Purification run:

- Wash frac tubing between samples
- Clean sample inlet tubing after sample loading
- Clean pump between samples


Post run (prepare for next run or for storage):

- Fill system/columns with ethanol
- CIP system, then fill with ethanol/buffer
- CIP columns, then fill with ethanol/re-equilibrate
- Re-equilibrate affinity or ion exchange columns
- Strip and recharge AC chelating columns
- Clean sample inlet tubing

1.7 Advanced zone

General

In the Advanced zone, all main parameter values can be viewed and modified if required. Examples of parameters that may be modified:



- Volumes
- Flow rates
- Gradient and slopes
- Watch parameters for peak detection

Note: Do not change default parameter values in a method plan **unless the result is clearly understood**.



Flow rates and pressure limits

No warning will be given if the flow rate or pressure limit is set higher than the values recommended for the columns used.

Getting help

To get help about parameter setting, click **Help** in the wizard boxes.



Tip

Use the Note field in the **Save As** page to specify changes made to the default settings.

Note
2005-12-19
Method

Specify path to import file

If an import file is going to be used for specification of sample information, the file path can be set in **Import File Location** in the **Miscellaneous Settings** page.

Tip

If the import files are saved in the default folder **MethodWizardImport**, there is no need to specify the path. For further information, see User Manual.

More information about Advanced zone

See ÅKTExpress User Manual.

2 Preparing solutions, system and columns

2.1 Print out a summary

A *Summary* can be printed from System Control and it is useful during the preparation of the system, columns and all solutions. The list includes:

- Buffers/solutions, volumes and inlets
- Sample ID, sample volume loaded and inlet position
- Usage of outlet and outlet position
- Type of column, column and column positions
- Loop positions
- Microplate type




a Select method plan

- In UNICORN™ **System Control**, click **Instant Run**.

Instant Run

- Select the required **Method plan** from the list.



b Select systems for the run

- Select on which **System(s)** the method plan will be run and the **Number of Samples** on each system.

c Print Summary

- Print out the Summary by clicking **Print**.

Note: It is also possible to copy the Summary information and paste into, for example, Microsoft® Excel®.

2.2 Prepare buffers, solutions and inlets

- Prepare buffers and solutions required for the run (see “A4 Buffers” on page 22).
- Immerse all inlet tubing in the appropriate liquid containers as described in the summary check list.

Inlets

A1-A4	See the Summary page for the specific run
A5	Water
A6	0.5 M NaOH
A7	See the Summary for the specific run
A8	20% EtOH
B1-B2	See the Summary for the specific run
S1-S4	See the Summary for the specific run

- Fill the pump and inlet tubing with liquid:

If...	then...
<i>small amounts of air need to be removed</i> or <i>buffer needs to be changed</i> <i>inlet tubing is empty</i>	perform Fill buffer inlet tubing in a Prepare and Maintain method: Standard System and Column Procedures
<i>the system has been unused for a week or longer</i> or <i>the pump has been run dry</i>	perform Purge pump with methanol in a Prepare and, Maintain method , see page 16.

- Immerse unused tubing in (e.g., 20% ethanol) to avoid air entering the tubing.
- If **Fill Sample inlet Tubings with Buffer** is not included in the method plan to use, make sure the sample inlet tubing is filled with affinity binding buffer or ion exchange buffer, see “5.3 System procedures Manual filling of the inlet tubing” on page 16.

2.3 Connect columns and superloop





For column positions, see the Summary and “A2 Column positions” on page 20.

a Remove air before connecting columns

Air remaining in the system may be removed by purging the pump and system manually by selecting **Pump Wash** and **System Wash**.

- Immerse A1 tubing in the buffer to be used.
- Select **System Control:Manual:PumpWash** or **SystemWash**.

b Connecting tubing to columns



① New models of HiTrap columns do not require the unions.

Tip


Replace the union's internal capillary with a 1 mm peek tubing (brown) to reduce backpressure, if required.

¹⁾ The standard system tubing (190 mm) should be replaced by the 350 mm capillary tubing (available in the accessory kit) when using a HiPrep 26/10 desalting column or gel filtration column in position 4 or 5.

c Column Attachment drop-to-drop

To avoid air bubbles, use the drop-to-drop procedure when connecting columns:

- In **Prepare and Maintain**, select **Column Attachment** to create a method plan which starts the pump with a low flow rate in the first column position.
- Click **Continue** to start filling Column position 1.
- Fix the tubing to the column drop-to-drop.
- Click **Next Breakpoint** on system to start filling the next column position.



Note: It is also possible to Attach Columns manually by starting a low flow and selecting **System Control:Manual:Flowpath:ColumnPosition**.

d Column preparations

- If preparations of the columns are not included in the method, prepare them according to the "5.2 Column procedures" on page 15.

Note: Uncharged IMAC columns must be loaded with metal ions in a Prepare and Maintain run before use to be functional.

e Superloop connection

Pre-filling the Superloop with buffer

To avoid pressing air into the system, pre-fill the Superloop with buffer manually before fitting it to the system.



- Remove the upper end pieces.
- Fill the upper chamber with buffer.
- Replace the end pieces and ensure there is no air trapped.
- Turn the Superloop upside down and fill the lower chamber the same way.
- Reinsert the inner end piece and turn it until the slot lines up with the countersink in the glass tube. Replace the other end pieces and ensure there is no air trapped.
- For further information, see Superloop Instruction.

Connecting the Superloop

- Fit the buffer-filled Superloop in the right hand holder.

Note: Make sure the Superloop is not hindering the ejection of the microplate.

- Connect the Superloop to the Loop Valve:
 - Top to **OUT 5**
 - Bottom to **IN 5**




Filling the Superloop

The Superloop is filled with protease/sample with assistance from a **Purify** method plan, if **Guided loading of Superloop** has been selected. To avoid filling before every run, it is also possible to include the filling in a **Prepare and Maintain** method (**Standard System**, and **Column Procedures**).

Guided loading of Superloop:

During **Guided loading of Superloop**, dialog boxes will give guidance to filling the Superloop at the beginning of the run. The filling is performed via a syringe connected to the injection valve.

- Follow the dialog boxes and press **Continue** between each step.
 - 1 Fill syringe with buffer and insert syringe in the manual inject port in injection valve.
 - 2 Inject buffer.
 - 3 Fill syringe with protease solution, or sample, and insert syringe in the manual inject port in injection valve.
 - 4 Inject protease solution, or sample.



Manual filling off-line:

The protease can be injected directly into the Superloop before starting the run.

- Connect a syringe to the lower port of the Superloop and fill it. Avoid introducing air bubbles.

f Filling capillary loop

The capillary loops are filled in the same way as Superloop.

g Calculating protease amount

- Calculate the volume and concentration of protease needed. Use, for example, the Microsoft Excel file available on the ÄKTExpress strategy CD.

What is the column volume (CV) of your HiTrap affinity column?	HiTrap 1ml
How many samples on this system will you be running?	1 Sample
How much protein will be applied to each HiTrap affinity column?	0 mg
How many CV of protease will be applied? (0.7 is default)	0.7 CV
How much protease do you need to cleave 1 mg of protein?	0 units or mg

h Optional loop volume


It is possible to increase the loop volume by prolonging each loop with an extra loop tubing. Each total loop volume will then be 20 ml. Following parameters have to be changed in Advanced Zone:

- Max Volume in Each Loop: 15 ml
- Loop Wash Volume: 40 ml
- Flush Volume Empty Loops: 20 ml
- Loop Wash Volume Empty Loops: 40 ml

3 Starting a run

3.1 Select method plan

- In UNICORN **System Control**, click **Instant Run**.
- Select the required **Method plan** from the list.




3.2 Select systems for the run

- Select on which **System(s)** the method plan will be run and the **Number of Samples** on each system.



3.3 Specify samples

- Enter identification names for the samples, either via the keyboard or using a bar code reader.



- For each sample, enter:
 - isoelectric point, **pI**,
 - extinction coefficient for the protein, **Ext Coeff**,
 - molecular weight of the protein, **MW**.

The data is automatically imported from an import file, if it has been prepared and placed in the specified folder and the location must be specified in **Import File Location** in Advanced zone settings pages, see p8. (See User Manual for further information.)

- Enter optional text, for example culture batch number.

3.4 Edit result file location and names

- If required, edit the folder path and file names of the result files to be created.

3.5 Print Summary

- Print out the Summary by clicking **Print**.


3.6 Preparations completed?

Ensure that the preparations according to “2 Preparing solutions, system and columns” on page 9 has been performed

3.7 Check the flow path

Check against the summary:


- Check that the tubing are fully immersed in the containers and that the tubing are filled with liquid.
- Immerse the waste and outlet tubing in appropriate containers.



- Make sure:
 - there is enough buffer available
 - the correct inlet is placed in each buffer
 - the outlets are placed in correct bottles
 - the columns are placed in correct positions.

3.8 Prepare the microplate

- Place an empty deep well microplate on the sled and check that the labelling **A** matches the labelling on the system.



3.9 Prepare the samples

- The samples should have been prepared and clarified using centrifugation and/or filtration through a 0.45 µm filter. If using crude columns, clarification is not needed.
- Place the sample tubes in the tube or (optional) flask holder or inject the sample into the Superloop.
- If **Fill sample inlets** is included in the method plan, the inlets should be placed in affinity binding buffer. During the run, the system will pause after initial buffer filling, and a message will appear requesting each sample inlet tubing to be gently moved to its sample.
- If **Fill sample inlets** is not included, gently move the already filled sample inlet tubing to each sample before starting the run.
- Make sure that no air enters the tubing. Place the tubing close to the bottom of the liquid container but

not too tight against the bottom.

- Secure the tubing with the tubing holder.

3.10 Final check

Run

- Perform a final check that tubing, columns and solutions are placed according to the *Summary*.

3.11 Start the run

- Click **Run** to start the run on the selected systems.


Be prepared for manual interactions (fill loops, guided loading of Superloop, fill sample inlets etc.).

4 Monitoring a run

Monitoring a run is performed in System Control.


View single or all systems

The run data of all systems can be viewed simultaneously by clicking **Expand >>**. Return to single system view by clicking **<< Collapse**.




Run data

Right-click the view pane and select **Properties** to select run data to be displayed.



Panes

Click Panes to select panes to display.




Zoom

Use the mouse drag-and-drop box to zoom in and out of the curve. Right-click to reset zoom.


Curves

Right-click to select which curves to be displayed.



Log book

Right-click to select log book properties to be displayed.



Run status

Each system's status indicator shows the run status:

Status.....	Color	Example of events
SYS_1	End.....	White
SYS_1	Run or Manual..	Green
SYS_1	Hold	Yellow Sample loading
SYS_1	Pause.....	Red Pump-wash
SYS_1	Error	Yellow, flashing Air detected (e.g., running out of buffer)

Tip


The chapter Troubleshooting in the User Manual can be useful if something unexpected happens.

5 Column and system procedures

5.1 Using a Prepare and Maintain method

Note: All preparations, except **Column Attachment, Fill Inlets Manually Using Syringe, Fill Buffer Inlets, Purge Pump with Methanol** and **Rinse All Outlets**, can also be included in a **Purify** method.

Create the method plan

- Click the **Method Wizard** icon in the **Method Editor** module. 
- Select **New** or open an existing method plan to edit.
- Select **Prepare and Maintain:Standard System and Column Procedures**.
- In **System Procedures**, select as appropriate.
- In the **Column Procedures** pages, select as appropriate. Each column type has its own page.
- On the last page, click **Finish** (or to enter the Advanced zone, click **Next**).
- In the **Save As** page, type the name of the method plan and click **OK**.

Perform the run

- Start the run using Instant Run. 

5.2 Column procedures

Equilibrate/re-equilibrate

Equilibration can be included in a Purification run. The procedure equilibrates the columns with buffer before the samples are loaded onto the columns. All four column types can be equilibrated.

An Equilibration method can also be created in **Prepare and Maintain:Standard System and Column Procedures**.

Customized equilibration can be performed in **Prepare and Maintain:Customized Column Procedures**. The procedure equilibrates up to 5 columns, of one column type, with up to 9 different solutions.

Affinity blank run



Affinity blank run can be included in a Purification run. If **Remove ethanol** is selected, EtOH will be washed out with water before equilibration.

It is recommended to include a blank run:

- Prior to a first time usage of an affinity column, or
- After metal ion charging.

Ion exchange blank run



Ion exchange blank run can be included in a Purification run. The procedure provides ion exchange columns with exchangeable counterions.

It is recommended to include a blank run:

- After long-term storage, or
- Prior to a first time usage of an ion exchange column.



Metal ion stripping (chelating columns)

To remove metal ions before regenerating the HiTrap chelating, HiTrap IMAC and HisTrap™ affinity columns, Metal Ion Stripping can be included in **Purify** or in a **Prepare and Maintain:Standard System and Column Procedures** method plan.

Strip of HisTrap, HiTrap IMAC and HiTrap Chelating (1 ml/5 ml)

Procedure	Step 1		Flow [ml/min]	
	Solution	Volume	RT	CR
Strip	His affinity binding buffer with 50mM EDTA	5 CV	1/5	0.8/4
	Step 2			
Strip	Deionized water	5–10 CV	1/5	0.8/4

Metal ion charging/re-charging

Metal ion charging/recharging of affinity columns can be performed in a Prepare and Maintain:Standard System and Column Procedures run or included in a Purify run as a post run (in combination with strip). The procedure is used to charge new or stripped HiTrap Chelating HP or HiTrap IMAC HP/FF, or stripped HisTrap HP columns with metal ions (e.g., Ni²⁺). In a Prepare and Maintain run several columns can be prepared simultaneously.

Note: Put the outlet tubing F11 in a liquid container for collection of any metal waste.

CIP Column (Cleaning-In-Place)

The purpose of CIP column is to clean contaminated columns.

CIP is normally included in a **Purify** method plan as a **Column Post Run**. However, **Customized column procedures** can be performed in a **Prepare and Maintain run**. The procedure cleans up to 5 columns of one type with up to 9 different solutions.

Note: Chelating columns have to be stripped before CIP with NaOH

Recommended solutions and inlets are listed in "A5 Recommended solutions and inlets for CIP in Customized column procedures (Prepare and Maintain)" on page 23. For detailed column cleaning procedures and column storage instructions, please refer to the instructions supplied with each column or to the web page: cytiva.com/protein-purification-labresearch


5.3 System procedures Manual filling of the inlet tubing

Either use the interactive procedure of a method plan **Fill Inlets Manually Using Syringe in Prepare and Maintain** or handle it manually in **System Control**.


Filling inlet tubing S1–S4, A1–A8, and B1–B2 using **Fill Inlets Manually Using Syringe** protocol

- Immerse tubing in the liquid containers.
- Start the **Fill Inlets Manually Using Syringe** protocol which will start with inlet **S1**.


- Connect a syringe to the purge valve.



- Turn the purge valve to open it.



- Slowly draw solution into the syringe. When fluid starts to enter the syringe, continue to draw a few milliliters before closing the purge valve. Check that there is no visible air left in the tubing.



- Press the **Next Breakpoint** button to turn the inlet valve to the next position, **S2**. Repeat the procedure for **S3–S4, A1–A8** and **B1–B2**.

Filling inlet tubing manually in **System Control**

To fill the inlets manually in **System Control**:

- Select which system to fill by left-clicking on the symbol for the system to the left in System control.
- Set the inlet valve to the appropriate position in **System Control:Manual:Flowpath**.
- Connect the syringe to the purge valve and fill the inlets as described above.
- To fill inlet **B1** and **B2**:
- Start a low flow in **System Control:Manual:Pump:Flow**.
- Set inlet valve to B1 or B2:
- In **System Control:Manual:Pump:Gradient**, select **Target 100% B** and **Mode A1/B1** to fill **B1** or **Mode A2/B2** to fill **B2**. Wait for the valve to turn (a clicking sound) before starting the purging procedure.
- When all inlets are filled, click **End**.

Purge pump with methanol

(Select on the **Prepare and Maintain** page.)

Both pump heads must be purged with methanol if:

- The system has been left unused for a week or longer, or
- The pump has been run dry.

Purging will maintain the pumping capacity and protect the pump piston seals.

Approximately 150 ml degassed 100% methanol from inlet A8 and 200 ml water from inlet A5 is needed.

Clean sample inlets

Cleaning of the sample inlet tubing can be included in a **Purify** run (**Clean Sample Inlets after Sample Loading**).

Note: The procedure is only available for protocols starting with AC/IEX.

After the samples have been loaded, the system will pause and a message will appear on the screen. Sample inlet tubing should then carefully be immersed into a separate bottle containing the appropriate wash solution, (e.g., water or buffer).

Clean pump

The procedure can be included in a **Purify** run (**System Procedures:Clean Pump with 0.5 M NaOH between Loading of Different Samples**)

After a sample has been loaded via the pump, the pump is cleaned with 0.5 M NaOH followed by water and AC/IEX binding buffer.

Wash frac tubing

The procedure can be included in a **Purify** run (**System Procedures:Wash Frac Tubing Between Sample**)

The fraction collector tubing, FracCollF2, will be washed with 2 ml buffer between fractionation of samples. One microplate well will be used for each wash.

Rinse all outlets

The procedure can be included in a **Prepare and Maintain** method by selecting **Standard System and Column Procedures:Rinse All Outlets**.

All outlet tubing F1-F11, including FracCollF2, will be washed with solution from inlet A1.

Note: A microplate must be inserted in the fraction collector to collect solution from FracCollF2.

Remove ethanol from the system

Removal of ethanol from the system can be included in a **Purify** run by selecting **Remove Ethanol From System** as **Preparation** on the **System Procedures** page.

In the beginning of the purification run, the system is washed with water. The loops to be used during purification will also be washed.

It is also possible to select the procedure in **Prepare and Maintain** by selecting **Standard System and Column Procedures:Remove Ethanol from System and Loops**.

Fill system with ethanol

The procedure can be included in a **Purify** run by selecting **Fill System with Ethanol** as **Post Run** on the **System Procedures** page.

After the purification has been completed, the system and the loops used during purification are filled with 20% ethanol.

It is also possible to select the procedure in **Prepare and Maintain** by selecting **Standard System and Column Procedures:Fill System+Loops+Outlets with Ethanol**.

Clean System

Cleaning of the system is usually named CIP System (Cleaning-In-Place). A CIP method can be created on different pages in the Wizard.

a Included in a purification run

A CIP System method can be included in a **Purify** run by selecting **CIP System then Fill with Ethanol** or **CIP System then Fill with Buffer** as **Post Run** on the **System Procedures** page.

CIP System then Fill with Ethanol: After the purification run, the system and the loops used

during purification are cleaned with 0.5 M NaOH followed by water, affinity binding buffer, water and finally 20% ethanol.

CIP System then Fill with Buffer: After the purification run, the system and the loops used during purification are cleaned with 0.5 M NaOH followed by water and affinity binding buffer.

b Standard method in Prepare and Maintain

The same type of CIP System methods can be created on the **System Procedures** page in **Prepare and Maintain** by selecting **Standard System and Column Procedures**.

c Customized method

Customized system procedures can be performed in a **Prepare and Maintain** run. **Customized System Procedures** is a procedure that cleans the system with up to 5 different cleaning solutions. The selected parts of the system will be washed with one cleaning solution at a time. The system will pause and a message will appear when the inlet tubing should be inserted into a new solution.

Note: Recommended cleaning solutions are 0.5 M NaOH, deionized water and buffer. When the system will be unused for a longer time, fill the system with water after recommended cleaning solutions and then 20% ethanol as the last step.

6 Evaluating the results

When a system has entered **End** status SYS_1 the result can be evaluated and the samples can be pooled.

6.1 Find and open the result files

- Use the **Recent Runs** or the **Find** tab, in the Evaluation module to locate the result file.

Note: Automatic pooling is not performed if the result file is opened via the **Files** tab.

- Click **+** to expand the list for the result file.
- Double-click a sample chromatogram to open it.

A result file may include chromatograms of:

- equilibration
- sample loading
- tag cleavage
- purification of sample(s)
- CIP/fill with ethanol etc. performed after the purification run

- If necessary, adjust the suggested pooling.
- If the extinction coefficient has been entered in the wizard, concentration and amount are automatically calculated. *Otherwise:* enter the extinction coefficient manually by marking a pool and then typing the value in the extinction coefficient field. The concentration is given in mg/ml, and amount is given in mg.
- The linearity of the UV monitor is limited. For more information, see Technical specification in the User Manual.
- Click the **Add to Pooling Protocol** button to add the adjusted pools to the Pooling protocol.
- Repeat the procedure for other chromatograms from the same, or other, result files.
- Click the **View Pooling Protocol** button.


6.3 Print or save the Pooling protocol

The Pooling protocol can be used as a help when making the physical pooling of the purified samples from the microplate.

- To print the Pooling protocol: Click the **Print** button to print the protocol on the default printer.
- To save the Pooling protocol as a file: Click **Export** and save the protocol in one of the following formats: text (.txt), Excel (.xls), HTML (.htm).

6.2 Adjust pooling and add to Pooling protocol

The chromatogram is displayed and UNICORN will automatically display a suggested pooling of the fractions. The pooled fractions are listed in a table below the chromatogram and the pooled peaks are numbered sequentially in the chromatogram.




Only adjacent fractions will be pooled. The fraction numbers for each pool are listed in the table as a range in retention order, (e.g., A6–A7 etc.)

If the pooling suggestion is not performed, choose **Operations:Pool**.


Pooling protocol

System	Result	Sample Id	Pool	Vol.	Conc.	Text	Target vol.
SYS_1							
	AC D59001						
		alpha-1					
			A2-A5	7.5620	2.9275		2.2138
		alpha-2					
			A7-B12	12.7663	3.7228		4.7527
		alpha-3					
			B10-B5	11.2592	3.5886		4.0405
		alpha-4					
			B3-C2	9.5732	0.8292		0.7938




6.4 Print the Report

- Select **File:Report** or click the **Report** icon.



- Select **xpress 1sample...**, **xpress 2samples...**, **xpress 3samples...**, or **xpress 4samples...** report format, depending on how many samples are included in the result file.
- Select **xpress 1sample** when the sample has been loaded with the sample pump. Select **xpress 1sample Loop** when the sample has been loaded from a loop or Superloop.
- The report format **xpressActiveChrom** can be selected if only the active zoomed window in the Evaluation module should be included in the report.
- Click **Edit** or **Preview** to inspect and modify the report format if needed.
- Click the **Print** button. Select the printer from the Main UNICORN window.
- Choose which pages and how many copies to print and click **OK**.




Getting help

Manuals and Cue Cards

- Open the **Help** menu and select


Manuals to access the :

- UNICORN User Manual
- ÄKTExpress User Manual
- ÄKTExpress Cue Cards




Context-specific Help

- Click the **Help** button in the dialog box, or press **F1**.



Example of Method Editor wizard help:



Appendix A

A1 Columns supported (examples)

Column	Volume [ml]
Any HiTrap AC ¹	1 or 5
GSTrap™ FF	1 or 5
GSTrap 4B	1 or 5
GSTrap HP	1 or 5
HiTrap HP	1 or 5
HiTrap FF	1 or 5
HiTrap FF Crude	1 or 5
HiTrap Chelating HP	1 or 5
HiTrap IMAC HP	1 or 5
HiTrap IMAC FF	1 or 5
2 × HiTrap Desalting	10
HiPrep 26/10 Desalting	53
Any HiTrap IEX ²	1 or 5
HiTrap Q HP	1 or 5
HiTrap SP HP	1 or 5
MonoQ 5/50 GL	1
MonoS 5/50 GL	1
RESOURCE Q	1 or 6
RESOURCE S	1 or 6
HiLoad 16/600 Superdex 75 or 200 prep grade	120
HiLoad 26/600 ³ Superdex 75 or 200 prep grade	320
HiPrep 16/60 ⁴ Sephacryl S-100 or S-200 or S-300 HR	120
HiPrep 26/60 Sephacryl S-100 or S-200 or S-300 HR	320

¹ Any HiTrap means that any HiTrap Affinity column can be used.

² Any HiTrap IEX means that any HiTrap ion exchange column can be used.

³ An optional column holder is available.

⁴ An optional column wrap is available.

A3 Typical run times

To estimate the total run time, see the Summary. Examples of estimated run times for some different types of runs are shown below.

A3.1 Typical run times for the different protocols

Protocol	Max. No. of samples	Typical run time ¹ (h)
1 Step		
AC/IEX	4	2.5
DS	4	0.5
GF	4	14
2 Step		
AC/IEX-DS	4	4.5
AC/IEX-GF	4	18
DS-AC/IEX	2	3
3 Step		
AC/IEX-(DS)-AC/IEX	3	7
DS-AC/IEX-DS	2	4
DS-AC/IEX-GF	2	11
4 Step		
AC/IEX-(DS)-AC/IEX-DS	2	5.5
AC/IEX-(DS)-AC/IEX-GF	2	12.5

¹ Run times are approximate and are valid at room temperature.

Sample loading time is not included. The largest columns have been used. AC: 5 ml, DS: HiPrep 26/10 Desalting, IEX: 6 ml, GF: HiPrep 26/60 Sephacryl,

A2 Column positions

Protocol	Column position					Max no of samples	On-column tag cleavage possible	
	1	2	3	4	5			
1-step	AC/IEX	AC/IEX	AC/IEX	AC/IEX	AC/IEX	-	4	Yes
	DS/GF	-	-	-	-	DS/GF	4 ¹	No
2-step	AC-DS/GF	AC	AC	AC	AC	DS/GF	4	Yes
	IEX-DS/GF	IEX	IEX	IEX	IEX	DS/GF	4	Yes ²
	DS-AC/IEX	AC/IEX	-	-	-	DS	2	No
3-step	AC-DS-AC/IEX	AC	AC	AC	AC(2)/IEX	DS	3	Yes
	IEX-DS-AC/IEX	IEX	IEX	IEX	AC/IEX(2)	DS	3	Yes ²
	DS-AC/IEX-DS/GF	AC/IEX	-	-	DS(1)	DS(2)/GF	2	No
4-step	AC-DS-IEX-DS/GF	AC	AC	IEX	DS(1)	DS(2)/GF	2	Yes
	AC-DS-AC-DS/GF	AC	AC	AC(2)	DS(1)	DS(2)/GF	2	Yes
	IEX-DS-AC-DS/GF	IEX	IEX	AC	DS(1)	DS(2)/GF	2	Yes ²
	IEX-DS-IEX-DS/GF	IEX	IEX	IEX(2)	DS(1)	DS(2)/GF	2	Yes ²

¹ Or one sample injected four times when using Superloop

² Note: This may not be a useful protocol for on-column tag-cleavage if a tagged protease is used.

A3.2 Typical run times for Included System procedures

SYSTEM PROCEDURES	Typical run time ¹ (min)
Preparation	
Remove EtOH	10
During Run	
Clean Sample inlets (4 samples)	5
Wash Frac Tubing (4 samples)	5
Post Run	
Fill with EtOH	10
CIP then fill with EtOH	35
CIP then fill with Buffer	25

¹ Run times are approximate and are valid at room temperature.

A3.3 Typical run times for Included Column procedures

COLUMN PROCEDURES	Typical run time ¹ /column (min)
Preparation	
AC Equilibration	5
AC Blank run	25
AC Remove EtOH before preparation	5
DS Equilibration	10
DS Remove EtOH before preparation	5
IEX Equilibration	5
IEX Blank run	25
IEX Remove EtOH before preparation	5
GF Equilibration	325
GF Remove EtOH before preparation	155

Typical run time¹
/column
(min)

COLUMN PROCEDURES	Typical run time ¹ /column (min)
Post run	
AC Fill with EtOH	5
AC CIP Fill with EtOH	25
AC CIP+ Re-equilibrate	20
AC Re-equilibrate	5
AC Strip+ Fill with EtOH	10
AC Strip+ Recharge+ Blank Run	40
AC Strip+ Recharge+ Fill with EtOH	20
DS Fill with EtOH	25
DS CIP Fill with EtOH	50
DS CIP+Re-equilibrate	25
IEX Fill with EtOH	5
IEX CIP Fill with EtOH	30
IEX CIP+Re-equilibrate	20
IEX Re-equilibrate	5
GF Fill with EtOH	485
GF CIP Fill with EtOH	1130
GF CIP+Re-equilibrate	650

¹ Run times are approximate and are valid at room temperature. Sample loading time is not included.

The largest columns have been used. AC: 5 ml, DS: HiPrep 26/10 Desalting, IEX: 6 ml, GF: HiPrep 26/60 Sephacryl.

A4 Buffers

Buffer types needed

for AC:	<ul style="list-style-type: none"> • binding buffer • cleavage buffer (optional) 	<ul style="list-style-type: none"> • elution buffer • extra wash buffer (optional)
for DS:	• one buffer type per desalting step	
for IEX:	• binding buffer	• elution buffer
for GF:	• one buffer type per run	

AC buffer suggestions for some His-tagged proteins

If performing...	suggested buffer
binding using HisTrap HP/FF, HisTrap FF crude, HiTrap IMAC HP/FF	20 mM sodium phosphate, pH 7.4, 0.5 M NaCl, 20–40 mM imidazole ¹⁾
binding using HiTrap Chelating	20 mM sodium phosphate, pH 7.4, 0.5 M NaCl, 10–20 mM imidazole ¹⁾
extra wash	20 mM sodium phosphate, pH 7.4, 0.5 M NaCl, 50 mM imidazole
cleavage using TEV protease on HisTrap HP or HiTrap Chelating HP ²⁾	20 mM sodium phosphate, pH 7.4, 0.5 M NaCl, 50 mM imidazole
step or gradient elution	20 mM sodium phosphate, pH 7.4, 0.5 M NaCl, 500 mM imidazole

AC buffer suggestions for GST-tagged proteins

If performing...	suggested buffer
binding using GSTrap HP/FF/4B	10 mM sodium phosphate, pH 7.4, 140 mM NaCl, 1–10 mM DTT ³⁾
cleavage using PreScission™ protease on GSTrap HP or FF ³⁾	50 mM Tris-HCl pH 7.5, 150 mM NaCl, 1 mM EDTA, 1 mM DTT
step or gradient	50 mM Tris-HCl,

DS buffer suggestions

If preparing for...	suggested buffer
AC	Use the AC binding buffer that will be used in the following step
IEX	Use the IEX binding buffer that will be used in the following step
protein storage	include e.g., 10% glycerol in a suitable buffer, (e.g., 50 mM Tris-HCl, pH 7.5, 150 mM NaCl).

IEX buffer suggestions

If for example...	suggested buffer (depends on protein pI)
binding to AIEX	50 mM Tris-HCl pH 8.0
binding to CIEX	20 mM MES pH 6.0
elution from AIEX	50 mM Tris-HCl pH 8.0, 1 M NaCl
elution from CIEX	20 mM MES pH 6.0, 1 M NaCl

GF buffer suggestions

If preparing for...	suggested buffer
further studies	a suitable buffer, e.g., 50 mM Tris-HCl pH 7.5, 150 mM NaCl
protein storage	include e.g., 10% glycerol in a suitable buffer (e.g., 50 mM Tris-HCl pH 7.5, 150 mM NaCl)

¹⁾ The imidazole concentration is protein dependent.

²⁾ In the affinity tag removal protocols, the affinity columns can be equilibrated with either AC-binding/wash buffer or an alternative cleavage buffer prior to protease injection.

³⁾ DTT enhance GST-binding.

⁴⁾ Higher amounts of reduced glutathions may give sharper elution peaks.

Tip

Buffer volumes required for a run

See the Summary for the specific run (see “2.1 Print out a summary” on page 9).

For other solutions (metal ion charging etc): see “5.2 Column procedures” on page 15.

Tip

Use high purity liquids

For best purification results, use high purity deionized water and chemicals. Filtering of liquids through a 0.45 µm filter and degassing the liquids is recommended.

A5 Recommended solutions and inlets for CIP in Customized column procedures (Prepare and Maintain)

AC DS IEX GF

Note: Ensure that the columns are filled with water prior to CIP. If required, select **Wash Columns with Water before CIP**.

If the columns will not be used for some days, fill with 5 CV 20% ethanol (HiTrap SP and RECOURCE S: 20% ethanol, 0.2 M sodium acetate), as step 9.

Column examples	Step 1					Step 2				
	Solution	Inlet	Volume	Flow [ml/min]		Solution	Inlet	Volume	Flow [ml/min]	
			RT	CR				RT	CR	
HiTrap HP/FF/FF crude, HiTrap Chelating/IMAC HP/IMAC FF (1 ml / 5 ml)	1 M NaOH	A7	1 CV ¹	1/5	0.8/4	Deionized water	A5	10 CV	1/5	0.8/4
GSTrap HP/FF/4B (1 ml / 5 ml)	6 M GuaHCl	A7	2 CV	1/5	0.8/4	GST A-buffer	A1	5 CV	1/5	0.8/4
HiTrap Desalting	1 mg pepsin ⁴ /ml in 0.1 M acetic acid, 0.5 M NaCl	A7 or	1 CV ²	5	4	0.2 M NaOH	A3	2 CV	5	4
HiPrep 26/10 Desalting	1 mg pepsin ⁴ /ml in 0.1 M acetic acid, 0.5 M NaCl	A7 or	1 CV ²	10	8	0.2M NaOH	A3	2 CV	10	8
HiTrap Q/SP HP (1ml / 5ml)	2 M NaCl	B2	4 CV	0.25/1.3	0.20/1	Deionized water	A5	2 CV	0.25/1.3	0.20/1
RESOURCE Q/S (1 ml / 6 ml)	1 M NaCl	A2	5 CV	4/6	3.2/4.8	Deionized water	A5	2 CV	4/6	3.2/4.8
Mono Q/S	1 M NaCl	A2	4 CV	0.5	0.4	Deionized water	A5	2 CV	0.5	0.4
HiLoad Superdex ⁵ pg (16/600 / 26/600)	1 M NaOH	A7	4 CV	0.5/1.3	0.4/1	Deionized water	A5	4 CV	0.5/1.3	0.4/1
HiPrep Sephacryl ⁶ HR (16/60 / 26/60)	0.5 M NaOH	A6	0.25 CV	0.3/0.8	0.2/0.6	Deionized water	A5	4 CV	0.3/0.8	0.2/0.6

Column examples	Step 3					Step 4				
	Solution	Inlet	Volume	Flow [ml/min]		Solution	Inlet	Volume	Flow [ml/min]	
			RT	CR				RT	CR	
HiTrap HP/FF/FF crude, HiTrap Chelating/IMAC HP/IMAC FF (1 ml / 5 ml)	30% isopropanol	A4	10 CV	0.5/2.5	0.4/2	Deionized water	A5	10 CV	1/5	0.8/4
GSTrap HP/FF/4B (1 ml / 5 ml)	1% Triton X-100 ³	A7	2 CV	1/5	0.8/4	GST A-buffer	A1	5 CV	1/5	0.8/4
HiTrap Desalting	Deionized water	A5	5 CV	7.5	6					
HiPrep 26/10 Desalting	Deionized water	A5	5 CV	15	12					
HiTrap Q/SP HP (1ml / 5ml)	1 M NaOH	A7	4 CV	0.25/1.3	0.20/1	Deionized water	A5	4 CV	0.25/1.3	0.20/1
RESOURCE Q/S (1 ml / 6 ml)	1 M NaOH	A7	5 CV	4/6	3.2/4.8	Deionized water	A5	2 CV	4/6	3.2/4.8
Mono Q/S	1 M NaOH	A7	4 CV	0.5	0.4	Deionized water	A5	2 CV	0.5	0.4
HiLoad Superdex ⁵ pg (16/600 / 26/600)	30% isopropanol	A4	0.5 CV	0.5/1.3	0.4/1	Deionized water	A5	2 CV	0.5/1.3	0.4/1
HiPrep Sephacryl ⁶ HR (16/60 / 26/60)	30% isopropanol	A4	0.5 CV	0.3/0.8	0.2/0.6	Deionized water	A5	4 CV	0.3/0.8	0.2/0.6

Column examples	Step 5					Step 6				
	Solution	Inlet	Volume	Flow [ml/min]		Solution	Inlet	Volume	Flow [ml/min]	
			RT	CR				RT	CR	
HiTrap HP/FF/FF crude, HiTrap Chelating/IMAC HP/IMAC FF (1 ml / 5 ml)	2 M NaCl	B2	1 CV	1/5	0.8/4	Deionized water	A5	10 CV	1/5	0.8/4
HiTrap Q/SP HP (1ml / 5ml)	30% isopropanol	A4	2 CV	0.25/1.3	0.20/1	Deionized water	A5	4 CV	0.25/1.3	0.20/1
RESOURCE Q/S (1 ml / 6 ml)	1 M HCl	A1	5 CV	4/6	3.2/4.8	Deionized water	A5	2 CV	4/6	3.2/4.8
Mono Q/S	1 M HCl	A1	4 CV	0.5	0.4	Deionized water	A5	2 CV	0.5	0.4

Column examples	Step 7					Step 8				
	Solution	Inlet	Volume	Flow [ml/min]		Solution	Inlet	Volume	Flow [ml/min]	
			RT	CR				RT	CR	
RESOURCE Q/S (1 ml / 6 ml)	1 M NaCl	A2	5 CV	4/6	3.2/4.8	Deionized water	A5	2 CV	4/6	3.2/4.8
Mono Q/S	1M NaCl	A2	4 CV	0.5	0.4	Deionized water	A5	2 CV	0.5	0.4

¹ Pause 1 hour.

² Pause over night at RT.

³ Alternatively 3-4 CV 70% ethanol.


⁴ If pepsin is used, ensure it has been completely eluted or has been denatured by NaOH.

⁵ HiLoad 16/600 and 26/600, Superdex 75 and 200 columns.

⁶ HiPrep 16/60 and 26/60, S-100, S-200 and S-300 columns.


Appendix B Flow charts

- Valve positions after reset (END), no flow




- Column equilibration
- Sample loading


- Wash out unbound sample
- Elution




- Intermediate peak collection in capillary loops
- System wash



- Sample loading from loops
- Loop wash





Ordering information

Product	Pack Size	Code Number
---------	-----------	-------------

Systems

ÅKTExpress	1	18-6645-01
USB/CAN device ³	1	28-9692-01
ÅKTExpress software & manuals ¹	1	28-9053-45

Accessories

Flask holder	1	18-1177-79
Large column holder ÅKTExpress	1	28-4007-37
Capillary loop	1 x 10 ml	11-0003-02
Red fingertights	8	28-4010-81
Tubeing 1 mm diameter	1 x 3 m	18-1142-38
Tubeing 1.6 mm diameter	1 x 3 m	18-1121-16
HiPrep 16/xx Column Wrap	1	28-9021-50
Union 1/16" female to 1/16" female	1	11-0003-39
Union M6 male to 1/16" female	1	18-1112-57
Union M6 female to 1/16" male	1	18-1112-58
Union M6 female to 1/16" female (PEEK)	1	18-1123-94
On-line filter	1	18-1112-44

Supported columns

Affinity chromatography

GSTrap FF	5 x 1 ml	17-5130-01
GSTrap FF	100 x 1 ml ²	17-5130-05
GSTrap FF	5 x 5 ml	17-5131-02
GSTrap FF	100 x 5 ml ²	17-5131-05
GSTrap 4B	5 x 1 ml	28-4017-45
GSTrap 4B	5 x 5 ml	28-4017-48
GSTrap 4B	100 x 1 ml ²	28-4017-46
GSTrap 4B	100 x 5 ml ²	28-4017-49
GSTrap HP	5 x 1 ml	17-5281-01
GSTrap HP	100 x 1 ml ²	17-5281-05
GSTrap HP	5 x 5 ml	17-5282-02
GSTrap HP	100 x 5 ml ²	17-5282-05
HisTrap excel	5 x 1 ml	17-3712-05
HisTrap excel	5 x 5 ml	17-3712-06
HisTrap HP	5 x 1 ml	17-5247-01
HisTrap HP	100 x 1 ml ²	17-5247-05
HisTrap HP	5 x 5 ml	17-5248-02
HisTrap HP	100 x 5 ml ²	17-5248-05
HisTrap FF	5 x 1 ml	17-5319-01
HisTrap FF	5 x 5 ml	17-5255-01
HisTrap FF	100 x 1 ml ²	17-5319-02
HisTrap FF	100 x 5 ml ²	17-5255-02

HisTrap FF Crude	5 x 1 ml	11-0004-58
HisTrap FF Crude	5 x 5 ml	17-5286-01
HisTrap FF Crude	100 x 1 ml ²	11-0004-59
HisTrap FF Crude	100 x 5 ml ²	17-5286-02
HiTrap Chelating	5 x 1 ml	17-0408-01
HiTrap Chelating	5 x 5 ml	17-0409-03
HiTrap Chelating	100 x 5 ml ²	17-0409-05
HiTrap IMAC HP	5 x 1 ml	17-0920-03
HiTrap IMAC HP	5 x 5 ml	17-0920-05
HiTrap IMAC FF	5 x 1 ml	17-0921-02
HiTrap IMAC FF	5 x 5 ml	17-0921-04
HiTrap TALON [®] crude	5 x 1 ml	28-9537-66
HiTrap TALON crude	5 x 5 ml	28-9537-67
HiTrap TALON crude	100 x 1 ml ²	28-9538-05
HiTrap TALON crude	100 x 1 ml ²	28-9538-09

Desalting

HiPrep 26/10 Desalting	1 x 53 ml	17-5087-01
HiPrep 26/10 Desalting	4 x 53 ml	17-5087-02
HiTrap Desalting	5 x 5 ml	17-1408-01
HiTrap Desalting	100 x 5 ml ²	11-0003-29

Ion Exchange chromatography

HiTrap Q HP	5 x 1 ml	17-1153-01
HiTrap Q HP	5 x 5 ml	17-1154-01
HiTrap SP HP	5 x 1 ml	17-1151-01
HiTrap SP HP	5 x 5 ml	17-1152-01
Mono Q 5/50 GL	1 x 1 ml	17-5166-01
Mono S 5/50 GL	1 x 1 ml	17-5168-01
RESOURCE Q	1 x 1 ml	17-1177-01
RESOURCE Q	1 x 6 ml	17-1179-01
RESOURCE S	1 x 1 ml	17-1178-01
RESOURCE S	1 x 6 ml	17-1180-01

Code numbers for more accessories and user replaceable spare parts can be found in the ÅKTExpress User Manual.

¹ Includes both ÅKTExpress and ÅKTExpress Mab software.

² 100-packs are special packs delivered on customer order. Includes connector package, domed nuts and instructions.

³ Needed between the first system and the computer.

Gel filtration

<i>HiLoad 16/600 Superdex 75</i>		
<i>prep grade</i>	1 x 120 ml	17-1068-01
<hr/>		
<i>HiLoad 16/600 Superdex 200</i>		
<i>prep grade</i>	1 x 120 ml	17-1069-01
<hr/>		
<i>HiLoad 26/600 Superdex 75</i>		
<i>prep grade</i>	1 x 320 ml	17-1070-01
<hr/>		
<i>HiLoad 26/600 Superdex 200</i>		
<i>prep grade</i>	1 x 320 ml	17-1071-01
<hr/>		
<i>HiPrep 16/60 Sephacryl</i>		
<i>S-100 HR</i>	1 x 120 ml	17-1165-01
<hr/>		
<i>HiPrep 16/60 Sephacryl</i>		
<i>S-200 HR</i>	1 x 120 ml	17-1166-01
<hr/>		
<i>HiPrep 16/60 Sephacryl</i>		
<i>S-300 HR</i>	1 x 120 ml	17-1167-01
<hr/>		
<i>HiPrep 26/60 Sephacryl</i>		
<i>S-100 HR</i>	1 x 320 ml	17-1194-01
<hr/>		
<i>HiPrep 26/60 Sephacryl</i>		
<i>S-200 HR</i>	1 x 320 ml	17-1195-01
<hr/>		
<i>HiPrep 26/60 Sephacryl</i>		
<i>S-300 HR</i>	1 x 320 ml	17-1196-01

Superloop

<i>Superloop</i>	1 x 10 ml	18-1113-81
<i>Superloop</i>	1 x 50 ml	18-1113-82

Proteases

<i>PreScission Protease</i>	500 units	27-0843-01
<i>Thrombin</i>	500 units	27-0846-01
<i>Factor Xa</i>	400 units	27-0849-01

Documents

<i>UNICORN user manual</i>	1	11-0003-68
<i>ÄKTExpress User Manual</i>	1	28-4090-22
<i>ÄKTExpress Operating instructions</i>	1	28-9579-08
<i>ÄKTExpress Installation Guide</i>	1	28-4090-29
<i>Superloop Instructions</i>	1	56-3015-99
<i>Affinity Chromatography Handbook</i>	1	18-1022-29
<i>Gel Filtration Handbook</i>	1	18-1022-18
<hr/>		
<i>Strategies for Protein Purification Handbook</i>	1	28-9833-31
<hr/>		
<i>Recombinant Protein Purification Handbook</i>	1	18-1142-75
<i>Antibody Purification Handbook</i>	1	18-1037-46
<hr/>		
<i>Ion Exchange Chromatography & Chromatofocusing</i>	1	11-0004-21

All documents can be viewed at:
cytiva.com



cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

GSTrap, HiLoad, HiPrep, HiTrap, HisTrap, Mono Q, Mono S, PreScission, RESOURCE, Sephacryl, Superdex, Superloop and UNICORN are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

TALON is a registered trademark of Clontech Laboratories Inc, Superflow is a trademark of Sterogene Bioseparations Inc, Microsoft and Excel are registered trademarks of Microsoft Corporation.

All other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

UNICORN © 2020–2021 Cytiva

Any use of UNICORN is subject to Cytiva Standard Software License End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact.

28409023 AC 04/2021