


Immobiline DryPlate

Instructions for Use



1 Introduction

Immobiline™ DryPlate offers a convenient and reliable way to obtain the utmost separation power of isoelectric focusing.

The Immobiline system has indefinitely stable pH gradients allowing high voltages for maximal separation and, when necessary, long focusing times (1, 2).

The rehydratable dry gels facilitate the use of additives such as urea, detergents, carrier ampholytes etc, for optimal performance, even for samples with poor solubility.

This manual gives general instructions on how to use Immobiline DryPlate for isoelectric focusing. Please consult the Application Notes and/or the articles in the reference list (3–21) for detailed instructions on specific applications.

Immobiline DryPlate is a polyacrylamide gel with an immobilized pH gradient. It is bound to plastic backing and is ready to use for isoelectric focusing after rehydration. The product is available with various pH gradients (see the Table below). The pH gradients are linear over the stated interval.

Table 1.1:

Product code	pH interval	Major application	Appl.Note
80112828	4 to 7	General purpose	-
80112829	4.2 to 4.9	α_1 -antitrypsin	470
80112830	4.5 to 5.4	Group specific component	471
80112831	5.0 to 6.0	Transferrin	472
80112832	5.6 to 6.6	Phosphoglucumutase	473

1.1 Package contents and technical data

Package contents

Each gel package contains 3 gels, filter papers, experimental result forms and instructions.

Designation	No. per pack.	Product code
Immobiline DryPlate	3	(See label)
Filter paper	50	
Experimental result form	3	
Instructions	1	71703001

Technical data

Gel dimensions:	Approx. 250 × 110 × 0.5 mm
Gel matrix:	Polyacrylamide T=4%, C=3%
Buffering capacity:	3 meqv/pH/L
Gel backing:	Polyester film
Storage:	-20°C
Shelf life:	18 months from manufacturing. Please observe "Expiry date" printed on each kit.

2 Preparing the gel

2.1 Preparing the rehydration solution

One of the advantages of the dry gel format is the opportunity to include different additives in the reswelling solution (10, 14–16, 18–19). The three options given in the Table below should therefore only be regarded as typical examples that will give good results for most applications. However, whenever required these recipes can easily be modified for further optimization. Consult the relevant Application Note ([Table 1.1, on page 2](#) and Reference list) or any of the cited references for instructions about specific applications. The rehydration process itself has also been investigated (20, 21).

Table 2.1:

	Alt.1	Alt.2	Alt.3
Application areas:	Water soluble proteins	Proteins with reduced solubility	Proteins with low solubility, e. g. Membrane proteins Lipoproteins
Composition:			
Distilled water	20.0 ml	19.5 ml	12.0 ml
Pharmalyte™ 3 to 10/			
Ampholine pH 3.5 to 9.5		0.5 ml	0.5 ml
Urea			9.6 g
Triton™ X-100			0.1 ml
DTT			60 mg
Total volume:	20.0 ml	20.0 ml	20.0 ml
Rehydration time:	1 to 2 h	1 to 2 h	15 to 18 h

Comments

Note: All chemicals should be of the highest purity. PlusOne chemicals are highly recommended. Double-distilled water should be used.

The presence of carrier ampholytes not only increases protein solubility but also their electrophoretic migration velocity resulting in shorter focusing times.

Mercaptoethanol (2%) or dithiothreitol (15 to 50 mmol/l) can be added to avoid oxidation of sensitive proteins.

Glycerol (20% to 25%) improves solubility of hydrophobic proteins and reduces the risk for urea crystallization. Lateral band spreading can be reduced by adding acetic acid (2 mmol/l) and applying the sample at the anodic side or adding Tris (2 mmol/l) and applying the sample at the cathodic side.

Triton X-100 can be replaced with other non-ionic or zwitterionic detergents, e.g. CHAPS. Other carrier ampholytes than Pharmalyte 3 to 10/Ampholine pH 3.5 to 9.5 may also be used.

Alt. 3 in [Table 2.1, on page 4](#) corresponds to what is used in the first dimension focusing in 2-D electrophoresis. This alternative can be regarded as a standard choice for focusing under denaturing conditions and will normally give high quality results with all kind of samples.

2.2 Opening the package

Note: *Wear clean gloves to avoid contamination of the gel surface, particularly when using silver stain.*

Note: *The gel is packed so that it is faced down to the aluminium foil backing of the package, and the gel support is uppermost.*

Note: *If only half of the gel is to be used, cut the package in half with sharp scissors, reseal the portion to be saved with tape, and store it at -20°C. Remember to identify the polarity of the remaining part.*

Open the gel package from the transparent side. Use scissors to cut around all four sides of the package, taking care not to cut either the gel or its transparent backing film.

To simplify gel handling later on, identify the polarity of the pH gradient. The support film has a precut corner which indicates the anodic side of the pH gradient.

2.3 Rehydrating Immobiline DryPlate

For this procedure the specially designed Reswelling Cassette is highly recommended. It allows fast, convenient, even and reproducible rehydration of the gel. It also facilitates including the additives necessary for optimal performance for each application. Proceed as follows: (See the instruction manual for the Reswelling Cassette or the Multiphor™ II Electrophoresis System user manual for detailed instructions).

Step	Action
1	To prevent the gel from adhering to the glass plate fitted with the U-frame, coat the plate with Repel-Silane.
2	Mark the cathodic side of the gel.

Step	Action
3	Wet a clean thick glass plate with a few drops of water and place the gel on the glass plate with the gel side up.
4	Roll the gel with a clean rubber roller (Product code 80110679) to remove all air bubbles from between the glass plate and the support film.
5	Mount the gel in the cassette taking particular care that the U-frame gasket seals also over the cut-off corner of the supporting plastic foil and that the clamps are mounted correctly to avoid leakage.
6	Fill the cassette with the desired rehydrating solution.
7	Leave the gel to rehydrate for the recommended time.
8	Open the cassette and check the gel surface. Remove excess liquid by placing a filter paper moistened in distilled water on top of the gel followed by a dry filter paper on top.
9	Blot the gel by gently rolling the rubber roller under slight pressure over the dry filter paper. Finally remove the filter papers carefully from the gel. (Since gels rehydrated in detergent containing solutions have less tendency to stick to dry filter paper, they can be dried with a simpler procedure: Place a piece of dry tissue paper (e.g. Kleenex) on the gel, press gently to ensure contact between tissue and gel, and remove the paper carefully).

3 Sample treatment

3.1 Sample preparation

Even if Immobiline DryPlate is exceptionally tolerant towards impure samples, best results are still obtained with samples that are free from precipitates. Should aggregation occur at the application point, this can often be overcome by diluting the sample or changing the sample application position.

Best results are generally obtained when the samples are solubilized in the rehydration solution. If this is not possible, the concentration of salt and buffer ions should still be kept at a minimum and, as a general rule, preferably below 50 mmol/l. Excess buffer and salt ions will cause local overheating due to high local currents, which can result in protein denaturation and/or prolonged running times.

Desalting and buffer exchange can be carried out by dialysis, or, more easily, by gel filtration using a prepacked Sephadex™ G-25 column (Choose NAP™-5 Column, NAP-10 Column or PD-10 Column depending on the sample size. See [Ordering information, on page 18](#)).

3.2 Sample concentration

In general, Immobiline gels can take much higher sample loads than corresponding gels with carrier ampholytes. Several factors will determine the optimal sample concentration and volume:

1. pH range.
2. Number and relative proportions of the components in the sample.
3. Sensitivity of the detection method used.

Guidelines: PhastGel™ Blue R stain detects proteins down to the µg range. Normally 15 to 20 µl of sample with a concentration of 0.5 to 10 mg/ml will give good results. Silver staining has about 20 times higher sensitivity. A suitable load in a narrow pH gradient is normally 2 to 3 times higher than the load in a pH 4 to 7 gradient.

4 Isoelectric focusing

4.1 Preparing the experiment

Setting the cooling temperature

Connect Multiphor II electrophoresis unit to MultiTemp™ II thermostatic circulator and set the desired running temperature. A running temperature of 10°C is often used except for gels containing urea, which are preferably run at somewhat higher temperatures (15°C or more) to avoid precipitation of the urea.

Switch on the thermostatic circulator 15 minutes before starting the run.

Positioning the gel on the cooling plate

Pipette a few milliliters of insulating fluid (kerosene or light paraffin oil) on the cooling plate of Multiphor II. Position the gel on the cooling plate so that the polarity of the gel corresponds with the polarity marked on the cooling plate. Ensure that no air bubbles are trapped between the gel and the cooling plate.

Electrode strips are used to ensure good electrical contact between the gel and the electrodes. This prevents sparking and allows salt ions from the gel to migrate into the electrode strips where they will stay and not interfere with the separation.

Soak the electrode strips evenly with approximately 3 ml distilled water and remove the excess by using a filter paper. The electrode strips should appear *very dry* before they are applied to the gel.

Lay the electrode strips along the long edges of the gel. Cut off the parts of the strips which protrude beyond the short ends of the gels using a pair of sharp scissors.

4.2 Sample application

There are basically three different methods for sample application. Which method to choose is determined primarily by the sample volume.

Step	Action
------	--------

- | | |
|---|---|
| 1 | For 5 to 20 µl sample volumes: Apply sample directly on the gel, using Immobiline applicator strip (Product code 18100276). This applicator strip makes it possible to use a multiple syringe which allows quick and simple sample loading, especially when working with large numbers of samples. Check that the contact between the gel and the applicator strip is uniform. Leave the applicator strip on the gel during focusing. |
|---|---|

Step	Action
2	For 15 to 20 μ l (and larger) sample volumes: Use sample application pieces (Product code 80112946). Apply the dry pieces to the gel surface at the desired position(s) in the gradient. Using a micropipette, apply 15 to 20 μ l volumes of sample solution on each piece. To apply larger volumes, use 2 or 3 pieces stacked or laid end-to-end. If you want to apply smaller volumes with by this method, trim the paper proportionally before applying it to the gel. Remove the pieces of paper after completing half the total focusing time.
3	For 2 to 10 μ l sample volumes: Apply the sample as droplets directly on the gel surface.

4.3 Running conditions

Place the electrode holder on the Multiphor II electrophoresis unit and align the electrodes with the center of the electrode strips. Finally, connect the two electrodes to the base unit and place the safety lid in position. Connect the electrode leads to the power supply.

Typical running conditions are listed in the Table below.

Table 4.1:

pH range	Voltage (V)	Current (mA)		Power (W)	Time (KVh)	Time (h)
4 to 7	3500	5.0	(1.0)	15.0	7 to 15	2 to 4
4.2 to 4.9	3500	5.0	(2.0)	15.0	15 to 25	4 to 7
4.5 to 5.4	3500	5.0	(2.0)	15.0	15 to 25	4 to 7
5.0 to 6.0	3500	5.0	(1.0)	15.0	15 to 25	4 to 7
5.6 to 6.6	3500	5.0	(1.0)	15.0	15 to 25	4 to 7

Comments

- Decrease the power and current settings proportionally if only part of the plate is being used.
- The settings above are only to be regarded as guidelines. Some proteins focus very slowly and may require as much as 50 to 60 KVh to give optimal sharp bands. This must be determined experimentally in each case: Run the sample for different times.
- Since there is no gradient drift in the Immobiline DryPlate there is no limitation in the electrophoresis system as such as to how long the experiment can be continued. The limits are set only by the risk of drying out the gel, oxidising or denaturing the sample.

These risks can be minimized by placing a plastic foil on top of the gel, running at low temperatures, flushing the unit with inert gas (N_2) and/or including a reducing agent (DTT or b-mercapto-ethanol) in the rehydration mixture. The surface can also be protected with DryStrip Cover Fluid (22).

5 Detection

All currently used detection methods can be applied on Immobiline DryPlate gels, including Coomassie™ Blue (4, 14), silver staining (23). Possible problems from extensive swelling of the gel can be reduced by adding ethanol (30%) to the washing solutions. Enzymatic- and immunologically-based staining procedures as well as blotting can also be used (10–12).

5.1 Silver Staining

This silver staining method is according to Heukeshoven and Dernick (23) with some modifications.

Silver staining solutions

Note: 250 ml of solutions are needed per gel. All steps should be done with gentle shaking of the tray.

Fixing solution:	Trichloroacetic acid	30.0 g
60 min	Ethanol	125 ml
	Dissolve in distilled water and make up to 250 ml.	
Wash:	Ethanol	150 ml
2 × 15 min	Acetic acid	50 ml
	Make up to 500 ml with distilled water	
Incubation solution:	Ethanol	75 ml
minimum 40 min	Sodiumacetate • 3 H ₂ O	17.0 g
	Glutardialdehyde (25%w/v) ¹	1.3 ml
	Sodium thiosulfate,	0.50 g
	Na ₂ S ₂ O ₃ • 5 H ₂ O	
	Make up to 250 ml with distilled water	
Wash:	Distilled water	
3 × 5 min		
Silver solution:	Silver nitrate	0.25 g
20 min	Formaldehyde ¹	50 µl
	Make up to 250 ml with distilled water	

¹ **Note:**
Add these components immediately before use.

Developing solution:	Sodium carbonate	6.25 g
5 to 15 min	Ethanol	75 ml
	Formaldehyde ²	25 µl
	Make up to 250 with distilled water	
Stop solution:	EDTA-Na ₂ • 2 H ₂ O	3.0 g
10 min	Ethanol	75 ml
	Make up to 250 ml with distilled water	
Wash:	Ethanol	150 ml
5 min	Make up to 500 ml with distilled water	
Preserving solution:	Glycerol (87% w/w)	25 ml
20 min	Ethanol	75 ml
	Make up to 250 ml with distilled water	

Step	Action
1	Fixation: remove the electrode strips by using forceps, thereafter immediately immerse the gel in the fixing solution for 60 minutes. This solution precipitates the proteins and allow detergents (if used) to diffuse out.
2	Washing: thereafter, wash the gel in washing solution for 2 × 15 minutes.
3	Incubation: place the gel in incubation solution for a minimum of 40 minutes. The gel can be left over night in the incubation solution.
4	Washing: thereafter, wash three times in distilled water, each time for 5 minutes.
5	Silver reaction: put the gel in silver solution for 20 minutes.
6	Developing: develop the gel in developing solution for 5 to 15 minutes. The protein bands should become intensively dark.
7	Stopping: stop the reaction by placing the gel in stop solution for 10 minutes.
8	Washing: wash in washing solution for 5 minutes.

² **Note:**
Add these components immediately before use.

Step	Action
9	Preserving: To preserve the silver stained gel, place the gel in preserving solution for 20 minutes. Then place the gel on a glass plate with the gel side up, and cover the gel with cellophane preserving sheet soaked in preserving solution. Allow the gel to dry in room temperature (Do not put the gel in a heating cabinet, the silver stain will bleach due to the increased temperature).

5.2 Coomassie staining

This is a general protein stain which detects protein concentrations down to the µg level. 250 ml of solution is used in each step.

Fixing solution:	Trichloroacetic acid	46 g
30 to 60 min	Sulphosalicylic acid	14 g
	Make up to 400 ml with distilled water	
Destaining solution:	1. Ethanol	500 ml
5 min	Make up to 1000 ml with distilled water	
	2. Acetic acid	160 ml
	Make up to 1000 ml with distilled water	
	Mix 1:1 before use	
Coomassie solution:	PhastGel Blue R	1 tablet
10 min	Dissolve 1 tablet in 400 ml destaining solution. Stir with a magnetic stirrer and heat the solution to 60°C. Filter before use. Use only once.	
	See above	
Destaining solution:	Glycerol	
Until background is clear		40 ml
Preserving solution:	Add 360 ml destaining solution and mix well.	

Step	Action
1	Fixation: remove the electrode strips with forceps. Immediately place DryPlate in Staining Kit containing fixing solution for 30–60 minutes. This solution precipitates the proteins.

Step	Action
2	Destaining: before staining, wash DryPlate in destaining solution for 5 minutes.
3	Staining: pour off the destaining solution and stain DryPlate for 10 minutes in staining solution which has been preheated to 60 °C.
4	Destaining: destain DryPlate with several changes of destaining solution until the background is clear.
5	Preserving: Place the destained DryPlate in the glycerol preserving solution for 30–60 minutes. Then place the gel on a glass plate with the gel side up, and cover the gel with a cellophane preserving sheet soaked in preserving solution. Allow it to dry at room temperature.

5.3 Electrophoretic transfer

Before electrophoretic blotting can take place, the support film must be removed to allow the current to pass through the gel. FilmRemover is ideal for this purpose. After the film and the gel have been separated, the proteins can be transferred to an immobilizing membrane by using the Multiphor II NovaBlot transfer kit. Complete information on running conditions is given in Multiphor II User Manual (Product code 18110343).

6 Evaluation

6.1 Determination of the isoelectric point

Because of the low ionic strength in the gel, the pH gradient cannot be directly measured with a standard surface pH electrode unless carrier ampholytes have been included in the rehydrating solution (21).

An alternative to direct pH measurement is to run pI calibration proteins in parallel with the experimental samples.

For details about the use of pI calibration proteins see the Instruction enclosed with each pI Calibration Kit.

6.2 Densitometric evaluation

ImageMaster 1D Software (Product code 80635037) is a powerful software package for protein quantitation and data analysis. By using ImageMaster 1D Software together with ImageScanner (Product code 18113445), you can capture, store, evaluate, and report all the information contained in your electrophoresis gels. ImageMaster 1D Software automatically selects your lanes, bands, subtracts the background, corrects the smiling, and integrates areas and band volume ($OD \times mm^2$). The software calculates relative amounts, percentages of totals, and amounts of proteins in real quantity units using a calibration curve.

ImageMaster 1D Software also calculates isoelectric points or molecular weights, compares bands across different lanes or gels, and produces hierarchical clustering information. Further lane comparison, databasing, and identification using an internal library can be done by using ImageMaster 1D Database (Product code 80635094) and an add-on package to ImageMaster 1D Software.

7 References

1. Isoelectric focusing in immobilized pH gradients. Principle, methodology and some applications. *J. Biochem. Biophys. Meth.* 6 (1981) 317–339, Bjellqvist B., Ek, K., Righetti, P.G., et al.
2. Righetti, P. G.. Immobilized pH gradients: Theory and Methodology. Vol. 20, Laboratory techniques in Biochemistry & Molecular Biology. Elsevier, 1990.
3. Analysis of alpha1-Antitrypsin phenotypes in Immobiline electrofocusing gels. Application Note 470, GE Healthcare AB.
4. Analysis of Gc phenotypes by IEF in Immobiline gels. Application Note 471, GE Healthcare AB.
5. Analysis of the Transferrin phenotypes in Immobiline electrofocusing gels. Application Note 472, GE Healthcare AB.
6. Analysis of Phosphoglucosyltransferase (PGM1) phenotypes in Immobiline isoelectric focusing gels. Application Note 473, GE Healthcare AB.
7. Subtyping of group specific component (Gc) in human semen, blood and vaginal fluid by isoelectric focusing in immobilized pH gradients. *Electrophoresis* 9 (1988) 602–605, Pötsch-Schneider, L. and Klein, H.
8. Carbohydrate analysis of transferrin subfractions isolated by preparative isoelectric focusing in immobilized pH gradients. *Electrophoresis* 4 (1992) 225–229, G. de Jong, W. L. van Noort and H. G. van Eijk.
9. Isoelectric focusing of apolipoproteins on immobilized pH gradients: Improved determination of apolipoprotein E phenotypes. *Electrophoresis* 9 (1988) 576–579, Baumstark, M. W., Berg, A., Halle, M. and Keul, J.
10. Microheterogeneity of apolipoprotein D as revealed by electroblotting following isoelectric focusing in Immobiline DryPlates. *Electrophoresis* 4 (1992) 262–264, Holmquist, L.
11. Phenotyping of Apolipoprotein-E. Immunoblotting After Isoelectric Focussing in Immobilized pH Gradients *Electrophoresis* 12 (1991) 59–63, Marz, W., Cezanne, S. and Gros, W.
12. Immobilized pH gradient isoelectric focusing and immunoblotting for investigations of anti-Borrelia burgdorferi IgG antibodies. *Electrophoresis* 4 (1992) 229–234, M. Cruz and Å. Sidén.
13. Isoelectric focusing in Immobilized pH Gradients – applications in Clinical Chemistry and Forensic Analysis (review). *J. Chromatogr.* 569 (1991) 197–228, Righetti, P. G., Gianazza, E., Bianchibosisio, A., et al.
14. The application and optimisation of solubilising additives in Immobiline gels. Application Note 345, GE Healthcare AB.

15. Membrane protein analysis by isoelectric focusing in immobilized pH gradients. Electrophoresis 6 (1985) 419–422, Rimpilainen, M. A. and Righetti, P. G.
16. Hybrid isoelectric focusing in rehydrated immobilized pH gradients with added carrier ampholytes: Demonstration of human globulins. Electrophoresis 6 (1985) 314–325, Altland, K. and Rossman, U.
17. Analysis of recombinant proteins by isoelectric focusing in immobilized pH gradients. Electrophoresis 4 (1992) 214–220, R. Bischoff, D. Roecklin and C. Roitsch.
18. Apolipoprotein E phenotyping by isoelectric focusing in immobilized pH gradients and silver staining. Electrophoresis 4 (1992) 252–258, R. Cartier and A. Sassolas.
19. Rapid and simple method for the identification of apolipoprotein E isomorphic phenotypes from whole serum. Electrophoresis 4 (1992) 258–262, M. Kohlmeier, H.-J. Drossel, P. Sinha and E. Köttgen.
20. Swelling kinetics of Immobiline gels for isoelectric focusing. Electrophoresis 5 (1984) 257–262, Gelfi, C and Righetti, P. G.
21. Improved rehydration procedure for polyacrylamide gels in presence of urea: Demonstration of inherited presence of prealbumin variants by isoelectric focusing in an immobilized pH gradient. Electrophoresis 5 (1984) 379–381, Altland, K., Bantzhoff, A., Hackler, R. and Rossman, U.
22. Immobiline DryStrip Kit Instructions (Product code 18103863).
23. Simplified method for silver staining of proteins in polyacrylamide gels and the mechanism of silver staining. Electrophoresis 6, (1985) 103–112, Heukeshoven, J. and Dernick, R.
24. pH measurements in ultranarrow immobilized pH gradients. J. Biol. Biophys. Meth. (1986) 113–124, Gelfi, C., Morelli, A., Rovida, E., and Righetti, P. G.

8 Ordering information

Designation	Product code
Immobiline DryPlate pH 4 to 7	80112828
Immobiline DryPlate pH 4.2 to 4.9	80112829
Immobiline DryPlate pH 4.5 to 5.4	80112830
Immobiline DryPlate pH 5.0 to 6.0	80112831
Immobiline DryPlate pH 5.6 to 6.6	80112832
Roller	80110679
IEF electrode strips (100)	18100440
Immobiline applicator strip (5)	18100276
IEF sample application pieces (200)	80112946
PD-10 column, Desalting samples $\leq$ 2.5 ml (30)	17085101
NAP-5 column, Desalting samples $\leq$ 0.5 ml (20/50)	17085301/02
NAP-10 column, Desalting samples $\leq$ 1.0 ml (20/50)	17085401/02
Multiphor II electrophoresis unit	18101806
EPS 3500 XL Power Supply	18113005
MultiTemp III Thermostatic Circulator, 115 V AC	18110277
MultiTemp III Thermostatic Circulator, 220 V AC	18110278
NovaBlot electrophoretic transfer kit	18101686
FilmRemover	18101375
Cellophane preserving sheets, 210×320 mm (50)	80112938
ImageMaster 1D Software	80635037
ImageScanner	18113445
PhastGel Blue R	17051801
Broad pI Calibration kit, pH 3 to 10	17047101
Low pI Calibration kit, pH 2.5 to 6.5	17047201
High pI Calibration kit, pH 5 to 10.5	17047301
<i>PlusOne chemicals</i>	
Urea	17131901
Glycerol 87% (w/w)	17132501
2-Mercaptoethanol	17131701

Designation	Product code
Dithiothreitol (DTT)	17131801



cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

Immobiline, MultiPhor, MultiTemp, NAP, Pharmalyte, PhastGel and Sephadex are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Coomassie is a trademark of Thermo Fisher Scientific. GE is a trademark of General Electric Company. TRITON is a trademark of Union Carbide Chemicals and Plastics Company Inc.

All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

71703001 AF V:6 12/2020