



Ready-to-use fluid management solutions for chromatography systems

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY13363-12May20-AN

Ready-to-use fluid management solutions for chromatography systems

A selection of fluid management components from GE Healthcare's ReadyToProcess platform of disposable and single-use equipment was used to arrange buffer and sample management solutions for four typical large-scale chromatography setups. These setups comprised different combinations of ReadyToProcess and AxiChrom™ columns and ÄKTA™ ready and ÄKTAprocess™ systems.

Ready-to-use plastic bags, plastic tubing, connectors and a mobile processing station provided flexible liquid-handling solutions that will simplify workloads and boost efficiency for all four configurations. The increasingly wide range of 'plug-and-play' bioprocessing components now available should thus find use in pilot and production facilities that utilize both disposable as well as fixed equipment.

Introduction

Interest in disposables and single-use bioprocessing equipment has increased over recent years. Today, new-build pilot and production facilities are often based on a flexible platform where unit operations occupy only a small footprint. To minimize space, for example, ready-to-use plastic bags and flexible plastic tubing increasingly replace fixed stainless steel buffer tanks and piping.

GE Healthcare's broad ReadyToProcess product platform brings such 'plug-and-play' options to biopharmaceutical manufacturing. Many disposable components are designed for single-use only while others can be re-used. All, however, are prepared for immediate use.

This Application note describes solutions for handling buffers and sample in a process setup for ÄKTA ready, a liquid chromatography system that operates with ready-to-use, disposable flow paths, as well as for ÄKTAprocess, a hard-piped chromatography system. Both systems operate



Fig 1. ReadyToProcess column connected to an ÄKTA ready chromatography system. Buffer and sample management is by a combination of ready-to-use liquid handling components such as ReadyCircuit™ bags and tubing, ReadyMate™ connectors, and ReadyKart mobile processing stations.

with disposable prepacked ReadyToProcess columns or conventional columns such as the AxiChrom range that offers Intelligent Packing with verified, preprogrammed packing methods.

Process using ÄKTA ready system System description

ÄKTA ready is a liquid chromatography system for process scale-up and production for early clinical phases. The system has a design, functionality and documentation to support its use in GLP and cGMP environments. As it operates with ready-to-use, disposable flow paths, the need for cleaning between unit operations/products/batches is eliminated and no development and validation of cleaning procedures is required. Using ÄKTA ready system with ReadyToProcess columns avoids the risk for cross-contamination.

Liquid-handling setup

The scenario illustrated in Figure 1 is a chromatographic purification process using a 10 L ReadyToProcess column. Sample and buffers are managed using ReadyCircuit, which are bags and tubing assemblies that quickly form sterile, self-contained purification modules.



Each bag is connected to the ÄKTA ready system inlet/outlets through a ReadyMate DAC 750 mini TC via a jumper with ReadyMate connectors at each end (Fig 2). Bag sizes for each solution and connection components are listed in Table 1. Figure 3 illustrates the complete operating setup.




Fig 2. ReadyCircuit 20 L bag and ReadyMate connection assembly. All four bag sizes (see Table 1) use a ReadyMate ½ inch inner diameter jumper with a 3 ft length.

Table 1. Bags and connection components used in the buffer and sample management setups

Process step	Bag description
Equilibration	100 L bag with 4 inlets
Sample	200 L bag with 4 inlets
Wash	50 L bag with 3 inlets
Elution	50 L bag with 3 inlets
Strip	20 L bag with 3 inlets
CIP	50 L bag with 3 inlets
Outlet waste	100 L bag with 4 inlets

Component	Connection
Jumper 7 pct	RMRM Jumper, 3/8 inch, C-FLEX, 3 ft
ReadyMate DAC 750 mini TC	To connect ReadyMate with mini TC
ReadyClamp	To secure ReadyMate connections




Fig 3. Complete operating setup with six bags connected to the ÄKTA ready flow kit via jumpers and ReadyMate connections.

Table 2. ReadyKart setup components used in buffer and sample management

Item/description
ReadyKart with three shelves
Top shelf assembly
ReadyKart shelf with tank hole
ReadyKart Mini frame for three shelves
ReadyKart Mini shelf no tank hole
ReadyKart 200 L Tank
ReadyKart 100 L Tank
ReadyKart 50 L Tray

As an alternative to the 10 L ReadyToProcess column shown in Figure 1, a conventional column can be used with ÄKTA ready system. In this case, the column is connected using a suitable length of tubing with mini TC connections of the same diameter as the installed system flow kit. Figure 4 shows how a conventional AxiChrom 140 column has replaced the ReadyToProcess column in the ÄKTA ready setup shown in Figure 1.



Fig 4. AxiChrom 140 column connected to an ÄKTA ready chromatography system. Buffer and sample management uses the same ready-to-use components as shown in Figure 1.

Process using ÄKTAprocess system

System description

ÄKTAprocess is an automated liquid chromatography system for process scale-up and large-scale biopharmaceutical manufacturing. Its proven stainless steel-based design can be configured to meet specific process demands. Three flow rate ranges extend up to 1800 L/h.

Liquid-handling setup

Buffer and sample management using ready-to-use bags can also be applied to a traditional hard-piped system like ÄKTAprocess. Figure 5 shows such an arrangement where the buffer assemblies for an AxiChrom 140 column are

managed in exactly the same way as for the ÄKTA ready setup shown in Figure 1. Note that ReadyKart and its component parts are also used in the same manner in the 'fully traditional' setup shown in Figure 5.



Fig 5. AxiChrom 140 column connected to an ÄKTAprocess chromatography system. Buffer and sample management use the same ready-to-use components as shown in Figure 1.

Replacing the AxiChrom 140 column with a 10 L ReadyToProcess column, thereby eliminating packing and repacking from the operating procedure, results in a buffer and sample management arrangement like the setup shown in Figure 6.



Fig 6. ReadyToProcess column connected to a traditional ÄKTAprocess chromatography system. Buffer and sample management use the same ready-to-use components as shown in Figure 1.

Conclusions

Disposables and single-use equipment are finding increasing use in downstream chromatographic processes, and are especially popular in new-build pilot and production facilities.

This Application note shows that ready-to-use plastic bags and flexible plastic tubing, plus supporting liquid handling systems, are sufficiently flexible to work effectively not only with the ReadyToProcess product platform but also with more traditional columns and 'hard-piped' systems for biopharmaceutical manufacturing.

Ordering information

Product	Code number
20 L bag with 3 inlets	12-4102-24
50 L bag with 3 inlets	12-4102-28
100 L bag with 4 inlets	12-4102-06
200 L bag with 4 inlets	12-4102-08
RMRM Jumper 0.375 in x 3 ft C-FLEX, 1 pct	12-4101-18
ReadyMate DAC 750 mini TC	28-9366-95
ReadyClamp	28-9366-90

Product	Code number
ReadyKart with three shelves, 2 pct	28-9778-94
ReadyKart top shelf assembly	28-9778-99
ReadyKart shelf with tank hole, 1 pct	28-9779-09
ReadyKart Mini frame for three shelves	28-9778-98
ReadyKart Mini shelf no tank hole	28-9779-10
ReadyKart 200 L Tank	28-9779-03
ReadyKart 100 L Tank	28-9779-04
ReadyKart 50 L Tray, 5 pct	28-9779-07

For local office contact information, visit
www.gelifesciences.com/contact

www.gelifesciences.com/readytoprocess

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden



imagination at work

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

AxiChrom, ReadyToProcess, ReadyCircuit, ReadyMate, ÅKTA, and ÅKTAprocess are trademarks of GE Healthcare companies.

ReadyMate is covered by US patent number 6,679,529 B2 owned by Johnson & Boleby Holdings, LLC and licensed to GE Healthcare companies.

The use of ÅKTA ready High Flow Kit and Low Flow Kit may require a license under the following US patents and their equivalents in other countries: US patent numbers 6,712,963 and 7,052,603 (assignee: SciLog, Inc.). A free, non-transferable license limited to using these products under above mentioned patents accompanies the purchase of the product from any GE Healthcare company and its licensed distributors.

© 2011 General Electric Company — All rights reserved.
First published Oct. 2011

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
UK

GE Healthcare Europe, GmbH
Munzinger Strasse 5
D-79111 Freiburg
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327
USA

GE Healthcare Japan Corporation
Sanken Bldg., 3-25-1, Hyakunincho
Shinjuku-ku, Tokyo 169-0073
Japan