

# VACU-GUARD vacuum protection filters

## FILTRATION

Protecting your equipment from harmful contaminants is more important than ever. For maximum protection of equipment, choose the Whatman™ VACU-GUARD family of products. These easy-to-use in-line filter devices help to remove particulate contamination from vacuum lines and protect equipment.

## VACU-GUARD

### Features and benefits

- **Excellent performance:** prevents fluid and aerosol contamination of vacuum pumps or aspiration suction systems while eliminating hazardous exhaust
- **Flexible:** designed for use with 6–10 or 10–12 mm ID tubing
- **Biosafe:** all materials pass USP Class VI Test for Plastics

### Applications

- Protects vacuum pumps and systems from aerosols and particulate contamination

## VACU-GUARD 150

### Features and benefits

- **Choice of media:** VACU-GUARD 150 capsule filters include all the features and benefits of standard VACU-GUARD disc filters, plus a range of media for specific applications
- **Added back-up protection:** use as a backup between a cold trap and pump to protect against moisture and organic vapors if cold trap fails

### Applications

- Activated carbon removes organic vapors from air
- Molecular sieve for removal of water and small organic and alkaline molecules from air streams
- Desiccant for use with high velocity acidic air



**Fig 1.** In-line disc filters protect vacuum systems from aqueous aerosols.

**Table 1.** Typical data VACU-GUARD In-Line Disc Filter — 50 and 60 mm

	50 mm	60 mm
Filtration area	16 cm <sup>2</sup>	25 cm <sup>2</sup>
Maximum operating pressure	1 bar (15 psi)	1 bar (15 psi)
Biosafety	All materials pass USP Class VI test for plastics	
Rated retention in air	99.99% particle retention for particles ≥ 0.1 μm	
Pore size (in liquid)	0.45 mm	0.45 mm
Housing	Polypropylene	Polypropylene
Filtration media	PTFE membrane	PTFE membrane
Connectors	1/4–3/8" (6–10 mm) SB (stepped barb) inlet and outlet	3/8–1/2" (10–12 mm) SB inlet and outlet
Flow rates (SLPM):		
2 psi (0.14 bar)*	15	27
4 psi (0.28 bar)*	27	57
6 psi (0.41 bar)*	38	83
10 psi (0.69 bar)*	53	139
Flow direction	Inlet to outlet	Inlet to outlet

\* Differential pressure

**Table 2.** Typical data VACU-GUARD 150 In-Line Capsule Filter

	<b>Activated carbon</b>	<b>Desiccant</b>	<b>Molecular sieve</b>
Chemical trap media	Activated carbon	Anhydrous calcium sulphate	Silico aluminate zeolite
Filter media	PTFE	PTFE	PTFE
Surface area or weight (nominal)	82 000 m <sup>2</sup> (carbon)	318 g (desiccant)	363 g (zeolite)
Flow rates (SLPM) (nominal):			
0.1 bar (1.45 psi)*	210	280	250
0.5 bar (7.25 psi)*	450	600	570
Maximum operating pressure:			
dry gas	4 bar (60 psi)	4 bar (60 psi)	4 bar (60 psi)
wet gas	1 bar (14 psi)	1 bar (14 psi)	1 bar (14 psi)
Connectors:			
inlet	Hose barb for 1/2" (12.7 mm) tube		
outlet	3/8-1/2" (10-12 mm) step barb		

\* Differential pressure

Note: As with any chemical reaction, care should be used to determine the safety and usefulness of VACU-GUARD 150 products prior to routine use. For example, the molecular sieve rapidly heats up when exposed to water.

## Ordering information

<b>Product</b>	<b>Quantity</b>	<b>Product code</b>
VACU-GUARD, 50 mm disc	10	6722-5000
VACU-GUARD, 60 mm disc	10	6722-5001
VACU-GUARD 150, activated carbon	1	6722-1001
VACU-GUARD 150, desiccant	1	6722-1002
VACU-GUARD 150, molecular sieve	1	6722-1003



**Fig 2.** In-line capsule filters trap chemicals in addition to aqueous aerosols.

[cytiva.com/whatman](http://cytiva.com/whatman)

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. Whatman is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](http://cytiva.com/contact)

CY17819-18Dec20-DF

