


Biacore™ sensor surface



Contents

01	Introduction	4			
	1.1 Principles of Biacore™ SPR systems	5			
	1.2 Biacore™ SPR terminology	5			
	1.3 Components of Biacore™ SPR systems	7			
	1.3.1 The SPR detection system	7			
	1.3.2 The sensor surface	8			
	1.3.3 The flow system	9			
	1.4 Scope of this handbook	9			
	1.5 Additional information	9			
02	Sensor chip overview	10			
	2.1 The sensor surface	11			
	2.1.1 The gold layer	11			
	2.1.2 The surface coating	12			
	2.2 Sensor chips and reagent kits	12			
	2.2.1 General purpose sensor chips	13			
	2.2.2 Sensor chips and reagent kits for antibody capture	13			
	2.2.3 Sensor chips for tagged ligand capture	14			
	2.2.4 Sensor chips for lipid attachment	14			
	2.2.5 Customizable sensor chips	14			
	2.3 Reference	14			
03	Surface preparation strategies	15			
	3.1 Covalent immobilization	16			
	3.2 Capture approaches	17			
	3.2.1 Streptavidin-biotin capture	19			
	3.2.2 NeutrAvidin-biotin capture	19			
	3.2.3 Reversible capture of biotinylated ligands	19			
	3.2.4 Antibody-based capture	20			
	3.2.5 Sensor chips and kits for antibody capture	21			
	3.2.6 Capturing tagged recombinant proteins	21			
	3.3 Hydrophobic attachment	22			
	3.3.1 Monolayer attachment	22			
	3.3.2 Bilayer attachment	22			
	3.4 Immobilizing small molecules	23			
	3.5 Techniques for stabilizing the surface	23			
	3.6 Strategies for surface preparation	24			
	3.6.1 How much ligand should be immobilized?	24			
	3.6.2 Choice of sensor surface	25			
	3.6.3 Choice of ligand attachment and capture method	25			
04	Covalent immobilization chemistry	26			
	4.1 Conditions for immobilization	27			
	4.1.1 Temperature	27			
	4.1.2 Buffer pH and ligand preconcentration	27			
	4.1.3 Other buffer conditions	28			
	4.1.4 Ligand concentration	28			
	4.1.5 Immobilization pH scouting	28			
	4.2 General immobilization procedures	29			
	4.2.1 Preparing solutions	29			
	4.2.2 Evaluating covalent immobilization	30			
	4.2.3 Testing the surface after immobilization	30			
	4.3 Activating the surface	31			
	4.4 Amine coupling	31			
	4.4.1 Chemistry	31			
	4.4.2 Required solutions	32			
	4.4.3 Typical immobilization protocol	32			
	4.4.4 Sensorgram	32			
	4.5 Reversed amine coupling	33			
	4.6 Micelle-mediated coupling of negatively charged ligands	33			
	4.7 Thiol coupling	34			
	4.7.1 Ligand thiol coupling	34			
	4.7.2 Surface thiol coupling	35			
	4.7.3 Maleimide coupling	38			
	4.8 Aldehyde coupling	41			
	4.8.1 Chemistry	41			
	4.8.2 Oxidizing the ligand	42			
	4.8.3 Required solutions	42			
	4.8.4 Typical immobilization protocol	42			
	4.8.5 Sensorgram	43			
	4.9 Ligand modification for capturing	43			
	4.9.1 Biotinylation for streptavidin-biotin capture	43			
	4.10 References	44			
05	Hydrophobic adsorption methods	45			
	5.1 Monolayer adsorption on Sensor Chip HPA	46			
	5.1.1 Preparation for use	46			
	5.1.2 Liposome preparation and adsorption	47			
	5.2 Liposome attachment on Sensor Chip L1	48			
	5.2.1 Preparations for use	48			
	5.2.2 Liposome preparation and adsorption	48			
	5.3 Liposome attachment on other sensor chips	49			
	5.4 On-surface reconstitution	49			
	5.4.1 Preparation of ligand and micelles	50			
	5.				

06	Regeneration	53
	6.1 Regeneration scouting strategy	54
	6.1.1 General considerations	54
	6.1.2 Choice of regeneration solution	55
	6.1.3 Factors affecting regeneration	55
	6.1.4 Regeneration scouting procedure	56
	6.2 Interpreting regeneration scouting	57
	6.2.1 Interpretation principles	57
	6.3 Verifying regeneration conditions	59
	6.4 Suggested regeneration solutions	59
	6.4.1 Protein surfaces	59
	6.4.2 Nucleic acids	60
	6.4.3 Serum and plasma	60
	6.4.4 Membrane-associated ligands	60
	6.4.5 Small molecules	61
	6.4.6 Other regeneration strategies	61
07	Sensor surface performance	62
	7.1 Testing the surface	63
	7.2 Troubleshooting surface preparation	64
	7.2.1 Low immobilization levels	64
	7.2.2 Low analyte responses	65
	7.2.3 When regeneration fails	65
	7.2.4 Problems with nonspecific binding	66
08	Storing and reusing sensor chips	68
	8.1 Recommended procedure for wet storage	69
	8.2 Recommended procedure for dry storage	70

01

Introduction




Fig 1.2. Schematic illustration of a sensorgram. The bars below the sensorgram curve indicate the solutions that pass over the sensor surface.

- **Regeneration** is the process of removing bound analyte from the surface after an analysis cycle without damaging the ligand, in preparation for a new cycle.
- Response is measured in **resonance units (RU)**. The response is directly proportional to the concentration of biomolecules on the surface.
- A **sensorgram** is a plot of response in RU against time (Fig 1.2), showing the progress of the interaction. This curve is displayed directly on the computer screen during an analysis. Sensorgrams may be analyzed to provide you with information on the rates of the interaction.
- In many assay situations, sample passes over two or more flow cells in series. The first flow cell in the series is usually blank and serves as a reference, while ligand is attached in the other flow cell(s). Surfaces with ligand are referred to as **active**. Blank surfaces used for reference purposes are referred to as **reference**.
- A particular sensorgram is referred to as a **curve** in several contexts in the software. This terminology is used to distinguish between different classes of sensorgram that recur within a run. For example, measurements on one active and one reference surface can generate separate curves for each of the two flow cells and a third **reference-subtracted curve** (active minus reference).
- A **report point** records the response on a sensorgram at a specific time averaged over a short time window, as well as the slope of the sensorgram over the window. The response may be absolute (above a fixed zero level determined by the detector) or relative to the response at another specified report point of your choice.

<

1.3.2 The sensor surface

The sensor chip is at the heart of Biacore™ SPR systems providing the physical conditions necessary to generate the SPR signal. The interaction that you want to study occurs on the surface of the sensor chip (Fig 1.4).




Fig 1.4. Illustration of a solution containing the analyte flowing over the sensor chip and interacting with the ligand.

The specificity of a Biacore™ analysis is determined through the nature and properties of the molecule attached to the sensor surface. You may attach biomolecules to the sensor chip surface using three different approaches (Fig 1.5):

- Covalent immobilization, where the molecule is attached to the surface through a covalent bond linkage.
- High affinity capture, where the molecule of interest is attached by noncovalent interaction with another molecule (which in turn is usually attached using covalent immobilization).
- Hydrophobic adsorption, which exploits hydrophobic interactions to attach either the molecule of interest or a hydrophobic carrier such as a lipid monolayer or bilayer to the sensor chip surface.

These approaches are discussed in more detail in Chapters 3 through 5.




Fig 1.5. Three approaches for attaching biomolecules to the sensor chip surface.

1.3.

02

Sensor chip overview

03


Surface preparation strategies

Molecules can be attached to the sensor chip surface using several techniques, including covalent immobilization using a variety of chemical methods, high affinity binding to a specific capturing molecule, and adsorption of lipid mono- or bilayers.

This chapter provides an overview of the methods available for attaching the ligand to the sensor chip surface. Detailed procedures are described in the next chapter.

3.

<



05

