


## Procedure

# GFP V<sub>H</sub>H capture surface for Biacore assays

The GFP V<sub>H</sub>H/nanobody is a single domain antibody derived from alpaca that selectively binds to green fluorescent protein (GFP)-tagged proteins. This nanobody is used for specific and stable capture of green fluorescent protein (GFP)-tagged proteins in Biacore™ assays (see Fig 1). GFP V<sub>H</sub>H also recognizes proteins tagged with enhanced GFP (eGFP) and yellow fluorescent protein derivatives (YFP, eYFP, and Venus).

The GFP V<sub>H</sub>H capture surface allows for easy and effective isolation and directed immobilization of fluorescent-tagged proteins from complex biological samples like prokaryotic or eukaryotic cell extracts. Immobilization of these proteins enables binding studies and kinetic analysis using Biacore systems.



**Fig 1.** (A) The principle behind immobilization of ligands using GFP V<sub>H</sub>H in a Biacore capture assay format. (B) A typical sensorgram of GFP V<sub>H</sub>H immobilization by amine coupling.

## Procedure

This procedure describes how to attach ligand, run the interaction analysis, and regenerate a GFP V<sub>H</sub>H surface.

### Attachment of GFP V<sub>H</sub>H

GFP V<sub>H</sub>H is attached on both active and reference surfaces of the following Biacore sensor chips: Sensor Chip CM5 (alternatively Sensor Chip CM4, Sensor Chip CM3, or Sensor Chip C1) and Series S sensor chips of the same series.

The conditions shown below should allow an attachment level of ~ 3000 RU on a Sensor Chip CM5 surface, which is close to saturation.

- Coupling chemistry: amine coupling
- Dilute GFP V<sub>H</sub>H to 50 µg/mL in 10 mM acetate pH 5.5 immobilization buffer
- Activation time: 7 min, ligand contact time: 7 min, deactivation time: 7 min
- Flow rate: 10 µL/min

### Capture of GFP-tagged protein

You can calculate how much of the GFP-tagged protein (ligand) you need to capture on the sensor surface to reach a theoretical R<sub>max</sub> (maximum binding capacity for analyte) of 50–100 RU with your respective analyte using the following equation:

$$R_L = \frac{\text{ligand MW}}{\text{analyte MW}} \times \frac{R_{\max}}{S_m}$$

R<sub>L</sub> = Attachment or capture level  
S<sub>m</sub> = Stoichiometric ratio  
R<sub>max</sub> (RU) = Maximal binding response

1. Test optimal conditions to reach the respective capture level for your GFP-tagged protein, for example, in a manual run addressing the active flow cell only.
2. Start with a picomolar (pM) dilution of the GFP-tagged protein, inject at 10 µL/min for 60 s.
3. If the concentration of the GFP-tagged protein in the extract is not known, set up a series of 10-fold dilutions in running buffer. Start testing with the lowest concentration.

## Investigate binding of analyte to GFP-tagged protein

1. Capture the desired amount of GFP-tagged protein on the active flow cell.
2. Inject dilutions of the analyte over both reference (GFP V<sub>H</sub>H) and active (GFP-tagged protein captured by GFP V<sub>H</sub>H) flow cells.
3. Start with a pM dilution of analyte in running buffer.
4. Set up a series of 10-fold dilutions in running buffer. Start testing with the lowest concentration.
5. Control reference flow cell for potential nonspecific binding of analyte to GFP V<sub>H</sub>H.

## Remove GFP-tagged protein from the GFP V<sub>H</sub>H surface

Applying the conditions below regenerates the GFP V<sub>H</sub>H surface and removes bound GFP-tagged protein. Apply regeneration over reference and active flow cells:

- Regeneration solution: 10 mM glycine pH 1.7
- Contact time: 2 × 30 s
- Flow rate: 30 µL/min

## References

1. Della Pia, E.A. and Martinez, K.L. Single domain antibodies as a powerful tool for high quality surface plasmon resonance studies. PLoS ONE **10(3)** (2015): e0124303. <https://doi.org/10.1371/journal.pone.0124303>

## cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. Biacore is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

CY17752-08Dec20-PD

