



Use of GD/XP syringe filters for accurate analysis of metals

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit cytiva.com/contact

CY14242-08Jun20-AN



Use of GD/XP syringe filters for accurate analysis of metals

Introduction

In order to be able to determine the inorganic ion content (e.g., metals) of a sample as accurately as possible, it is important to avoid contamination during sample preparation such as a filtration step. Whatman™ GD/XP syringe filters (Fig 1) have been specifically designed to minimize inorganic extractables.



Fig 1. GD/XP syringe filter.

GD/XP syringe filters comprise a two-layer prefilter stack consisting of 20 µm and 5 µm polypropylene filters and a final membrane filtration layer positioned below the prefilter stack (Fig 2). This multilayer construction increases the volume of liquid that can be filtered compared to a filter with a single membrane layer.

Whatman GD/XP 0.45 µm PP filters (cat. number 6993-2504) were analyzed to determine the quantity of metals that might be released from the filters. The method involved flushing the filters with an acid solution and then analyzing the acid solution for metals by ICP-OES. These conditions should provide a realistic view of the quantity of metals leached during filtration prior to metal analysis.




Fig 2. GD/XP syringe filters contain three filtration layers and this leads to an increase in the volume of liquid that can be filtered. Layer 1 is a 20 µm polypropylene filter, layer 2 is a 5 µm polypropylene filter, and layer 3 is a 0.45 µm filter (such as polypropylene, PTFE etc.).

Materials and methods

Three lots of GD/XP syringe filters (catalog number 6993-2504: 0.45 µm Polypropylene) were used in this study. Lot numbers used were O684, L925, and W466. The tests were performed at GE China Technology Center in Shanghai, China.

Whatman GD/XP syringe filters were washed with 10 ml of HNO₃/HCl (1:2). This solution was evaporated at 90°C until almost dry. The residue was digested in 1 ml of HNO₃/HCl (1:2) until the residue completely dissolved. Each sample was diluted to 15 ml with deionized water. Blanks were prepared in an identical manner to the samples. Each lot of GD/XP filters were tested in duplicate. ICP-OES (Inductively Coupled Plasma Optical Emission) was used for determining the metal content of the samples.

ICP-OES conditions

RF power:	1400 W
Coolant flow:	13.00 L/min
Auxiliary flow:	0.80 L/min
Nebulizer flow:	0.80 L/min
Nebulizer:	Cross flow
Spray chamber:	Quartz scott type
Plasma torch:	Quartz standard

Discussion and conclusion

Critical metals such as Cd, As, and Pb were below the detection limits. Zn, Fe, Ba, Al, Cu, and Sr could be detected but they were present at low levels. Consider the following scenario: assume a GD/XP filter is used to filter a 10 ml solution. The filter will introduce 2.12 µg of Zn. Thus the back ground concentration of Zn due to the filter would be 0.21 ppm. In fact, Zn had the highest concentration. Hence the background levels for other metals that were detected would be lower. Therefore we may conclude that GD/XP, 0.45 µm PP syringe filters are suitable for the filtration of samples prior to metal analysis.

Results

Table 1. GD/XP metal content

GD/XP Lot #	GD/XP metal content (µg per filter)						Ave.	Std dev.
	L925 (Filter 1)	L925 (Filter 2)	O684 (Filter 1)	O684 (Filter 2)	W466 (Filter 1)	W466 (Filter 2)		
Metal								
As	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	-
Cd	< 0.005	< 0.005	< 0.005	< 0.005	< 0.005	< 0.005	< 0.005	-
Cr	< 0.01	< 0.01	< 0.01	0.03	< 0.01	< 0.01	< 0.01	-
Cu	0.04	0.04	0.04	0.05	0.03	0.05	0.04	0.01
Mo	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	-
Ni	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	-
Pb	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	-
V	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	< 0.01	-
Zn	1.33	1.58	1.56	1.85	3.51	2.94	2.12	0.88
Fe	4.35	1.27	0.68	2.29	1.12	0.71	1.74	1.41
Ba	< 0.01	< 0.01	< 0.01	< 0.01	0.03	0.02	0.02	-
Al	0.95	0.53	0.3	0.79	0.45	0.36	0.56	0.26
Sr	0.06	0.03	0.02	0.04	0.02	0.02	0.03	0.02

< = below detection limit

For local office contact information, visit

www.gelifesciences.com/contact

www.gelifesciences.com/LaboratoryFiltration

GE Healthcare UK Limited

Amersham Place

Little Chalfont

Buckinghamshire, HP7 9NA

UK

GE and GE monogram are trademarks of General Electric Company.

Whatman is a trademark of General Electric Company or one of its subsidiaries.

© 2014 General Electric Company – All rights reserved.

First published Sep. 2014.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

GE Healthcare Europe, GmbH
Munzinger Strasse 5
D-79111 Freiburg
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327
USA

GE Healthcare Japan Corporation
Sanken Bldg., 3-25-1, Hyakunincho
Shinjuku-ku, Tokyo 169-0073
Japan