

# UNICORN 6.3.2 Service Pack 2

## Installation instruction

**Note:** *UNICORN 6.3.2 or UNICORN 6.3.2 SP1 must already be installed prior to installing Service Pack 2. For comprehensive instructions on how to install UNICORN 6.3.2 please refer to the UNICORN 6.3 Administration and Technical Manual.*

Follow the instructions below to download and install UNICORN 6.3.2 Service Pack 2.

Step	Action
1	Go to the website <a href="http://www.gelifesciences.com/unicorn">http://www.gelifesciences.com/unicorn</a> .
2	At the bottom of the main panel, click <b>Previous UNICORN Versions</b> . <i>Result:</i> A new main panel opens.
3	Click <b>UNICORN 6</b> , then click <b>UNICORN 6.3.2</b> .
4	Click <b>Related Documents</b> tab in the middle of the page.
5	Scroll down to <b>Installation Guide</b> .
6	Click the link <b>UNICORN 6.3.2 Service Pack 2 software and Installation Instruction</b> . Download the zip file to a suitable folder on your computer.
7	Make sure UNICORN is not running.
8	Locate the downloaded zip file on your computer.
9	On the <b>File</b> menu, click <b>Extract All....</b> Choose a suitable folder to save the extracted files. <i>Result:</i> The software file <b>UNICORN 6.3.2 SP2.exe</b> and UNICORN 6.3.2 Service Pack 2 installation instruction <b>29138154AA.pdf</b> are saved on your computer.
10	Double-click the <b>UNICORN 6.3.2 SP2.exe</b> to start the installation. Follow the instructions on the screen to install the service pack.



Step	Action
11	When the installation has finished, verify that UNICORN 6.3.2 Service Pack 2 has been correctly installed.
a	Start UNICORN. Make sure the <b>Administration</b> module is selected.
b	Click the <b>Administration</b> module icon on the taskbar.
c	On the <b>Administration</b> module header toolbar, select <b>Help:About UNICORN</b> . A pop-up window opens.
d	Verify that the software version is shown as <b>UNICORN 6.3.2 SP2</b> .

## Implement changes

Follow the instructions below to implement the UNICORN Windows authentication corrections on your computer.

Step	Action
1	Make sure you are logged in as a domain user with local administrator rights on the computer running the UNICORN application.
2	Remove all Windows® users, who are not able to log in, from the UNICORN access group.
3	Add the removed Windows users back to the UNICORN access group.

## More information

See *UNICORN 6.3 Administration and Technical Manual* for more information about software installations.

For local office contact information,  
visit

[www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

GE Healthcare Bio-Sciences AB

Björkgatan 30  
751 84 Uppsala  
Sweden

[www.gelifesciences.com/unicorn](http://www.gelifesciences.com/unicorn)

GE and GE monogram are trademarks of General Electric Company.  
UNICORN is a trademark of General Electric Company or one of its subsidiaries.

Windows is a registered trademark of Microsoft Corporation.

All other third party trademarks are the property of their respective owner.

Any use of UNICORN is subject to GE Healthcare Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

UNICORN 6 © 2009-2014 General Electric Company.

© 2014 General Electric Company – All rights reserved.

First published Dec. 2014

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Limited  
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Japan Corporation  
Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan

