

HyClone media and supplements

Reconstitution protocol

ActiCHO P Powder media

This protocol provides instructions for reconstituting ActiCHO™ P from powder. The instructions apply for all filling sizes.

Content

HyClone™ ActiCHO P Powder is a chemically-defined, powdered cell culture medium containing poloxamer 188, without insulin or L-glutamine. It is hydrolysate-free and animal-derived component-free (ADCF). There are no proteins, peptides or growth factors in ActiCHO P.

Storage and shelf life

Store ActiCHO P Powder dry and protected from light at 2°C to 8°C.

Refer to the product label for the expiry date.

Required materials

- Mixing vessel
- Stirrer

Note: *A magnetic stirrer bar can be used for small-scale reconstitutions (up to 5 L). An overhead or bottom-mounted impeller is recommended for larger volumes.*

- Calibrated pH meter
- Calibrated osmometer
- Highly purified water, such as WFI (water for injection)

Preparation of liquid medium

Step	Action				
1	<p>Fill a clean mixing vessel to 96% of the final volume with high quality purified water, such as WFI (water for injection) at ambient temperature (18°C to 25°C) and start stirring. For example, to prepare a 1 L solution of ActiCHO P take approximately 960 mL of water.</p> <p>Note:</p> <p><i>The stirring should be vigorous enough to quickly draw the powder below the liquid surface, but not too fast to cause formation of air bubbles and excessive foaming.</i></p>				
2	<p>Add 22.36 g/L ActiCHO P Powder slowly to the vessel, avoiding formation of clumps. Stir for 30 minutes.</p> <p>Result:</p> <p>The solution will remain cloudy in this step, but should be clear of any clumps or dry powder residues.</p>				
3	<p>Slowly add 6.5 mL/L of a 5 N NaOH solution or 3.25 mL/L of a 10 N NaOH solution and continue stirring for 60 minutes.</p> <p>Result:</p> <p>After this step the solution will be clear.</p>				
4	<p>Add 1.8 g/L sodium bicarbonate (NaHCO₃) into the vessel and continue stirring for 10 minutes.</p>				
5	<p>Adjust the pH to between 6.90 and 7.35 by dropwise addition of 5 N or 10 N NaOH. After adjusting, continue stirring for an additional 60 minutes to ensure that all components are completely dissolved.</p> <p>Note:</p> <p><i>The pH will gradually increase with longer mixing times.</i></p>				
6	<p>Adjust to the final volume with high quality purified water, such as WFI (water for injection) and stir for an additional 10 minutes.</p>				
7	<p>Measure and record the final pH and osmolality.</p> <p>Result:</p> <table><tr><td>pH</td><td>6.90 to 7.55</td></tr><tr><td>Osmolality</td><td>285 to 330 mOsmol/kg</td></tr></table>	pH	6.90 to 7.55	Osmolality	285 to 330 mOsmol/kg
pH	6.90 to 7.55				
Osmolality	285 to 330 mOsmol/kg				

Step	Action
8	<p>Sterilize immediately by membrane filtration. Use a low-binding filter membrane type, such as PVDF, PES or cellulose acetate.</p> <p><i>Result:</i></p> <p>The medium is a clear yellow/brown liquid.</p>
9	<p>Store the medium protected from light at 2°C to 8°C until use.</p>



cytiva.com/bioprocess

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

ActiCHO and HyClone are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

All other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

29136548 AC V:6 02/2021