Amersham™ QuickStain

ELECTROPHORESIS REAGENTS

The Amersham™ QuickStain Kit contains a $Cy^{TM}5$ fluorophore and labeling buffer for easy detection of proteins in SDS-PAGE gels and on Western blotting membranes. The ready-to-use $Cy^{TM}5$ N-hydroxysuccinimide (NHS) ester and Tris labeling buffer ensure robust and consistent labeling for detection of proteins in diverse samples.

Key benefits:

- Ease of use: Ready-to-use reagents for fast and convenient analysis of proteins following one-dimensional (1D) electrophoresis and/or Western blotting
- Fast protocols: Choose between a quick five minute protocol and a standard 30 min quantitative protocol
- **Wide sample range:** Stain complex and purified samples with protein concentrations ranging from 1 µg/mL to 20 mg/mL, without the need to pre-measure protein concentration
- Sensitivity: As good as silver staining, sub-ng sensitivity
- Broad dynamic range: With a fluorescence imager or scanner, it is possible to simultaneously detect main bands and impurity levels down to parts per thousand (0.1%)
- Robustness: A 10-fold dilution of samples with the labeling buffer, which has a high buffering capacity, ensures reproducible labeling at optimal pH for a wide variety of samples
- Versatility: Most materials exhibit low Cy[™]5 fluorescence background, which makes the kit compatible with many common gel formats and membranes

With the easy Cy™5 prelabeling of Amersham™ QuickStain, there is no need to do post-staining (e.g., Coomassie™ or silver staining) of gels or Ponceau S staining of membranes (see Fig 1). The kit is compatible with common gel-buffer electrophoresis systems, including Tris-acetate and Tris-glycine buffer systems. Many plastic materials exhibit a low fluorescence background in the Cy™5 spectral region, which makes Cy™5 staining well suited for precise quantitation of proteins using a wide variety





Fig 1. Amersham[™] QuickStain is designed for easy detection of proteins after SDS-PAGE and/or after Western blotting. Simply add the provided labeling buffer and Cy[™]5 to the sample prior to electrophoresis.

of gel formats (i.e., precast, homecast, cassettes, backed, and nonbacked). Using CCD imagers and laser scanners from Cytiva, proteins can be detected with high sensitivity and a broad linear dynamic range. As a result, Amersham™ QuickStain can be used for a wide variety of protein analysis applications, from quick visualization of protein content in different samples to accurate quantitation of target proteins using calibration curves.

Traditional staining methods often limit the dynamic range, typically to one or two orders of magnitude for silver and Coomassie™ staining. In contrast, fluorescence detection has a much broader detection window. The protocols used with the Amersham™ QuickStain kit have been optimized to take advantage of this broad dynamic range. As shown in Figure 2, Amersham™ QuickStain gives substantially better signal-to-noise and signal-to-background ratios compared with Coomassie™ staining.








Fig 2. Comparison of gel images and lane profiles using (A) Amersham™ QuickStain and (B) Coomassie™ staining. The two gels were loaded with the same quantities of six proteins per well, as indicated above each lane. Novex™Tris-glycine gels were used, and gel images were generated using Typhoon™ FLA 9500 scanner (QuickStain) and ImageQuant™ LAS 500 imager (Coomassie™). The lane profiles in (C) show the intensity of bands in lane 3 for each image.

Labeling protocols for SDS-PAGE

The Amersham™ QuickStain kit is based on Cy™5 NHS ester labeling of reactive amino groups in proteins. To ensure minimal interference from sample components and to maintain an optimal pH for labeling, the prelabeling protocol starts with a 10-fold dilution of the sample using the Amersham™ QuickStain labeling buffer. There is no need to remove excess dye, because free dye migrates faster than the proteins of interest and therefore ends up at the bottom of the gel. The compositions of the two kit components are provided in Table 1.

Table 1. Amersham™ QuickStain components

| Су™5 | Four vials of Cy™5 NHS ester in anhydrous dimethyl sulfoxide (DMSO) each vial contains 35 microliters |
|-----------------|---|
| Labeling Buffer | Three vials of Tris-HCl labeling buffer with 0.1% sodium dodecyl each vial contains 700 microliters |

Quick protocol - for qualitative analysis

- 1. Dilute sample 10-fold with labeling buffer
- 2. Add Cy™5 to the sample
- 3. Heat sample at 95°C for 3-5 min
- 4. Add sample loading buffer containing lysine*
- 5. Heat sample at 95°C for 3 min
- 6. Load sample on gel

Total time less than 15 min

Quantitative protocol

- 1. Dilute sample 10-fold with labeling buffer
- Add Cy[™]5 to the sample
- 3. Incubate sample for 30 min at room temperature
- 4. Add sample loading buffer containing lysine*
- 5. Heat sample at 95°C for 3 min
- 6. Load sample on gel


Total time less than 40 min

* Lysine quenches the labeling reaction. Add DTT to loading buffer if running reducing SDS-PAGE.

Note: In Western blotting applications we recommend performing the labeling of cell lysates in the sample lysis buffer. Doing so will lead to a lower degree of labeling and a lower signal that is suitable for detection on membranes.

Fixed amount of dye in the labeling reaction

Keeping the amount of dye fixed in the labeling reaction and varying the protein concentration (Fig 3) leads to linear dose-response curves and well-resolved protein bands on the gel (1). For accurate quantitation, label calibrants under conditions identical to those used for the samples, and make a calibration curve.



Samples with different protein concentration (c)

Fig 3. The Amersham[™] QuickStain protocols are based on adding a fixed amount of dye to samples with different protein concentrations. For a 20 μ L reaction volume, add 1 μ L of Cy[™]5.

Cytiva imagers and scanners for detection of Amersham QuickStain

Both CCD imagers and laser scanners can be used to detect and quantitate signals from Cy™5 labeled proteins. Amersham™ Imager 600RGB (Fig 4, Fig 5) is an advanced CCD imager that can be used to detect three different fluorophores in multiplexing experiments. Amersham™ brand Typhoon™ RGB and 5 systems (Fig 4) are high-end variable-mode laser scanners for the most demanding applications. Both Amersham Typhoon systems have a large working area and many laser-filter combinations. All three instruments are well suited for detection of proteins prelabeled using Amersham™ QuickStain in gels and on membranes. The high-end detectors (CCDs and photomultiplier tubes) and the 16bit depth file format enable a linear dynamic range of up to five orders of magnitude. This wide linear range makes it possible to do quantitative measurements of both strong and weak signals in the same image, which is a key benefit for many applications (e.g., purity analysis).



Fig 4. Cytiva's Amersham[™] Imager 600RGB (left) and Amersham[™] brand Typhoon[™] scanner (right) can be used in Cy[™]5 imaging mode to detect proteins prelabeled using Amersham[™] QuickStain.




Fig 5. Linearity, dynamic range, and limit of detection for proteins using Amersham[™] QuickStain with Amersham[™] Imager 600RGB. A 2-fold dilution series of Cy[™]5 labeled phosphorylase b sample shows a linear dynamic range of 3.3 orders of magnitude and a detection limit of 100 pg.


Application examples

Amersham™ QuickStain can be used to detect proteins in different types of SDS-PAGE analysis, including:

- Monitoring of bioprocessing steps, including protein quantitation with calibrants and standard curves
- Purity analysis and fraction screening in protein purification processes
- Detection of proteins on Western blot membranes for comparison to antibody-based detection of target protein

Quantitation of antibodies

Antibodies play a key role in biotechnology and therapeutics. It is important to measure the amount of antibody in a sample and to analyze its purity. Amersham™ QuickStain was used to detect antibody expressed in Chinese hamster ovary (CHO) cells over time. Samples were labeled with Cy™5 prior to SDS-PAGE. Both non-reduced and reduced samples were run on the gels to check the integrity of the antibodies. It was possible to accurately measure the amount of antibody in the samples, using a calibration curve, and to evaluate the purity (Fig 6). The quantities of antibody determined from the SDS-PAGE analysis agreed well with enzyme-linked immunosorbent assay (ELISA) measurements. In this way the production process of an antibody was monitored (2).







Fig. 6. Detection of antibody in CHO cells; protein was prelabeled using Amersham™ QuickStain. (A) Images of SDS-PAGE gels. Calibrant samples, samples collected at different times from the cell culture, and markers were loaded onto each gel. The non-reduced gel shows bands of intact antibody, as expected. Reduced samples exhibit the heavy chain band at approximately M_r 50 000 and the light chain band at approximately M_r 23 100. (B) Example of a dose-response curve for a therapeutic monoclonal antibody from a "reduced" gel image similar to the one in the top panel. This curve is for the heavy chain of trastuzumab.

Purity analysis

It is standard practice to use SDS-PAGE to measure the purity of different fractions during protein purifications. The results are used to evaluate the purification process and form the basis for fraction picking and selection of a purification strategy. Thus, it is important to assess both purity and total amount of protein in each fraction. The excellent sensitivity of Amersham™ QuickStain reveals purity with details that are not seen with Coomassie™ staining. In fact, signals from impurities as low as 0.1% of the main band's signal intensity can be detected (see Fig. 7). The detection limit is similar to the limit for silver staining, but the wide dynamic range allows for quantitation of both weak and strong signals and a more stringent determination of purity. In addition, with SDS-PAGE slab gel analysis of multiple samples loaded in adjacent lanes, it is easy to assign well-resolved bands to protein impurities. In contrast, lane profiles in chip- and capillary-based electrophoresis are more difficult to evaluate.






Fig 7. Gel images showing Amersham™ QuickStain results of different fractions from purification of an MBP-His dual-tagged protein using different immobilized affinity chromatography (IMAC) media. Images (A) and (B) are the same image with different contrast settings. Marker samples are loaded in the outer lanes. The signals from highlighted impurity bands (red box) fall within 0.1%–0.7% of the main band's signal intensity.





Fig 8. SDS-PAGE analysis of the recovery and purity of his-tagged GFP in the elution pool after a Capto™ Chelating IMAC step using different imidazole concentrations in the column equilibration and wash buffers. Red arrow indicates position of full-length his-GFP. Buffer exchange prior to labeling was performed using the Amersham™ WB MiniTrap™ Kit. Prior to electrophoresis, Cy™5 labeling was performed using Amersham™ QuickStain.


In the example in Figure 8, the purity of histidine-tagged (his-tagged) GFP in different eluates was analyzed to optimize immobilized metal ion affinity chromatography (IMAC) elution conditions. Amersham $^{\text{TM}}$ QuickStain was used to analyze both the purity and the amount of protein in the different fractions.

Total protein detection on Western blotting membranes

Amersham™ QuickStain can also replace Ponceau S staining of proteins on membranes. Multiplexing is possible using the multichannel detection capabilities of many imagers and scanners. For example, a Cy™3 labeled secondary antibody can be used to detect the target protein of interest, and Amersham™ QuickStain can be used to measure the total protein content of a sample (Fig 9). This multiplexing capability provides a way to evaluate the transfer efficiency and to adjust for loading errors using normalization (3). It is also possible to perform chemiluminescence Western blotting experiments and compare these results with the total protein content measured using Amersham™ QuickStain and a suitable CCD imager or scanner.




Amersham™ QuickStain total protein (Cy™5)



Western Blot DHFR (Cy™3)

Overlay



1 2 3 4 5 6 7 8 9

Fig 9. Multiplex detection of total protein and target protein with Amersham™ QuickStain (Cy™5), Amersham ECL Plex™ (Ab-Cy™3), and Amersham™ Imager 600RGB. DHFR (Cy™3, green) was detected in nine different samples from a growth optimization of *E coli*. Total protein in the samples was prelabeled with Cy™5 (red). The overlay image shows the DHFR band in yellow. Membrane blocking was performed with 3% BSA in PBS-Tween™; rabbit anti-DHFR C-terminal primary antibody was diluted 1:1000; ECL Plex™ Goat anti-rabbit-Cy™3 secondary antibody was diluted 1:2500.

Data storage and image analysis

Data are stored either in linear 16-bit grayscale TIFF file format when using a CCD camera (.tif files) or in square root-encoded 16-bit TIFF file formats when using a laser scanner (.gel files). For image analysis we recommend either Amersham™ Imager 600 analysis software or ImageQuant™ TL software. ImageQuant™ TL is general-purpose software that can be used to analyze 1D electrophoresis gels, membranes, dot blots, and arrays.

References

- Bjerneld, E. et al. Pre-labeling of diverse protein samples with a fixed amount of Cy™5 for sodium dodecyl sulfate-polyacrylamide gel electrophoresis analysis. Anal. Biochem. 484, 51-57 (2015).
- 2. Application note: Analysis of therapeutic antibodies using Amersham™ WB system, Cytiva, 29114027, Edition AA (2014).
- 3. Hagner-McWhirter, Å. et al. Cy™5 total protein normalization in Western blot analysis. Anal. Biochem. **486**, 54–61 (2015).

Ordering information

| Quantity | Product code | |
|---------------|--|--|
| 1 | RPN4000 | |
| | | |
| | | |
| 50 columr | ns 28918007 | |
| er 500 mL | 28941279 | |
| 1 × 500 μL | _ 28444008 | |
| 2 × 1 mL | 28985738 | |
| | | |
| 150 µg | PA43009 | |
| oat-α- 150 μg | 28901106 | |
| | | |
| escent 120 µL | RPN850E | |
| escent 500 μL | RPN851E | |
| | 50 column er 500 mL 1 × 500 μL 2 × 1 mL 150 μg oat-α- 150 μg | |

100 sheets

3030-861

Related products

| Blotting membranes | Quantity | Product code |
|---|----------------------|-----------------|
| Amersham™ brand Hybond™ LFP | 1 roll/PK | 10600022 |
| Amersham [™] brand Hybond [™] LFP, 30×60 cm | 5 sheets /PK | 10600040 |
| Amersham™ brand Hybond™ LFP, 20 × 20 cm | 10 sheets/PK | 10600060 |
| Amersham™ brand Hybond™ LFP, 10 × 10 cm | 10 sheets/PK | 10600091 |
| Amersham™ brand Hybond™ LFP Sandwich 0.2 PVDF, 8 × 9 cm | For 10 minigel blots | 10600123 |
| Amersham™ brand Protran™ Premium 0.2 NC | 1 roll/PK | 10600004 |
| Amersham [™] brand Protran [™] Premium 0.2 NC 30×60 cm | ,5 sheets/PK | 10600035 |
| Amersham [™] brand Protran [™] Premium 0.2 NC 20×20 cm | ,10 sheets/PK | 10600050 |
| Amersham™ brand Protran™ Premium 0.2 NC 10 × 10 cm | ,10 sheets/PK | 10600081 |
| Amersham™ brand Protran™ Premium Sandwich 0.2 NC, 8 × 9 cm | For 10 minigel blots | 10600118 |
| Blocking agent | | |
| Amersham ECL™ Prime Blocking Reagent | 40 g | RPN418 |
| Imagers | | |
| Amersham™ brand Typhoon™ 5 | 1 | 29187191 |
| Amersham™ brand Typhoon™ RGB | 1 | 29187193 |
| Amersham™ brand Typhoon™ 600RGB | 1 | 29083467 |

Related literature

| Noiatoa iitorataro | |
|---|----------|
| Data file: Amersham™ Imager 600 | 29098107 |
| Data file: Amersham™ brand Typhoon™ Biomolecular Imager | 29229537 |
| Handbook: Imaging, principles and methods | 29020301 |
| Handbook: Western blotting, principles and methods | 28999897 |
| Application note: Accurate comparability assessment of a biosimilar interferon in process development | 29115478 |
| Application note: Analysis of therapeutic antibodies using Amersham™ WB system | 29114027 |

cytiva.com/amersham

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. Amersham, Capto, Cy, CyDye, ECL, ECL Plex, Hybond, ImageQuant, MiniTrap, Protran, Sepharose, and Typhoon are trademarks of Global Life Sciences Solutions USA LLC or an affiliate.

Coomassie and Novex are trademarks of Thermo Fisher Scientific LLC. Tween is a trademark of Croda Group of Companies. All other third-party trademarks are the property of their respective owners.

© 2021 Cytiv

3MM Chr

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

