



Flow accuracy and gradient performance of ÄKTA pilot 600

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY14004-02Jun20-AN

Flow accuracy and gradient performance of ÄKTA™ pilot 600

Accurate flow is important for your chromatography step to provide reproducible results at different scales. For ÄKTA pilot 600 chromatography system, flow and gradient formation is controlled by the UNICORN™ system control software. In this work, flow accuracy was investigated at a combination of different backpressures, flow rates, temperatures, and liquid viscosities. Gradient accuracy and range were calculated from the flow accuracy. Gradient accuracy was also evaluated by running step and linear gradients at different flow rates. Within the buffer B gradient range (%B) from 1%B to 99%B at a system flow rate of 5–600 mL/min, gradient accuracy was found to be $\pm 1\%$ B.

Introduction

The ability to scale up your bioprocess in an efficient way is important to increase speed to market. A comparable flow accuracy between scales can facilitate scaling of your methods. When scaling, it is important that the results (as visualized in the chromatogram) are comparable between scales. Gradient formation should be fast and reproducible in its full length at all available flow rates, backpressures, temperatures, and liquid viscosities. Examples of applications where maintained flow accuracy and gradient performance are especially important when scaling up your process are separation of mAb charge variants, of closely related entities, and of proteins of similar charge and size.

ÄKTA pilot 600 is an automated bench-top chromatography system equipped with two high-performance, low-pulsation piston pump modules to ensure accurate and reproducible results (Fig 1). The piston pump modules enable creation of binary gradients with efficient mixing at all flow rates and pressures allowed with the system. For high pump accuracy, even at low flow rates, and to prevent siphoning when the pumps are not running, a restrictor module is situated directly after the pump modules. The maximum flow rate of the system is 1200 mL/min when operated in isocratic mode. In gradient mode, the maximum flow rate is 600 mL/min.



Fig 1. ÄKTA pilot 600 high-accuracy chromatography system designed to facilitate sanitization (1). 1= pumps, 2 = restrictor, and 3 = column valve.

Maximum operating pressure is 20 bar. This work demonstrates flow accuracy and gradient performance of ÄKTA pilot 600.

Designation of flow accuracy

Flow accuracy can be specified either as percentage of full scale (% FS) or as percentage of reading (% RD), or a combination of both. If an instrument has a flow accuracy specified as % FS, the error will be of a fixed value irrespective of flow rate. For example, if the full-scale flow rate is 1000 mL/min and the system has an accuracy of 2% FS, the deviation will be ± 20 mL/min at all flow rates. This means that at a flow rate of 100 mL/min, the deviation will also be ± 20 mL/min (or $\pm 20\%$ deviation from the read value).

If a system instead has a flow accuracy specified as % RD, the error will always be the same percentage of the read value. In such a case, with a system accuracy of 2% RD, the deviation would be ± 20 mL/min at a flow rate of 1000 mL/min, but only ± 2 mL/min at 100 mL/min.

Designation of gradient accuracy

Gradient accuracy is often presented as a gradient acceptance range, with a defined accuracy within the acceptance range or as a chromatogram over a programmed and a performed gradient (Fig 2). The gradient acceptance range is given as % of buffer B (%B) for the flow rate range of the system, and the gradient accuracy is defined as the deviation from %B. The chromatogram is usually presented as the programmed %B and the performed %B at a single flow rate. The higher the flow rate; the lower the system dead volume; and the longer the gradient formation time, the better the gradient performance of the system.




Fig 2. (A) Example of gradient acceptance range for diaphragm pumps. (B) Example of how diaphragm pump gradient accuracy can be illustrated.

Materials and methods

Flow accuracy of ÄKTA pilot 600 was investigated at a combination of different backpressures (0 to 20 bar), flow rates (0.2 to 600 mL/min), temperatures (4°C to 35°C), and liquid viscosities (0.7 to 5 cP). The gradient accuracy and acceptance range were calculated from the flow accuracy. Gradient accuracy was also investigated by running step and linear gradients at different flow rates, using H₂O as buffer A and 1% acetone in H₂O as buffer B.

Results

Flow accuracy of ÄKTA pilot 600 was found to be $\pm 2\%$ RD. Gradient deviation, within the acceptance range 1%B–99%B, was less than 1%B, at system flow rates of 5–600 mL/min. At dual pump flow rates between 600–1200 mL/min, the gradient accuracy was the same, but the gradient range was limited by the maximum flow rate of 600 mL/min of each pump. Figure 3 depicts the gradient range at flow rates between 0–1200 mL/min. The gradient formation at low, intermediate, and high flow rates is illustrated in Figure 4.




Fig 3. Gradient range of ÄKTA pilot 600. Blue = normal flow, red = dual flow.




Fig 4. Chromatograms from ÄKTA pilot 600, showing the programmed and the resulting step and linear gradients at 100, 300, and 600 mL/min, respectively.

Conclusion

This work demonstrates the flow and gradient accuracy of ÄKTA pilot 600. The flow accuracy of ÄKTA pilot 600 was found to be $\pm 2\%$ RD. Gradient deviation was found to be less than 1%B within the gradient range from 1%B to 99%B at a system flow rate of 5–600 mL/min. The high accuracy over a multitude of parameters such as flow rate, viscosity, and pressure indicates that the system has a highly stable performance, omitting the need of a separate flowmeter to feedback control the pumps. The high accuracy, together with small hold-up volumes, strong sanitary performance (1), and a similar gradient performance as the ÄKTA avant and ÄKTA pure chromatography systems, makes ÄKTA pilot 600 well-suited for scaling, while minimizing the uncertainties associated with process transfer between GMP and non-GMP environments.

Reference

1. Application note: Sanitization of ÄKTA pilot 600 using sodium hydroxide. GE Healthcare, KA453DDMM17AN.

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

gelifesciences.com/bioprocess

GE, the GE Monogram, ÅKTA, and UNICORN are trademarks of General Electric Company.
Any use of UNICORN software is subject to GE Healthcare Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.
© 2018 General Electric Company
All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.
A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.
GE Healthcare UK Ltd., Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK
GE Healthcare Europe GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany
GE Healthcare Bio-Sciences Corp., 100 Results Way, Marlborough, MA 01752, USA
HyClone Laboratories Inc., 925 W 1800 S, Logan, UT 84321, USA
GE Healthcare Japan Corp., Sanken Bldg., 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan
For local office contact information, visit gelifesciences.com/contact.
KA454070318AN