



What do cells need from a bioreactor?

Intellectual Property Notice: The Biopharma business of GE Healthcare was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company. Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. All other third-party trademarks are the property of their respective owners. © 2020 Cytiva All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information. For local office contact information, visit cytiva.com/contact

CY15165-13Jul20-WP



What do cells need from a bioreactor?

Bioreactors are widely used to cultivate cells during the development and manufacturing of modern biopharmaceuticals. Cells are very sensitive to changes in the culture environmental conditions, such as aeration, agitation, nutrients, and pH. This paper discusses the importance of aeration and available options to control the oxygen mass transfer coefficient ($k_L a$) within a bioreactor.

Delivering oxygen to cells

In cell culture, oxygen is a key substrate for growth, production, and maintenance activities. Cells obtain their oxygen in free and noncompound forms, called dissolved oxygen (DO). One of the most important functions of bioreactors is providing dissolved oxygen to cells continuously through a process called aeration. Aeration in the bioreactor typically occurs when:

1. Oxygen diffuses through overlay to the cell culture medium interface
2. Oxygen from the spargers dissolves in the cell culture through convection with the help of agitation

Agitation disperses the oxygen bubbles and promotes mass transfer of the gas bubbles through the gas-liquid (cell culture medium) interface (Fig 1). The rate of oxygen transfer (OTR) from gas to liquid interface is a function of physicochemical properties of the cell culture medium, the geometrical parameters of the bioreactor, and presence of cells.




Fig 1. Diagram of a gas bubble in liquid, showing how the bubble is released, solubilized, and transferred to a cell.

Oxygen utilization rate (OUR) is often cell line-dependent. Table 1 lists the rates for common industrial cell lines (1, 2).

Table 1. OUR of cell lines typically used in biomanufacturing

Cell line	OUR [10^{-13} mol/cell/h]
DG44 (CHO)	2
CHO	5.0–8.04
NSO (myeloma)	2.19–4.06
MAK (hybridoma)	4.6
FS-4 (human diploid cells)	0.5
HFN7.1 (hybridoma)	2

Due to its low solubility in liquid phase (Fig 1) and increasing metabolic consumption by the cells with time, oxygen is supplied continuously to the cell culture. Oxygen supply is carefully controlled for optimal cell growth by manipulating bioreactor parameters.

During batch cell culture, OUR (or OTR) is initially low during the lag phase (Fig 2), where cells are self-synthesizing and there is little gain of cell density. As cell density increases during the exponential phase, OUR increases until OTR becomes a limiting rate, as determined by the mass transfer of oxygen into the bulk liquid.




Fig 2. Phases of cell growth.

The OTR and OUR rates are correlated by the oxygen mass transfer coefficient, $k_L a$. Therefore, the OTR, through its correlation to $k_L a$, defines a theoretical maximum cell density that could be achieved in cell culture.

Utility of $k_L a$ values

Because of this association with cell density, $k_L a$ values are particularly useful in the following scenarios:

Scenario 1: Evaluating scalability within the same bioreactor platform

The conventional scale-up of bioprocesses is based on physicochemical and geometric similarity. $k_L a$ is kept constant for this scenario. The OTR should remain constant for a bioreactor platform with geometric similarity. Bioreactor physical characteristics at the different scales are altered to provide the necessary OTR at controlled temperature, pH, and DO to achieve the target cell density.

Scenario 2: Technical transfer across different bioreactor designs

During the comparison, $k_L a$ is utilized as a target performance metric when a process is transferred from one bioreactor platform to another design. Bioreactor hardware design (e.g., stirrer geometry and aeration-sparger option) and running parameters (e.g., gas flow rate or power input) are altered to achieve a similar $k_L a$, providing a similar cell density.

Equation for $k_L a$

Figure 1 is a diagram of a gas bubble in liquid. For this discussion imagine that the gas bubble contains oxygen, and the liquid is the liquid in a bioreactor. $k_L a$ can be represented by the following equation:

$$k_L a = k_L \times a$$

- Where $k_L a$ is the mass transfer coefficient from the gas to liquid phase, given in sec^{-1}
- k_L = liquid side mass transfer coefficient (resistance in gas side film can be neglected)
- a = bubble surface (available for diffusion)

Key variables that impact k_La values

Any change to process and engineering parameters or to physical characteristics will have an impact on k_La and should be considered when evaluating bioreactor platforms and performing scaling calculations. Here are four key variables that can affect k_La values (Fig 3):




Fig 3. Diagram of a bioreactor process showing key factors that can influence k_La values.

Gas bubble size

When gas bubble size decreases, surface area and gas residency time increases, causing bubbles to stay in the culture longer. Thus, there is a greater opportunity for oxygen to release mass transfer into the cell culture medium. An increase in this oxygen residence time improves k_La .

Mixing

In a bioreactor, mixing is used to eliminate gradients of concentration (cell, gas, medium, and nutrient), temperature, and other properties. Mixing time is widely used to characterize mixing efficiency in a bioreactor. Mixing efficiency is one of the most significant factors affecting both performance and scale-up in a bioreactor.

Gas bubble size and residency time are highly dependent upon three mixing conditions: impeller type, speed, and location(s).

k_La values generally increase as tip speed increases. However, tip speed is proportional to shear forces that can lead to cell death. Bioreactors, therefore, are designed with different impeller types, combinations, and locations to achieve target k_La values without creating these shear forces. Generally, k_La values are closely associated with impeller design, with Rushton typically higher than paddle, which is typically higher than marine and pitched impeller.

Air flow rate

Higher oxygen availability drives k_La increases. Increasing oxygen supply to a bioreactor drives this availability and can be controlled by modifying concentration (air vs O₂ enrichment) and volumetric flow. Although high k_La values are desirable, it is important to consider the actual operating conditions and implications to cell viability and associated process costs. For example, high air flow rates can cause cell damage due to shear forces. Excessive foam might also be generated, requiring a high concentration of antifoam that could hinder downstream processing. Additionally, higher air flow rates require a larger exhaust filter area, driving consumable cost increases.

Properties of the liquid or medium

During cell culture, small bubbles collide and coalesce to form larger bubbles, decreasing surface area (a) and subsequently k_La . Be aware of reported k_La values in which high salt concentrations are used, because this can prevent bubble coalescing.

Antifoaming agents are used to influence surface tension, resulting in reduced bubble coalescence and foaming. However, this principle does not always lead to increases in OTR wherein antifoam also reduces bubble mobility, which subsequently reduces the k_La (3).

Other factors that affect cell culture k_La

Measurement method

Several different methods are used. Most commonly the nitrogen stripping (i.e., gassing-out) method is employed.

When scaling a process within the same platform, it is important to use an identical method for measuring k_La . k_La , when combined with process engineering parameters (i.e., tip speed, power input), can be used to experimentally determine the cell density in a larger bioreactor compared with a smaller bioreactor.

Temperature

Increasing temperatures inversely affects both the volumetric mass transfer coefficient and oxygen solubility in culture medium. Oxygen solubility in pure water falls with increasing temperature (i.e., $-0.5 \times 10^{-3} \text{ kg/m}^3$ between 35°C and 30°C; 3). Therefore, it is important to note the temperature conditions from vendor-supplied characterization data.

Sparger characteristics

k_La values will vary widely with sparger characteristics, including number, pore size, and surface area, because these factors affect bubble size, gas velocity, and flow rates.

Conclusion

As discussed here, k_La is impacted by multiple factors. A thorough understanding of a bioreactor platform's physical design, mixing mechanism, sparging options, as well as cell line characteristics will inform decision-making when scaling cell culture processes up and out.

References

1. Ruffieux P. A. *et al.* Measurement of volumetric (OUR) and determination of specific (qO_2) oxygen uptake rates in animal cell cultures. *Journal of Biotechnology* **63**: 85–95 (1998).
2. Xiu Z. L. *et al.* Estimation of rates of oxygen uptake and carbon dioxide evolution of animal cell culture using material and energy balances. *Cytotechnology* **29**: 159–166 (1999).
3. Doran P. M. *Bioprocess engineering principles*, **Second edition**, Academic Press, Massachusetts (2012).

gelifesciences.com

GE and the GE Monogram are trademarks of General Electric Company.

© 2018 General Electric Company

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.

A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Ltd., Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Europe GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp., 100 Results Way, Marlborough, MA 01752, USA

HyClone Laboratories Inc., 925 W 1800 S, Logan, UT 84321, USA

GE Healthcare Japan Corp., Sanken Bldg., 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan

GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden

For local office contact information, visit gelifesciences.com/contact.

KA1217020318WP