




Preparative purification of proteins with size exclusion chromatography columns



Content

For direct access, click on your size exclusion chromatography (SEC) topic of interest

When and why should you use SEC

1. SEC fundamentals [>>](#)
2. Why use SEC in a protein purification protocol? [>>](#)

Cytiva SEC columns

3. How to select Cytiva SEC columns [>>](#)
4. HiPrep™ Sephacryl™ columns [>>](#)
5. HiLoad™ columns [>>](#)
6. Superdex™ Increase and Superose™ Increase for volumes < 0.5 mL [>>](#)
7. Superdex™ Increase and Superose™ Increase for volumes > 0.5 mL [>>](#)

SEC tips, tools, and summary

8. Tips for successful SEC [>>](#)
9. Useful SEC tools [>>](#)
10. Summary [>>](#)

Appendix

11. Appendix 1 — Help for column selection [>>](#)
12. Appendix 2 — Application examples [>>](#)


1

SEC fundamentals


Chromatography techniques commonly used for protein purification

Technique:


**Ion exchange
Chromatography (IEX)**




**Hydrophobic interaction
Chromatography (HIC)**



**Size exclusion
Chromatography (SEC)**



**Affinity
Chromatography (AC)**



Separation principle:

Charge

Hydrophobicity



Size

Biorecognition

**Chromatography techniques enable separation of proteins
based on differences in specific properties**

Principles of SEC

Separates molecules based on size



Largest molecules elute first

Characteristics


- Nonbinding technique — the separation takes place in only 1 column volume (CV)
- Mild conditions — good for sensitive biomolecules
- Any buffer can be used
- Limited in sample volume
- By nature a slow technique

A well-packed SEC column is critical for high-resolution separations

More on SEC animation and Cytiva handbook


Watch the video on YouTube

The principle of gel filtration (size exclusion chromatography)



Download the handbook

Size Exclusion Chromatography Principles and Methods




2

Why use SEC in a protein purification protocol?

The CiPP model for simplified protein purification planning

Capture, intermediate Purification, and Polishing (CiPP)




The initial capture stage isolates, concentrates, and stabilizes the target protein.

Intermediate purification removes bulk contaminants.

The final polishing step removes the most difficult impurities, such as aggregates or isoforms of the target protein.

SEC is widely used for polishing

Capture, intermediate Purification, and Polishing (CiPP)




SEC is widely used for polishing:

- It effectively removes dimers and aggregates of the target protein
- The target protein will be size homogeneous


SEC also simultaneously enables the transfer of the target protein to the buffer of choice.

How many chromatography steps should be used in a purification protocol?




The number of steps to be included will depend on the purity requirements and intended use of the protein.

Addition of chromatography steps will increase purity at the cost of decreased yield of active protein.



Use of SEC to improve antibody purity



B = buffer exchange to neutralize low pH Ab elution buffer

C = concentration for sample volume reduction. May also be performed before SEC.


After the first affinity chromatography step that isolates the antibody from initial sample, a SEC step will remove antibody aggregates and/or fragments to obtain monomeric antibodies.

After SEC, the purity is very high: 95% to 99%.

Use of SEC to improve antibody purity

The addition of SEC removed antibody aggregates to obtain monomeric antibodies

Step 1 = AC




Column: [HiTrap™ MabSelect™ Prisma 1 mL](#)
Binding buffer: 20 mM phosphate, 150 mM NaCl pH 7.4,
Elution buffer: 50 mM sodium acetate pH 3.5
Sample: 6 mL of supernatant containing polyclonal human IgG
Flow rate: 0.5 mL/min


¹ These columns are produced on-demand.



Step 2 = SEC



Use of SEC to improve antibody purity



B1 = buffer exchange to remove imidazole or salts

B2 = buffer exchange to prepare for IEX

C = concentration for sample volume reduction. May also be performed before SEC.

HCP = host-cell proteins

A single IMAC step delivers a moderate protein purity (> 80%).


Whether you choose a 2-step or 3-step protocol, SEC will be used as a last step for removal of remaining impurities.

In the 3-step purification protocol, IEX enables removal of impurities such as HCP.

Use of SEC to improve purity of his-tagged proteins

The addition of a SEC step removes impurities such as truncated forms and aggregates of your target protein


Step 1 = IMAC



Column: HisTrap™ FF 1 mL

Sample: 50 mL of (his)₁₀-Trx-P 450 in *E. coli* lysate


Step 2 = SEC



Column: HiLoad™ 16/60 Superdex™ 200 pg

Sample: 5.2 mL of eluted pool from HisTrap FF

SDS-PAGE analysis



Step 1


Step 2

3

**How to select the most
suitable Cytiva SEC column
for your protein purification**

Cytiva offers prepacked SEC columns for user convenience and reproducible results

For a broad range of biomolecules and resolutions




Colored bars indicate the fractionation range for each resin



How to select the best SEC column for a specific application


- Choose the resin that has a fractionation range where the target molecule falls in the middle of the range
- If contaminants are close in size to the target molecule, choose a resin with higher resolution
- Choose column type depending on the sample volume that should be applied



The image shows the cover of a Cytiva selection guide. On the left, the text 'Selection guide' is in small green font, followed by 'Size exclusion chromatography columns and resins' in large, bold, green font. Below this is the Cytiva logo. On the right, there is a photograph of a person in a lab coat and blue gloves holding a test tube. At the bottom right, a green button contains the text 'Download our selection guide >>'.

Selection guide


Size exclusion chromatography columns and resins

 **cytiva**

Download our selection guide >>

General guidelines for SEC column selection

For more information, click on the column of your interest



¹These columns are produced on-demand. [Contact your Cytiva representative](#)


For more detailed guidance, see [Appendix](#)



For volumes > 0.5 mL, the choice of column depends on the resolution needed


For good resolution: HiPrep™ Sephacryl™

HiPrep 16/60 Sephacryl S-300 HR
M_r range ~ 10 000–1 500 000




For high resolution: HiLoad™ Superdex™/Superose™ pg

HiLoad 16/600 Superdex 200 pg
M_r range ~ 10 000–600 000



For highest resolution: HiScale™ with “Increase” resin

Superdex 200 Increase HiScale 16/40
M_r range ~ 10 000–600 000



Highest resolution →

Protein mix (HiPrep Sephacryl + HiLoad Superdex)

- | | |
|---------------------------------------|---|
| 1. Ferritin (M _r 440 000) | 4. Ovalbumin (M _r 44 000) |
| 2. Aldolase (M _r 158 000) | 5. Carbonic anhydrase (M _r 29 000) |
| 3. Conalbumin (M _r 75 000) | 6. Ribonuclease (M _r 13 700) |

Running conditions	HiPrep Sephacryl and HiLoad Superdex pg	Superdex Increase HiScale
Sample volume	0.5 mL	1.6 mL
Flow rate	1 mL/min	2 mL/min

Protein mix (HiScale Superdex Increase)

1. Thyroglobulin (M_r 669 000)
2. Aldolase (M_r 158 000)
3. Ovalbumin (M_r 44 000)
4. Ribonuclease (M_r 13 700)

Data files available for download: sample volume > 0.5 mL

Click on the document of interest

HiLoad™ Superdex™ columns

Superdex prep grade and prepacked HiLoad columns

SIZE EXCLUSION CHROMATOGRAPHY

Superdex™ prep grade (pgl) is a high resolution size exclusion chromatography resin (Fig 1). It is composed of cross-linked agarose and dextran. The steep selectivity of the dextran component and the high chemical and physical stability of the agarose give high-resolution separations at flow velocities up to 50 cm/h. Three resin types are available in laboratory and larger pack sizes: Superdex 50 prep grade, Superdex 75 prep grade, and Superdex 200 prep grade. Their main features are:

- Steep selectivity provides high resolution
- High chemical stability

The resins are also available in prepacked, high-performance HiLoad™ columns offered in two different column sizes, 16 and 26 mm diameter, both with 600 mm bed height (Fig 2). The HiLoad 16/600 and 26/600 columns provide a number of significant advantages for high resolution work:

- Prepared for convenience and reproducibility
- High-resolution separation of biomolecules
- High chemical stability and easy scale-up
- Easy connection to, for example, AKTA™ chromatography systems

Each column is expertly packed and individually tested. This combination of prepacked convenience and reproducibility makes HiLoad Superdex pg columns a confident choice for fast, high-resolution size exclusion chromatography at preparative laboratory scale.

The columns run with a wide variety of equipment: AKTA systems or simple pump-based configurations.




Fig 1. Superdex size exclusion chromatography resins.





Fig 2. HiLoad Superdex 50, 75, and 200 pg columns bring convenience and high resolution to size exclusion chromatography. Each is available in two column sizes: HiLoad 16/600 and HiLoad 26/600.

Download



HiLoad 16/600 Superose™ 6 pg

HiLoad 16/600 Superose 6 prep grade column

HiLoad™ 16/600 Superose™ 6 prep grade (pgl) prepacked column (Fig 1) is designed for purification of large proteins and protein complexes by size exclusion chromatography (SEC) in sample volumes up to 5 mL. The column is packed with Superose 6 prep grade chromatography resin having a fractionation range suitable to separate biomolecules with molecular weights (M_r) from ~5000 to 5 000 000 by size exclusion chromatography (SEC). This very wide fractionation range also makes the column suitable for purification of membrane proteins and other macromolecules.

HiLoad 16/600 Superose 6 pg column offers:

- Convenient prepacked column format that removes the need for self-packing and enables highly reproducible protein separations
- High-resolution separation of protein complexes and other large biomolecules for high protein purity
- High chemical stability enabling thorough cleaning for long column lifetime and minimal carry-over
- Easy connection to AKTA™ and other liquid chromatography systems

HiLoad 16/600 Superose 6 pg column is an empty KX column that has been expertly packed and individually tested. The combination of prepacked convenience and reproducibility makes HiLoad 16/600 Superose 6 pg column a confident choice for fast, high-resolution SEC at preparative scale. At maximum recommended sample load (5 mL), up to about 250 mg of proteins can be purified.

The column operates with a wide variety of equipment, including simple pump-based configurations and AKTA chromatography systems.

Resin characteristics

The chromatography resin in the column is Superose 6 prep grade and this is based on cross-linked agarose. The size and the distribution of the particles allow high flow, high efficiency, and good capacity. Further characteristics of Superose 6 prep grade resin are shown in Table 1.





Fig 1. HiLoad 16/600 Superose 6 pg column brings convenience and high resolution to size exclusion chromatography (also known as gel filtration).

Matrix	Composite of cross-linked agarose
Particle size, μm	50.1 to 50
Particle size, nm	50.1 to 50
Fractionation range	6×10^3 to 5×10^6 (globular proteins)
Exclusion limit (M _r)	$\sim 4 \times 10^6$ (globular proteins)
pH stability	3 to 12
operational ^a	cleaning-in-place (CIP) ^b
Temperature	operational 4°C to 40°C storage 4°C to 30°C

^a Media particle size of the sorbent is not a function of flow rate.
^b pH range where resin can be subjected to cleaning- or sanitization-in-place without significant change in function.

Download



HiPrep™ Sephacryl™ columns

Sephacryl High Resolution resins HiPrep Sephacryl HR columns

SIZE EXCLUSION CHROMATOGRAPHY

Sephacryl™ High Resolution (HR) chromatography resins allow fast and reproducible purification of proteins, polysaccharides, and other macromolecules by size exclusion chromatography (SEC) at laboratory and industrial scale. Five Sephacryl HR chromatography resins are available as prepacked columns and in laboratory and larger pack sizes: Sephacryl S-100 HR, S-200 HR, S-300 HR, S-400 HR, and S-500 HR. Characteristics of the resins include:

- Purification over a wide molecular weight range
- High reproducibility due to high stability
- High recoveries
- Well-suited to industrial-scale use

All five Sephacryl HR resins are available in prepacked HiPrep™ Sephacryl HR SEC columns. Each chromatography resin is available in two different prepacked column sizes: 16/60 (120 mL) and 26/100 (320 mL). HiPrep Sephacryl HR SEC columns provide the excellent purification properties of Sephacryl HR resins in a convenient, ready-to-use format.

Characteristics of the columns include:

- Convenient, easy-to-use, prepacked SEC columns in two different column sizes
- Choice of five selectivities covering a wide molecular weight range
- Reliable and reproducible preparative purification
- Easy connection to AKTA™ chromatography systems




Fig 1. Sephacryl HR chromatography resins and HiPrep Sephacryl HR prepacked columns offer the user a wide range of choice and reliable purification by SEC.


Characteristics of resins

The matrix of Sephacryl HR resins is a cross-linked copolymer of allyl dextran and N,N'-methylene bisacrylamide. This crosslinking gives good rigidity and chemical stability. The narrow particle size distribution, together with steep selectivity curves, results in good preparative characteristics with maintained resolution. The hydrophilic nature of the resins minimizes nonspecific adsorption and maximizes recovery.

The high resolution and flow characteristics, long-term physical and chemical stability, and ease of handling make Sephacryl HR the resin of choice for routine purification.

Sephacryl HR SEC resins fulfill process chromatography requirements in terms of stability, scalability, and bulk availability. As members of the BioProcess™ family of chromatography resins, they carry full technical and regulatory support for production-scale operations.

Download



Data files available for download: sample volume < 0.5 mL

Click on the document of interest

Superdex™ 30 Increase >>

Superdex 30 Increase columns

SIZE EXCLUSION CHROMATOGRAPHY

Superdex™ 30 Increase prepacked columns (Fig 1) are designed for high-resolution, small-scale purification and analysis of peptides and other biomolecules with molecular weights (M) from ~100 to 7000. These new generation size exclusion chromatography (SEC) columns replace their predecessors, Superdex Peptide columns, delivering improved performance in preparative and analytical purification.

Superdex 30 Increase columns provide:

- Small-scale preparative purification (sample volumes of 4 to 500 µL) and characterization of peptides and small biomolecules
- Higher resolution compared with Superdex Peptide for improved purity and analysis results
- Three times faster separations than Superdex Peptide with the same resolution
- Tolerance to repeated harsh cleaning procedures at high pH (1 M NaOH), giving long column life and minimal carry-over

Resin characteristics

Superdex 30 Increase resin is based on a high-flow agarose base matrix with good pressure-flow properties. The small bead size with a narrow bead size distribution allows for high-resolution separations. In addition, low nonspecific interaction permits high recovery of biological materials. The characteristics of the Superdex 30 Increase resin are shown in Table 1.



Fig 1. Superdex 30 Increase 10/300 GL and Superdex 3.2/300 GL columns for high-resolution separation and analysis of peptides and small biomolecules.

Matrix	Composite of cross-linked agarose and dextran
Particle size, d _{av} ¹	9 µm
Fractionation range	M _w ~100 to 7000
pH stability	3 to 12
operational ²	cleaning-in-place (CIP) ³ 1 to 14
temperature	4°C to 40°C
storage	4°C to 30°C

¹ Median particle size of the cumulative volume distribution.
² pH range where resin can be operated without significant change in function.
³ pH range where resin can be subjected to cleaning or sanitation in place without significant change in function.

cytiva.com

cytiva

DOWNLOAD

Superdex 75 Increase >>

Superdex 75 Increase columns

SIZE EXCLUSION CHROMATOGRAPHY

Superdex™ 75 Increase prepacked columns are designed for rapid separation and analysis of proteins with molecular weights ranging from M_w 3000 to 70 000 by size exclusion chromatography (SEC), also called gel filtration (Fig 1). These versatile columns offer separations with high resolution for a variety of applications including preparative protein purification (µg to mg quantities), aggregate analysis, studies of complex formation, and screening of samples. Superdex 75 Increase belongs to the new generation of SEC columns that replaces the well-known Superdex 75 columns.

Superdex 75 Increase columns offer:

- Increased resolution compared with Superdex 75, for higher purity and improved analysis results
- Reduced runtime compared with Superdex 75, for obtaining results faster
- Versatile use in both preparative and analytical applications
- Tolerance to repeated harsh cleaning procedures at high pH, giving long column lifetime and minimal carry-over

Resin characteristics

Superdex 75 Increase resin is based on a high-flow agarose base matrix with good pressure-flow properties. The small bead size with a narrow size distribution allows for high-resolution separations. In addition, low nonspecific interactions permit high recovery of biological materials. The characteristics of the Superdex 75 Increase resin are shown in Table 1.




Fig 1. Superdex 75 Increase 10/300 GL and Superdex 3.2/300 GL columns.

Matrix	Composite of cross-linked agarose and dextran
Fractionation range	M _w 3000 to 70 000 (globular proteins)
Exclusion limit	M _w ~100 to 30 000 (dextran)
Particle size, d _{av} ¹	9 µm (wet)
pH stability	3 to 12
operational ²	cleaning-in-place (CIP) ³ 1 to 14
Storage temperature	4°C to 40°C
Storage temperature	4°C to 30°C

¹ Median particle size of the cumulative volume distribution.
² pH range where resin can be operated without significant change in function.
³ pH range where resin can be subjected to cleaning or sanitation in place without significant change in function.

cytiva.com

cytiva

DOWNLOAD

Superdex 200 Increase

Superdex 200 Increase columns

SIZE EXCLUSION CHROMATOGRAPHY

Superdex™ 200 Increase prepacked columns (Fig 1) are designed for size exclusion chromatography (SEC)/high-resolution gel filtration in small-scale (µg to mg), preparative purification, as well as for characterization and analysis of proteins with molecular weights (M) from 10 000 to 600 000. These versatile columns offer rapid separations with high resolution for a variety of applications including protein purification, aggregate analysis, studies of complex formation, and screening of samples.

Superdex 200 Increase columns offer:

- Versatile use for both preparative and analytical purposes, for a large variety of proteins
- Increased resolution compared with Superdex 200, for higher purity and improved analysis results
- Reduced runtime compared with Superdex 200, to get results faster

Size exclusion chromatography

SEC separates molecules according to their differences in size as they pass through a chromatography resin packed in a column. Unlike ion exchange or affinity chromatography, molecules do not bind to the chromatography resin and the buffer composition has no direct effect on resolution (the degree of separation between peaks). Consequently, a significant advantage of SEC is that conditions can be varied to suit the sample type or the requirements for further purification, analysis, or storage, without altering the separation.

SEC is an excellent technique for discriminating between monomer, oligomer, and aggregated forms of a target protein.

Resin characteristics

The chromatography resin in Superdex 200 Increase columns is based on a high-flow agarose base matrix with good pressure-flow properties and a small particle (bead) size (median particle size, 8.6 µm). In addition, the low nonspecific interaction




Fig 1. Superdex 200 Increase 10/300 GL and Superdex 3.2/300 GL columns.

Matrix	Composite of cross-linked agarose and dextran
Fractionation range	M _w 10 000 to 600 000 (globular proteins)
Exclusion limit	M _w ~1000 to 100 000 (dextran)
Particle size, d _{av} ¹	8.6 µm
pH stability	Operational ² 3 to 12 Cleaning-in-place (CIP) ³ 1 to 14
Operating temperature	4°C to 40°C
Storage temperature	4°C to 30°C

¹ Median particle size of the cumulative volume distribution.
² pH range where resin can be operated without significant change in function.
³ pH range where resin can be subjected to cleaning or sanitation in place without significant change in function.

cytiva.com

cytiva

DOWNLOAD

Superose™ 6 Increase >>

Superose 6 Increase columns

SIZE EXCLUSION CHROMATOGRAPHY

Superose™ 6 Increase prepacked columns (Fig 1) are designed for rapid separation and analysis of proteins and other biomolecules by size exclusion chromatography (SEC, also known as gel filtration), in the molecular weight (M_w) range for globular proteins between 5000 and 5 000 000. This very wide fractionation range of the chromatographic resin makes it suitable for purification of protein complexes, membrane proteins, and other macromolecules. The columns are also useful as a screening tool to explore the molecular-weight distribution of unknown samples. This new generation SEC resin replaces the well-known Superose 6 columns.

Superose 6 Increase columns offer:

- Increased resolution compared with Superose 6, for higher purity and improved analysis results
- Reduced runtime compared with Superose 6, for rapid results
- Versatile use for both preparative and analytical applications, especially for large proteins and complexes
- Tolerance to repeated harsh cleaning procedures at high pH, giving long column lifetime and minimal carry-over

Resin characteristics

The chromatography resin in Superose 6 Increase columns is based on a high-flow agarose base matrix with good pressure-flow properties and small bead size. The small beads of the resin allow high-resolution analytical separations. The low nonspecific interaction of the resin permits high recovery of biological materials. Further characteristics of the Superose 6 Increase resin are shown in Table 1.




Fig 1. Superose 6 Increase 10/300 GL and Superose 3.2/300 GL columns.

Fractionation range	M _w 5000 to 5 × 10 ⁶ (globular proteins)
Exclusion limit (M _w)	Approx. 4 × 10 ⁶ (globular proteins)
pH stability	3 to 12 (long term)
Operating temperature	4°C to 40°C
Storage temperature	4°C to 30°C
Matrix	Composite of cross-linked agarose
Average bead size	8.6 µm

¹ Median particle size.

cytiva.com

cytiva


DOWNLOAD

4

HiPrep Sephacryl SEC columns


HiPrep Sephacryl columns deliver good resolution over a broad fractionation range for routine SEC

- **Five different fractionation ranges**
- **Two column dimensions**
- **Sample volumes up to 5 mL and 13 mL**



Good price/resolution compromise!

Sephacryl resins technical specifications



- **Matrix:** Cross-linked copolymer of allyl dextran and N,N-Methylene bisacrylamide
- **Particle size, d_{50V}^1 :** ~ 50 μm
- **pH stability, operational²:** 3 to 11
- **pH stability, CIP³:** 2 to 13


Fractionation range/Resin	Sephacryl™ S-100 HR	Sephacryl S-200 HR	Sephacryl S-300 HR	Sephacryl S-400 HR	Sephacryl S-500 HR
Fractionation range (M_r) Globular proteins	1000 to 100 000	5000 to 250 000	10 000 to 1 500 000	20 000 to 8 000 000	No data
Fractionation range (M_p) Dextrans	No data	1000 to 80 000	2000 to 400 000	10 000 to 2 000 000	40 000 to 20 000 000

¹ Median particle size of the cumulative volume distribution.

² pH range where resin can be operated without significant change in function.

³ pH range where resin can be subjected to cleaning-in-place (CIP) without significant change in function.

Technical specifications — HiPrep prepacked columns



- **Max. pressure over the packed bed during operation:**
0.15 MPa, 1.5 bar, 22 psi
- **Column hardware pressure limit:**
0.5 MPa, 5 bar, 73 psi


Parameter	HiPrep™ 16/60	HiPrep 26/60
Bed dimensions	16 mm × 600 mm	26 mm × 600 mm
Approximate bed volume	120 mL	320 mL
Recommended sample volume	Up to 5 mL	Up to 13 mL
Recommended operating flow rate	0.5 mL/min	1.3 mL/min
Max. operating flow rate	1.0 mL/min	2.7 mL/min

All different Sephacryl™ resins are available in both column sizes

Comparing separations on the different Sephacryl resins


Sephacryl™ S-100 HR

HiPrep 16/60 Sephacryl S-100 HR




Sephacryl S-200 HR

HiPrep 16/60 Sephacryl S-200 HR



Sephacryl S-300 HR

HiPrep 16/60 Sephacryl S-300 HR



Sample:

Standard proteins

500 µL of a mixture comprising :

IgG	(M _r 160 000)
BSA	(M _r 67 000)
β-lactoglobulin	(M _r 35 000)
cytochrome C	(M _r 12 400)
cytidine	(M _r 240)

Buffer:

50 mM sodium phosphate, 150 mM NaCl, pH 7.0

Flow rate:


0.8 mL/min

Detection:

A₂₈₀

Sephacryl S-400 HR

HiPrep 16/60 Sephacryl S-400 HR



Sample: Dextrans

1 mL of a mixture containing
Dextran > 1 × 10⁷
Dextran 410 (M_p 276 500)
Dextran 12 (M_p 9890)

Buffer:

0.25 M NaCl

Flow rate:


0.5 mL/min

Detection:

Refractive Index (RI)

Sephacryl S-500 HR

HiPrep 16/60 Sephacryl S-500 HR



Sample: Dextrans

1 mL of a mixture containing
Dextran > 1 × 10⁷
Dextran DXT1185K (M_p 1.05 × 10⁶)
Dextran 50 (M_p 43 500)


Check chromatograms to see which resin best fits with your sample.

5

HiLoad SEC columns

HiLoad Superdex prep grade and Superose prep grade¹ SEC columns deliver high resolution

- Four different fractionation ranges
- Two column dimensions
- Sample volumes up to 5 mL and 13 mL




¹ Superose 6 prep grade is from Nov. 2018 available in prepacked format



For high resolution and high recovery needs

Superdex prep grade and Superose prep grade resins technical specifications




- **Matrix:** Composite of cross-linked agarose and dextran (Superdex pg)
Composite of cross-linked agarose (Superose pg)
- **Particle size, d_{50v} ¹:** ~ 34 μm (Superdex pg); ~ 30 to 40 μm (Superose pg)
- **pH stability, operational²:** 3 to 12
- **pH stability, CIP³:** 1 to 14

Fractionation range/Resin	Superdex™ 30 prep grade	Superdex 75 prep grade	Superdex 200 prep grade	Superose™ 6 prep grade
Fractionation range (M _r) Globular proteins	< 1000	~ 3000 to 70 000	~ 10 000 to 60 000	~ 5000 to 5 000 000

¹ Median particle size of the cumulative volume distribution.
² pH range where resin can be operated without significant change in function.
³ pH range where resin can be subjected to cleaning-in-place (CIP) without significant change in function.

Technical specifications — HiLoad prepacked columns




- **Max pressure over the packed bed during operation:**
0.3 MPa, 3 bar, 42 psi
- **Column hardware pressure limit:**
0.5 MPa, 5 bar, 73 psi

Parameter	HiLoad™ 16/600 ¹	HiLoad 26/600 ¹
Bed dimensions	16 mm × 600 mm	26 mm × 600 mm
Approximate bed volume	120 mL	320 mL
Recommended sample volume	Up to 5 mL	Up to 13 mL
Recommended operating flow rate	1.0 mL/min	2.6 mL/min
Max. operating flow rate	1.7 mL/min	4.4 mL/min

¹ HiLoad columns are called XK when sold as empty columns. Superose™ 6 prep grade in XK 26/70 is available as a custom column from CDP with code no. 90100043.


HiLoad 16/600 Superose 6 pg column — what results to expect?

Small sample volume and low flow




Sample volume: 0.5 mL
Flow rate: 0.8 mL/min

Small sample volume and high flow



Sample volume: 0.5 mL
Flow rate: 1.6 mL/min

Large sample volume and low flow



Sample volume: 5 mL
Flow rate: 0.8 mL/min

Protein mix (same in all three runs):

1. Thyroglobulin (M_r 669 000)
2. Ovalbumin (M_r 44 000)
3. Ribonuclease (M_r 13 700)


General rule: small sample volume and low flow rate gives the best resolution

6

**Superdex Increase and
Superose Increase columns
for sample volumes < 0.5 mL**


Superdex Increase and Superose Increase for highest resolution

- Four different fractionation ranges
- Three standard column dimensions + dimensions on demand
- Sample volumes from 4 μL up to 0.5 mL



For highest resolution and speed

Superdex Increase and Superose Increased resins technical specifications



- **Matrix:** Composite of cross-linked agarose and dextran (Superdex Increase)
Composite of cross-linked agarose (Superose Increase)
- **Particle size, d_{50v}^1 :** ~ 9 μm
- **pH stability, operational²:** 3 to 12
- **pH stability, CIP³:** 1 to 14

Fractionation range/Resin	Superdex™ 30 Increase	Superdex 75 Increase	Superdex 200 Increase	Superose™ 6 Increase
Fractionation range (M_r) Globular proteins	~ 100 to 7000	~ 3000 to 70 000	~ 10 000 to 600 000	~ 5000 to 5 000 000

¹ Median particle size of the cumulative volume distribution.
² pH range where resin can be operated without significant change in function.
³ pH range where resin can be subjected to cleaning-in-place (CIP) without significant change in function.

Technical specifications — Superdex Increase and Superose Increase standard columns



- **Typical pressure over the packed bed during operation:**
3 MPa, 50 bar, 435 psi (10/300 and 5/150)
2 MPa, 20 bar, 290 psi (3.2/300)
- **Column hardware pressure limit:**
5 MPa, 50 bar, 725 psi (10/300 and 3.2/300)
10 MPa, 100 bar, 1450 psi (5/150)

Parameter	Tricorn™ 10/300 GL	Tricorn 5/150 GL	3.2/300
Bed dimensions	10 mm × 300 mm	5 mm × 150 mm	3.2 mm × 300 mm
Approximate bed volume	24 mL	3 mL	2.4 mL
Recommended sample volume	25 to 500 µL	4 to 50 µL	4 to 50 µL
Recommended flow rate	0.8, 0.75, or 0.5 mL/min ¹	0.45 or 0.3 mL/min ²	0.075 or 0.04 mL/min ³
Max. operating flow rate	1.8, 1.6, 1.5, or 1.2 mL/min ⁴	0.75 mL/min	0.15 mL/min

¹ 0.8 for Superdex™ 30 and 75 Increase, 0.75 for Superdex 200 Increase, 0.5 for Superose™ 6 Increase

² 0.45 for Superdex 75 and 200 Increase, 0.3 for Superose 6 Increase

³ 0.075 for all Superdex Increase, 0.04 for Superose 6 Increase

⁴ 1.8 for Superdex 200 Increase, 1.6 for Superdex 75 Increase, 1.5 for Superose 6 Increase, 1.2 for Superdex 30 Increase

[Download data files >>](#)

7

**Superdex Increase and
Superose Increase columns for
larger volumes (0.5 to 8.5 mL)**

Superdex Increase and Superose Increase resins in HiScale columns for larger sample volumes

**We listened
to you!**

**Now available
on demand!**



If you like our new generation Superdex™ Increase and Superose™ Increase columns, but need to purify larger volumes of up to 8.5 mL, we can pack the resin in HiScale™ columns.

On-demand service¹. Contact your Cytiva representative.

¹ Custom Designed Products offer this as nonstandard products. This means that delivery could be longer than catalog products and we do not provide Instruction for use for these columns.

Technical specifications: HiScale prepacked columns



- **Max pressure over the packed bed during operation:**
2 MPa, 20 bar, 290 psi
- **Column hardware pressure limit:**
2 MPa, 20 bar, 290 psi

Parameter	HiScale™ 16/40	HiScale 26/40
Bed dimensions	16 mm × 400 mm	26 mm × 400 mm
Approximate bed volume	80 mL	212 mL
Recommended sample volume	Up to 3.2 mL	Up to 8.5 mL
Max. operating flow rate	~ 2 mL/min [†]	~ 4 mL/min [†]

* Max pressure over the packed bed is limited by the pressure limit of the column hardware.


† This value could differ between different resins and resin lots. Max. operating flow rate for the individual column is stated in the documentation included with each column.

High resolution is maintained when scaling up Superdex 200

Increase with different columns formats

Standard: 10/300 GL


Sample: 0.5 mL
Flow rate: 0.75 mL/min
Run time: ~ 32 min



Resolution	
Thyro.	NA
Aldo.	2.25
Ovalb.	1.55
Rnase	2.17

HiScale™ 16/40


Sample: 1.6 mL
Flow rate: 2 mL/min
Run time: ~ 40 min



Resolution	
Thyro.	NA
Aldo.	2.64
Ovalb.	1.95
Rnase	2.79

HiScale 16/40


Sample: 3.2 mL
Flow rate: 2 mL/min
Run time: ~ 40 min



Resolution	
Thyro.	NA
Aldo.	1.87
Ovalb.	1.08
Rnase	1.47

HiScale 26/40


Sample: 4.25 mL
Flow rate: 4 mL/min
Run time: ~ 53 min



Resolution	
Thyro.	NA
Aldo.	2.71
Ovalb.	2.01
Rnase	2.81

HiScale 26/40

Sample: 8.5 mL
Flow rate: 4 mL/min
Run time: ~ 53 min




Resolution	
Thyro.	NA
Aldo.	1.97
Ovalb.	1.17
Rnase	1.62

Sample: Thyroglobulin, aldolase, ovalbumin, and ribonuclease A in PBS buffer, pH 7.4

Antibody purification with Superdex 200 Increase in different columns with maintained high resolution

Standard: 10/300 GL


Sample: 0.5 mL IgG (2 mg/mL)



Flow rate: 0.75 mL/min

HiScale™ 16/40


Sample: 1.6 mL IgG (2 mg/mL)



Flow rate: 2 mL/min

HiScale 26/40

Sample: 4.25 mL IgG (2 mg/mL)



Flow rate: 4 mL/min


Sample: IgG 2 mg/mL in PBS-buffer, pH 7.4

8

Tips for successful size exclusion chromatography


Protect the packed bed in the column — decrease flow rate when working in cold room or with viscous liquids

Low temperature increases pressure



Pressure over the column at different flow rate and temperature on Superdex™ 200 Increase 10/300 GL column in water

High viscosity increases pressure



Pressure over the column when increasing the amount of ethanol (viscosity increases up to ~ 50% ethanol) on Superdex 200 Increase 5/150 GL column

Recommendation

When working with

- viscous liquids
- or at low temperature

Our recommendation is to lower the flow rate to avoid damaging the packed bed in the column.

Save your column by cleaning with NaOH every 10 to 20 SEC cycles


What to do?

Sodium hydroxide is a very efficient cleaning solution.

See chromatograms to the right for example of a column that could be refreshed.


[Learn more >>](#)

Fresh column



Superose™ 6 Increase 10/300 GL


Dirty column



Column used for a very long time without cleaning

- Resolution gradually decreases
- A gap was formed resulting in this very poor resolution

Refreshed column




The column performance could be restored

after adjusting the adapter and NaOH cleaning.

Tips for taking care of your precious SEC columns

More on [cytiva.com/ProteinResearch](https://www.cytiva.com/ProteinResearch)

Tips for maximizing your SEC column lifetime



Tips and hints for protecting your SEC column investment

READ MORE

Inspect your column to ensure you get it right



Have you ever had to deal with bubbles in the resin bed of your chromatography column? Or even a gap between the resin bed and the adapter?

READ MORE

Prevent SEC columns drying out



How to connect storage device to your SEC column to protect your column during storage

WITH A SYRINGE

WITH AN ÄKTA™ SYSTEM

Excluding air from the column



Tutorial showing how to connect SEC column drop-to-drop to ÄKTA systems

WATCH NOW

9

Summary

Adding a SEC polishing step will improve your proteins' purity



SEC gives highly size homogeneous samples.

Cytiva offers a variety of SEC resins and column formats to meet your recovery, speed, and purity needs.

It is critical to regularly maintain your SEC column to ensure high performance.

10

Useful tools to ensure successful SEC runs




Cytiva expertise made available for you at

[cytiva.com/SEC](https://www.cytiva.com/SEC)

Principles and methods
Cytiva Handbook >>

Help for column selection
Interactive selector >>

Ordering information
Cytiva online shop >>



Ordering information

Columns with SEC "Increase" resins

Product code	Product name
28990944	Superdex™ 200 Increase 10/300 GL
29148721	Superdex 75 Increase 10/300 GL
29219757	Superdex 30 Increase 10/300 GL
29091596	Superose™ 6 Increase 10/300 GL
29321903*	Superose 6 Increase HiScale™ 16/40
29321904*	Superose 6 Increase HiScale 26/40
29321905*	Superdex 200 Increase HiScale 16/40
29321906*	Superdex 200 Increase HiScale 26/40
29321907*	Superdex 75 Increase HiScale 16/40
29321908*	Superdex 75 Increase HiScale 26/40

HiLoad™ columns

Product code	Product name
28989335	HiLoad 16/600 Superdex 200 pg
28989333	HiLoad 16/600 Superdex 75 pg
28989331	HiLoad 16/600 Superdex 30 pg
29323952	HiLoad 16/600 Superose 6 pg†
28989336	HiLoad 26/600 Superdex 200 pg
28989334	HiLoad 26/600 Superdex 75 pg
28989332	HiLoad 26/600 Superdex 30 pg

HiPrep™ columns

Product code	Product name
17116501	HiPrep 16/60 Sephacryl™ S-100 HR
17116601	HiPrep 16/60 Sephacryl S-200 HR
17116701	HiPrep 16/60 Sephacryl S-300 HR
28935604	HiPrep 16/60 Sephacryl S-400 HR
28935606	HiPrep 16/60 Sephacryl S-500 HR
17119401	HiPrep 26/60 Sephacryl S-100 HR
17119501	HiPrep 26/60 Sephacryl S-200 HR
17119601	HiPrep 26/60 Sephacryl S-300 HR
28935605	HiPrep 26/60 Sephacryl S-400 HR
28935607	HiPrep 26/60 Sephacryl S-500 HR

† New prepacked column, launched Nov. 2018

* These products are available on-demand. [Contact your Cytiva representative](#)

11

Appendix 1 — Column selection

We help you to select the right SEC column!

What molecule do you need to purify?
Click on your choice.

**Peptides or
other small
biomolecules >>**

**Recombinant
tagged
proteins
>>**

**Antibodies
or
similar
>>**


**Macro-molecules,
viruses, protein
complexes
>>**

**Smaller
molecules**

Larger molecules

SEC columns for peptides or other small biomolecules

For more information, click on the column of your interest




¹These columns are produced on-demand.

[Contact your Cytiva representative](#)




SEC columns for recombinant tagged proteins

For more information, click on the column of your interest




SEC columns for monoclonal antibodies or other antibodies

For more information, click on the column of your interest



Macromolecules, viruses, large proteins, and protein complexes

For more information, click on the column of your interest



12

Appendix 2 — Application examples

Screening for optimal sample load of virus-like particles on HiPrep 16/60 Sephacryl S-500 HR column

Purpose of the study


Virus-like particles (VLP)s are used as vaccines.

To increase productivity, it is important to determine the maximum amount of feed per milliliter of chromatography resin that can be loaded to give an acceptable level of purification.

The effect of increased sample load was evaluated using HiPrep™ 16/60 Sephacryl™ S-500 HR column.


Resolution decreased with increased sample volume as expected. Product purity was analyzed by SDS-PAGE.

Sample load on SEC column



(A) Purification of a virus-like particle (VLP) by SEC using HiPrep 16/60 Sephacryl S-500 HR. Various sample volumes, previously purified on a Capto™ Q column, were loaded on the column.
(B) Enlargement of peaks presented in (A).

Purity check (SDS-PAGE)



SDS-PAGE analysis (reducing conditions, 4% to 12% polyacrylamide gel, Coomassie stained) of eluted pools where the arrows indicate surface proteins of the VLP (Mr 69 000, 54 000, and 27 000).

Purification of insulin using HiPrep 26/60 Sephacryl S-100 HR


Background

Insulin consists of two chains (A and B) held together by -S-S bonds.

When these links have been broken, the two chains can, despite small differences in molecular weight, be separated by SEC.


Purification of insulin

Column: HiPrep 26/60 Sephacryl S-100 HR
Sample: 1 mL of a mixture comprising bovine insulin chain A (M_r 2532) and chain B (M_r 3496), 0.5 mg/mL of each
Buffer: 50 mM sodium phosphate, 150 mM NaCl, pH 7.0
Flow rate: 2.0 mL/min




Purification of a tagged protein using HiLoad Superdex 200 pg column as polishing step

1. Capture: AC




Column: MBPTrap™ HP 5 mL
Sample: 15 mL of MBP-MCAD in *E. coli* lysate,
 $M_r \sim 85\,500$

2. Polishing: SEC




Column: HiLoad™ 16/600 Superdex™ 200 pg
Sample: 2 mL eluted fraction from AC

Purity check (SDS-PAGE)



Purification of an untagged protein using HiLoad Superdex 75 pg column as polishing step


1. Capture: IEX



Column: HiPrep™ Q XL 16/10

Sample: 40 mL of clarified *E. coli* extract with DAOCS


2. Intermediate purification: HIC



Column: SOURCE™ 15ISO, packed in HR column 16/10

Sample: 40 mL of DAOCS pool from IEX


3. Polishing: SEC



Column: HiLoad™ 16/600 Superdex™ 75 pg

Sample: 3 mL of concentrated DAOCS pool from HIC

Purity check (SDS-PAGE)



Improved purity of a his-tagged protein using Superdex 75 Increase 10/300 GL column

Purpose

A sample of a purified his-tagged protein that had oligomerized during storage and freeze-thawing was run on Superdex™ 75 Increase 10/300 GL column to remove aggregates.

Columns: Superdex 75 Increase 10/300 GL for preparative SEC (A)
Superdex 75 Increase 5/150 GL for analysis (B)

Samples: Concentrated, partially purified histidine-tagged protein (A)
Samples of SEC fractions and from pooled peaks (B)


Sample volumes: 500 μ L (A)
25 μ L (B)

Buffer: PBS

Flow rates: 0.2 mL/min (A)
0.45 mL/min (B)


System: ÄKTA pure 25 with 2 mm UV cell to extend the linear absorbance range

Purification



The final pool (fractions marked with red bar), with target protein, contained approximately 5 mg of target protein.

Purity check



Analysis of pool of both peaks from the preparative SEC purification.

SDS-PAGE separation of fluorescent prestained samples on an 8% to 18% gel.

Thank you



Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. ÄKTA, Capto, HiLoad, HiPrep, HisTrap, HiScale, HiTrap, MabSelect, MBPTrap, Sephacryl, SOURCE, Superdex, Superose, and Tricorn trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)