

BioProcess™ IC system

INLINE BUFFER PREPARATION

The BioProcess™ inline conditioning (IC) system (Fig 1) is a buffer management system for inline buffer preparation using concentrated single component stock solutions of acid, base, salt, and water for injection (WFI). The use of concentrates will significantly reduce buffer volumes, saving both floor space and tank volumes. In addition, many different buffers can be prepared from the same set of concentrates, streamlining buffer preparation even further. For accuracy in formulation and consistency between preparations it is possible to select different feedback modes featuring the built-in dynamic control functionality.

Intensifying large scale buffer management:

- Automated inline buffer preparation
- Quality rigor to ensure buffers meet specifications
- Integrated buffer preparation and chromatography operations

Freeing up resources

Traditional buffer preparation is a highly manual and resource-intensive activity in biomanufacturing. Large volumes of different buffers, formulated according to specific recipes, are required to produce a biopharmaceutical and due to the large quantities used, buffer management can easily become a bottleneck in production. In addition to high labor costs, there is also a risk of losing an entire buffer batch if it does not pass the release criteria.

As an extension of the concept of inline dilution (ILD), IC is a buffer management strategy that allows for buffer production on demand from concentrated single component stock solutions. Implementing IC can help to streamline the entire buffer preparation process and reduce manual handling by automating several steps, not only reducing the risk of human errors, but also making it possible to reassign personnel to other tasks.




Fig 1. BioProcess IC system for inline buffer preparation.

The BioProcess IC system has dedicated pump lines for WFI, acid, base and salt. It is also possible to include additional pumps for additives and/or feed (Fig 2). Each pump line has several inlet ports to allow different buffer families to be formulated with the same setup. This makes it possible to produce exact amounts of different buffers sequentially during a run (Fig 3). The process is controlled by a range of sensors; flow meters on each pump, and pH and conductivity sensors for accurate buffer preparation. The pH and conductivity sensors give feedback to the pumps to adjust the flow of each component. The Buffer Verification Function (BVF) in the UNICORN™ software monitors the pH and/or the conductivity of the mixed buffer in real time. The system will automatically bypass the column and redirect from outlets to waste if the detected values are outside set limits. This ensures that only the correct buffer reaches the column or collecting vessel, creating a robust and reliable production process.




Fig 2. Schematic flow diagram of a BioProcess IC system.




Fig 3. Preparation of buffers as well as strip and cleaning in place (CIP) solutions required for a mAb capture step. The arrows indicate preparation of formulations within specifications.

Preparing several buffers from the same stock solutions

While ILD requires one concentrate for each final buffer, a wide range of buffers can be prepared from a few single component stock solutions using IC (Fig 4). This will save time and space, and drastically reduce documentation and QC testing requirements.

Over 100 unique buffers of five different buffer systems have been formulated in an automated manner using the BioProcess IC system. Some examples are listed in Table 1 and 2, showing buffers with a great variety of pH values, buffer concentrations, salt concentrations, and additives have been formulated. Glycerol mixtures with acetate and dilution gradients with potassium phosphate have also been successfully prepared using the BioProcess IC system. Salt gradients are also accurately created by mixing the relevant stock solutions (Fig 5).

Monitoring real time data for full process control

The UNICORN control software is used to create recipes, control the system, and analyze each process run. All buffer characteristics are recorded, allowing operators to monitor the trending of those properties by looking at chromatography results in real-time.

There is also a possibility to integrate into a DeltaV automation system or to use a standard for open platform communication (OPC) for compatibility with other third-party system control softwares. This allows for real-time data to be shared effectively, for run data to be stored within the automation system, and for appropriate start and stop signals to be properly sequenced.




Fig 4. Different combinations of stock solutions and WFI generates different buffers in an automated manner.

Table 1. Example of phosphate buffer range that can be prepared using only three stock solutions

Input Stock solutions	Output		
	Buffer conc. (mM)	Salt conc. range (mM)	pH range
0.3 M NaH ₂ PO ₄	20	0–500	6.8–7.3
0.3 M Na ₂ HPO ₄	30	0–500	6.8–7.4
3.5 M NaCl	50	0–500	6.8–7.4

Table 2. Experimental data on buffers typically formulated with the BioProcess IC system from Cytiva

Buffer conc. (mM)	pH	NaCl conc. range (mM)
Phosphate buffers		
10*	6–8	0–6
20 [†]	6.5–7.4	0–500
25 [‡]	6–8	0–1000
30	6–7	0–1000
35	7.2	50
50 [§]	6–7.4	0–150
200	6.8	-
Sodium acetate buffers		
1.8	3.6	0–100
10	3.5	0–107
20	5.2	100
25	3.5–5.5	0–500
30	3.6	100
38	5.1	100
42.3	5.3	-
50	3.5–7.5	0–500
150	4	0–500
Sodium citrate buffers		
10 [¶]	4.5–5	0–300
25	3.2–4.2	0
100	3.5–5	-
Formic acid buffers		
390	2	2000
15	3.5	50–250
Tris buffers		
10	8.2	-
16	8	13
16.2	9	0–1000
20	7.5	20–500
25	7–9	0–1000
50**	8–9	50–1000
80	8	-
100 ^{††}	8	-

* Prepared with 0–2 M urea and 0–3 mM MES.

[†] Prepared with 0–3 M urea and 0%–0.04% Tween.

[‡] Prepared with 0–0.1 M (NH₄)₂SO₄.

[§] Prepared with 0–1 M (NH₄)₂SO₄ and 0%–0.1% Tween.

[¶] Prepared with 0–7 M urea.

** Prepared with 0–0.1 M (NH₄)₂SO₄.

^{††} Prepared with 0.1% Tween.




Fig 5. Gradient. Sodium phosphate buffer with NaCl—conductivity increases while total flow and pH remains unchanged. Acid and base are adjusted to maintain pH at the set point, while conductivity changes during the gradient run (sodium phosphate pH flow with gradient 30 mM sodium phosphate pH 6). Stock concentrations: 500 mM NaH₂PO₄, 350 mM Na₂HPO₄. Stock solution NaCl: 3.5M = gradient 0.1 to 1 M.).

Fulfilling quality rigor


Because buffers are used to maintain purification conditions, as well as to stabilize the product and preserve its functional characteristics, correct buffer formulation is crucial for success in bioproduction. The BioProcess IC system features built-in dynamic control to ensure accuracy in buffer formulation. Three modes of dynamic feedback control can be used: flow; pH and flow; and pH and conductivity.

- **Flow feedback:** A known buffer formulation is entered in the system control software. The software adjusts the flow rates of the specified stock solutions to achieve the desired formulation. This control mode is useful when the temperature is constant, and the stock solutions are accurate.
- **pH and flow feedback:** The method specifies target pH and the software adjusts the flow rates of the acid and base stock solutions to achieve the desired pH and buffer concentration in the final formulation.
- **pH and conductivity feedback:** The method specifies the target pH and conductivity, and the dynamic control functionality of the software uses feedback from flow, conductivity, and pH sensors to adjust flow rates of the stock solutions to achieve the desired conductivity and pH. In this control mode both the temperature and the concentrations of the stock solution can vary without affecting the pH and conductivity accuracy of the final buffer formulation.

Feedback control minimizes buffer variability and compensates for variability in the stock solutions, and thus improves process robustness and quality of products (Figs 6 and 7).

The BioProcess IC system is equipped with two dual pH modules, controlling and monitoring. Each module includes both CPS11E and CPS61E pH probes. With dual flow paths, the system can effectively run buffers where salt is present and those where it's not without the pH sensor suffering from salt memory.

CPS11E is used for measuring buffers without salt; CPS61E requires some conductivity to measure properly and is used when salt is present. CPS61E is more resistant to sodium hydroxide and can be left in the flow path during column CIP.



Buffer	pH	Conductivity (mS/cm)
20 mM sodium phosphate, 150 mM NaCl, pH 7.4 buffer	7.41 ± 0.05	17.3 ± 0.3
20 mM sodium phosphate, 500 mM NaCl, pH 7 buffer	7.44 ± 0.05	17.3 ± 0.2
20 mM sodium phosphate, 150 mM NaCl, pH 7.4 buffer	7.40 ± 0.04	17.3 ± 0.2
20 mM sodium phosphate, 500 mM NaCl, pH 7 buffer	7.00 ± 0.03	47.2 ± 0.2
20 mM sodium phosphate, 500 mM NaCl, pH 7 buffer	7.01 ± 0.04	47.3 ± 0.3
20 mM sodium phosphate, 500 mM NaCl, pH 7 buffer	6.99 ± 0.02	47.3 ± 0.3

Fig 6. Robustness, reliability and reproducibility: The image shows three runs of preparation 20 mM sodium phosphate, 150 mM NaCl, pH 7.4 followed by preparation of 20 mM sodium phosphate, 500 mM NaCl, pH 7. The switch between the two buffers takes 2 minutes.




Fig 7. Starting point: Three different concentrations of base and acid stock solutions and seven different combinations were tested. All combinations reached target values of pH and conductivity.

Several sensors throughout the flow path

When using pH feedback the system uses two sensors, one for controlling and one for monitoring. The accuracy of each pH sensor in the system is ± 0.10 pH units, which is also the accuracy of the buffer. By adding the variances, the monitoring pH can at the most be ± 0.15 pH units from the target pH.

The accuracy of the conductivity monitors between 0.1 and 300 mS/cm can be found in the technical specifications table.

When selecting pH and conductivity feedback control, the flow rates will be adjusted to meet the pH and conductivity target values.

The accuracy of the flow rate for each pump is $\pm 1\%$ of the reading or 0.1 L/h for small pumps and 1.0 L/h for larger pumps (whichever is greater), meaning that a smaller pump has a smaller error. To optimize flow accuracy the smaller pumps are used for the highly concentrated acid and base stock solutions.

Integrating buffer production and chromatography

The BioProcess IC system has the functionality required to operate as a chromatography unit, allowing direct connection to a chromatography column to make it possible to deliver buffers directly on to the column without the need for storage in bags or tanks (Fig 8).

Since a column can be connected directly to the BioProcess IC system, the buffer mixing and the chromatography run can be done in an integrated manner. With all steps controlled by the same software, all relevant data for buffer characteristics and for the run itself are recorded at the same time, resulting in one batch record for both operations, simplifying documentation.




Fig 8. Inline conditioning system can be used as (A) central buffer preparation station or (B) as an integral part of a chromatography or filtration unit operation.

In combination with the column option it is possible to add a titration pump, placed after the column connection. This pump is designed to enable pH adjustment of the product eluted from the column, using flow or pH feedback independently of the other system pumps.

Technical specifications

	BioProcess IC system 10 mm	BioProcess IC system 1 in.	BioProcess IC system 1½ in.
Flow rate range	60–600 L/h	200–2000 L/h	1000–5000 L/h
Dimensions (W × D × H)*	1450 × 2119 × 1911 mm	1712 × 2437 × 2196 mm	1720 × 2900 × 2100 mm
Operating max. pressure	6 barg (max 3 barg at 40°C–60°C)	6 barg (max 3 barg at 40°C–60°C)	5 barg (max 3 barg at 40°C–60°C)
Pressure sensor			
Range	0.1–6 barg	0.1–6 barg	0.1–5 barg
Accuracy	± 0.12 barg	± 0.12 barg	± 0.12 barg
Total flow accuracy	± 2% or 0.5 L/h	± 2% or 2 L/h	± 2% or 4 L/h
Temperature sensor			
Range	2°C–60°C		
Accuracy	± 2°C		
Conductivity sensor			
Range	0.1 to 49.99 mS/cm, precision: ± 2% or 0.05 mS/cm 50 to 99.99 mS/cm, precision: ± 2.5% 100–300 mS/cm, precision: ± 4%		
UV sensor - fixed or variable/ multiple wavelength			
Range	0–2 AU		
Accuracy	linearity ± 5%		
pH sensors			
Range	2–12		
Accuracy within IC application range [†]	± 0.10 pH units		
Lowest conductivity limit	0.6 mS/cm (CPS11E) and 0.5 mS/cm (CPS61E)		
Wetted parts			
Main piping size	10 mm	1 in.	1½ in.
Piping	Polypropylene (PP)	PP	PP pre-air trap and stainless steel (SS) post-air trap
Flange type	Hygienic tube clamp connection (TC)		
Gaskets	EPDM		
Flow meter	UNS N06022 (Hastelloy)		
Air trap	Ethylene propylene rubber (EPDM), polyamid (PA), PP	EPDM, PA, PP	EPDM, PA, SS
Air sensor	PP	PP	PP pre-air trap and SS post-air trap
pH sensor	Glass, Polytetrafluoroethylene (PTFE)		
Conductivity sensor	Polyetheretherketone (PEEK)		
Temperature sensor			
Integrated in conductivity sensor	PEEK		
Safety	UNS N06022		
UV sensor	PEEK, EPDM, quartz, titanium		
Pressure sensor	UNS N06022		
Pump chamber	PP, PTFE, PEEK, oxide ceramics, EPDM	PP, PTFE, PEEK, oxide ceramics, EPDM	ASTM 316L (1.4435), EPDM, PP, Santoprene
Material requirements for wetted parts	USP class VI, 21 CFR 177, animal origin free or in compliance with EMA/410/01		
Wetted and pressure retaining materials	PP, 316L (EN 1.4435), glass, quartz glass, EPDM, PEEK, FPM, FKM, PTFE, oxide ceramics, titanium grade 2, Santoprene		

AxiChrom™ column compatibility

AxiChrom column range (inner diameter in mm)	300–600	600–1200	1000–1600
--	---------	----------	-----------

System control

Software	UNICORN 7.10 or higher
----------	------------------------

*Dimensions excluding HMI (rotatable), stack light (foldable) and feed pump.

[†] Buffer capacity of at least 5 mmol NaOH per pH unit for the controlling pH electrode.



Ordering information

For ordering information and further technical details and buffer list evaluation, please contact your regional Cytiva sales office.

The system can be optimized to fit the process and facility by either integrating it into the chromatography unit operation or by using it as a standalone buffer preparing unit. To further intensify the buffer preparation process, ready-made HyClone™ stock solutions that can be used, eliminating all raw material handling and manual preparation of stock solutions. Advantages of using single component stock solutions as opposed to concentrated buffers for dilution, are that the common ion effect will not limit the maximum concentration to the same extent and there is no need to handle pH and conductivity shifts due subsequent dilution.

For more information, please visit: [cytiva.com/buffer-management](https://www.cytiva.com/buffer-management)

[cytiva.com/buffer-management](https://www.cytiva.com/buffer-management)

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva. AxiChrom, BioProcess, HyClone, and UNICORN are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Any other trademarks are the property of their respective owners.

Any use of UNICORN is subject to Cytiva Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

© 2026 Cytiva

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY11780-12May26-DF

