

Install Injection Valve V9-J

Instructions

Scope

This document describes how to install Injection Valve **V9-J** on ÄKTA™ systems.

Instruction

Follow the steps below to replace or install a module or module panel. Module panels must be used in positions not occupied by modules.

Note: *The illustrations show the principle of how to replace and install a module. The position of the module on the instrument and the type of module used will depend on the module being installed.*



CAUTION

Disconnect power. Always switch off power to the ÄKTA instrument before replacing any of its components, unless otherwise stated in the user documentation.

Step	Action
1	Disconnect power from the instrument by using the instrument power button.
2	If a module is to be replaced, loosen the tubing connectors and remove the tubing from the existing module.

Step	Action
------	--------

- | | |
|---|--|
| 3 | Loosen the screw in the module or module panel with a Torx™ T20 screwdriver. |
|---|--|



- | | |
|---|---|
| 4 | Remove the module or module panel and disconnect the cable at the back. |
|---|---|




Step	Action
------	--------

- | | |
|---|--|
| 5 | Connect the cable to the new module or module panel. |
|---|--|




- | | |
|---|--|
| 6 | Insert the module or module panel and fasten it with a Torx T20 screwdriver. |
|---|--|




Ports and flow paths

The tables below describe the ports of and different flow paths through Injection Valve **V9-J**.


	Port	Description
	Col	Outlet to the Column Valve or column.
	LoopF	Port for connection of a loop. Used to fill a loop.
	Syr	Port for connection of a syringe.
	W2	Loop waste
	W1	Pump waste
	LoopE	Port for connection of a loop. Used to empty a loop into the flow path.
	Pump	Inlet from the Pump.

The image below shows the flow paths in Injection Valve **V9-J** when (A) loading a sample to a loop, (B) injecting a sample onto the column, and (C) washing the system.



Connect tubing

The illustration shows the tubing connecting Injection Valve **V9-J**.



Part	Function
1	Sample loop or Superloop, connected to ports LoopF and LoopE
2	Tubing 4 , connects port Pump to the mixer
3	Tubing W1 , directs pump waste to a waste bottle through port W1
4	Tubing W2 , directs loop waste to a waste bottle through port W2
5	Tubing 5 , connects port Col to a column or column valve



CAUTION

Fasten the waste tubing. Make sure that tubing is securely fastened to the waste ports **W**, **W1**, and **W2**.



CAUTION

Make sure there is a sample loop or a stop plug in place in the loop positions in the injection valve to avoid liquid spurting out during valve turns.

Node ID

All modules have a pre-configured Node ID according to their function. The Node ID is used by the instrument to distinguish between several units of the same type.

The Node ID for the Injection Valve **V9-J** is 25. Check the Node ID of the module when troubleshooting. See the *User Manual* of your instrument for more information.

Note: *The function of a module is defined by the module type and the Node ID, not by its physical position.*

Page intentionally left blank



cytiva.com/akta

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

ÄKTA and UNICORN are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Torx is a trademark of Acument Intellectual Properties, LLC.

All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

Any use of UNICORN is subject to Cytiva Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

29325017 AD V:4 11/2020