


Introduction to protein purification

Part 2: Protein purification




Chromatography is the most common technique for purifying proteins

What is chromatography?

Separation of components using differential affinities of the components for a mobile phase and for a stationary phase, through which they pass

"Drawing with colour"
(Tswett, inventor of chromatography, 1903)




Liquid chromatography

- Buffer (mobile phase) pass through a medium (resin) containing chromatography beads (stationary phase)
- Different bead properties can be used to purify different target proteins
 - For example, a bead with small pores can be used to purify a small target from larger molecules in the same sample



Chromatography principle




Protein properties used in chromatography

Protein property	Technique
Biospecific affinity (tags)	Affinity chromatography (AC)
Hydrophobicity	Hydrophobic interaction chromatography (HIC)
Net charge	Ion exchange chromatography (IEX)
Size/volume	Size exclusion chromatography (SEC) (also called gel filtration [GF])

We offer handbooks that are a great reference for practical tips and in-depth information about common chromatography techniques. To download, visit our [handbook page](#).


Schematics of common chromatography techniques

AC



[Watch AC principle video](#)

IEX




[Watch IEX principle video](#)

SEC



[Watch SEC \(or GF\) principle video](#)

HIC



[Watch HIC principle video](#)

Comparison: Manual vs automated purification

Manual purification

Little training or start-up time required

Easy to do parallel runs for increased throughput



Automated purification

Automation — convenience


Gradient elution — high resolution

Reproducibility

Documentation



Summary protein purification



Efficient purification of proteins is required both in research and in production of biopharmaceuticals.

Chromatography is a technique to purify proteins based on properties like size and net charge.

Protein purification can be performed manually or automatically.

Life science online courses

Browse our learning courses and grow your skills

Self-paced learning

No hidden cost

Up-to-date courses

Catalog

Search...

CATEGORIES

- ☐ Bioprocess manufacturing
- ☐ Buffer preparation and management
- ☐ Cell culture
- ☐ Cell processing for cell and gene therapy
- ☐ Chromatography software and systems
- ☐ Chromatography techniques
- ☐ Cryopreservation/cord blood banking
- ☐ Genomics
- ☐ Point of care cell and gene therapy
- ☐ Process development in bioprocessing
- ☒ Protein purification
- ☐ Surface plasmon resonance
- ☐ Viral vector and vaccine manufacturing
- ☐ Western blotting



Load More

To browse our learning courses visit
<https://www.cytivalifesciences.com/education/online-learning/courses>

Thank you



Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. ÄKTA is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)