

Liquid-Liquid Mixing at 500 L Scale

Mixing system: Pall Magnetic Mixer system

Mixing biocontainer: 500 L Round mixing biocontainer (off-center impeller)

Application mixing type: Liquid-liquid

The Pall Magnetic Mixer platform is a scalable single-use mixing system. The heart of this system is a mixing biocontainer incorporating an innovative bottom-mounted, magnetically-driven impeller capable of providing efficient, high speed, high-torque mixing in powder-liquid and liquid-liquid applications.

Introduction

Liquid-liquid mixing is a common requirement in biopharmaceutical processing. Biotech applications involving stress-sensitive materials (e.g. some proteins) often require the use of reduced impeller mixing speeds to minimize risk of yield loss. In this experiment, mixing performance was evaluated at a power input of approximately 1 W/m^3 which is apparently very soft as it presents no significant visible movement of the liquid surface; sodium chloride (NaCl) aqueous solution was added to water at the 250 L and 500 L scales, with the time to reach homogeneity determined by conductivity measurement.

Experimental

A Magnetic Mixer 500 L biocontainer was installed in a 500 L round mixing tank and then filled with 250 L of room-temperature water via a top-mounted liquid inlet port. Mixing was then initiated with an impeller speed of 45 rpm (equivalent to $\sim 1 \text{ W/m}^3$ power input) and then 500 mL of 20% (w/w) NaCl solution was added to the top of the biocontainer over a duration of ~ 1 second on the opposite side from the impeller just above the liquid surface. The test was repeated five times. Solution homogeneity was monitored using two conductivity sensors located at the top and bottom of the biocontainer; homogenization was deemed complete when both sensors had reached a stable reading within $\pm 5\%$ of the conductivity signal shift. Figure 2 shows the test setup. The test was repeated at the 500 L scale with an impeller speed of 60 rpm.





Figure 1
500 L tank and drive

Figure 2

Test setup



Results

The homogenization times are reported in Table 1. As an example, the conductivity vs time chart for the 250 L tests is shown in Figure 3. With the impeller speeds tested, at the 250 L scale (50% of tank capacity) it took in average 116 seconds to reach homogeneity, and at the 500 L scale (100% of tank capacity) it took in average 198 seconds to reach homogeneity. The longest homogenization time (at the 500 L scale) was 238 seconds.


Table 1

Homogenization time in seconds (s) in a 500 L round tank. Each test was repeated 5 times

Filling Level	Speed (rpm)	Test 1 (s)	Test 2 (s)	Test 3 (s)	Test 4 (s)	Test 5 (s)
50%	45	100	136	113	113	118
100%	60	222	174	170	187	238

Figure 3

Conductivity evolution during the 250 L test. The blue lines indicate the addition times and the red lines indicate the calculated mixing times



Conclusion

The mixing capabilities of the Magnetic Mixer system facilitate rapid liquid-liquid mixing of water-like solutions at the 250 L and 500 L scales, even when low impeller speeds are required due to product stress-sensitivity concerns. Under the conditions evaluated, homogenization was typically completed within 2 and 4 minutes respectively.



Visit us on the Web at www.pall.com/biotech
Contact us at www.pall.com/contact

Corporate Headquarters

Port Washington, NY, USA
+1 800 717 7255 toll free (USA)
+1 516 484 5400 phone

European Headquarters

Fribourg, Switzerland
+41 (0)26 350 53 00 phone


Asia-Pacific Headquarters

Singapore
+65 6389 6500 phone

Pall Corporation has offices and plants throughout the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

IF APPLICABLE Please contact Pall Corporation to verify that the product conforms to your national legislation and/or regional regulatory requirements for water and food contact use.

© Copyright 2021, Pall Corporation. Pall  are trademarks of Pall Corporation. ® Indicates a trademark registered in the USA.

USD 3482 GN21.07867
May 2021