



Biotech

USD3285

Pall Kleenpak™ Nova Capsules with Supor® EX Grade ECV Filter Membrane



High Capacity, High Flow Rate Filter Capsules

Supor EX grade ECV filters are high capacity, high flow rate 0.2 µm sterilizing grade filters validated for the retention of *Brevundimonas diminuta* (ATCC 19146) at a challenge level of 10⁷ colony forming units (CFU) per cm² membrane. A unique polyethersulfone (PES) membrane pairing ensures Supor EX grade ECV filters deliver highly efficient filtration of cell harvest material, process intermediates, growth media, buffers and final bulk biological process fluids. Supor EX grade ECV filters also deliver unrivalled efficiency with more viscous process streams.

Available in capsules with filter areas of 0.1 m² to 3 m², in in-line and T-style configurations, and with a variety of inlet and outlet options, Kleenpak Nova capsules with Supor ECV membranes have a size and format to meet your process needs today and can scale as your process grows to meet your future needs without a change of membrane or materials.

Filtration. Separation. Solution.SM

Key Features and Benefits

- ▶ Asymmetric PES membrane arranged in laid over pleat geometry for excellent flow rates and high capacity per unit membrane
- ▶ Low protein and preservative binding
- ▶ Available as 0.1 m², 0.5 m², 1 m², 2 m² and 3 m² filter capsules for coverage of a broad range of fluid volumes
- ▶ Encapsulated format for minimized cleaning and low installation costs
- ▶ Easy integrity testing through Staubli vent valve
- ▶ Fluid compatibility over a broad pH range for use in multiple applications

Quality Standards

A Certificate of Test provided with every filter documents the following:

- ▶ Filters are manufactured in a controlled environment
- ▶ Manufactured for use in conformance with cGMP
- ▶ Materials of construction meet USP<88> Biological Reactivity Test *In Vivo* for Class VI-121°C plastics
- ▶ Q.C lot release tests
 - Fabrication Integrity – each filter has successfully passed a Forward Flow integrity test correlated to bacterial retention per modified ASTM F838-05
 - Lot sample bacterial retention tests with *Brevundimonas diminuta* (ATCC 19146) at a challenge level of 10⁷ CFU/cm²
 - Effluent quality for cleanliness, TOC and water conductivity, pH and pyrogens per USP standard methods

Specifications

Materials of Construction

Filter Membrane	Dual-layer hydrophilic polyethersulfone
Support and Drainage Layers	Polypropylene
Core, End Caps, Fin, Adapter	Polypropylene
Cage	Polypropylene with TiO ₂ whitener ¹
O-rings	Silicone elastomer
Sealing Technology	Thermal bonding without adhesives
Housing Bowl	Polypropylene
Housing Head	Polypropylene with TiO ₂ whitener ¹

¹ TiO₂ is an insoluble inorganic mineral filler that does not contribute to organic extractables

Effective Filtration Area

Format	NP1L	NP5L	N(P/T)6	N(P/T)7	N(P/T)8
Nominal Filter Length	1 in.	5 in.	10 in.	20 in.	30 in.
Effective Filtration Area	0.1 m ²	0.52 m ²	1.04 m ²	2.08 m ²	3.12 m ²

Operating Parameters²

Format	NP1L	NP5L	N(P/T)6	N(P/T)7	N(P/T)8
Maximum Differential Pressure	4.1 bar (60 psi) at 20 °C	5.0 bar (72.5 psi) at 40 °C	3.0 bar (43.5 psi) at 40 °C	3.0 bar (43.5 psi) at 40 °C	3.0 bar (43.5 psi) at 40 °C
Maximum Operating Pressure	5.2 bar (75.4 psi) at 20 °C 4.0 bar (58 psi) at 40 °C	5.0 bar (72.5 psi) at 40 °C	3.0 bar (43.5 psi) at 40 °C	3.0 bar (43.5 psi) at 40 °C	3.0 bar (43.5 psi) at 40 °C

² In compatible fluids which do not soften, swell or adversely affect the filter or its materials of construction

Sterilization³

Format	NP1L	NP5L	N(P/T)6	N(P/T)7	N(P/T)8
Autoclave (G option only)	5 x 60 minute cycles at 125 °C 1 x 60 minute cycle at 135 °C	5 x 60 minute cycles at 125 °C	3 x 60 minute cycles at 125 °C	3 x 60 minute cycles at 125 °C	3 x 60 minute cycles at 125 °C
Gamma Irradiation (G option only)	Maximum of 50 kGy	Maximum of 50 kGy	Maximum of 50 kGy	Maximum of 50 kGy	Maximum of 50 kGy

³ Pre-sterilized Kleenpak Nova capsules must not be re-sterilized. Kleenpak Nova capsules must not be sterilized in-situ by passing steam under pressure

Typical Non Volatile Residue (NVR) Extractables in Water at 20 °C⁴

Format	NP1L	NP5L	N(P/T)6	N(P/T)7	N(P/T)8
Post 1 x 60 minute autoclave cycle at 125°C	< 15 mg per capsule	< 75 mg per capsule	< 150 mg per 10 in. capsule	< 150 mg per 10 in. capsule	< 150 mg per 10 in. capsule
Post gamma irradiation at 50 kGy	< 2 mg per capsule	< 15 mg per capsule	< 25 mg per 10 in. capsule	< 25 mg per 10 in. capsule	< 25 mg per 10 in. capsule

⁴ Tested on elements after 24 hr extraction

Integrity Test Values (Air test gas, water wet)*

Format	NP1L	NP5L	N(P/T)6	N(P/T)7	N(P/T)8
Maximum Allowable Forward Flow at 2760 mbar (40 psi)	2.4 mL/min	11 mL/min	21 mL/min	37 mL/min	54 mL/min

* Values correct at 20 °C. Contact Pall for multi-element integrity test values or other fluid values and recommended test procedures

Clean Water Flow Rates

Format	NP1L	NP5L	N(P/T)6	N(P/T)7	N(P/T)8
Clean Water Flow at 100 mbar Differential Pressure	1.4 L/min	8.5 L/min	17 L/min	34 L/min	51 L/min

Nominal Dimensions

In-Line	NP1L	NP5L	NP6	NP7	NP8
Maximum Diameter Including Valves	154 mm (6.1 in.)	154 mm (6.1 in.)	154 mm (6.1 in.)	154 mm (6.1 in.)	154 mm (6.1 in.)
Length with Hose Barb Inlet/ Outlet	–	275 mm (10.8 in.)	397 mm (15.6 in.)	644 mm (25.4 in.)	895 mm (35.2 in.)
Length with Sanitary Inlet / Outlet	128 mm (5.0 in.)	213 mm (8.4 in.)	332 mm (13.1 in.)	584 mm (23.0 in.)	834 mm (32.8 in.)
T-Style	–	–	NT6	NT7	NT8
Maximum Diameter Including Valves	–	–	240 mm (9.5 in.)	240 mm (9.5 in.)	240 mm (9.5 in.)
Length	–	–	349 mm (13.7 in.)	598 mm (23.5 in.)	848 mm (33.4 in.)

Ordering Information

N []		[]		UECVP		[]		[]		[]	
Code	Description	Code	Nominal Filter Length	Nominal Filter Area	Code	Inlet / Outlet Configurations	Code	Format	Code	Vent and Drain Options	
P	In-line	1L*	1 in.	0.1 m ²	1	1 – 1½ in. sanitary flange inlet and outlet	Blank	Non-sterile, autoclavable	Blank	Staubli vent and stepped hose barb drain	
T	T-style	5L*	5 in.	0.5 m ²	9	1 in. (25 mm) single barb hose barb inlet and outlet	G	Non-sterile, autoclavable / irradiatable	A	Staubli vent and drain	
		6	10 in.	1 m ²	19	1 – 1½ in. sanitary flange inlet and 1 in. (25 mm) single barb hose barb outlet	S	Presterilized by irradiation			
		7	20 in.	2 m ²	6*	½ in. (13 mm) single barb hose barb inlet and outlet					
		8	30 in.	3 m ²	16*	1 – 1½ in. sanitary flange inlet and ½ in. (13 mm) single barb hose barb outlet					
					1H**	1 – 1½ in. sanitary flange inlet and outlet, with ½ in. sanitary port on inlet					
					1H9**	1 – 1½ in. sanitary flange inlet and 1 in. (25 mm) single barb hose barb outlet, with ½ in. sanitary port on inlet					

*In-line only

** T-Style only



Biotech

Corporate Headquarters

Port Washington, NY, USA
+1.800.717.7255 toll free (USA)
+1.516.484.5400 phone

European Headquarters

Fribourg, Switzerland
+41 (0)26 350 53 00 phone

**Asia-Pacific Headquarters
Singapore**

+65 6389 6500 phone

Visit us on the Web at www.pall.com/biotech

Contact us at www.pall.com/contact

International Offices

Pall Corporation has offices and plants throughout the world in: Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, New Zealand, Norway, Philippines, Poland, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, United Kingdom, and Vietnam. Distributors in all major industrial areas of the world. To locate the Pall office or distributor nearest you, visit www.pall.com/contact.

The information provided in this literature was reviewed for accuracy at the time of publication. Product data may be subject to change without notice. For current information consult your local Pall distributor or contact Pall directly.

© 2018, Pall Corporation. Pall, , Kleenpak, and Supor are trademarks of Pall Corporation. ® indicates a trademark registered in the USA. TM indicates a common law trademark. *Filtration.Separation.Solution.* is a service mark of Pall Corporation.