

# Mustang™ Q

## CHROMATOGRAPHY CAPSULES

### High-throughput single-use capsules for process scale anion exchange chromatography

Mustang™ Q single-use products are designed to capture negatively charged biomolecules in downstream bioprocessing by anion exchange. These products will be used either for the initial capture of the target biomolecule in a purification process or for the removal of negatively charged contaminants in subsequent purification steps.

Mustang Q membrane is an anion exchange support with pendant quaternary amine functional groups in a cross-linked polymeric coating. The 0.8 µm pores in Mustang Q membranes are large enough to allow large biomolecules to access all the binding sites by direct fluid convection. This produces a very high dynamic capacity for very high molecular weight molecules (e.g., DNA, plasmid), or even particles as large as viruses. This is in contrast to conventional chromatographic beads with diffusive pores that typically have lower throughputs and dynamic binding capacities.

Combining this high-capacity Q membrane with a pleated design has resulted in the Mustang Q single-use capsule filters. Each Mustang Q element contains 16 layers of pleated membrane that will bind negatively charged biomolecules. These units are specifically designed to be single use to eliminate cleaning and cleaning validation. They are available in capsule or cartridge format to accommodate the entire range of volumes used in the biopharmaceutical processes. The 16-layer construction has been held constant from the laboratory product (0.35 mL) to the industrial scale. The ease of linear scale-up enables a shortened development time.

### DNA removal efficiency

Mustang Q units will typically reduce the DNA level in a solution of sonicated calf thymus DNA from 10 µg/mL to less than 10 pg/mL in one pass. This is a removal efficiency of 6 logs. The presence of host cell protein or other contaminants, and the flow rates used, may influence test sensitivity, efficiency of binding and capacity. Performance capability in a specific solution must be determined and validated by the user. For more information, contact us.



Fig 1. The range of Mustang Q single-use capsule filters.

### Features and benefits

- Binding efficiency – negatively charged biomolecules are readily bound in a single pass
- Speed – high flow rates enable the processing of large volumes in less than one shift
- Scalability – a range of sizes accommodate different volumes and capacities needed in biopharmaceutical processing
- Convenience – ready to use and autoclavable
- Flexibility – available in capsule or cartridge format
- Cost – lower operating costs and capital investment compared to conventional columns that need validated packing and cleaning

## High quality standards

- Manufactured to high quality assurance standards in accordance with BS EN ISO9002 1997
- Membrane lots tested for dynamic protein binding capacity, and peak position using standard proteins
- Identified by a lot number and a specific serial number for traceability of manufacturing history, satisfying QC/QA requirements
- Supplied with certificate of analysis to confirm quality and quality control
- Meets USP biological reactivity tests *in vivo* in accordance with USP class VI 50°C and all materials listed in the drug master file submitted to the FDA

## Comprehensive validation

- Extensive validation to ensure consistent and reliable performance
- Comprehensive validation guide available on request

## Technical specifications

### Material of construction

|  | <b>Novasip™ capsule<br/>CLxMSTGQP1</b> | <b>Kleenpak™<br/>Nova capsule<br/>NPxMSTGQP1</b> |
|--|--|--|
| Membrane                                     | Modified hydrophilic polyethersulfone  | Modified hydrophilic polyethersulfone            |
| Membrane support and drainage layer assembly | Polypropylene                          | Polypropylene                                    |
| Endcaps, core and cage                       | Polypropylene                          | Polypropylene                                    |
| Housing bowl                                 | Polyetherimide                         | Polypropylene                                    |
| Housing head                                 | Polyetherimide with TiO <sub>2</sub>   | Polypropylene with TiO <sub>2</sub>              |
| O-rings                                      | Silicone elastomer                     | Silicone elastomer                               |

Note: For CLxMSTGQP1, x = M05 for a unit with a bed volume of 10 mL, 3 for a unit with a bed volume of 60 mL. For NPxMSTGQP1, x = 6 for a unit with a bed volume of 260 mL, 7 for a unit with a bed volume of 520 mL, 8 for a unit with a bed volume of 780 mL.

### Capsules operating characteristics <sup>(1)</sup>

|                                 |   |
|---------------------------------|---|
| Maximum operating pressure      | 3.0 bar (43.5 psi) at 40°C (104°F)          |
| Maximum differential pressure   | 3.0 bar (43.5 psi) at 40°C (104°F)          |
| Maximum sanitization conditions | 1 M NaOH for 30 min for one cycle only      |
| Maximum autoclave conditions    | 121°C (250°F) for 30 min for one cycle only |


<sup>(1)</sup> With fully compatible fluids that do not soften, swell or adversely affect the capsule or its materials of construction  
1 bar = 100 kPa

## Available sizes and typical dynamic binding capacity

| Bed volume (mL) | Typical BSA binding capacity (mg) | Flow rate (L/min/bar Delta P) <sup>(2)</sup> | Product code                |
|-----------------|-----------------------------------|--|-----------------------------|
| 0.35            | 17                                | 6.400 mL/min (0.430 mL/min/psi)              | MSTG18Q16                   |
| 10              | 500                               | 0.245 (17 mL/min/psi)                        | CLM05MQTGSP1                |
| 60              | 3000                              | 1.250 (86 mL/min/psi)                        | CL3MSTGQP1                  |
| 260             | 13 000                            | 3.900 (269 mL/min/psi)                       | NP6MSTGQP1/<br>AB1MSTGQ7PH4 |
| 520             | 26 000                            | 7.800 (538 mL/min/psi)                       | NP7MSTGQP1/<br>AB2MSTGQ7PH4 |
| 780             | 39 000                            | 11.700 (807 mL/min/psi)                      | NP8MSTGQP1/<br>AB3MSTGQ7PH4 |

<sup>(2)</sup> For liquid of 1 cP. For other liquids, divide the flow rate by the viscosity expressed in cP.

## DNA breakthrough curves on Mustang Q capsules



**Fig 2.** DNA breakthrough curves on Mustang Q capsules.

The DNA breakthrough curves above demonstrate the high breakthrough capacity of Mustang Q capsules. After preconditioning and equilibration with 25 mM Tris-HCl pH 8.0, a solution of 140 µg/mL, 210 µg/mL and 338 µg/mL of herring sperm DNA in 25 mM Tris-HCl pH 8.0 was pumped respectively through the 10 mL, 60 mL and 260 mL capsules at a flow rate of 10-20 CV/min. The absorbance at 260 nm was continuously monitored downstream of the capsule.

Scaling down unit characteristics

| Materials of Construction     |  |
|-------------------------------|--|
| Membrane                      | Modified hydrophilic polyethersulfone                                    |
| Housing                       | 316L stainless steel   |
| Support frit                  | 316B stainless steel   |
| Housing clamp                 | 300 series stainless steel housing                                       |
| O-rings                       | Silicone elastomer   |
| Nominal Dimensions            |  |
| Coin device                   | 18 mm outer diameter   |
| Housing                       | 51 mm height   |
| Fittings                      | M6 male to 1⁄8 in. OD tubing and<br>M6 male to luer lock female fittings |
| Operating Characteristics     |  |
| Maximum temperature           | 40°C (104°F)   |
| Maximum differential pressure | 17 bar (247 psi)   |

Nominal dimensions for capsules

| Capsule type                               | CLM05MSTGEP1     | CL3MSTGEP1       | NP6MPSTGEP1       | NP7MSTGEP1        | NP8MSTGEP1        |
|--|------------------|------------------|-------------------|-------------------|-------------------|
| Maximum diameter including valves          | 123 mm (4.8 in.) | 123 mm (4.8 in.) | 154 mm (6.1 in.)  | 154 mm (6.1 in.)  | 154 mm (6.1 in.)  |
| Length with 38 mm (1½ in.) sanitary flange | 84 mm (3.3 in.)  | 157 mm (6.2 in.) | 333 mm (13.1 in.) | 581 mm (22.9 in.) | 831 mm (32.7 in.) |
| Bed volume (mL)                            | 10               | 60               | 260               | 520               | 780               |

Ordering information

| Product   | Packaging     | Product code |
|---|---------------|--------------|
| Mustang Q coin  | 10/box        | MSTG18Q16    |
| Stainless steel coin housing assembly   | 1/box         | MSTG18H16    |
| Mustang coin housing frits  | 2 frits/box   | MSTG18H16F   |
| Mustang coin housing gaskets  | 2 gaskets/box | MSTG18H16G   |
| Mustang Q capsule, 10 mL, connections 1 to 1½ in. compatible sanitary connection  | 1/box         | CLM05MSTGQP1 |
| Mustang Q capsule, 60 mL, connections 1 to 1½ in. compatible sanitary connection  | 1/box         | CL3MSTGQP1   |
| Mustang Q capsule, 260 mL, connections 1 to 1½ in. compatible sanitary connection | 1/box         | NP6MSTGQP1   |
| Mustang Q capsule, 520 mL, connections 1 to 1½ in. compatible sanitary connection | 1/box         | NP7MSTGQP1   |
| Mustang Q capsule, 780 mL, connections 1 to 1½ in. compatible sanitary connection | 1/box         | NP8MSTGQP1   |
| Mustang Q cartridge, 260 mL, code 7   | 1/box         | AB1MSTGQ7PH4 |
| Mustang Q cartridge, 520 mL, code 7   | 1/box         | AB2MSTGQ7PH4 |
| Mustang Q cartridge, 780 mL, code 7   | 1/box         | AB3MSTGQ7PH4 |

For more information, contact us.

cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.  
Kleenpak, Mustang and Novasip are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.  
Any other third-party trademarks are the property of their respective owners.  
© 2024 Cytiva  
For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

