MabSelect™ SuRe 70

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Matrix	Highly cross-linked agarose
Ligand	Alkaline-stabilized Protein A-derived (<i>E. coli</i>)
Ligand coupling	Single-point attachment
Coupling chemistry	Epoxy
Particle size d_{50v} *	50 μ m
DBC QB _{10%} [†]	> 70 g trastuzumab per liter resin at 6 min RT
Recommended maximum operating flow velocity [‡]	220 cm/h (20 cm bed height), 400 cm/h (10 cm bed height)
Chemical stability	Stable in aqueous buffers commonly used in Protein A chromatography
pH stability, operational [§]	3 to 12
pH stability, CIP [¶]	2 to 13.7
CIP stability [¶]	0.1 to 0.5 M NaOH
Delivery conditions	20% ethanol

* Median particle size of the cumulative volume distribution

[†] DBC at 10% breakthrough by frontal analysis at a RT of 6 min (100 cm/h) in a HiScreen™ column at 10 cm bed height

[‡] Base matrix packed in an AxiChrom™ 1000 column with 100 cm i.d. at 20 cm bed height, operating pressure up to 3 bars, using buffers with the same viscosity as water at 20°C

[§] pH range where resin can be operated without significant change in function

[¶] pH range where resin can be subjected to cleaning-in-place (CIP), 0.1 M NaOH without significant change in function, 0.5 M NaOH can be used with some loss in capacity

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
그림 2. DBC, QB_{10%} for MabSelect SuRe 70, MabSelect SuRe LX, and MabSelect Prisma resins in a HiScreen column at 10 cm bed height using trastuzumab at different RTs.

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
그림 3. A study of accelerated alkaline stability for MabSelect SuRe 70 resin showed that the resin can be used for 200 cycles at 0.1 M NaOH and for 100 cycles at 0.5 M NaOH with more than 80% relative remaining capacity.

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표 3. Typical protein A capture protocol

Step	Column volume (CV)	RT (min)/linear flow velocity (cm/h) at 20 cm bed height	Buffer
Equilibration	3	6/200	20 mM sodium phosphate, pH 7.4 + 150 mM NaCl*
Load	70% to 80% of QB ^{10%}	6/200	As required
Wash 1	5	6/200	20 mM sodium phosphate, pH 7 + 500 mM NaCl
Wash 2	1	6/200	50 mM sodium acetate, pH 6
Elution step	3	6/200	50 mM sodium acetate, pH 3.5
Strip	2	6/200	100 M acetic acid, pH 2.9
Cleaning in place (CIP)	3 (15 min contact time in upflow)	5.2/220	0.1 M NaOH to 0.5 M NaOH
Re-equilibration	3 to 5 (in upflow)	5.2/220	20 mM sodium phosphate, pH 7.4 + 150 mM NaCl*
Only after last run/for storage	4	7.5/160	20% ethanol

* For lab-scale use, 20 mM sodium phosphate, pH 7.4 + 500 mM NaCl can be used also for equilibration and reequilibration to reduce the number of different buffers.

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
그림 4. Normalized protein A consumable cost using the protein MabSelect SuRe 70 protein A resin for 5 batches and utilizing the full lifetime, respectively, at 6 min and 2 min RT. Resin and column packing cost per gram of mAb produced is calculated based on the number of cycles used. Using short RT will reduce resin cost at low resin lifetime utilization and increase costs at high lifetime utilization.


그림 5. Duration of the MabSelect SuRe 70 Protein A chromatography step decreases with shorter RT, despite more cycles being run. Resin utilization does not impact the times required to run a batch.

표 2. Process parameters

Parameter	6 min RT		2 min RT	
	5 batches	Full lifetime	5 batches	Full lifetime
Resin utilization	5 batches	Full lifetime	5 batches	Full lifetime
Volume clarified cell culture feed	2000 L			
Titer	4.25 g/L			
Column diameter	45 cm			
Column bed height	20 cm		10 cm	

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Product	Size	Product code
HiTrap MabSelect SuRe 70	1 × 1 mL	17542411
HiTrap MabSelect SuRe 70	5 × 1 mL	17542412
HiTrap MabSelect SuRe 70	1 × 5 mL	17542413
HiTrap MabSelect SuRe 70	5 × 5 mL	17542414
HiScreen MabSelect SuRe 70	1 × 4.7 mL	17542410
MabSelect SuRe 70	25 mL	17542401
MabSelect SuRe 70	200 mL	17542402
MabSelect SuRe 70	1 L	17542403
MabSelect SuRe 70	5 L	17542404
MabSelect SuRe 70	10 L	17542405
ReadyToProcess MabSelect SuRe 70 columns	Contact us for sizes and product codes	

Please [contact us](#) to request samples of MabSelect SuRe 70 resins.

cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva. AxiChrom, BioProcess, Capto, Fast Trak, HiScreen, HiTrap, MabSelect, MabSelect Prisma, MabSelect SuRe, and ReadyToProcess are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Any other trademarks are the property of their respective owners. The Danaher trademark is a proprietary mark of Danaher Corporation. © 2025 Cytiva Any use of software may be subject to one or more end-user license agreements, a copy of, or notice of which, are available on request. For local office contact information, visit cytiva.com/contact