

# ReadyToProcess™ hollow fiber cartridges

## BIOPROCESS™ FILTRATION

ReadyToProcess™ hollow fiber cartridges are designed to help you simplify and speed up drug development and production by minimizing device preparation steps such as wetting and sterilizing. Because of the single-use nature of the cartridges, you can eliminate the need for cleaning protocols and cleaning validation. ReadyToProcess™ hollow fiber cartridges incorporate ReadyMate™ Disposable Aseptic Connectors (DAC), which allow you to install and aseptically connect the cartridges to other components within Cytiva's ReadyToProcess™ product platform, including BioProcess™ containers and tubing assemblies.

For your peace of mind, ReadyToProcess™ hollow fiber cartridges are supplied preconditioned for reduced conductivity and total organic carbon (TOC) level and processed to minimize endotoxin and bioburden levels. The cartridges are packaged and sterilized by gamma irradiation.

### ReadyToProcess™ hollow fiber cartridges offer the following benefits:

- Preconditioned for reduced conductivity and TOC level
- Gamma irradiated for sterility to an assurance level of  $10^{-6}$
- Processed to minimize endotoxin and bioburden levels
- Supplied fully wet for quick start-up
- Incorporate ReadyMate™ DACs for easy installation and aseptic connections to ReadyCircuit™ products
- Enable optimized process throughput

ReadyToProcess™ hollow fiber cartridges are available in a wide variety of ultrafiltration molecular weight cut-offs and microfiltration pore sizes, to allow scalability from laboratory to production scale.

### Compatible systems

ReadyToProcess™ hollow fiber cartridges are compatible with the following filtration systems:

- Xcellerex™ Automated Perfusion System (APS)
- ÄKTA™ brand readyflux™ automated filtration system



Fig 1. ReadyToProcess™ hollow fiber cartridges.

### ReadyToProcess™ product platform

ReadyToProcess™ hollow fiber cartridges are part of the ReadyToProcess™ product platform, a suite of products that brings a comprehensive offering of ready-to-use cell culture, purification, and fluid management products such as:

- ReadyToProcess™ columns and ÄKTA™ ready automated chromatography system
- WAVE Bioreactor™ systems and Cellbag™ bioreactor chambers
- ReadyCircuit™ bioprocessing containers and tubing assemblies
- ReadyKart hardware and portability solutions
- ReadyMate™ Disposable Aseptic Connectors

# Applications

Typical applications for ReadyToProcess™ hollow fiber cartridges are:

- Aseptic clarification and purification of vaccines, monoclonal antibodies, recombinant proteins, and plasmids
- Aseptic cell processing
- Operations in environments where terminal sterilization is not feasible
- Across all stages of production, from preclinical to commercial stages
- Processing in areas where the risk of cross-contamination is high

# Specifications

Specifications for ReadyToProcess™ hollow fiber cartridges are listed in Tables 1 to 3.

**Table 1.** Operating specifications

<b>Ultrafiltration</b>	
Temperature range	4°C to 80°C
Maximum feed pressure	< 10°C: 5.2 barg (75 psig)
	10°C to 25°C: 4.5 barg (65 psig)
	25°C to 80°C: 3.4 barg (50 psig)
Maximum transmembrane pressure	10 000 to 30 000 MWCO:
	< 10°C: 4.1 barg (60 psig)
	10°C to 25°C: 3.4 barg (50 psig)
	25°C to 50°C: 3.1 barg (45 psig)
	50°C to 80°C: 2.4 barg (35 psig)
	50 000 to 750 000 MWCO:
< 10°C: 3.4 barg (50 psig)	
	10°C to 25°C: 3.1 barg (45 psig)
	25°C to 50°C: 2.4 barg (35 psig)
	50°C to 80°C: 1.7 barg (25 psig)
<b>Microfiltration</b>	
Temperature range	0.1 to 0.2 µm: up to 80°C
	0.45 to 0.65 µm: up to 50°C
Maximum feed pressure	At 25°C:
	0.1 µm: 2.1 barg (30 psig)
	0.2 µm: 1.7 barg (25 psig)
	0.45 µm: 1 barg (15 psig)
	0.65 µm: 1 barg (15 psig)
Maximum transmembrane pressure	At 25°C:
	0.1 µm: 1.4 barg (20 psig)
	0.2 µm: 1 barg (15 psig)
	0.45 µm: 0.7 barg (10 psig)
	0.65 µm: 0.7 barg (10 psig)

**Table 2.** Technical specifications

<b>Size</b>	<b>Cartridge diameter</b>
3, 3X2	0.9 cm (0.375 in)
4, 4X2	1.9 cm (0.75 in)
5, 6	3.2 cm (1.25 in)
8, 9	5.1 cm (2 in)
35, 55	7.6 cm (3 in)
45, 65	10.8 cm (4.25 in)
152	16.8 cm (6.6 in)
<b>Size</b>	<b>Path length</b>
3, 4, 5, 8, 35, 45	30 cm (12 in)
3X2, 4X2, 6, 9, 55, 65, 152	60 cm (24 in)
<b>Fluid connections (excludes size 152)</b>	
Feed/retentate	ReadyMate™ connector
Permeate	ReadyMate™ connector
<b>Materials of construction</b>	
Cartridge housing, fittings, fiber	Polysulfone
Cartridge sleeve	Polypropylene
Fiber net	Polyethylene
Fiber potting	Epoxy
Permeate tubing	C-Flex™ (3, 3X2) AdvantaPure™ (4, 4X2, 5, 6, 8, 9, 35, 55, 45, 65)
ReadyMate™	Polycarbonate, silicone
Gaskets	Silicone
<b>Regulatory conformance</b>	
EMA/410/01	Compliant
Biological Reactivity Test, In Vivo per USP <88> Class VI	Compliant
21CFR177 Indirect Food Additives	Compliant (except Size 3 and 3X2)
L929 MEM Elution Test – ISO 10993-5 (Cytotoxicity)	Compliant
Hemolysis – rabbit blood (direct contact) – ISO 10993-4	Compliant

**Table 3.** Nominal feed stream flow rates in mL/min\*

Nominal lumen i.d. (mm)	Shear rate ~ 2000 s <sup>-1</sup>	Shear rate ~ 4000 s <sup>-1</sup>	Shear rate ~ 8000 s <sup>-1</sup>	Shear rate ~ 16 000 s <sup>-1</sup>
<b>Size 3 and 3X2</b>				
0.5	60	120	250	500
0.75	100	200	400	800
1	150	300	600	1200
<b>Size 4 and 4X2</b>				
0.5	300	600	1200	2400
0.75	400	800	1500	3000
1	600	1200	2500	5000
<b>Size 5 and 6</b>				
0.5	1100	2100	4300	8600
0.75	1400	2800	5600	11 200
1	2000	4000	8000	16 000
<b>Size 8 and 9</b>				
0.5	2700	5400	10 600	21 500
0.75	4400	8800	18 000	35 000
1	6100	12 200	24 500	49 000
<b>Size 35 and 55</b>				
0.5	6600	13 200	26 000	52 000
0.75	10 000	20 000	40 000	80 000
1	15 000	30 000	60 000	120 000
<b>Size 45 and 65</b>				
0.5	14 000	28 000	55 000	111 000
0.75	19 000	39 000	77 000	154 000
1	31 000	61 000	122 000	245 000
<b>Size 152</b>				
1	70 000	140 000	280 000	560 000

\*Based on a fluid viscosity of 1 cP.

# Ordering information

## ReadyToProcess™ hollow fiber ultrafiltration cartridges

Model number	Pore size (MWCO)	Fiber i.d. (mm)	Membrane surface area (m <sup>2</sup> )	Membrane surface area (ft <sup>2</sup> )	Nominal flowpath length (cm)
RTPUFP-10-C-3X2MS	10 000	0.5	0.029	0.31	60
RTPUFP-30-C-3X2MS	30 000	0.5	0.029	0.31	60
RTPUFP-50-C-3X2MS	50 000	0.5	0.029	0.31	60
RTPUFP-100-C-3X2MS	100 000	0.5	0.029	0.31	60
RTPUFP-300-C-3X2MS	300 000	0.5	0.029	0.31	60
RTPUFP-500-C-3X2MS	500 000	0.5	0.029	0.31	60
RTPUFP-750-C-3X2MS	750 000	0.5	0.029	0.31	60
RTPUFP-750-E-3X2MS	750 000	1	0.023	0.24	60
RTPUFP-10-C-4MS	10 000	0.5	0.065	0.7	30
RTPUFP-30-C-4MS	30 000	0.5	0.065	0.7	30
RTPUFP-50-C-4MS	50 000	0.5	0.065	0.7	30
RTPUFP-100-C-4MS	100 000	0.5	0.065	0.7	30
RTPUFP-300-C-4MS	300 000	0.5	0.065	0.7	30
RTPUFP-500-C-4MS	500 000	0.5	0.065	0.7	30
RTPUFP-750-C-4MS	750 000	0.5	0.065	0.7	30
RTPUFP-750-E-4MS	750 000	1	0.042	0.45	30
RTPUFP-10-C-4X2MS	10 000	0.5	0.14	1.5	60
RTPUFP-30-C-4X2MS	30 000	0.5	0.14	1.5	60
RTPUFP-50-C-4X2MS	50 000	0.5	0.14	1.5	60
RTPUFP-100-C-4X2MS	100 000	0.5	0.14	1.5	60
RTPUFP-300-C-4X2MS	300 000	0.5	0.14	1.5	60
RTPUFP-500-C-4X2MS	500 000	0.5	0.14	1.5	60
RTPUFP-750-C-4X2MS	750 000	0.5	0.14	1.5	60
RTPUFP-750-E-4X2MS	750 000	1	0.085	0.9	60
RTPUFP-10-C-5S	10 000	0.5	0.2	2.1	30
RTPUFP-30-C-5S	30 000	0.5	0.2	2.1	30
RTPUFP-50-C-5S	50 000	0.5	0.2	2.1	30
RTPUFP-100-C-5S	100 000	0.5	0.2	2.1	30
RTPUFP-300-C-5S	300 000	0.5	0.2	2.1	30
RTPUFP-500-C-5S	500 000	0.5	0.2	2.1	30
RTPUFP-750-C-5S	750 000	0.5	0.2	2.1	30
RTPUFP-750-E-5S	750 000	1	0.12	1.3	30
RTPUFP-10-C-6S	10 000	0.5	0.48	5.2	60
RTPUFP-30-C-6S	30 000	0.5	0.48	5.2	60
RTPUFP-50-C-6S	50 000	0.5	0.48	5.2	60
RTPUFP-100-C-6S	100 000	0.5	0.48	5.2	60
RTPUFP-300-C-6S	300 000	0.5	0.48	5.2	60
RTPUFP-500-C-6S	500 000	0.5	0.48	5.2	60
RTPUFP-750-C-6S	750 000	0.5	0.48	5.2	60
RTPUFP-750-E-6S	750 000	1	0.28	3	60
RTPUFP-10-C-8S	10 000	0.5	0.53	5.7	30
RTPUFP-30-C-8S	30 000	0.5	0.53	5.7	30
RTPUFP-50-C-8S	50 000	0.5	0.53	5.7	30
RTPUFP-100-C-8S	100 000	0.5	0.53	5.7	30
RTPUFP-300-C-8S	300 000	0.5	0.53	5.7	30
RTPUFP-500-C-8S	500 000	0.5	0.53	5.7	30
RTPUFP-750-C-8S	750 000	0.5	0.53	5.7	30
RTPUFP-750-E-8S	750 000	1	0.36	3.9	30

## ReadyToProcess™ hollow fiber ultrafiltration cartridges continued

Model number	Pore size (MWCO)	Fiber i.d. (mm)	Membrane surface area (m <sup>2</sup> )	Membrane surface area (ft <sup>2</sup> )	Nominal flowpath length (cm)
RTPUFP-10-C-9S	10 000	0.5	1.15	12.5	60
RTPUFP-30-C-9S	30 000	0.5	1.15	12.5	60
RTPUFP-50-C-9S	50 000	0.5	1.15	12.5	60
RTPUFP-100-C-9S	100 000	0.5	1.15	12.5	60
RTPUFP-300-C-9S	300 000	0.5	1.15	12.5	60
RTPUFP-500-C-9S	500 000	0.5	1.15	12.5	60
RTPUFP-750-C-9S	750 000	0.5	1.15	12.5	60
RTPUFP-750-E-9S	750 000	1	0.84	9	60
RTPUFP-500-C-35S	500 000	0.5	1.35	14.5	30
RTPUFP-750-C-35S	750 000	0.5	1.35	14.5	30
RTPUFP-500-C-55S	500 000	0.5	3.25	35	60
RTPUFP-750-C-55S	750 000	0.5	3.25	35	60
RTPUFP-10-C-45S	10 000	0.5	3.5	37	30
RTPUFP-750-C-45S	750 000	0.5	3.5	37	30
RTPUFP-10-C-65S	10 000	0.5	6.1	66	60
RTPUFP-30-C-65S	30 000	0.5	6.1	66	60
RTPUFP-100-C-65S	100 000	0.5	6.1	66	60
RTPUFP-500-C-65S	500 000	0.5	6.1	66	60
RTPUFP-750-C-65S	750 000	0.5	6.1	66	60
RTPUFP-10-E-65S	10 000	1	4.4	47	60
RTPUFP-30-E-65S	30 000	1	4.4	47	60
RTPUFP-750-E-65S	750 000	1	4.4	47	60
RTPUFP-10-E-152 <sup>†</sup>	10 000	1	8.3	89.3	60
RTPUFP-30-E-152 <sup>†</sup>	30 000	1	8.3	89.3	60
RTPUFP-750-E-152 <sup>†</sup>	750 000	1	8.3	89.3	60

<sup>†</sup>Size 152's do not carry a sterile claim due to their non-self-contained design but are exposed to the same radiation dose as the remainder of the product line.

## ReadyToProcess™ hollow fiber microfiltration cartridges

Model number	Pore size (µm)	Fiber i.d. (mm)	Membrane surface area (m <sup>2</sup> )	Membrane surface area (ft <sup>2</sup> )	Nominal flowpath length (cm)
RTPCFP-1-E-3MS	0.1	1	0.011	0.12	30
RTPCFP-2-E-3MS	0.2	1	0.011	0.12	30
RTPCFP-4-E-3MS	0.45	1	0.011	0.12	30
RTPCFP-6-D-3MS	0.65	0.75	0.012	0.13	30
RTPCFP-1-E-3X2MS	0.1	1	0.023	0.24	60
RTPCFP-2-E-3X2MS	0.2	1	0.023	0.24	60
RTPCFP-4-E-3X2MS	0.45	1	0.023	0.24	60
RTPCFP-6-D-3X2MS	0.65	0.75	0.027	0.29	60
RTPCFP-1-E-4MS	0.1	1	0.042	0.45	30
RTPCFP-2-E-4MS	0.2	1	0.042	0.45	30
RTPCFP-4-E-4MS	0.45	1	0.042	0.45	30
RTPCFP-6-D-4MS	0.65	0.75	0.046	0.5	30
RTPCFP-1-E-4X2MS	0.1	1	0.085	0.91	60
RTPCFP-2-E-4X2MS	0.2	1	0.085	0.91	60
RTPCFP-4-E-4X2MS	0.45	1	0.085	0.91	60
RTPCFP-6-D-4X2MS	0.65	0.75	0.095	1.02	60

## ReadyToProcess™ hollow fiber microfiltration cartridges continued

Model number	Pore size (µm)	Fiber i.d. (mm)	Membrane surface area (m <sup>2</sup> )	Membrane surface area (ft <sup>2</sup> )	Nominal flowpath length (cm)
RTPCFP-1-E-5S	0.1	1	0.12	1.3	30
RTPCFP-2-E-5S	0.2	1	0.12	1.3	30
RTPCFP-4-E-5S	0.45	1	0.12	1.3	30
RTPCFP-6-D-5S	0.65	0.75	0.16	1.7	30
RTPCFP-1-E-6S	0.1	1	0.28	3	60
RTPCFP-2-E-6S	0.2	1	0.28	3	60
RTPCFP-4-E-6S	0.45	1	0.28	3	60
RTPCFP-6-D-6S	0.65	0.75	0.37	4	60
RTPCFP-1-E-8S	0.1	1	0.36	3.9	30
RTPCFP-2-E-8S	0.2	1	0.36	3.9	30
RTPCFP-4-E-8S	0.45	1	0.36	3.9	30
RTPCFP-6-D-8S	0.65	0.75	0.41	4.4	30
RTPCFP-1-E-9S	0.1	1	0.84	9	60
RTPCFP-2-E-9S	0.2	1	0.84	9	60
RTPCFP-4-E-9S	0.45	1	0.84	9	60
RTPCFP-6-D-9S	0.65	0.75	0.93	10	60
RTPCFP-6-D-35S	0.65	0.75	1	10.8	30
RTPCFP-6-D-55S	0.65	0.75	2.5	27	60
RTPCFP-1-E-45S	0.1	1	2.5	27	30
RTPCFP-2-E-45S	0.2	1	2.5	27	30
RTPCFP-4-E-45S	0.45	1	2.5	27	30
RTPCFP-6-D-45S	0.65	0.75	2.65	28.5	30
RTPCFP-1-E-65S	0.1	1	4.4	47	60
RTPCFP-2-E-65S	0.2	1	4.4	47	60
RTPCFP-4-E-65S	0.45	1	4.4	47	60
RTPCFP-6-D-65S	0.65	0.75	4.9	53	60
RTPCFP-1-E-152 <sup>†</sup>	0.1	1	8.2	89.3	89.3
RTPCFP-2-E-152 <sup>†</sup>	0.2	1	8.2	89.3	89.3
RTPCFP-4-E-152 <sup>†</sup>	0.45	1	8.2	89.3	89.3

<sup>†</sup>Size 152's do not carry a sterile claim due to their non-self-contained design but are exposed to the same radiation dose as the remainder of the product line.

## cytiva.com/bioprocess

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate. ÅKTA, BioProcess, Cellbag, ReadyCircuit, readyflux, ReadyMate, ReadyToProcess, WAVE Bioreactor, and Xcellerex are trademarks of Global Life Sciences Solutions USA LLC or an affiliate.

ReadyMate is covered by US patent number 6,679,529 B2 owned by Johnson & Boley Holdings, LLC and licensed to Cytiva.

AdvantaPure is a trademark of NewAge Industries, Inc. C-Flex is a trademark of Saint-Gobain Performance Plastics. Any other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

CY11887-29Apr21-DF

