

Envirochek™ and Envirochek HV Sampling Capsules and Accessories

For the concentration and recovery of *Cryptosporidium* and *Giardia* cysts from source or finished water

- Typically greater than 70% recovery of target organisms. Capsules replace traditional string wound filter method and increase capture and recovery of organisms by 15 times.
- Capsule method results in a smaller pellet after elution, requiring less analysis time and cost.
- Disposable design eliminates false positives from cross-contamination.
- The 1 µm pore size membrane eliminates false negatives for 100% collection of *Cryptosporidium* and *Giardia*.
- Simple and safe capsules eliminate contact with organisms. No filter holder assembly or cleaning required, and the self-contained capsule means that the potentially contaminated filter element does not need to be handled or cut apart.
- Saves time by allowing the processing of multiple samples at the same time. Process up to eight samples at one time with Pall's Laboratory Shaker.



- Capsule method does not affect the viability of the organisms like harsh chemicals used in other techniques.
- Serialized for traceability.
- Approved by regulatory agencies worldwide, including U.S. EPA methods 1622 and 1623, United Kingdom DWI standard operating protocols, and ISO/DIS 15553.

Applications

Envirochek capsules are used for collection and recovery of *Cryptosporidium* oocysts and *Giardia* cysts in both finished and source water, including surface water, municipal water supplies, samples in containers, or wells.

- **Envirochek capsules** are validated and listed in U.S. EPA Methods 1622 and 1623, and used for sampling source water for *Cryptosporidium* and *Giardia*. Together, the EPA method and Envirochek capsule present a major improvement over the previous string wound cartridge method. The Envirochek capsule is also listed in ISO/DIS 15553.
- **Envirochek HV capsules** are designed for sampling up to 1,000 liters or more of treated water and validated for up to 50 liters of source water. The Envirochek HV incorporates a track etched membrane designed to process high volumes of treated water while maintaining high recovery characteristics (see Tables 1 and 2) and meeting U.S. EPA requirements. Envirochek HV capsules are also approved by the United Kingdom DWI standard operating protocols for monitoring drinking water for *Cryptosporidium*.

In the Envirochek sampling capsule method, the water sample to be tested is passed through the capsule either by sampling at the source or by taking a "grab" sample back to the laboratory and filtering it on the bench top. The capsule is then filled with an elution solution, placed on a laboratory shaker, and vigorously shaken to elute any captured oocysts and cysts. Up to eight capsules can be processed at once using Pall's Laboratory Shaker. The elution solution is decanted and centrifuged to a pellet for further examination by the user's method of choice.

Table 3 shows the time savings from using the patented Envirochek sampling capsule and method. This simple process saves time and money by reducing steps and allowing the technician to attend to other tasks. The disposable capsule eliminates the arduous cleaning required by other products that use reusable filter housings.

Specifications

Materials of Construction

Envirochek™ Capsule (PN 12110, 12107)
 Filter Media: Supor® membrane, hydrophilic polyethersulfone
 Housing: Polycarbonate
 Filter Support Material: Polypropylene
 End Caps: Green vinyl
 Adhesive: Urethane

Envirochek HV Capsule (PN 12096, 12097, 12098, 12099)
 Filter Media: Polyester, hydrophilic membrane
 Housing: Polycarbonate
 Filter Support Material: Polypropylene
 End Caps: Blue vinyl
 Adhesive: Urethane

Effective Filtration Area

1,300 cm²

Dimensions

Length: 21.6 cm (8.5 in.)
 Diameter: 6.1 cm (2.4 in.)

Inlet/Outlet Connections

12.7 mm (1/2 in.) straight hose barb

Elution Capacity

Minimum of 127 mL

Maximum Operating Pressure

Envirochek Capsule (PN 12110, 12107)
 2.1 bar (30 psid) at 21 °C (70 °F)

Envirochek HV Capsule (PN 12099, 12097, 12096)
 4.2 bar (60 psid) at 21 °C (70 °F)

Sterilization

Provided non-sterile

Performance

Table 1

IPR and MS/MSD Validation Data for *Cryptosporidium* Recovery Using Method 1622 and the Envirochek HV Capsule with 50 Liters Source Water Samples

Sample Description	Lab	Sample Turbidity (ntu)	Spike Dose (#)	Mean % Rec.	RSD or RPD
IPR	1	< 0.1	95.9	44.9	9.5%
IPR	2	0.05	95.9	67.6	17.0
IPR	3	0.03	95.9	60.5	18.8
IPR Acceptable Range				24-100	< 55%
Matrix Blank	1	10.6	0		
MS1	1	9.4	95.9		
MS2	1	10.1	95.9	29.2	0.0%
Matrix Blank	2	2.1	0		
MS1	2	2.1	95.9		
MS2	2	2.1	95.9	72.0	20.3
Matrix Blank	3	1.8	0		
MS1	3	1.8	95.9		
MS2	3	1.8	95.9	56.3	11.1
IPR Acceptable Range				13-111	< 61%

IPR = Initial Precision and Recovery

MS/MSD = Matrix Spike/Matrix Spike Duplicate

RSD = Relative Standard Deviation

RPD = Relative Percent Difference

Table 2

IPR Tier 1 Validation Data for *Cryptosporidium* Recovery Using Method 1622 and the Envirochek HV Capsule for 1,000 Liters Finished Drinking Water Samples

Sample Description	Sample Turbidity (ntu)	Spike Dose (#)	% Rec.	Mean % Rec.	RSD or RPD
Reagent Blank	< 0.1	0			
IPR2	< 0.1	99.3	68.5		
IPR2	< 0.1	99.3	58.4		
IPR3	< 0.1	99.3	60.4		
IPR4	< 0.1	99.3	61.4	62.2	7.0
IPR Acceptable Range				13-143	< 67%

Table 3


Time Savings

Step	Approximate Time to Process Eight Samples (Minutes)	
	Envirochek Capsule	Other Major Protocol
Set-up	5	50
Elution	55	240
Concentration	65	240
Cleaning Equipment	0	320
Total Time	125	850
Time/Test (minutes)	16	106

Figure 1

Envirochek™ HV Capsule Produces a Smaller Pellet Size When Sampling Turbid Water, Reducing Time and Cost for Analysis

A smaller pellet following concentration of the sample can reduce cost by eliminating the need for additional ISM kits and fewer slides. Testing time is reduced because additional slides do not have to be prepared and visually analyzed. The Envirochek HV capsule is designed to trap 100% of *Cryptosporidium* and *Giardia* while allowing finer particulate to pass through the filter that would add to pellet size. Each capsule is integrity tested to assure no bypass, providing highly accurate recovery rates as well as time and cost savings.



Laboratory Shaker

For effective elution of *Cryptosporidium* oocysts and *Giardia* cysts off Envirochek sampling capsules

- Meets U.S. EPA requirements for Methods 1622 and 1623 for recovery of organisms.
- Processes up to eight Envirochek sampling capsules at one time.
- Accurate and consistent shaking speed with a digital read-out.
- Simulates human wrist action to save time, reduce fatigue, and provide consistent results.
- CE marked for European community use.
- UL-listed power cord also included.



Applications

- Specifically designed for use with Pall's Envirochek sampling capsules for the elution of *Cryptosporidium* oocysts and *Giardia* cysts. Add an elution buffer to the capsule and the shaking action elutes the captured

organisms off the membrane surface for further concentration and identification.

- Shaker accommodates up to eight Envirochek sampling capsules at one time.

Specifications

Timer

Includes 30-minute timer, calibrated at five minute intervals, or may be run continuously.

Shaking Speed

Variable speed regulation from 100 to 700 osc/minute (1,400 vibrations). Speed is monitored within 1 rpm of setpoint by digital tachometer.

Stroke

8 mm

Set-up

8 clamps included for attachment of Envirochek sampling capsules.

Type of Mounting

Benchtop. Suction cup feet prevent creeping.

Load Carrying Capacity

Accepts loads up to 4 kg (9 lb.)

Power Requirements

PN 4821: 115V, 50/60 Hz
PN 4822: 230V, 50/60 Hz
3-wire UL listed power cord

Power Consumption

80 watts

Dimensions

Height: 23 cm (9.1 in.)
Width: 77 cm (30.3 in.)
Depth: 27 cm (10.6 in.)
Weight: 11 kg (24.2 lb.)

Laureth-12

Laureth-12 is a component required in the elution solution. Available in paste form, one 50-gram bottle per package. One bottle has sufficient volume for approximately 200 tests. For the complete formulation, refer to the Envirochek™ HV sampling capsule protocol, PN 33210, or the Pall web site at www.pall.com/lab. The formulation can also be found in U.S. EPA Methods 1622 and 1623.



Ordering Information

Product No.	Description	Packaging	Product No.	Description	Packaging
12110	Envirochek sampling capsule	1/pkg	4820	Laureth-12 Paste, 50 g bottle	1/pkg
12107	Envirochek sampling capsule, bulk pack, individually bagged	25/pkg	4821	Laboratory shaker, 115 V, 50/60 Hz	1/pkg
12099	Envirochek HV sampling capsule	1/pkg	4822	Laboratory shaker, 230 V, 50/60 Hz	1/pkg
12098	Envirochek HV sampling capsule, bulk pack, individually bagged	25/pkg			
12097	Envirochek HV sampling capsule, for U.K. DWI	1/pkg			
12096	Envirochek HV sampling capsule, for U.K. DWI	25/pkg			

Reference Literature

- The Envirochek HV Capsule: Recovering *Cryptosporidium* from High Volume Source and Finished Drinking Water Samples, Technical Article, PN 33249
- Envirochek HV Sampling Capsule Protocol, PN 33210



Life Sciences

Pall Life Sciences
600 South Wagner Road
Ann Arbor, MI 48103-9019 USA

800.521.1520 toll free in USA
734.665.0651 phone
734.913.6114 fax

Australia – Lane Cove, NSW
Tel: 02 9428-2333
1800 635-082 (in Australia)
Fax: 02 9428-5610

Austria – Wien
Tel: 043-1-49 192-0
Fax: 0043-1-49 192-400

Canada – Ontario
Tel: 905-542-0330
800-263-5910 (in Canada)
Fax: 905-542-0331

Canada – Québec
Tel: 514-332-7255
800-435-6268 (in Canada)
Fax: 514-332-0996

800-808-6268 (in Canada)

China – P. R., Beijing
Tel: 86-10-8458 4010
Fax: 86-10-8458 4001

France – St. Germain-en-Laye
Tel: 01 30 61 39 92
Fax: 01 30 61 58 01
Lab-FR@pall.com

Germany – Dreieich
Tel: 06103-307 333
Fax: 06103-307 399
Lab-DE@pall.com

India – Mumbai
Tel: 91-22-55995555
Fax: 91-22-55995556

Italy – Milano
Tel: 02-477961
Fax: 02-47796394

or 02-4122985

Japan – Tokyo
Tel: 03-6901-5800
Fax: 03-5322-2134

Korea – Seoul
Tel: 2-569-9161
Fax: 2-569-9092

Poland – Warszawa
Tel: (+ 48 22) 510 21 00
Fax: (+ 48 22) 510 21 01

Russia – Moscow
Tel: 00 75 01 787 76 14
Fax: 00 75 01 787 76 15

Singapore
Tel: (65) 389-6500
Fax: (65) 389-6501

South Africa – Johannesburg
Tel: +27-11-2662300
Fax: +27-11-3253243

Spain – Madrid
Tel: 91-657-9876
Fax: 91-657-9836

Sweden – Lund
Tel: +46 (0)46 158400
Fax: +46 (0)46 320781

Switzerland – Basel
Tel: 061-638 39 00
Fax: 061-638 39 40

Taiwan – Taipei
Tel: 2-2545-5991
Fax: 2-2545-5990

United Kingdom – Farlington
Tel: 02392 302600
Fax: 02392 302601

Lab-UK@europe.pall.com

Visit us on the Web at www.pall.com/lab

E-mail us at Lab@pall.com

Filtration. Separation. Solution.SM

© 2005 Pall Corporation. Pall, PALL, Envirochek, and Supor are trademarks of Pall Corporation. Filtration. Separation. Solution. is a service mark of Pall Corporation. ® indicates a registered trademark in the U.S.