


Skip to the good bit with 4-in-1 Mini-UniPrep™ Syringeless Filters

Avoid the unnecessary steps of analytical sample prep



We're over the old, wasteful way to filter samples?

Our 4-in-1 Syringeless Filters help you speed through analytical sample prep

They're preassembled with fewer consumables and less contamination risk. So, you can simply skip to analysis faster and with less waste.

The benefits of Mini-UniPrep Syringeless Filters include:



4-in-1 syringeless filters, syringe, membrane filters, vial, and septum



Enhances HPLC column life



Faster and effortless workflows



Significantly reduce your plastic footprint




Follow the QR code to learn more and sign up for a free sample

Smarter filtration just makes sense

Using a preassembled Mini-UniPrep Syringeless Filter means improved clarity of your sample.

Reasons for smarter filtration include:


- Protect your UHPLC or HPLC autosampler's pump, valves, and tubing from damage or clogging.
- Prevent particulates from building up in the column that can affect data quality and reduce column life.



Condensing workflows the more sustainable way

Mini-UniPrep Syringeless Filters are fast and easy to use. They help avoid repetitive sample prep tasks that cause muscle fatigue.

Mini-UniPrep Syringeless Filters also generate less waste than traditional methods. So your team can work towards your lab's sustainability goal.



[cytiva.com](https://www.cytiva.com)

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva.

Mini-UniPrep is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Any other third-party trademarks are the property of their respective owners.

© 2024 Cytiva

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY46152-27Aug2024-FL