


eQualification

Your digital qualification experience



Forget binders and paper documentation. Digital documentation is the power of now. Here are **four compelling reasons to move to digital qualification today.**



1 Improve security


With eQualification, you have the benefit of increased accuracy without the risk of information loss. Avoid the challenge of updating documents with handwritten changes. Maintain data and documentation integrity with GDP, GAMP5, and ALCOA+ compliance standards.



2 Be more flexible

eQualification allows you to:

- Select flexible qualification offerings including physical or electronic, which help you meet your process and budget needs.
- Replace paper documentation with electronic storage.
- Follow the process live or remotely.



3 Save time

With eQualification your digital qualification binder is available to you as soon as your content is ready. Reviewing, approving, and sharing documents is easier and faster compared to paper binders. Your final qualification reports include hyperlinks for improved navigation during your audits.




4 Contribute to sustainability

eQualification eliminates shipping and paper waste helping to reduce CO₂ emissions.

How does eQualification work?

To eQualification



Transition from conventional binder

1

Your qualification service content is created and sent to you for review

2

You approve content by signing it digitally and one of our engineers is scheduled to complete it

3

Efficient completion of IQ/OQ, RQ, or CCP by our field service engineers

4

Once approved, you can download the complete report in PDF format

Additional benefits

- eQualification is available for installation and operational qualification (IQ/OQ), requalification (RQ), and change control procedure (CCP).
- The Cytiva eQualification solution is validated to 21 CFR Part 11 requirements. Non-product software is validated to Cytiva's quality management system process for computer system validation (CSV).
- You can access the full audit trail.
- Your final report includes all appendices and attachments in one document.

cytiva.com/iqoq

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corporation or an affiliate doing business as Cytiva. Any other trademarks are the property of their respective owners. Any use of software may be subject to one or more end user license agreements, a copy of, or notice of which, are available on request.

© 2025 Cytiva
For local office contact information, visit cytiva.com/contact

CY44776-04Dec25-IG

Learn more about eQualification, your digital qualification experience

Learn more