

# Instructions: Staining Tray Set

## 2-D electrophoresis

The 80-6468-17 Staining Tray Set provides a convenient means of staining up to four large format (26 × 20 cm) gels.

The set includes two stainless steel trays for staining solutions; one perforated stainless steel tray that seats within the staining trays; and a transparent plastic cover. The perforated insert supports and restrains gels for transfer between staining trays while allowing staining solutions to drain rapidly.

### Recommended stain solution volumes

Number of gels	Stain/destain solution (litres)
1	1.0
4	1.5

**Amersham Biosciences UK Limited**  
Amersham Place Little Chalfont  
Buckinghamshire England HP7 9NA

**Amersham Biosciences Biotech AB**  
SE-751 84 Uppsala Sweden

**Amersham Biosciences Biotech Inc.**  
800 Centennial Avenue PO Box 1327  
Piscataway NJ 08855 USA

**Amersham Biosciences Biotech Europe GmbH**  
Munzinger Strasse 9 D-79111 Freiburg Germany

[www.amershambiosciences.com](http://www.amershambiosciences.com)

Amersham Biosciences is a trademark of Amersham plc.

© 2000 Amersham Biosciences Inc. All rights reserved.

All goods and services are sold subject to the terms and conditions of sale of the company within the Amersham Biosciences group that supplies them. A copy of these terms and conditions is available on request.

Printed in the USA.