

UNICORN 5.11

Software Change Description 28-4094-64 AA



This document describes the implemented changes and improvements in UNICORN™ v5.11 compared to v5.10.

- Maintenance module now reads serial numbers for all components in ÄKTAcrossflow™
- Reports from Maintenance module now print out complete list of ÄKTApilot™ components
- Profibus system protocol extended
- Improved Method wizard for ÄKTAcrossflow and ÄKTExpress™ allows faster start-up in the UNICORN network
- Filter area removed from ÄKTAcrossflow report
- Method wizard for ÄKTAcrossflow compatible with Microsoft™ Windows™ XP SP2



www.gehealthcare.com/unicorn

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

ÅKTAcrossflow, ÅKTApilot, ÅKTAexpress, Drop design, and UNICORN are trademarks of GE Healthcare companies. GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Windows and Microsoft are trademarks of Microsoft Corp.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare that supplies them. General Electric Company reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your local GE Healthcare representative for the most current information.

Purification and preparation of fusion proteins and affinity peptides comprising at least two adjacent histidine residues may require a license under US patent 5,284,933 and US patent 5,310,663, including corresponding foreign patents (assignee: Hoffmann-La Roche Inc).

© 2006 General Electric Company – All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Münzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327, Piscataway,
NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo 169-0073, Japan

Asia Pacific Tel +65 6275 1830 Fax +65 6275 1829 • Australasia Tel +61 2 9899 0999 Fax +61 2 9899 7511 • Austria Tel 01/57606-1619 Fax 01/57606-1627 • Belgium Tel 0800 73 888 Fax 02 416 82 06 • Canada Tel 800 463 5800 Fax 800 567 1008 • Central, East, & South East Europe Tel +43 1 972720 Fax +43 1 97272 2750 • Denmark Tel 45 16 2400 Fax 45 16 2424 • Finland & Baltics Tel +358 (0)9 512 39 40 Fax +358 (0)9 512 39 439 • France Tel 01 69 35 67 00 Fax 01 69 41 96 77 Germany Tel: 089 96281 660 Fax 089 96281 620 • Greater China Tel +852 2100 6300 Fax +852 2100 6338 • Italy Tel 02 27322 1 Fax 02 27302 212 • Japan Tel +81 3 5331 9336 Fax +81 3 5331 9370 • Latin America Tel +55 11 3933 7300 Fax +55 11 3933 7304 • Middle East & Africa Tel +30 210 9600 687 Fax +30 210 9600 693 • Netherlands Tel 0800 82 82 82 1 Fax 0800 82 82 82 4 • Norway Tel 815 65 555 Fax 815 65 666 • Portugal Tel 21 417 7035 Fax 21 417 3184 • Russia & other C.I.S. & N.I.S Tel +7 (495) 956 5177 Fax +7 (495) 956 5176 • Spain Tel 93 594 49 50 Fax 93 594 49 55 • Sweden Tel 018 612 1900 Fax 018 612 1910 • Switzerland Tel 0848 8028 12 Fax 0848 8028 13 • UK Tel 0800 616928 Fax 0800 616927 • USA Tel 800 526 3593 Fax 877 295 8102



imagination at work