

## Addendum

---

### 1. 1.3 Product markings:

- GE Healthcare Code No: 18-1163-74 (Pack 50)  
18-1162-08 (Pack 100)  
Custom order (Pack 200)
- Serial number
- Year of manufacture
- Operating temperature: +4°C to +30°C
- Supply air pressure: 6 to 7 bar g
- Air consumption: 500 NI/min\* (Pack 50)  
1000 NI/min\* (Pack 100)  
1250 NI/min\* (Pack 200)
- Maximum outlet pressure: 7 bar g
- Pump capacity: 50 l/min (Pack 50)  
100 l/min (Pack 100)  
200 l/min (Pack 200)
- EX category: II 3G

### 2. Table 1-1: Dimensions and design specifications for Chromaflow Packing stations.

Description	
Footprint Pack 50 Pack 100 Pack 200	810 x 715 mm 810 x 715 mm 830 x 910 mm
Height Pack 50 Pack 100 Pack 200	1175 mm 1175 mm 1290 mm
Weight Pack 50 Pack 100 Pack 200	115 kg 130 kg 170 kg
Operating temperature	+4°C to +30°C
Supply air pressure	6 to 7 bar g
Noise level (maximum)	85 dB (A)



Description	
Air consumption Pack 50 Pack 100 Pack 200	500 NI/minute 1000 NI/minute 1250 NI/minute
Max. outlet pressure	7 bar g
Pump capacity* Pack 50 Pack 100 Pack 200	50 l/min 100 l/min 200 l/min

\*: Flow valid for water at 20°C with 7 bar g supply air pressure, 1 bar g pressure resistance. Sufficient air supply needed.

### 3. Table 1-3: Interface connection dimensions and sizes for the Chromaflow Packing station.

Interface	Dimension, internal diameter	Flange type and size
Packing station inlet B	Pack 50: 1" (22.1 mm) Pack 100: 1.5" (34.8 mm) Pack 200: 1.5" (34.8 mm)	Pack 50: TC 1" Pack 100: TC 1.5" Pack 200: TC 1.5"
Packing station inlet C	Pack 50: 1" (22.1 mm) Pack 100: 1.5" (34.8 mm) Pack 200: 1.5" (34.8 mm)	Pack 50: TC 1" Pack 100: TC 1.5" Pack 200: TC 1.5"
Packing station inlet SOB	Pack 50: 1" (22.1 mm) Pack 100: 1" (22.1 mm) Pack 200: 1.5" (34.8 mm)	Pack 50: TC 1" Pack 100: TC 1" Pack 200: TC 1.5"
Packing station outlet D	Pack 50: 1" (22.1 mm) Pack 100: 1" (22.1 mm) Pack 200: 1.5" (34.8 mm)	Pack 50: TC 1" Pack 100: TC 1" Pack 200: TC 1.5"
Packing station outlet SIT	Pack 50: 1" (22.1 mm) Pack 100: 1" (22.1 mm) Pack 200: 1.5" (34.8 mm)	Pack 50: TC 1" Pack 100: TC 1" Pack 200: TC 1.5"
Packing station outlet SIB	Pack 50: 1" (22.1 mm) Pack 100: 1" (22.1 mm) Pack 200: 1.5" (34.8 mm)	Pack 50: TC 1" Pack 100: TC 1" Pack 200: TC 1.5"