

Attach Bottle holder and connect Adapter for air sensor

Instructions

Scope

Instructions on how to use Bottle holder (28956327) and Adapter for air sensor (28956342).

Note: *The bottle holder is designed for the ÄKTA™ avant and ÄKTA pure instruments and may not be possible to attach to other ÄKTA instruments.*

Attach Bottle holder


The Bottle holder is used for holding bottles and can be attached to the holder rails.



Connect Adapter for air sensor

The air sensor adapter is used to hold an external air sensor.

The air sensor with adapter is connected to the Bottle holder, see illustrations below.



Part	Description
1	Bottle holder
2	Air sensor
3	Adapter for air sensor

Page intentionally left blank



cytiva.com/AKTA

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

AKTA is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

28968642 AC V:4 10/2020