

# ReadyToProcess™ columns

Overview and ordering information

## Prepacked chromatography columns for scale-up and GMP manufacturing

- Save time and maximize throughput by eliminating in-house packing.
- Choose from a wide range of high-performance resins and column dimensions to suit your process needs.

See available dimensions on [page 2](#).

## Predictable production planning

- Manage timelines and minimize delays with a single supplier for resins and columns.
- Get consistent performance with verified packing methods.
- Get support for purchase, applications, and troubleshooting from one supplier.

## Short lead times

- Select configurations can ship 1 or 4 weeks after order acceptance.\*

### Look for the lead time tags in the column list



### Click the chromatography mode to jump to column list:

Affinity (mAbs and antibody variants) →

Affinity (others) →

Anion exchange →

Cation exchange →

Hydrophobic interaction →

Multimodal →

Size exclusion →

#### \*Details

*Lead time* refers to the time between when we accept your purchase order to the time your columns can be shipped from our facility in Sweden.

#### Maximum order quantities

4-week lead time: 10 columns


1-week lead time: 1 column

#### Exceptions

Shorter lead times do not apply to orders placed in late June, July, late November, or December. For help planning your orders or for more details, reach out to your sales representative.




## Inner diameters and bed heights



\*Technical feasibility of the combination of diameters, bed heights, and resin will be evaluated by our team to ensure our quality standards can be fulfilled.

## Resin bed volumes from 0.4 L to 71 L



If you can't find the configuration you need in this document, reach out to us about customization.

**MabSelect Prisma™ resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess™ MabSelect Prisma column	80	80	0.4	29475859	*
		100	0.5	29453291	*
		150	0.8	29420989	4
		200	1	17549861	4
	126	80	1	29944817	*
		100	1.2	29475856	*
		120	1.5	29512381	*
	178	200	2.5	17549862	4
		100	2.5	29475852	*
		150	3.7	29804293	4
	251	200	5.0	17549863	4
		100	4.9	29475842	*
	359	200	10	17549864	4
		100	10	29802518	*
	450	150	15	29601750	4
		200	20	17549865	1
	600	100	16	29703720	*
		200	32	17549866	4
	600	140	40	29698529	*
		150	42	29639034	*
		200	57	17549867	4
		250	71	29697135	*

\*Please check lead time with your Cytiva representative.



If you can't find the configuration you need, reach out to us about customization.

**MabSelect Prisma X resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess MabSelect Prisma X column	80	150	0.8	30042317	*
		200	1	30042315	*
	126	150	1.9	30042314	*
		200	2.5	30042313	*
	178	150	3.7	30042312	*
		200	5	29944818	4
	251	150	7.4	30042310	*
		200	10	29947596	*
	359	150	15	30042309	*
		200	20	30042308	*
	450	150	24	30042302	*
		200	32	30042307	*
	600	150	42	30042301	*
		200	57	30042306	*

\*Please check lead time with your Cytiva representative.

**MabSelect™ SuRe 70 resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess MabSelect SuRe 70 column	80	150	0.8	30042324	*
		200	1	30042323	*
	126	150	1.9	30042322	*
		200	2.5	30012494	*
	178	150	3.7	30042321	*
		200	5	29940012	*
	251	150	7.4	30042319	*
		200	10	29940011	4
	359	150	15	30042318	*
		200	20	29938388	4
	450	150	24	30042305	*
		200	32	29941295	*
	600	150	42	30042303	*
		200	57	29941294	*

\*Please check lead time with your Cytiva representative.

**MabSelect SuRe LX resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29026927	4
MabSelect SuRe LX column	126	200	2.5	29027644	4
	178	200	5	29145983	4
	251	200	10	29172004	4
	359	200	20	29172006	1
		250	25	29610439	4
	450	150	24	29378061	*
		200	32	29256260	4
		150	42	29742730	*
	600	200	57	29376119	4
		250	71	29698518	*

\*Please check lead time with your Cytiva representative.

 If you can't find the configuration you need, reach out to us about customization.

**MabSelect SuRe resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		150	0.8	29287558	4
MabSelect SuRe column	80	200	1	28951110	4
		250	1.3	29287559	4
		80	1	29357019	*
	126	100	1.2	29351070	*
		150	1.9	29287560	4
		200	2.5	28901717	4
		250	3.1	29287561	4
		150	3.7	29287555	4
	178	200	5	29145980	4
		250	6.2	29287556	4
		100	4.9	29357048	*
		120	5.9	29587576	*
	251	150	7.4	29287562	4
		160	7.9	29452019	*
		200	10	28901718	4
		250	12.4	29287557	4
		100	10	29660006	*
	359	120	12	29647158	*
		130	13	29705618	*
		150	15	29304086	4
		200	20	28901719	4
		250	25	29304087	4
	450	150	24	29354399	*
		200	32	29212474	4
		160	45	29498040	*
	600	200	57	29376118	4
		250	71	29738873	*

\*Please check lead time with your Cytiva representative.

**MabSelect Xtra™ resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	359	200	20	29097172	*
MabSelect Xtra NS column					

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

**MabSelect resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	28951128	4
MabSelect NS column	126	200	2.5	28941522	4
	251	200	10	28941523	4
	359	200	20	28941524	4
	450	200	32	29256261	4
		250	40	29237040	*
	600	250	71	29703498	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

 If you can't find the configuration you need, reach out to us about customization.

**MabSelect VH3 resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29804777	4
MabSelect VH3 column	359	150	15	29939903	4

**MabSelect VL resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	251	200	10	29732180	4
MabSelect VL NS column	359	150	15	29801782	4
		200	20	29731571	4
	450	200	32	17542066	4

NS = non-sanitized.

**Capto™ L resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		90	1.1	29277419	*
Capto L NS column	126	200	2.5	29015989	*
		178	200	5	29146142
	251	200	10	29015990	*
	359	200	20	29015991	*
	450	200	32	29256262	*
	600	200	57	29719608	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

**KappaSelect resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	126	150	1.9	29588695	*
KappaSelect NS column	359	200	20	29106344	*
	600	150	42	29516452	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

**LambdaFabSelect resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	126	200	2.5	29168997	*
LambdaFabSelect NS column	359	200	20	29339637	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.



## Capto AVB resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	126	200	2.5	29501102	*
Capto AVB NS column	178	100	2.5	29747481	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

 If you can't find the configuration you need, reach out to us about customization.

## Capto Blue resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	359	200	20	29695369	*
Capto Blue column					

\*Please check lead time with your Cytiva representative.

## Capto Chelating resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29719177	*
Capto Chelating NS column	126	200	2.5	29719178	*
	251	200	10	29597509	*
	359	250	25	29940154	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

## Capto DeVirS resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29606061	*
Capto DeVirS column		150	3.7	29801814	*
	178	200	5	29266282	*
		250	6.2	29803363	*
	251	200	10	29304512	*
	600	250	71	29691384	*

\*Please check lead time with your Cytiva representative.

## Capto Lentil Lectin resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	150	0.8	29697174	*
Capto Lentil Lectin NS column	450	130	21	29624834	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

## Alpha-1 Antitrypsin Select resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Alpha-1	80	200	1	29512525	*
Antitrypsin Select NS column	359	200	20	29300743	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

## AVB Sepharose High Performance resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		100	0.5	29282169	*
AVB Sepharose HP NS column	80	180	0.9	29351777	*
	126	200	2.5	29304511	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.



## Blue Sepharose 6 Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29620541	*
Blue Sepharose 6	251	200	10	29034314	*
Fast Flow column	359	200	20	29027418	*

\*Please check lead time with your Cytiva representative.

 If you can't find the configuration you need, reach out to us about customization.

## Chelating Sepharose Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess					
Chelating Sepharose	450	200	32	29700559	4
Fast Flow NS column					

NS = non-sanitized.

## Heparin Sepharose 6 Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess					
Heparin Sepharose 6	178	200	5	29717295	4
Fast Flow NS column					

NS = non-sanitized.

## IMAC Sepharose 6 Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29318226	4
IMAC Sepharose 6	178	200	5	29336847	4
Fast Flow NS column	251	200	10	29624097	4

NS = non-sanitized.

## VIISelect resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess					
VIISelect NS column	600	100	28	29658522	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

## VIIISelect resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29175834	*
VIIISelect NS column	126	200	2.5	29037897	*
	178	200	5	29296446	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.



## Capto DEAE resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Capto DEAE column	80	150	0.8	29657118	4
		200	1	29719080	4
	126	150	1.9	29802805	4
		200	2.5	29731668	4
		251	10	28929127	4

 If you can't find the configuration you need, reach out to us about customization.

## Capto Q resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Capto Q column	80	100	0.5	29804778	*
		200	1	28951090	4
	126	150	1.9	29287566	4
		200	2.5	28901723	4
		150	3.7	29287567	4
	178	200	5	29146143	4
		250	6.2	29287568	4
		150	7.4	29287569	4
	251	200	10	28901724	4
		250	12.4	29287570	4
		150	15	29304088	4
	359	200	20	28901725	4
		250	25	29304089	4
	450	200	32	29212473	4
		250	40	29620561	*
	600	150	42	29661597	*
		200	57	29376129	4

\*Please check lead time with your Cytiva representative.

## Capto Q ImpRes resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Capto Q ImpRes column	80	200	1	29101654	4
		126	200	2.5	29101655
	178	200	5	29146150	4
		150	7.4	29610087	4
	251	200	10	29101657	4
		200	20	29101658	4
	359	220	22	29611696	*
		200	32	29256251	4
		220	35	29402830	*
	450	250	40	29939673	*
		200	57	29376127	4
	600	250	71	29691383	*

\*Please check lead time with your Cytiva representative.

## Capto Q XP resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Capto Q XP column	359	200	20	29714881	*

\*Please check lead time with your Cytiva representative.

## CM Sepharose Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess CM Sepharose Fast Flow column	126	200	2.5	28929117	4
	251	200	10	28929118	4
	359	200	20	28929119	4
	450	200	32	29409615	4
	600	110	31	29511633	*

\*Please check lead time with your Cytiva representative.

**DEAE Sepharose Fast Flow resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	28951126	4
DEAE Sepharose	126	200	2.5	28929114	4
Fast Flow column	178	150	3.7	29704025	4
		150	7.4	29704024	4
	251	200	10	28929115	4
		150	15	29703329	4
	359	200	20	28929116	4
		130	37	29511635	*
600	200	57	29697781	4	

\*Please check lead time with your Cytiva representative.

 If you can't find the configuration you need, reach out to us about customization.

**Q Sepharose Fast Flow resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	28951125	4
Q Sepharose		250	1.3	29287576	4
Fast Flow column	126	150	1.9	29588338	4
		200	2.5	28929076	4
		100	2.5	29660005	*
	178	150	3.7	29711780	4
		200	5	29146155	4
	251	200	10	28929079	4
		250	12.4	29358340	4
	359	150	15	29317031	4
		200	20	28929082	4
		250	25	29304095	4
	450	200	32	29212479	4
		250	40	29366934	*
600	150	42	29516451	*	
	200	57	29376120	4	

\*Please check lead time with your Cytiva representative.

**Q Sepharose High Performance resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)	
ReadyToProcess	80	100	0.5	29417695	*	
Q Sepharose HP		150	0.8	29805177	*	
NS column		200	1	29101715	*	
		126	200	2.5	29101716	*
		178	200	5	29204338	*
		251	200	10	29101717	*
			150	24	29716840	*
		450	200	32	29212478	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

**Q Sepharose XL resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29032670	4
Q Sepharose XL column	126	200	2.5	29000545	4
		178	200	5	29296929
	251	200	10	29015616	4
		100	10	29353387	*
	359	200	20	29695976	4
		100	28	29704278	*
	600	200	57	29623200	4

\*Please check lead time with your Cytiva representative.



## Capto S resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	28951093	4
Capto S column	126	200	2.5	28901729	4
	178	200	5	29482665	4
	251	200	10	28901730	4
	359	200	20	28901731	4
	450	200	32	29615551	4
	600	200	57	29485444	4

 If you can't find the configuration you need, reach out to us about customization.

## Capto S ImpAct resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		150	0.8	29287577	4
Capto S ImpAct column	80	200	1	17371761	4
		250	1.3	29717186	4
		126	200	2.5	17371762
	178	200	5	29146146	4
	251	200	10	17371763	4
	359	200	20	17371764	4
		190	30	29871571	*
		200	32	29212475	4
	450	220	35	29875368	*
		250	40	29691087	*
		600	200	57	29376130

\*Please check lead time with your Cytiva representative.

## Capto SP ImpRes resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29101661	4
Capto SP ImpRes column	126	200	2.5	29101662	4
		250	3.1	29715992	4
	178	200	5	29146147	4
	251	100	4.9	29728648	*
		200	10	29101663	4
	359	200	20	29101665	4
	450	200	32	29251611	4
	600	200	57	29376128	4
		250	71	29715925	*

\*Please check lead time with your Cytiva representative.

## SP Sepharose Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		100	0.5	29697210	*
SP Sepharose Fast Flow column	80	200	1	28951097	4
		126	200	2.5	28929105
	178	200	5	29146156	4
	251	150	7.4	29655825	4
		200	10	28929106	4
		150	15	29372451	4
	359	200	20	28929107	4
		250	25	29375266	4
	450	200	32	29256256	4
		250	40	29276957	*
	600	200	57	29376121	4

\*Please check lead time with your Cytiva representative.

## SP Sepharose High Performance resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29101711	*
SP Sepharose HP NS column	126	200	2.5	29101712	*
	178	200	5	29204339	*
	251	200	10	29101713	*
	450	150	24	29611698	*
		200	32	29256257	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.



## Capto Butyl resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	126	200	2.5	29011998	4
Capto Butyl column	251	200	10	29696789	4
	359	200	20	29594551	4
		250	25	29660869	*
	450	250	40	29873456	*

\*Please check lead time with your Cytiva representative.

 If you can't find the configuration you need, reach out to us about customization.

## Capto Butyl ImpRes resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29713752	4
Capto Butyl ImpRes column	126	200	2.5	29655954	4
	178	200	5	29647159	4
		250	6.2	29746199	4
	251	200	10	29138138	4
	359	200	20	29229399	4
	450	200	32	29256254	4
		170	48	29871335	*
	600	200	57	29474655	4

\*Please check lead time with your Cytiva representative.

## Capto Butyl-S resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	251	200	10	29718930	*
Capto Butyl-S column					

\*Please check lead time with your Cytiva representative.

## Capto Octyl resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	178	200	5	29436142	*
Capto Octyl column	600	200	57	29437660	*

\*Please check lead time with your Cytiva representative.

## Capto Phenyl (high sub) resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29316634	4
Capto Phenyl (high sub) column	126	200	2.5	29207816	4
		150	3.7	29287578	4
	178	200	5	29225988	4
		250	6.2	29287579	4
		150	7.4	29287580	4
	251	200	10	29194484	4
		250	12.4	29287581	4
		150	15	29304091	4
	359	200	20	29119051	4
		250	25	29304092	4
		180	29	29722538	*
	450	200	32	29321929	4
		200	32	29321929	4
	600	200	57	29376124	4

\*Please check lead time with your Cytiva representative.



## Capto Phenyl ImpRes resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	180	0.9	29730658	*
Capto Phenyl ImpRes column	80	200	1	29101697	4
		150	1.9	29609021	4
	126	200	2.5	29101698	4
		178	200	5	29642661
	251	150	7.4	29696391	4
		200	10	29101700	4
	359	210	10.4	29660003	*
		150	15	29943675	4
		200	20	29101702	4
	450	200	32	29256253	4
		250	40	29729894	*
	600	200	57	29649594	4

\*Please check lead time with your Cytiva representative.

 If you can't find the configuration you need, reach out to us about customization.

## Capto PlasmidSelect resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		100	0.5	29702352	*
Capto PlasmidSelect column	80	150	0.8	29716910	4
		200	1	29415793	4
	126	150	1.9	29512474	4
		100	2.5	29700170	*
		130	3.2	29605795	*
	178	150	3.7	29435027	4
		200	5	29513793	4
	251	150	7.4	29435206	4
	359	150	15	29397750	4
	450	150	24	29397609	*

\*Please check lead time with your Cytiva representative.

## Butyl Sepharose 4 Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29318210	4
Butyl Sepharose 4 Fast Flow NS column	126	200	2.5	29698799	4
		100	4.9	29690006	*
	251	200	10	29698651	4
		359	200	20	29717004
	600	150	42	29722865	*
160		45	29731499	*	

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

## Butyl Sepharose High Performance resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29713754	*
Butyl Sepharose HP NS column	126	200	2.5	29450716	*
		178	200	5	29204341
	251	200	10	29713753	*
		150	24	29402837	*
	450	200	32	29256259	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

**Phenyl Sepharose 6 Fast Flow (low sub) resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Phenyl Sepharose 6 Fast Flow (low sub) column	80	200	1	28951111	4
	126	200	2.5	28901735	4
	251	200	10	28901736	4
	359	200	20	28901737	4

 If you can't find the configuration you need, reach out to us about customization.

**Phenyl Sepharose 6 Fast Flow (high sub) resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Phenyl Sepharose 6 Fast Flow (high sub) column	80	200	1	28951197	4
	126	200	2.5	28929108	4
	178	200	5	29453083	4
	251	200	10	28929109	4
	359	200	20	28929110	4
	450	100	16	29689470	*
		200	32	29398527	4
	600	200	57	29489504	4

\*Please check lead time with your Cytiva representative.

**Phenyl Sepharose High Performance resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Phenyl Sepharose HP NS column	80	200	1	29101719	*
	126	200	2.5	29101720	*
	178	200	5	29204340	*
	251	200	10	29101721	*
		150	24	29657449	*
	450	200	32	29256258	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

**PlasmidSelect Xtra resin**

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess PlasmidSelect Xtra NS column	80	200	1	29244205	*
	178	200	5	29400484	*
		150	7.4	29397204	*
	251	200	10	29600356	*
		250	12.4	29400485	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.



## Capto adhere resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	28951109	4
Capto adhere column	80	80	1	29271532	*
		150	1.9	29294361	4
126		200	2.5	28901714	4
		250	3.1	29287571	4
		150	3.7	29615116	4
		200	5	29146144	4
178		250	6.2	29714372	4
		150	7.4	29714373	4
		200	10	28901715	1
251		150	15	29304090	4
		200	20	28901716	1
359		250	25	29317029	4
		150	24	29378062	*
		160	25	29659080	*
		200	32	29212477	4
450		210	33	29707318	*
		200	57	29376122	4
		250	71	29711668	*

\*Please check lead time with your Cytiva representative.

 If you can't find the configuration you need, reach out to us about customization.

## Capto adhere ImpRes resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29101703	4
Capto adhere ImpRes column	126	100	1.2	29730498	*
		200	2.5	29101704	4
178		150	3.7	29287573	4
		200	5	29146151	4
		250	6.2	29287574	4
251		150	7.4	29287572	4
		200	10	29101705	4
		250	12.4	29287575	4
		150	15	29304093	4
359		200	20	29101706	4
		250	25	29304094	4
450	200	32	29212476	4	
600	200	57	29376125	4	

\*Please check lead time with your Cytiva representative.

## Capto Core 400 resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		100	0.5	29517132	*
Capto Core 400 column	80	160	0.5	29724312	*
		200	1	29697211	4
		126	200	2.5	29639038
178	250	6.2	29641373	4	
251	150	7.4	29751769	4	
359	200	20	29453290	4	
450	200	32	29364879	4	

\*Please check lead time with your Cytiva representative.

## Capto Core 700 resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess		100	0.5	29263633	*
Capto Core 700 column	80	150	0.8	29704354	4
		200	1	29020257	4
		150	1.9	29464142	4
	126	200	2.5	29119812	4
		250	3.1	29874192	4
		150	3.7	29714758	4
	178	200	5	29410473	4
		250	6.2	29610456	4
		251	10	29175700	4
	359	200	20	29255522	4
	450	200	32	29445827	*
	600	200	57	29631005	4

\*Please check lead time with your Cytiva representative.

If you can't find the configuration you need, reach out to us about customization.

## Capto MMC resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	28951118	4
Capto MMC column	126	200	2.5	28929120	4
		178	200	5	29146145
	251	200	10	28929121	4
		250	12.4	29589044	4
	359	150	15	29372600	4
		200	20	28929122	4
		150	24	29716843	*
	450	200	32	29256255	4
		250	40	29943933	*
		200	57	29428374	4
	600	250	71	29661048	*

\*Please check lead time with your Cytiva representative.

## Capto MMC ImpRes resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess	80	200	1	29101707	4
Capto MMC ImpRes column	126	200	2.5	29101708	4
		178	200	5	29146153
	251	200	10	29101709	4
		150	15	29661505	4
	359	200	20	29101710	4
		450	200	32	29256252
	600	200	57	29505248	4

## Sephacryl™ S-500 HR resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Sephacryl S-500 HR NS column	178	200	5	29274322	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.

If you can't find the configuration you need, reach out to us about customization.

## Sepharose 6 Fast Flow resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Sephacryl S-500 HR NS column	80	250	1.3	29339640	4
Sephacryl S-500 HR NS column	126	200	2.5	28935148	4
	251	200	10	28935149	4
	450	200	32	29512774	4

## Superdex™ 200 Prep Grade resin

Column	Inner diameter (mm)	Bed height (mm)	Bed volume (L)	Product code	Lead time (week)
ReadyToProcess Superdex 200 Prep Grade NS column	126	250	3.1	29357044	*
Superdex 200 Prep Grade NS column	178	250	6.2	29432991	*
	251	200	10	29946367	*

\*Please check lead time with your Cytiva representative. NS = non-sanitized.



## [cytiva.com/ReadyToProcess-columns](https://www.cytiva.com/ReadyToProcess-columns)

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva. Capto, MabSelect, MabSelect PrismA, MabSelect Xtra, ReadyToProcess, Sephacryl, Sepharose, and Superdex are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. Any other trademarks are the property of their respective owners.

© 2026 Cytiva

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

CY1797-30Jan26-FL