


# HyClone SFM4Insect



**HyClone™ SFM4Insect is a versatile animal-derived component-free cell culture medium developed to support the growth of multiple insect cell lines and production of a variety of recombinant proteins without compromising performance. HyClone SFM4Insect has been successfully tested in a variety of culture systems ranging from T-flasks, shaker flasks, and various traditional and disposable bioreactors for production of a variety of recombinant proteins using the baculovirus expression vector system (BEVS).**

## HyClone SFM4Insect medium



Growth kinetics of batch cultures of Sf9 cells comparing SFM4Insect to our traditional serum-free SFX-Insect medium.

SFM4Insect provides excellent growth of many key insect cell lines, including Sf9, Sf21, High Five™ (T.ni), and D.mel-2 cells, while requiring minimal adaptation.

### Key features


- Supports multiple insect cell lines
- High-performance growth in multiple culture systems including high-density bioreactors
- Complexed lipids for enhanced stability
- Available in bottles, single-use containers, or custom packaging
- Manufactured from traceable components according to cGMP and ISO guidelines

To learn more about HyClone SFM4Insect, visit [www.cytiva.com/hyclone](http://www.cytiva.com/hyclone)

Product	Size	Code number
HyClone SFM4Insect liquid medium with L-glutamine	500 mL bottle*	SH30913.01
	1000 mL bottle*	SH30913.02
	6 × 1000 mL bottle*	SH30913.LS
	5 L bag	SH30913.04
	10 L bag	SH30913.05
	20 L bag	SH30913.06
	50 L bag	SH30913.07
	100 L bag	SH30913.08
	200 L bag	SH30913.09
HyClone SFM4Insect powder medium with L-glutamine	5 L bottle†	SH30912.01
	10 L bottle†	SH30912.02
	50 L bottle†	SH30912.03
	100 L polybag/pail	SH30912.04
	500 L polybag/pail	SH30912.05
	1000 L polybag/pail	SH30912.06

\* Polyethylene terephthalate (PETE) bottle

† High-density polyethylene (HDPE) bottle



Doubling times from batch cultures of Sf9 cells comparing SFM4Insect to our traditional serum-free SFX-Insect medium.

## cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. HyClone is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

High Five is a trademark of Invitrogen Corporation. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

CY16146-10Dec20-FL

