

Benchkote sheet for ÄKTA systems

The superior surface protection and convenience of Benchkote™ sheet is now available for ÄKTA™ systems. Benchkote sheet for ÄKTA systems is designed to cover the top buffer tray on ÄKTA protein purification systems.

Benchkote surface protection

Benchkote sheets are manufactured from absorbent, hydrophobic material designed to protect lab surfaces against hazardous spills. The high-quality, smooth, and absorbent Whatman™ paper quickly absorbs liquid spills, and the laminated polyethylene layer prevents flow-through to the working surface.

For more details, [click here](#)

Benchkote surface protectors

Benchkote sheet for ÄKTA systems

Lab hygiene made easy

- Protects the top buffer tray from spillages and salt deposits
- Simplifies daily and weekly preventive maintenance of the top buffer tray
- Comes in three convenient pack sizes: 10, 25, and 50 sheets
- Fits ÄKTA start, ÄKTA pure, ÄKTA avant, and ÄKTA go protein purification systems

Lab scale ÄKTA systems

ÄKTA lab-scale protein purification systems are designed for purification of biomolecules, providing speed, performance, and flexibility in research and process development. UNICORN™ software is used to plan, run and evaluate chromatography runs on ÄKTA systems and combines simplicity with power in protein purification. ÄKTA systems, backed by Cytiva columns and chromatography resins, are designed for a broad range of applications and purification tasks.

For more details, click here

Chromatography systems

Ordering information

Product	Pack size	Material	Product code
Benchkote sheet for ÄKTA pure	10 sheets	Absorbent paper and laminated polyethylene	2300-10060
Benchkote sheet for ÄKTA pure	25 sheets	Absorbent paper and laminated polyethylene	2300-10061
Benchkote sheet for ÄKTA pure	50 sheets	Absorbent paper and laminated polyethylene	2300-10062
Benchkote sheet for ÄKTA start	10 sheets	Absorbent paper and laminated polyethylene	2300-10063
Benchkote sheet for ÄKTA start	25 sheets	Absorbent paper and laminated polyethylene	2300-10064
Benchkote sheet for ÄKTA start	50 sheets	Absorbent paper and laminated polyethylene	2300-10065
Benchkote sheet for ÄKTA avant	10 sheets	Absorbent paper and laminated polyethylene	2300-10072
Benchkote sheet for ÄKTA avant	25 sheets	Absorbent paper and laminated polyethylene	2300-10073
Benchkote sheet for ÄKTA avant	50 sheets	Absorbent paper and laminated polyethylene	2300-10074
Benchkote sheet for ÄKTA go	10 sheets	Absorbent paper and laminated polyethylene	2300-10092
Benchkote sheet for ÄKTA go	25 sheets	Absorbent paper and laminated polyethylene	2300-10093
Benchkote sheet for ÄKTA go	50 sheets	Absorbent paper and laminated polyethylene	2300-10094

Benchkote sheet for ÄKTA start, ÄKTA pure, ÄKTA avant, and ÄKTA go systems



For more details, click here

Benchkote sheet for ÄKTA systems

Bench protector sheet



Benchkote surface protectors

cytiva.com/AKTA

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. ÄKTA, Benchkote, UNICORN, and Whatman are a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2020 Cytiva

For local office contact information, visit cytiva.com/contact

CY12814-20Aug20-FL

