

Shelf Life Statement

Regarding: Hollow fiber cartridges and adaptor kits

Hollow fiber cartridges and adaptor kits are made from materials that are considered chemically inert and thus have no assigned expiration date. Product degradation is not expected when stored as recommended: room temperature and away from direct sunlight.

Please note the statement above does not include ReadyToProcess™ Hollow Fiber cartridges.

For shelf life of ReadyToProcess Hollow Fiber cartridges and accessories, please refer to the product specific Certificate of Quality (CoQ). CoQs are available for download by production lot at www.cytiva.com/certificates.

Signature: <i>Rebecca Moore</i>	Date: 26 October 2023
Name: Rebecca Moore	Title Regulatory Support Manager

General statements, including this one, are available on the Customer Regulatory Support website at www.cytiva.com/rsf. See "Download without a subscription" then "Statements."

The information herein is provided to the best of our knowledge and belief, based on the information we have been given and the information available to us at the time that this document was issued. The information may change as additional facts come to light, or as the circumstances warrant. Your use hereof is hence at your own risk and does not relieve you of your own responsibility for checking the suitability of our products for your intended use.

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate.
For additional information on trademarks see cytiva.com/trademarks.
© 2023 Cytiva