

19 May 2023

RE: Discontinuation Notice Filtration Devices containing Polypropylene Membrane

Dear Valued Customer,

Cytiva has been notified by our supplier that a critical raw material used in Filtration Devices will no longer be available. The raw material affects Filtration Devices containing **Polypropylene Membrane**. The below tables reference all part numbers that are affected. Cytiva has limited material available and will produce as many filters as possible. Once supply is consumed, the parts will be discontinued. We expect most of the supply to be consumed during 2023.

Polypropylene membranes exhibit a mildly hydrophobic filter surface that offers broad chemical resistance to many organic solvents. The filters were especially suited to filtration of aqueous /organic solvent suspensions.

Parts to be discontinued once supply is depleted:

6788-1302	PURADISC 13/0.2 PP 100/PK
6785-1302	PURADISC 13/0.2 PP 500/PK
6788-1304	PURADISC 13/0.45 PP 100/PK
6785-1304	PURADISC 13/0.45 PP 500/PK
6786-2502	PURADISC 25/0.2 PP 50/PK
6788-2502	PURADISC 25/0.2 PP 200/PK
6790-2502	PURADISC 25/0.2 PP 1000/PK
6878-1302	Whatman GD/X 13/0.2 PP 150/PK
6878-1304	Whatman GD/X 13/0.45 PP 150/PK
6878-2502	Whatman GD/X 25/0.2 PP 150/PK
6878-2504	Whatman GD/X 25/0.45 PP 150/PK
6879-2504	Whatman GD/X 25/0.45 PP 1500/PK
6978-2504	Whatman GD/XP 25/0.45 PP 150/PK
UN203NPEPP	MUP 0.2uM PP 100/PK
UN503NPEPP	MUP 0.2uM PP 1000/PK
UN203NPUPP	MUP 0.45uM PP 100/PK
UN503NPUPP	MUP 0.45uM PP 1000/PK
US203NPEPP	MUP SS 0.2uM PP 100/PK
US503NPEPP	MUP SS 0.2uM PP 1000/PK
US203NPUPP	MUP SS 0.45uM PP 100/PK
US503NPUPP	MUP SS 0.45uM PP 1000/PK
AV125EPP	Autovial 12 0.2uM PP 50/PK
AV125UPP	Autovial 12 0.45uM PP 50/PK

Suggested replacements are available: select the following for a hydrophilic PTFE membrane that offers broad chemical resistance, low protein binding and low extractables and can be used for both aqueous and aggressive organic solvents.

Discontinued Part		Recommended Alternative	
Part Number	Item Description	Part Number	Item Description
6788-1302	PURADISC 13/0.2 PP 100/PK	6772-1302	PURADISC 13/0.2 H-PTFE 100/PK
6785-1302	PURADISC 13/0.2 PP 500/PK	6773-1302	PURADISC 13/0.2 H-PTFE 500/PK
6788-1304	PURADISC 13/0.45 PP 100/PK	6772-1304	PURADISC 13/0.45 H-PTFE 100/PK
6785-1304	PURADISC 13/0.45 PP 500/PK	6773-1304	PURADISC 13/0.45 H-PTFE 500/PK
6786-2502	PURADISC 25/0.2 PP 50/PK	6772-2502	PURADISC 25/0.2 H-PTFE 50/PK
6788-2502	PURADISC 25/0.2 PP 200/PK	6773-2502	PURADISC 25/0.2 H-PTFE 200/PK
6790-2502	PURADISC 25/0.2 PP 1000/PK	6774-2502	PURADISC 25/0.2 H-PTFE 1000/PK

Suggested alternative replacement: select the following for a hydrophilic RC or PES membrane that can be used when low protein binding, low extractables is required. RC can be used with both aqueous / organic solvents while PES is primarily chosen for aqueous buffers.

Discontinued Part		Recommended Alternative	
Part Number	Item Description	Part Number	Item Description
6788-1302	PURADISC 13/0.2 PP 100/PK	6756-1302	Puradisc 13/0.2 RC 100/pk
6785-1302	PURADISC 13/0.2 PP 500/PK	6757-1302	Puradisc 13/0.2 RC 500/pk
6788-1304	PURADISC 13/0.45 PP 100/PK	6756-1304	Puradisc 13/0.45 RC 100/pk
6785-1304	PURADISC 13/0.45 PP 500/PK	6757-1304	Puradisc 13/0.45 RC 500/pk
6786-2502	PURADISC 25/0.2 PP 50/PK	6756-2502	Puradisc 25/0.2 RC 50/pk
6788-2502	PURADISC 25/0.2 PP 200/PK	6757-2502	Puradisc 25/0.2 RC 200/pk
6790-2502	PURADISC 25/0.2 PP 1000/PK	6758-2502	Puradisc 25/0.2 RC 1000/pk
6878-1302	Whatman GD/X 13/0.2 PP 150/PK	6887-2502	Whatman GD/X 25 RC 0.2UM 150/PK
6878-1304	Whatman GD/X 13/0.45 PP 150/PK	6882-2504	Whatman GD/X 25 RC 0.45UM 150/PK
6878-2502	Whatman GD/X 25/0.2 PP 150/PK	6887-2502	Whatman GD/X 25 RC 0.2UM 150/PK
6878-2504	Whatman GD/X 25/0.45 PP 150/PK	6882-2504	Whatman GD/X 25 RC 0.45UM 150/PK
6879-2504	Whatman GD/X 25/0.45 PP 1500/PK	6883-2504	Whatman GD/X 25 RC 0.45UM 1500/PK
6978-2504	Whatman GD/XP 25/0.45 PP 150/PK	6994-2504	Whatman GD/XP 25/0.45 PES 150/PK
UN203NPEPP	MUP 0.2uM PP 100/PK	UN203NPERC	MUP 0.2uM RC 100/PK
UN503NPEPP	MUP 0.2uM PP 1000/PK	UN503NPERC	MUP 0.2uM RC 1000/PK

UN203NPUPP	MUP 0.45uM PP 100/PK	UN203NPURC	MUP 0.45uM RC 100/PK
UN503NPUPP	MUP 0.45uM PP 1000/PK	UN503NPURC	MUP 0.45uM RC 1000/PK
US203NPEPP	MUP SS 0.2uM PP 100/PK	US203NPEPES	MUP SS 0.2uM PES 100/PK
US503NPEPP	MUP SS 0.2uM PP 1000/PK	US503NPEPES	MUP SS 0.2uM PES 1000/PK
US203NPUPP	MUP SS 0.45uM PP 100/PK	US203NPUPES	MUP SS 0.45uM PES 100/PK
US503NPUPP	MUP SS 0.45uM PP 1000/PK	US203NPUPES	MUP SS 0.45uM PES 100/PK

Suggested alternative replacements: select the following for a hydrophobic PTFE membrane that can be used with most organic solvents.

Discontinued Part

Part Number Item Description

6788-1302	PURADISC 13/0.2 PP 100/PK
6785-1302	PURADISC 13/0.2 PP 500/PK
6788-1304	PURADISC 13/0.45 PP 100/PK
6785-1304	PURADISC 13/0.45 PP 500/PK
6786-2502	PURADISC 25/0.2 PP 50/PK
6788-2502	PURADISC 25/0.2 PP 200/PK
6790-2502	PURADISC 25/0.2 PP 1000/PK
6878-1302	Whatman GD/X 13/0.2 PP 150/PK
6878-1304	Whatman GD/X 13/0.45 PP 150/PK
6878-2502	Whatman GD/X 25/0.2 PP 150/PK
6878-2504	Whatman GD/X 25/0.45 PP 150/PK
6879-2504	Whatman GD/X 25/0.45 PP 1500/PK
6978-2504	Whatman GD/XP 25/0.45 PP 150/PK
UN203NPEPP	MUP 0.2uM PP 100/PK
UN503NPEPP	MUP 0.2uM PP 1000/PK
UN203NPUPP	MUP 0.45uM PP 100/PK
UN503NPUPP	MUP 0.45uM PP 1000/PK
US203NPEPP	MUP SS 0.2uM PP 100/PK
US503NPEPP	MUP SS 0.2uM PP 1000/PK
US203NPUPP	MUP SS 0.45uM PP 100/PK
US503NPUPP	MUP SS 0.45uM PP 1000/PK

Recommended Alternative

Part Number Item Description

6784-1302	PURADISC 13/0.2 PTFE 100/PK
6783-1302	PURADISC 13/0.2 PTFE 500/PK
6784-1304	PURADISC 13/0.45 PTFE 100/PK
6783-1304	PURADISC 13/0.45 PTFE 500/PK
6784-2502	PURADISC 25/0.2 PTFE 50/PK
6785-2502	PURADISC 25/0.2 PTFE 200/PK
6798-2502	PURADISC 25/0.2 PTFE 1000/PK
6874-1302	Whatman GD/X 13/0.2 PTFE 150/PK
6874-1304	GD/X 13/0.45 PTFE 150/PK
6874-2502	Whatman GD/X 25/0.2 PTFE 150/PK
6874-2504	Whatman GD/X 25/0.45 PTFE 150/PK
6875-2504	Whatman GD/X 25/0.45 PTFE 1500/PK
6974-2504	Whatman GD/XP 25/0.45 PTFE 150/PK
UN203NPEORG	MUP 0.2uM PTFE 100/PK
UN503NPEORG	MUP 0.2uM PTFE 1000/PK
UN203NPUORG	MUP 0.45uM PTFE 100/PK
UN503NPUORG	MUP 0.45uM PTFE 1000/PK
US203NPEORG	MUP SS 0.2uM PTFE 100/PK
US503NPEORG	MUP SS 0.2uM PTFE 1000/PK
US203NPUORG	MUP SS 0.45uM PTFE 100/PK
US503NPUORG	MUP SS 0.45uM PTFE 1000/PK

Parts to be discontinued once supply is depleted:

Discontinued Part		Recommended Alternative	
Part Number	Item Description	Part Number	Item Description
AV125EPP	Autovial 12 0.45uM PP 50/PK	No suggested alternative	-
AV125UPP	Autovial 12 0.45uM PP 50/PK	No suggested alternative	-

Please contact your local Cytiva sales team to request pricing or samples of the suggested alternatives. I regret any inconvenience these events may cause. Thank you for your continued support and business with Cytiva.

Sincerely

Daniel McElroy

Daniel McElroy
Global Product Leader - Analytical
Cytiva Lab Filtration
Daniel.mcelroy@cytiva.com
www.cytiva.com