



His 标签蛋白纯化层析组合精选



Cytiva 用于组氨酸标签蛋白纯化的产品具有以下优点:

- 未螯合或预螯合镍/钴填料的灵活选择
- 高蛋白质结合能力，低的镍或钴离子脱落
- 产品范围涵盖了来自原核和真核表达系统的标签蛋白的纯化
- 不同的规格可以进行筛选、捕获纯化和放大





Ni Sepharose™ excel 革新您的 His 蛋白纯化

产品特点:

- 耐受高达 100 mM EDTA 和高浓度还原剂，简化样品前处理，并最大限度保护蛋白活性
- 无需脱 Ni 直接进行 NaOH 彻底清洗，有效避免交叉污染，延长使用寿命
- Ni 脱落率显著降低，无需反复再生，使用保存更轻松
- 可适用于真核分泌表达的 His 标签蛋白的捕获和纯化，例如昆虫细胞或 CHO 细胞

两步法或三步法纯化策略



图(A). 在 ÄKTExpress™ 系统上，使用 HisTrap™ crude 和 HiLoad™ 16/60 Superdex™ 75 pg 对未过滤的大肠杆菌裂解液进行两步纯化

图(B). SDS-PAGE 结果

订购信息

实验阶段

实验阶段	产品类型	货号	产品名称	产品规格		
STEP 1	His 标签蛋白纯化预装柱	★ 17524701	HisTrap HP	5x1 mL		
		17524802	HisTrap HP	5x5 mL		
		17531901	HisTrap FF	5x1 mL		
		17525501	HisTrap FF	5x5 mL		
		11000458	HisTrap FF crude	5x1 mL		
		17528601	HisTrap FF crude	5x5 mL		
		17092003	HiTrap IMAC HP	5x1 mL		
		17092005	HiTrap IMAC HP	5x5 mL		
		17092102	HiTrap IMAC FF	5x1 mL		
		17092104	HiTrap IMAC FF	5x5 mL		
		17040801	HiTrap Chelating HP	5x1 mL		
		17040903	HiTrap Chelating HP	5x5 mL		
		★ 17371205	HisTrap excel	5x1 mL		
		17371206	HisTrap excel	5x5 mL		
		28953766	HiTrap TALON crude	5x1 mL		
		28953767	HiTrap TALON crude	5x5 mL		
		His 标签蛋白手动纯化	28401353	His SpinTrap	50x100 µL	
			28932171	His SpinTrap Kit	1 套	
			★ 28400989	His MultiTrap HP	4x96 孔板	
			28400990	His MultiTrap FF	4x96 孔板	
			11003399	His GraviTrap	10x1 mL	
			29000593	His SpinTrap TALON	50x100 µL	
			29000594	His GraviTrap TALON	10x1 mL	
			28401351	His GraviTrap KIT	1 套	
			29000596	His MultiTrap TALON	4x96 孔板	
			His 标签蛋白纯化填料	17526802	Ni Sepharose HP	100 mL
		17531802		Ni Sepharose 6FF	100 mL	
	★ 17371202	Ni Sepharose excel		100 mL		
	★ 28934388	HiTrap Capto IEX Selection Kit		5 根不同 1 mL 预装柱		
	STEP 2	Capto 离子交换系列	11001303	HiTrap Capto Q	5x5 mL	
			17544123	HiTrap Capto S	5x5 mL	
			28916540	HiTrap Capto DEAE	5x5 mL	
			11003275	HiTrap Capto MMC	5x5 mL	
			17371755	HiTrap Capto S ImpAct	5x5 mL	
			28405846	HiTrap Capto adhere	5x5 mL	
			17531602	Capto Q	100 mL	
			17546802	Capto SP ImpRes	100 mL	
17531702			Capto MMC	100 mL		
17546855			HiTrap Capto SP ImpRes	5x5 mL		
17547055			HiTrap Capto Q ImpRes	5x5 mL		
STEP 3			凝胶过滤预装柱	★ 29148721	Superdex 75 Increase 10/300 GL	1x24 mL
				28990944	Superdex 200 Increase 10/300 GL	1x24 mL
				29091596	Superose 6 Increase 10/300 GL	1x24 mL
	28989333	HiLoad 16/600 Superdex 75 pg		1x120 mL		
	28989335	HiLoad 16/600 Superdex 200 pg		1x120 mL		
	29323952	HiLoad 16/600 Superose 6 pg		1x120 mL		
	28989334	HiLoad 26/600 Superdex 75 pg		1x320 mL		
	28989336	HiLoad 26/600 Superdex 200 pg		1x320 mL		

快速订购通道

- 1 扫描右侧 Cytiva 微信公众号，进入 **Cytiva 在线商城**，即刻了解消耗品最新促销活动；
- 2 如您有产品应用及技术问题，欢迎拨打 **Cytiva 智荟专线**：400 810 9118



cytiva.com.cn/purify

Cytiva 和 Drop 标识是 Global Life Sciences IP Holdco LLC 或其附属公司的注册商标。Cytiva 版权所有商品和服务的销售需遵守在 Cytiva 企业中运营之供应商公司的销售条款与条件。可应要求提供这些条款与条件的副本。如需了解最新信息，请联系您当地的 Cytiva 代表。如需查看当地办公室的联系信息，请访问 cytiva.com.cn/contact。

CY20757-25Aug21-FL

