

PharmAssure™ 25 mm Hydrophobic Vent Filter

Product Code: HP2002

Description

PharmAssure 25 mm Hydrophobic Vent Filter has a 0.2 µm hydrophobic acrylic copolymer membrane.

This product is CE marked under the Medical Device Regulation (EU) 2017/745 and FDA 510k registered.

Configuration

HP2002: PharmAssure 25mm Hydrophobic Vent filter device with female inlet Luer lock port and male outlet Luer lock port (50 units per case).

Indications

The PharmAssure 25 mm Hydrophobic Vent Filter is designed for sterile venting of vial air during the preparation of injectable drug solutions. A supported 0.2 µm hydrophobic acrylic copolymer membrane protects the user from exposure to potentially toxic aerosols when attached to a needle and inserted through the rubber vial septum.

Contraindications

None.

Precautions

FOLLOW INSTRUCTIONS FOR USE CAREFULLY



Material of Construction

Refer to the Cytiva website (www.cytiva.com) for the Product Safety Data Information in the PSDI_Syringe filters w_Versapor RC_Family.

Performance and Specifications

- **Maximum Recommended Working Pressure:** 15 psi (1.0 bar)
- **Maximum Operating Temperature:** 55°C
- **Pore Size:** 0.2 µm
- **Filter Housing Colour (Upper/Lower):** Clear/clear
- **Filter Housing Diameter:** 28 mm
- **Effective Filtration Area:** 2.8 cm²

Sterility

- Sterile and non-pyrogenic
- **Sterilization method:** Gamma irradiation
- **Endotoxin activity:** < 0.25 EU/mL using LAL test


Shelf Life

- 3 years

Quality

- All materials in the fluid pathway meet relevant sections of ISO 10993 series of standards
- Male and female Luer connectors conform to ISO 80369-7
- Sterilised in accordance with ISO 11137-1, ISO 11137-2 and ISO 11137-3
- Designed and manufactured using quality systems approved to ISO 9001 and ISO 13485
- Manufacturing Environment: ISO 14644 Class 8

Filter Components



Visit us on the Web at www.cytiva.com/pallmedical

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

PharmAssure is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Pall is a trademarks of Pall Corporation. All other third-party trademarks are the property of their respective owners.

© 2025 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit cytiva.com/contact

250322.11GL
04/2025

