

Selection guide


Protein and Nucleic Acid Sample Preparation



Nucleic Acid and Protein Workflows

Nucleic Acid Workflow
Genome and Gene Analysis

Protein Workflow
Differential Expression Analysis Protein Mapping




Nucleic Acid and Protein Workflows

Nucleic Acid Workflow

Construction of Expression-ready DNA Clones

Protein Workflow

Overexpression and Purification of Specific Proteins: Structure and Function Studies



Protein Sample Preparation

	Code no.	2-D Electrophoresis	1-D Electrophoresis	Liquid Chromatography	From 2-D Electrophoresis to MS	From 1-D Electrophoresis to MS	From LC to MS	Protein Expression	Mass Spectrometry	
Untagged Protein Enrichment										
NHS HP SpinTrap™	28-9031-28		•						•	Bulk NHS HP Sepharose™ and empty spin columns for coupling of an affinity molecule with a primary amine
NHS Mag Sepharose	28-9440-09 28-9513-80		•						•	Magnetic beads for coupling of an affinity molecule with a primary amine, pull-down applications
NHS HP SpinTrap Buffer Kit	28-9135-69		•						•	Buffers for immuno-affinity enrichment using NHS HP SpinTrap
Streptavidin HP SpinTrap	28-9031-30		•						•	Coupling of a biotinylated affinity molecule, spin column format
Streptavidin HP SpinTrap Buffer Kit	28-9135-68		•						•	Buffers for immuno-affinity enrichment using Streptavidin HP SpinTrap
Streptavidin HP MultiTrap™	28-9031-31		•						•	Coupling of a biotinylated affinity molecule in 96-well filter plates by vacuum or centrifuge, larger sample series or automation
Streptavidin Mag Sepharose	28-9857-38 28-9857-99		•						•	Magnetic beads for coupling of a biotinylated affinity molecule for IP experiment or direct capture of biotinylated biomolecules
Protein A HP SpinTrap	28-9031-32		•						•	Coupling of an IgG affinity molecule, spin column format
Protein A HP MultiTrap	28-9031-33		•						•	Coupling of an IgG affinity molecule in 96-well filter plates by vacuum or centrifuge, large sample series or automation
Protein A Mag Sepharose	28-9440-06 28-9513-78		•						•	Magnetic beads for coupling of an IgG affinity molecule, immunoprecipitation applications
Protein G HP SpinTrap	28-9031-34		•						•	Coupling of an IgG affinity molecule, spin column format
Protein A/G HP SpinTrap Buffer Kit	28-9135-67		•						•	Buffers for immuno-affinity enrichment using the Protein A HP SpinTrap, Protein G HP SpinTrap, and Ab SpinTrap
Protein G HP MultiTrap	28-9031-35		•						•	Coupling of an IgG affinity molecule in 96-well filter plates by vacuum or centrifuge, larger sample series or automation
Protein G Mag Sepharose	28-9440-08 28-9513-79		•						•	Magnetic beads for coupling of an IgG affinity molecule, immunoprecipitation applications
Ab SpinTrap	28-4083-47		•						•	Coupling of an IgG affinity molecule, spin column format, larger pack size
Immunoprecipitation Starter Pack	17-6002-35		•						•	Coupling of an IgG affinity molecule, batch application
Phos SpinTrap Fe	28-9298-81								•	Enrichment of phosphorylated peptides, spin column format, IMAC precharged with Fe
TiO ₂ Mag Sepharose	28-9440-10 28-9513-77		•						•	Magnetic beads, titanium dioxide chromatography for enrichment of phosphorylated peptides
MagRack 6	28-9489-64									Magnetic rack designed for Mag Sepharose products, enables preparation of up to six samples in 1.5 ml microcentrifuge tubes
MagRack Maxi	28-9864-41									Magnetic rack designed for Mag Sepharose, enables preparation of sample volumes up to 50 ml

For more information on protein sample preparation, please visit www.cytiva.com/sampleprep

Protein Sample Preparation (cont.)

	Code no.	2-D Electrophoresis	1-D Electrophoresis	Liquid Chromatography	From 2-D Electrophoresis to MS	From 1-D Electrophoresis to MS	From LC to MS	Protein Expression	Mass Spectrometry	
Small-scale Antibody Purification										
Protein A HP SpinTrap	28-9031-32							•		Purification of IgG, spin column format
Protein A HP MultiTrap	28-9031-33							•		Purification of IgG in 96-well filter plates by vacuum or centrifuge, larger sample series or automation
Protein A Mag Sepharose Xtra	28-9670-56 28-9670-62							•		Magnetic beads for small-scale purification/screening of monoclonal and polyclonal antibodies from various species
Protein G HP SpinTrap	28-9031-34							•		Purification of IgG, spin column format
Protein G HP MultiTrap	28-9031-35							•		Purification of IgG in 96-well filter plates by vacuum or centrifuge, larger sample series or automation
Protein G Mag Sepharose Xtra	28-9670-66 28-9670-70							•		Magnetic beads for small-scale purification/screening of monoclonal and polyclonal antibodies from various species
Ab SpinTrap	28-4083-47							•		Purification of IgG, spin column format, larger pack size
Ab Buffer Kit	28-9030-59							•		Buffers for antibody purification using the Protein A HP SpinTrap, Protein G HP SpinTrap, and Ab SpinTrap
rProtein A GraviTrap™	28-9852-54		•					•		Purification of antibodies in milligram scale, gravity-flow column
Protein G GraviTrap	28-9852-55		•					•		Purification of antibodies in milligram scale, gravity-flow column
rProtein A/Protein G GraviTrap	28-9852-56		•					•		Purification of antibodies in milligram scale, gravity-flow column
Desalting/Buffer Exchange/Clean-up										
Disposable PD-10 Desalting	17-0851-01	•	•	•			•	•	•	Clean-up of proteins/oligosaccharides, sample volumes up to 2.5 ml, gravity-flow/centrifugation
LabMate™ PD-10 Buffer	18-3216-03	•	•	•			•	•	•	Buffer reservoir for easy equilibration
PD MidiTrap™ G-25	28-9180-08	•	•	•			•	•	•	Clean-up of proteins/oligosaccharides, sample volumes up to 1.0 ml, gravity-flow/centrifugation
PD MiniTrap™ G-25	28-9180-07	•	•	•			•	•	•	Clean-up of proteins/oligosaccharides, sample volumes up to 0.5 ml, gravity-flow/centrifugation
PD SpinTrap G-25	28-9180-04	•	•	•			•	•	•	Clean-up of proteins/oligosaccharides, sample volumes up to 130 µl, microcentrifugation
PD MultiTrap G-25	28-9180-06	•	•	•			•	•	•	Clean-up of proteins/oligosaccharides in 96-well plates by centrifuge, sample volumes up to 130, larger sample series or automation
PD MidiTrap G-10	28-9180-11	•	•	•			•	•	•	Clean-up of peptides/small proteins, sample volumes up to 800 µl, gravity flow
PD MiniTrap G-10	28-9180-10	•	•	•			•	•	•	Clean-up of peptides/small proteins, sample volumes up to 300 µl, gravity flow
Mini Dialysis Kit – 1 kDa cut-off	80-6483-75 80-6483-94	•	•	•			•	•	•	Dialysis of small sample volumes, 250 µl and 2 ml Dialysis of small sample volumes, max. sample volume 2 ml
Mini Dialysis Kit – 8 kDa cut-off	80-6484-13 80-6484-32	•	•	•			•	•	•	Dialysis of small sample volumes, 250 µl and 2 ml Dialysis of small sample volumes, max. sample volume 2 ml
2-D Clean-Up Kit	80-6484-51	•							•	Removal of interfering contaminants and concentration of total protein
SDS-PAGE Clean-Up Kit	80-6484-70		•						•	Removal of interfering contaminants

For more information on protein sample preparation, please visit www.cytiva.com/sampleprep

Protein Sample Preparation (cont.)

	Code no.	2-D Electrophoresis	1-D Electrophoresis	Liquid Chromatography	From 2-D Electrophoresis to MS	From 1-D Electrophoresis to MS	From LC to MS	Protein Expression	Mass Spectrometry	
Enzyme Regulation										
Protease Inhibitor Mix	80-6501-23	•	•	•				•	•	Protease inhibitor cocktail
Nuclease Mix	80-6501-42	•	•	•				•	•	Cocktail of nucleases for the removal of DNA and RNA from the protein sample
Fractionation										
2-D Fractionation Kit	80-6501-04	•	•	•					•	Fractionation of the total protein in the sample into six discrete fractions
Total Protein Quantitation										
2-D Quant Kit	80-6483-56	•	•	•	•	•	•	•	•	Quantitation of protein amount using a colorimetric assay
NanoVue™ Plus	28-9569-65	•	•	•			•	•	•	UV/Visible spectrophotometer for quick and accurate quantitation of nucleic acids and proteins, 0.5 to 2 µl sample volumes
Histidine-tagged Protein Capture										
His GraviTrap	11-0033-99							•		Purification of histidine-tagged proteins, with high binding capacity, gravity flow column
His GraviTrap TALON	29-0005-94									Gravity flow columns for purification of histidine-tagged proteins to high purity
His Mag Sepharose Ni	28-9673-88 28-9673-90 28-9799-17							•		Magnetic beads for small-scale purification/screening of histidine-tagged proteins
His Mag Sepharose excel	17-3712-20 17-3712-21 17-3712-22							•		Magnetic beads for small-scale purification/screening of histidine-tagged proteins from expression in mammalian and insect cell systems
His GraviTrap Kit	28-4013-51							•		Gravity flow columns for histidine-tagged protein purification and buffers
His MultiTrap HP	28-4009-89							•		Purification of histidine-tagged proteins with high binding capacity in 96-well filter plates by vacuum or centrifuge, for larger sample series or automation
His MultiTrap FF	28-4009-90							•		Purification of histidine-tagged proteins with high binding capacity in 96-well filter plates by vacuum or centrifuge, for larger sample series or automation
His MultiTrap TALON	29-0005-96							•		Purification of histidine-tagged proteins with high purity in 96-well filter plates by vacuum or centrifuge, for larger sample series or automation
His SpinTrap	28-4013-53							•		Purification of histidine-tagged proteins with high binding capacity, spin column format
His SpinTrap TALON	29-0005-93							•		Purification of histidine-tagged proteins with high purity, spin column format
His SpinTrap Kit	28-9321-71							•		Microspin™ columns and premade buffers for histidine-tagged protein purification
His Buffer Kit	11-0034-00							•		Buffers for the purification of histidine-tagged proteins using His GraviTrap, HisTrap, or His SpinTrap
Anti-His Antibody	27-4710-01							•		Unconjugated antibody with affinity for the histidine tag
HisTrap™ HP columns	17-5247-01							•		Simple purification with a syringe or chromatography system such as ÄKTA™
HisTrap FF columns	17-5319-01							•		Simple purification with a syringe or chromatography system such as ÄKTA
HisPrep™ FF 16/10 column	17-5256-01							•		Easy scale-up purification of histidine-tagged proteins
HisTrap FF crude Kit	28-4014-77							•		Purification with a syringe or chromatography system such as ÄKTA, no filtration of the sample
Maltose Binding Protein (MBP)-tagged Protein Capture										
MBPTrap™ HP	28-9187-78							•		Simple purification with a syringe or chromatography system such as ÄKTA
Strep-tag™ II-tagged Protein Capture										
StrepTrap™ HP	28-9075-46							•		Simple purification with a syringe or chromatography system such as ÄKTA

For more information on protein sample preparation, please visit www.cytiva.com/sampleprep

Protein Sample Preparation (cont.)

Code no.	2-D Electrophoresis	1-D Electrophoresis	Liquid Chromatography	From 2-D Electrophoresis to MS	From 1-D Electrophoresis to MS	From LC to MS	Protein Expression	Mass Spectrometry	
GST-tagged Protein Capture									
GST Detection Module	27-4590-01						•		Identification of GST-tagged proteins by biochemical or immunological assay
GST 96-well Detection Module	27-4592-01						•		Reagents and plates with precoated wells for the detection of GST-tagged proteins
GSTPrep™ FF 16/10 column	28-9365-50						•		Simple purification with a syringe or chromatography system such as ÄKTA
GSTrap™ HP columns	17-5281-01						•		Simple purification with a syringe or chromatography system such as ÄKTA
GSTrap FF columns	17-5130-02						•		Simple purification with a syringe or chromatography system such as ÄKTA
GSTrap 4B columns	28-4017-45						•		Simple purification with a syringe or chromatography system such as ÄKTA
GST GraviTrap	28-9523-60						•		Purification of GST-tagged proteins in gravity-flow format
GST MultiTrap FF	28-4055-01						•		Purification of GST-tagged proteins in 96-well filter plates by vacuum or centrifuge, for larger sample series or automation
GST MultiTrap 4B	28-4055-00						•		Purification of GST-tagged proteins in 96-well filter plates by vacuum or centrifuge, for larger sample series or automation
GST SpinTrap	28-9523-59						•		Purification of GST-tagged proteins, spin column format
GST Buffer Kit	28-9523-61						•		Buffers for purification of GST-tagged proteins using GSTrap, GST GraviTrap, or GST SpinTrap
pGEX vectors	multiple						•		Expression vectors for GST fusion proteins
GST Vector primers for sequencing	multiple						•		Sequencing primers for the pGEX vectors
<i>E. coli</i> BL21	27-1542-01						•		Bacterial strain for the expression of GST fusion protein
M13K07 Helper Phage	27-1524-01						•		M13 phage for the GST expression system
GST Bulk Kit	27-4570-01						•		Bulk Glutathine Sepharose 4B and empty gravity columns
PreScission™ Protease	27-0843-01						•		Purified enzyme for site-specific cleavage of the GST tag from fusion proteins expressed using pGEX-P vectors
Thrombin	27-0846-01						•		Purified enzyme for site-specific cleavage of the GST tag from fusion proteins expressed using pGEX-T vectors
Factor Xa	27-0849-01						•		Purified enzyme for site-specific cleavage of the GST tag from fusion proteins expressed using pGEX-X vectors
Anti-GST Antibody	27-4577-01						•		Unconjugated antibody with affinity for the GST tag
Protein Depletion									
HiTrap™ Albumin & IgG Depletion	28-9466-03	•	•	•					Prepacked 1 ml column for the depletion of albumin and IgG from larger sample volumes (~150 µl) of human plasma and serum, for use with a syringe or chromatography system such as ÄKTA
Albumin & IgG Depletion SpinTrap	28-9480-20	•	•	•					Depletion of albumin and IgG from smaller sample volumes (~50 µl) of human plasma and serum, spin column format
Lysis/Protein Extraction									
Sample Grinding Kit	80-6483-37	•	•	•			•	•	The mechanical breaking up of cells or tissue releasing the total protein content
Yeast Protein Extraction Buffer Kit	28-9440-45	•	•	•			•	•	Extraction of proteins from yeast cells
Mammalian Protein Extraction Buffer	28-9412-79	•	•	•			•	•	Extraction of proteins from mammalian cells
2-D Protein Extraction Buffer Trial Kit	28-9435-22	•	•	•			•	•	Protein extraction from tissues and cell lines, trial kit with six different buffers for optimization of the method
2-D Protein Extraction Buffer-I	28-9435-23	•	•	•			•	•	Protein extraction buffer - larger volume of Buffer-I from the trial kit

For more information on protein sample preparation, please visit www.cytiva.com/sampleprep

Protein Sample Preparation (cont.)

Code no.	2-D Electrophoresis	1-D Electrophoresis	Liquid Chromatography	From 2-D Electrophoresis to MS	From 1-D Electrophoresis to MS	From LC to MS	Protein Expression	Mass Spectrometry	
Lysis/Protein Extraction									
2-D Protein Extraction Buffer-II	28-9435-24	•	•	•			•	•	Protein extraction buffer - larger volume of Buffer-II from the trial kit
2-D Protein Extraction Buffer-III	28-9435-25	•	•	•			•	•	Protein extraction buffer - larger volume of Buffer-III from the trial kit
2-D Protein Extraction Buffer-IV	28-9435-26	•	•	•			•	•	Protein extraction buffer - larger volume of Buffer-IV from the trial kit
2-D Protein Extraction Buffer-V	28-9435-27	•	•	•			•	•	Protein extraction buffer - larger volume of Buffer-V from the trial kit
2-D Protein Extraction Buffer-VI	28-9435-28	•	•	•			•	•	Protein extraction buffer - larger volume of Buffer-VI from the trial kit
triplePrep Kit	28-9425-44	•	•	•				•	Extraction of DNA, RNA, and protein from the same undivided sample
Protein Concentration									
Vivaspin 500, 3 kDa MWCO PES	28-9322-18	•	•	•	•	•	•	•	Concentration of up to 500 µl sample through ultrafiltration, 3 kDa cut-off
Vivaspin 500, 5 kDa MWCO PES	28-9322-23	•	•	•	•	•	•	•	Concentration of up to 500 µl sample through ultrafiltration, 5 kDa cut-off
Vivaspin 500, 10 kDa MWCO PES	28-9322-25	•	•	•	•	•	•	•	Concentration of up to 500 µl sample through ultrafiltration, 10 kDa cut-off
Vivaspin 500, 30 kDa MWCO PES	28-9322-35	•	•	•	•	•	•	•	Concentration of up to 500 µl sample through ultrafiltration, 30 kDa cut-off
Vivaspin 500, 50 kDa MWCO PES	28-9322-36	•	•	•	•	•	•	•	Concentration of up to 500 µl sample through ultrafiltration, 50 kDa cut-off
Vivaspin 500, 100 kDa MWCO PES	28-9322-37	•	•	•	•	•	•	•	Concentration of up to 500 µl sample through ultrafiltration, 100 kDa cut-off
Vivaspin 2, 3 kDa MWCO PES	28-9322-40	•	•	•	•	•	•	•	Concentration of up to 2 ml sample through ultrafiltration, 3 kDa cut-off
Vivaspin 2, 5 kDa MWCO PES	28-9322-45	•	•	•	•	•	•	•	Concentration of up to 2 ml sample through ultrafiltration, 5 kDa cut-off
Vivaspin 2, 10 kDa MWCO PES	28-9322-47	•	•	•	•	•	•	•	Concentration of up to 2 ml sample through ultrafiltration, 10 kDa cut-off
Vivaspin 2, 30 kDa MWCO PES	28-9322-48	•	•	•	•	•	•	•	Concentration of up to 2 ml sample through ultrafiltration, 30 kDa cut-off
Vivaspin 2, 50 kDa MWCO PES	28-9322-57	•	•	•	•	•	•	•	Concentration of up to 2 ml sample through ultrafiltration, 50 kDa cut-off
Vivaspin 2, 100 kDa MWCO PES	28-9322-58	•	•	•	•	•	•	•	Concentration of up to 2 ml sample through ultrafiltration, 100 kDa cut-off
Vivaspin 6, 3 kDa MWCO PES	28-9322-93	•	•	•	•	•	•	•	Concentration of up to 6 ml sample through ultrafiltration, 3 kDa cut-off
Vivaspin 6, 5 kDa MWCO PES	28-9322-94	•	•	•	•	•	•	•	Concentration of up to 6 ml sample through ultrafiltration, 5 kDa cut-off
Vivaspin 6, 10 kDa MWCO PES	28-9322-96	•	•	•	•	•	•	•	Concentration of up to 6 ml sample through ultrafiltration, 10 kDa cut-off
Vivaspin 6, 30 kDa MWCO PES	28-9323-17	•	•	•	•	•	•	•	Concentration of up to 6 ml sample through ultrafiltration, 30 kDa cut-off
Vivaspin 6, 50 kDa MWCO PES	28-9323-18	•	•	•	•	•	•	•	Concentration of up to 6 ml sample through ultrafiltration, 50 kDa cut-off
Vivaspin 6, 100 kDa MWCO PES	28-9323-19	•	•	•	•	•	•	•	Concentration of up to 6 ml sample through ultrafiltration, 100 kDa cut-off
Vivaspin 20, 3 kDa MWCO PES	28-9323-58	•	•	•	•	•	•	•	Concentration of up to 20 ml sample through ultrafiltration, 3 kDa cut-off
Vivaspin 20, 5 kDa MWCO PES	28-9323-59	•	•	•	•	•	•	•	Concentration of up to 20 ml sample through ultrafiltration, 5 kDa cut-off
Vivaspin 20, 10 kDa MWCO PES	28-9323-60	•	•	•	•	•	•	•	Concentration of up to 20 ml sample through ultrafiltration, 10 kDa cut-off
Vivaspin 20, 30 kDa MWCO PES	28-9323-61	•	•	•	•	•	•	•	Concentration of up to 20 ml sample through ultrafiltration, 30 kDa cut-off
Vivaspin 20, 50 kDa MWCO PES	28-9323-62	•	•	•	•	•	•	•	Concentration of up to 20 ml sample through ultrafiltration, 50 kDa cut-off
Vivaspin 20, 100 kDa MWCO PES	28-9323-63	•	•	•	•	•	•	•	Concentration of up to 20 ml sample through ultrafiltration, 100 kDa cut-off

For more information on protein sample preparation, please visit www.cytiva.com/sampleprep

Buffer and Sample Filtration, Sample Collection

	Pack Size	Code No.
Klari-Flex™ Bottle Top Filtration System for Sterile Filtration		
KF 250 ml 0.22 PES Funl Sterile	12/pk	6515-2502
KF 500 ml 0.22 PES Funl Sterile	12/pk	6515-5002
KF 1000 ml 0.22 PES Funl Sterile	12/pk	6515-1002
Klari-Flex Bottle Top System Hardware for Sterile Filtration		
KF Cradle Ring	1/pk	6517-0001
KF Pedestal Stand	1/pk	6517-0002
Polycap PES Filter Capsules for Sterile Filtration		
Polycap 36 0.2/0.2 PES with filling bell/Sterile	1/pk	6715-3602
Polycap 75 0.2/0.2 PES with filling bell/Sterile	1/pk	6715-7502
Polycap 150 0.2/0.2 PES with filling bell/Sterile	1/pk	6718-9502
SPARTAN 13 mm Syringe Filters		
13 mm Spartan™, 0.2 µm RC	100/pk	10463100
13 mm Spartan, 0.45 µm RC	100/pk	10463110
SPARTAN 30 mm Syringe Filters		
30 mm Spartan, 0.2 µm RC	100/pk	10463060
30 mm Spartan, 0.45 µm RC	100/pk	10463050
GD/X Syringe Filter for Lysate Clarification		
25 mm GD/X 0.2 µm PES (Sterile)	50/pk	6896-2502
25 mm GD/X 0.45 µm PES (Sterile)	50/pk	6896-2504
Syringe filters for Sterile Filtration		
4 mm Puradisc 0.2 µm PVDF (Sterile)	50/pk	6791-0402
25 mm Puradisc 0.2 µm PES (Sterile)	50/pk	6780-2502
25 mm Puradisc 0.45 µm PES (Sterile)	50/pk	6780-2504
Mini-UniPrep™ Syringeless Filters for LC and LC-MS		
Mini-UniPrep 0.2 µm PES (with Slit Septum Cap)	100/pk	US203NPEPES
Mini-UniPrep 0.2 µm RC (with Standard Cap)	100/pk	UN203NPERC
Sample Collection Card		
ProteinSaver US	100/pk	10534612
ProteinSaver EU	100/pk	10531018

For more information on Buffer & Sample Filtration, please visit www.cytiva.com/whatman

Nucleic Acid Sample Preparation

PCR and RT-PCR

Amplification method	Starting quantity	Hot Start Mix Ready-To-Go*	Hot Start Master Mix	PuReTaq Ready-To-Go* PCR	RTG RT-PCR*
PCR	Basic PCR (amplicon up to 3 kb)	•	•	•	
	Hot start PCR (amplicon up to 3 kb)	•	•		
	Long-range PCR (amplicon up to 20 kb)				
RT-PCR	Basic RT-PCR (amplicon up to 3 kb)				•
	Long-range RT-PCR (amplicon up to 6 kb)				
Downstream applications	Genotyping	•	•	•	•
	Cloning and sequencing	•	•	•	•
	Gene expression	•	•	•	•

* Ready-To-Go (RTG™) is a single-dose, room temperature stable bead format that contains all PCR or RT-PCR components. You only need to add DNA template, primer, and water for a reaction.

DNA Clean-up

Sample types	Starting quantity	Exo-ProStar™	GFX PCR and Gel Band	GFX 96 PCR	MicroSpin G-25	AutoSeq G-50	ProbeQuant G-50	NICK columns	NAP-5 columns	NAP-10 columns	NAP-25 columns	CyScribe GFX	MicroSpin G-50	MicroSpin S-200	MicroSpin S-300	MicroSpin S-400
PCR products	Any size, for sequencing only															
	50 bp to 10 kb, few samples	•	•													
	100 bp to 10 kb, many samples	•		•												
	10 to 50 bp	•			•											
Agarose gel slices or enzyme removal	50 kbp to 10 bp		•													
Sequencing reactions	12 to 25 µl	•				•										
Labeled DNA fragments	> 20 bases	•					•									
	> 10 bases	•			•											
	100 µg sample in 100 µl or no microcentrifuge							•								
Oligonucleotides	100 to 150 µl				•											
	0.1 to 0.5 ml								•							
	0.5 to 1 ml									•						
	1 to 2.5 ml										•					
cDNA	CyDye™ labeled probes										•					
A range of sample types	25 to 50 µl for labeled DNA	•											•			
	PCR, sequencing, and labeling reactions												•			
	10 to 100 µl for buffer exchange or desalting												•			

For more information on protein sample preparation, please visit www.cytiva.com/sampleprep

Nucleic Acid Sample Preparation (cont.)

DNA Clean-up

Sample types	Starting quantity	Exo-ProStar™	GFX PCR and Gel Band	GFX 96 PCR	MicroSpin G-25	AutoSeq G-50	ProbeQuant G-50	NICK columns	NAP-5 columns	NAP-10 columns	NAP-25 columns	CyScribe GFX	MicroSpin G-50	MicroSpin S-200	MicroSpin S-300	MicroSpin S-400
A range of volumes and sample sizes	PCR (25 to 50 µl) or labeling reaction (25 to 50 µl) < 100 bp	•												•		
	PCR (25 to 50 µl) or labeling reaction (50 to 75 µl) or fragment < 200 bp	•													•	
	PCR (50 to 100 µl) or labeling reaction (75 to 100 µl) or remove primers > 24 bases	•														•
Downstream applications	Cloning and sequencing	•	•	•	•	•			•	•	•		•	•	•	•
	Gene expression		•	•	•		•	•	•	•	•	•	•	•	•	•
	Genotyping		•	•	•								•	•	•	•

For more information on protein and nucleic acid sample preparation, please visit www.cytiva.com/sampleprep

Plasmid DNA Preparation

Sample types	Starting quantity	plasmid-Prep Mini	plasmid-Prep Midi	TempliPhi 100/500*	TempliPhi Large Construct*	TempliPhi Sequence Resolver Kit*
Bacterial culture	1 µl (high-throughput)					
	1 µl (low- to medium-throughput)			•		
	1 to 3 ml	•				
	25 to 500 ml		•			
Bacterial colony	< 1 colony (high-throughput)					
	< 1 colony (low- to medium-throughput)			•		
Bacteria glycerol stock	< 1 µl (high-throughput)					
	< 1 µl (low- to medium-throughput)			•		
Purified DNA (small vectors)	> 1 ng (high-throughput)					
	> 1 ng (low- to medium-throughput)			•		
Purified BAC/fosmid DNA	> 1 ng				•	
Difficult templates (GC rich, secondary structures)	0.1 to 1 ng					•
M13 phage plaque	< 1 plaque			•		
M13 phage glycerol stock	< 1 µl			•		
	Cloning	•	•			
Downstream applications	Sequencing	•	•	•	•	•

* Can be used directly for PCR and subcloning. Additional steps will be needed for transformation or transfection.

Nucleic Acid Sample Preparation (cont.)

RNA Preparation

Sample type and format	Starting quantity	RNAspin Mini	RNAspin Midi	RNAspin 96	QuickPrep mRNA	QuickPrep Micro mRNA	mRNA Purification Kit
Total RNA							
Cultured cells and tissue	Up to 200 mg tissue or 5×10^7 cells		•				
Tube format	10 to 30 mg or up to 2×10^6 cells	•					
96-well plate vacuum	10 to 30 mg or up to 2×10^6 cells			•			
96-well plate centrifuge	30 mg or up to 1×10^7 cells			•			
mRNA purification	1 to 5×10^7 cells or 0.5 g tissue				•		
Eukaryotic cells or tissue	1 to 1×10^7 cells or 100 mg tissue					•	
Eukaryotic total RNA	Total RNA or 25 mg to 1 g tissue						•
Downstream applications	Gene expression	•	•	•	•	•	•

DNA, RNA, and Protein Preparation

Sample type	Starting quantity	triplePrep kit
Cultured cells	0.3 to 5×10^7 cells	•
Animal tissue	1 to 20 mg tissue	•
Downstream applications	Cloning and sequencing	•
	Gene expression	•
	Genotyping	•
	Protein analysis (1-D, 2-D, LCMS)	•

For more information on nucleic acid sample preparation, please visit www.cytiva.com/sampleprep

Nucleic Acid Sample Preparation (cont.)

Genomic RNA Preparation

Sample types	Starting quantity (DNA yield)	tissue & cells Mini Spin	tissue & cells Midi Flow	blood Mini Spin	blood Midi Flow	bacteria Mini Spin	BACC1	BACC2 & BACC3	Nucleon HT	Phytopure	GenomiPhi V2	GenomiPhi HY
Animal tissue	> 300 cells (4 to 7 µg DNA yield)										•	•
	> 300 cells (40 to 50 µg DNA yield)											•
	5 to 50 mg	•										
	Up to 200 mg		•									
Paraffin-embedded tissue or difficult samples	Up to 25 mg								•			
	20 to 30 µm thick paraffin section								•			
Cultured cells	> 300 cells (4 to 7 µg DNA yield)										•	
	> 300 cells (40 to 50 µg DNA yield)											•
	Up to 5.0 × 10 ⁵ cells	•										
	1 to 3 × 10 ⁵						•					
	Up to 2.0 × 10 ⁷ cells		•									
Blood	3 × 10 ⁶ to 1 × 10 ⁷ cells							•				
	5 to 10 µl (4 to 7 µg DNA yield)										•	
	5 to 10 µl (40 to 50 µg DNA yield)											•
	50 to 300 µl			•								
	1 ml						•					
	1 to 8 ml				•							
Buffy coat	10 ml							•				
	5 to 10 µl (4 to 7 µg DNA yield)										•	
	5 to 10 µl (40 to 50 µg DNA yield)											•
Nucleated red blood cells	50 to 300 µl			•								
	10 µl (4 to 7 µg DNA yield)										•	
	10 µl (40 to 50 µg DNA yield)											•
	10 µl			•								
Bone marrow (suspended cells)	25 to 200 µl (tested)				•							
	5 to 10 µl (4 to 7 µg DNA yield)										•	
	5 to 10 µl (40 to 50 µg DNA yield)											•
Bacteria (Gram - and +)	200 µl			•								
	> 300 cells (4 to 7 µg DNA yield)										•	
	> 300 cells (40 to 50 µg DNA yield)											•
Plants	Up to 4.0 × 10 ⁹					•						
	1 cm leaf (4 to 7 µg DNA yield)										•	
	1 cm leaf (40 to 50 µg DNA yield)											•
	1 seed (4 to 7 µg DNA yield)										•	
	1 seed (40 to 50 µg DNA yield)											•
Up to 1.0 g									•			

For more information on nucleic acid sample preparation, please visit www.cytiva.com/sampleprep

Nucleic Acid Sample Preparation (cont.)

Genomic RNA Preparation

Sample types	Starting quantity (DNA yield)	tissue & cells Mini Spin	tissue & cells Midi Flow	blood Mini Spin	blood Midi Flow	bacteria Mini Spin	BACC1	BACC2 & BACC3	Nucleon HT	Phytopure	GenomiPhi V2	GenomiPhi HY
FTA paper/Guthrie card	3 × 3 mm piece (4 to 7 µg DNA yield)										•	
	3 × 3 mm piece (40 to 50 µg DNA yield)											•
Buccal swab	Single swab			•								
	Single swab (4 to 7 µg DNA yield)										•	
	Single swab (40 to 50 µg DNA yield)											•
Purified genomic DNA	> 10 ng (4 to 7 µg DNA yield)										•	
	> 10 ng (40 to 50 µg DNA yield)											•
Downstream applications	Genotyping	•	•	•	•	•	•	•	•	•	•	•
	Cloning and sequencing	•	•	•	•	•	•	•	•	•	•	•

For more information on nucleic acid sample preparation, please visit www.cytiva.com/sampleprep

Ordering information

Columns	Pack size	Code No.
Genomic DNA Preparation		
Plasmid DNA Purification		
Amersham tissue and cells genomicPrep™ Mini Spin Kit	50	28-9042-75
	250	28-9042-76
Amersham tissue and cells genomicPrep Midi Flow Kit	25	28-9042-73
Amersham blood genomicPrep Mini Spin Kit	50	28-9042-64
	250	28-9042-65
Amersham blood genomicPrep Midi Flow Kit	25	28-9042-61
	100	28-9042-62
Amersham bacteria genomicPrep Mini Spin Kit	50	28-9042-58
	250	28-9042-59
Nucleon™ BACC1 Genomic DNA Kit	50	RPN8501
Nucleon BACC2 Genomic DNA Kit	50	RPN8502
Nucleon BACC3 Genomic DNA Kit	50	RPN8512
Nucleon PhytoPure	50 × 0.1 g	RPN8510
	50 × 1.0 g	RPN8511
Nucleon HT (hard tissue)	50	RPN8509
Whole Genome Amplification		
GenomiPhi™ V2 DNA Amplification Kit	25	25-6600-30
	100	25-6600-31
	500	25-6600-32
GenomiPhi HY DNA Amplification Kit	25	25-6600-22
	100	25-6600-20
	1000	25-6600-25

Columns	Pack size	Code No.
Plasmid DNA Preparation		
Plasmid DNA Purification		
Amersham PlasmidPrep Mini Spin Kit	50	28-9042-69
	250	28-9042-70
Amersham PlasmidPrep Midi Flow Kit	25	28-9042-67
	100	28-9042-68
Sephacryl™ S-1000 SF	750 ml	17-0476-01
Yeast Plasmid Isolation Kit	50	US79220-50RXNS
Rolling Circle Amplification – Template Amplification		
TempliPhi™ 100 DNA Amplification Kit	100	25-6400-10
TempliPhi 500 DNA Amplification Kit	500	25-6400-50
TempliPhi Large Construct Kit	1000	25-6400-80
TempliPhi Sequence Resolver Kit	20 rxns	28-9035-29
	50 rxns	28-9035-30
	200 rxns	28-9035-31
RNA Preparation		
Total RNA		
Amersham RNAspin Mini Kit	20	25-0500-70
	50	25-0500-71
	250	25-0500-72
Amersham RNAspin Midi Kit	20	25-0500-73
Amersham RNAspin 96 Kit	4 × 96	25-0500-75
Amersham RNAspin 96 Filter Plate	1	25-0500-88
Amersham CsTFA (solution)	100 ml	17-0847-02
mRNA		
QuickPrep™ Micro mRNA Purification Kit	1	27-9255-01
QuickPrep mRNA Purification Kit	4	27-9254-01
mRNA Purification Kit	2	27-9258-01
	4	27-9258-02

Columns	Pack size	Code No.
PCR		
PCR Kits & Components		
Amersham Hot Start Mix RTG 0.5 ml	100 rxn	28-9006-46
Amersham Hot Start Mix RTG 0.2 ml	96 rxn	28-9006-53
	480 rxn	28-9006-54
Amersham Hot Start Master Mix	100	25-1500-01
PuReTaq Ready-To-Go™ PCR Beads 0.2 ml	96 rxn plate	27-9557-01
	5 × 96 rxn plate	27-9557-02
PuReTaq Ready-To-Go PCR Beads 0.5 ml	100 rxn tubes	27-9558-01
PuReTaq Ready-To-Go PCR Beads 0.2 ml	96 rxn h-tube	27-9559-01
Amersham Ready-To-Go RAPD Analysis Beads	100	27-9500-01
Amersham Ready-To-Go RAPD Analysis Kit	100 + 6 Primers	27-9502-01
RT-PCR		
RT-PCR Kits & Components		
Amersham Ready-To-Go RT-PCR Beads 0.5 ml	100 rxn tubes	27-9266-01
Amersham Ready-To-Go RT-PCR Beads 0.2 ml	96 rxn tube	27-9267-01
	96 rxn h-tube	27-9259-01

Ordering information

Columns	Pack size	Code No.
DNA Clean-up		
Purification of PCR Products and Restriction Fragments		
GFX™ PCR DNA and Gel Band Purification Kit	100	28-9034-70
	250	28-9034-71
GFX 96 PCR Purification Kit	10 × 96	28-9034-45
Amersham MicroSpin S-400 HR columns	50	27-5140-01
Amersham MicroSpin S-300 HR columns	50	27-5130-01
Amersham MicroSpin G-25 columns	50	27-5325-01
ExoProStar™	20	US78220
	100	US78210
	500	US78211
	2000	US78212
	5000	US78225
ExoProStar 1-Step*	20	US77701
	100	US77702
	500	US77705
	2000	US77720
	5000	US77750
Dye Terminator Removal		
AutoSeq G-50	50	27-5340-01
	250	27-5340-02
	1000	27-5340-03
Amersham MicroSpin G-50 columns	50	27-5330-01
	250	27-5330-02
Amersham MicroSpin S-200 HR columns	50	27-5120-01

Columns	Pack size	Code No.
DNA Clean-up		
Labeling and Clean-up		
CyScribe GFX Purification Kit	25	27-9606-01
	50	27-9606-02
Amersham NICK columns	20	17-0855-01
	50	17-0855-02
ProbeQuant™ G-50 Micro columns	50	28-9034-08
Desalting and Nucleotide Removal		
Amersham MicroSpin G-25 columns	50	27-5325-01
Amersham NAP™-5	20	17-0853-01
	50	17-0853-02
Amersham NAP-10	20	17-0854-01
	50	17-0854-02
Amersham NAP-25	20	17-0852-01
	50	17-0852-02
Sephadex™ G-25 DNA Grade SF	100 g	17-0572-02
Sephadex G-50 DNA Grade F	25 g	17-0573-01
	100 g	17-0573-02
Sephadex G-100 DNA Grade SF	100 g	17-0574-02

Columns	Pack size	Code No.
Nucleotides		
dATP, 100 mM	25 µmol	28-4065-01
	100 µmol	28-4065-02
	500 µmol	28-4065-03
dCTP, 100 mM	25 µmol	28-4065-11
	100 µmol	28-4065-12
	500 µmol	28-4065-13
dGTP, 100 mM	25 µmol	28-4065-21
	100 µmol	28-4065-22
	500 µmol	28-4065-23
dTTP, 100 mM	25 µmol	28-4065-31
	100 µmol	28-4065-32
	500 µmol	28-4065-33
dUTP, 100 mM	25 µmol	28-4065-41
	100 µmol	28-4065-42
dNTP Set (100 mM each A,C,G,T)	4 × 25 µmol	28-4065-51
	4 × 100 µmol	28-4065-52
	4 × 500 µmol	28-4065-53
Amersham DNA Polymerization Mix dNTP Set (20 mM each A,C,G,T)	10 µmol	28-4065-57
	40 µmol (4 × 10 µl)	28-4065-58
Amersham PCR Nucleotide Mix dNTP Set (25 mM each A,C,G,T)	500 µl	28-4065-60
Amersham PCR Nucleotide Mix dNTP Set (2 mM each A,C,G,T)	1 ml	28-4065-62
Amersham PCR Nucleotide Mix dNTP Set (10 mM each A,C,G,T)	500 µl	28-4065-64
DNA, RNA and protein preparation		
triplePrep Kit	50	28-9425-44

* Not available in North America or Japan

cytiva.com/sampleprep

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate. ÄKTA, CyDye, ExoProStar, genomicPrep, GenomiPhi, GFX, GraviTrap, GSTprep, GSTrap, HiTrap, HisTrap, HisPrep, Klari-Flex, LabMate, MBPTrap, MicroSpin, MidiTrap, MiniTrap, MultiTrap, Mini-UniPrep, NanoVue, NAP, PreScission, ProbeQuant, QuickPrep, Ready-To-Go, RTG, Sephacryl, Sephadex, Sepharose, Spartan, SpinTrap, StrepTrap, TempliPhi, and Whatman are trademarks of Cytiva companies.

CyDye: This product, or portions thereof, is manufactured under an exclusive license from Carnegie Mellon University under US patent number 5,268,486 and equivalent patents in the US and other countries.

IMAC Sepharose products and Ni Sepharose products (Histidine- tagged protein purification, Lab products): Purification and preparation of fusion proteins and affinity peptides comprising at least two adjacent histidine residues may require a license under US patent numbers 5,284,933 and 5,310,663, and equivalent patents and patent applications in other countries (assignee: Hoffman La Roche, Inc).

pGEX Vectors are to be used for scientific investigation and research and for no other purpose whatsoever and a license for commercial use of the licensed products and the processes claimed in US patent 5,654,176 and equivalent patents and patent applications in other countries must be negotiated directly with Millipore Corp (formerly Chemicon International Inc) by the purchaser prior to such use.

Phi 29 DNA polymerase and its use for DNA synthesis is covered by US patent numbers 5,854,033, 5,198,543 and 5,576,204.

Ready to Go RT-PCR Beads: Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, and 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

StrepTrap HP: This product is covered by US patent number 6,103,493 and equivalent patents and patent applications in other countries. The purchase of StrepTrap HP includes a license under such patents for non-profit and in-house research only. Please contact IBA (info@iba-go.com) for further information on licenses for commercial use of StrepTactin.

2-D Protein Extraction Buffers and Mammalian Extraction Buffer and Yeast Extraction Buffer Kit are manufactured by g-Biosciences, St Louis, MO, USA. All other third-party trademarks are the property of their respective owners.

© 2020 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information. For local office contact information, visit cytiva.com/contact

CY13985-24Jul20-SG

